

# FORMFIT® OA EASE

FORMFIT® OA EASE



The Formfit OA Ease is a single upright brace with a Dynamic Force System™ strap that may help relieve pain from mild-to-moderate unicompartmental knee osteoarthritis (OA). Its wraparound design is easy to apply and remove using simple rubber pull-tabs, colour-coded guide dots, and strap tension indicator lines. This low-profile knee brace is comfortable, easy to use, and adapts to a variety of leg anatomies. The Formfit OA Ease can also be machine washed for effortless care.

## INDICATIONS

Unicompartmental osteoarthritis of the knee

Other conditions that may benefit from unicompartmental unloading, such as degenerative meniscus

## FEATURES & BENEFITS

Easy Wraparound Design - Simple rubber pull-tabs, colour-coded guide dots and a convenient hinge alignment guide make application and removal simple

3-Point Leverage System with Dynamic Force System™ strap and shells unload the affected knee compartment

Stretchy, lightweight, and breathable fabric feels great against the skin and can be machine washed for effortless care

Tool-free polycentric hinge with available ROM stop options for range-of-motion restriction

Fine-tune strap tension for custom pain relief with marked dosing indicator lines

Extension kits are available to further enhance fit and comfort

KNEE OA SOLUTIONS

# FORMFIT® OA EASE


## ORDERING INFORMATION

Size	Calf circ. 15 cm below mid patella	Thigh circ. 15 cm above mid patella	Right Medial	Left Medial	Right Lateral	Left Lateral
XS	28 - 32 cm	36 - 41 cm	FOA0518711	FOA0618711	FOA0528711	FOA0628711
S	32 - 36 cm	41 - 46 cm	FOA0518712	FOA0618712	FOA0528712	FOA0628712
M	36 - 41 cm	46 - 53 cm	FOA0518713	FOA0618713	FOA0528713	FOA0628713
L	41 - 46 cm	53 - 60 cm	FOA0518714	FOA0618714	FOA0528714	FOA0628714
XL	46 - 51 cm	60 - 67 cm	FOA0518715	FOA0618715	FOA0528715	FOA0628715
XXL	51 - 56 cm	67 - 75 cm	FOA0518716	FOA0618716	FOA0528716	FOA0628716
XXXL	56 - 61 cm	75 - 81 cm	FOA0518717	FOA0618717	FOA0528717	FOA0628717

## ACCESSORIES

Description	Part #
Small Extension Kit 2 Pieces	FOA000301
Small Extension Kit 10 Pieces	FOA000302
Large Extension Kit XS-SM (Right Medial/Left Lateral)	FOA000401
Large Extension Kit MD-LG (Right Medial/Left Lateral)	FOA000402
Large Extension Kit XL-XXXL (Right Medial/wLeft Lateral)	FOA000403
Large Extension Kit XS-SM (Right Lateral/Left Medial)	FOA000404
Large Extension Kit MD-LG (Right Lateral/Left Medial)	FOA000405
Large Extension Kit XL-XXXL (Right Lateral/Left Medial)	FOA000406
Wash Bag	WB000001
Suspension Strap	OAW00007
Flexion/Extension Stop Kit	B-301040047

## REQUIRED MEASUREMENTS



Measure calf circumference 15 cm below mid-patella. Leg in extension.



Measure thigh circumference 15 cm above mid-patella. Leg in extension.