

ÖSSUR FORMFIT® PRO KNEE

Össur Formfit® Pro Knee is a compressive knee sleeve that features a SmartPad™ buttress for targeted support and cushioning. Our patented design realigns the patella, suitable for conditions such as Patella Tendonitis (Jumper’s Knee) and Chondromalacia Patella. The dynamic, 3-D MotionTech™ knit is supportive, meets certified medical grade compression requirements and has been proven to be superior to the leading brands at moisture wicking and breathability, keeping your knee cool during wear. Lightweight, removable spiral stays ensure optimum shape and fit, while distinct patella graphics make donning intuitive.



MotionTech™



CoolVent™



StayBilizer™



SmartPad™

INDICATIONS FOR USE

For knee conditions that may benefit from compression, such as;

- Patella tendonitis
- Chondromalacia Patella
- Mild sprains and strains
- Knee pain
- Feeling of instability

CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit with superior moisture wicking and breathability, keeping your knee cool during wear
- Certified medical grade compression (20-30 mmHg)
- Lightweight CoolVent™ knit in popliteal area for added breathability and comfort
- Patented SmartPad™ for targeted support and cushioning
- Patella graphics for intuitive placement
- Patented FoldLock™ system for easier and more secure application
- Lightweight, removable spiral steel stays

REQUIRED MEASUREMENTS



Measure calf circumference 15 cm (6") below mid-patella. Leg in extension.



Measure thigh circumference 15 cm (6") above mid-patella. Leg in extension.

Size	Calf circumference 15 cm (6") below mid patella		Thigh circumference 15 cm (6") above mid patella	
	cm	inch	cm	inch
XS	26 - 31	10¼ - 12¼	36 - 41	14⅛ - 16⅛
S	31 - 36	12¼ - 14⅛	41 - 46	16⅛ - 18⅛
M	36 - 41	14⅛ - 16⅛	46 - 51	18⅛ - 20⅛
L	41 - 46	16⅛ - 18⅛	51 - 56	20⅛ - 22
XL	46 - 50	18⅛ - 19¾	56 - 60	22 - 23⅝

ÖSSUR FORMFIT® PRO KNEE

ORDERING INFORMATION

BLACK

BLUE

Size	Part#	Size	Part#
XS	FP100001	XS	FP101001
S	FP100002	S	FP101002
M	FP100003	M	FP101003
L	FP100004	L	FP101004
XL	FP100005	XL	FP101005

