

FORMFIT PRO ANKLE

Össur Formfit® Pro Ankle is a compressive ankle sleeve with 2 integrated malleoli pads, designed to worn after mild ankle injuries. The dynamic 3-D MotionTech™ knit is supportive and meets certified medical grade compression requirements. Proven superior to leading brands for moisture wicking and breathability, your ankle stays cool during wear. CoolVent™ technology in the area in front of the ankle further enhances breathability and comfort.



MotionTech™



CoolVent™



SmartPad™



FoldLock™

INDICATIONS FOR USE

For Ankle conditions that may benefit from compression, such as:

- Ankle injury rehabilitation
- Mild sprains and strains
- Mild instability

CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated

FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit with superior moisture wicking and breathability, keeping your ankle cool during wear
- Certified medical grade compression (20-30 mmHg)
- CoolVent™ technology in the area in front of the ankle for added breathability and comfort
- 2 integrated malleoli SmartPads™ for targeted support and cushioning
- Patented FoldLock™ system for easier and more secure application

REQUIRED MEASUREMENTS



Measure ankle circumference.



BLUE

Size	Ankle circumference		Blue		Black	
	cm	inch	Left	Right	Left	Right
XS	18 - 20	7 - 7¾	FP201201	FP201101	FP200201	FP200101
S	20 - 23	7¾ - 9	FP201202	FP201102	FP200202	FP200102
M	23 - 26	9 - 10¼	FP201203	FP201103	FP200203	FP200103
L	26 - 29	10¼ - 11½	FP201204	FP201104	FP200204	FP200104
XL	29 - 32	11½ - 12½	FP201205	FP201105	FP200205	FP200105