

FORMFIT® PRO ELBOW

Formfit® Pro Elbow provides relief from elbow pain. This 3D compressive elbow sleeve is manufactured using our MotionTech™ technology and features two integrated epicondyle SmartPads™ for more targeted support and cushioning. The RadialStrap™ provides additional compression and support.



MotionTech™



CoolVent™



SmartPad™



FoldLock™

INDICATIONS FOR USE

- For elbow conditions that may benefit from compression, such as:
- Medial and lateral epicondylitis
- Elbow pain and swelling


CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

FEATURES AND BENEFITS

- 3D MotionTech™ for compression and dynamic fit keep
- Epicondyle SmartPads™ for targeted support and cushioning
- Patented FoldLock™ system for easier and more secure application
- RadialStrap™ for increased compression and support

REQUIRED MEASUREMENTS



Measure arm circumference exactly at the elbow. Elbow in full extension.

ORDERING INFORMATION

Size	Elbow circumference		Black	Blue
	cm	inch		
XS	17 - 20	6¾ - 7⅞	FP300001	FP301001
S	20 - 23	7⅞ - 9	FP300002	FP301002
M	23 - 26	9 - 10¼	FP300003	FP301003
L	26 - 29	10¼ - 11¾	FP300004	FP301004
XL	29 - 32	11¾ - 12½	FP300005	FP301005

