

FORMFIT® PRO KNEE

Össur Formfit® Pro Knee is a compressive knee sleeve that features a SmartPad™ buttress for targeted support and cushioning. Our patented design realigns the patella, suitable for conditions such as Patella Tendinitis (Jumper's Knee) and Chondromalacia Patella. The dynamic, 3D MotionTech™ knit is supportive, meets certified medical grade compression requirements and has been proven to be superior to the leading brands at moisture wicking and breathability, keeping your knee cool during wear. Lightweight, removable spiral stays ensure optimum shape and fit, while distinct patella graphics make donning intuitive.



MotionTech™



CoolVent™



StayBilizer™



SmartPad™

INDICATIONS FOR USE

For knee conditions that may benefit from compression, such as:

- Patella tendinitis
- Chondromalacia Patella
- Mild sprains and strains
- Knee pain
- Feeling of instability


CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit with superior moisture wicking and breathability, keeping your knee cool during wear
- Certified medical grade compression (20-30 mmHg)
- Lightweight CoolVent™ knit in popliteal area for added breathability and comfort
- Patented SmartPad™ for targeted support and cushioning
- Patella graphics for intuitive placement
- Patented FoldLock™ system for easier and more secure application
- Lightweight, removable spiral steel stays

REQUIRED MEASUREMENTS



Measure calf circumference 15 cm below mid-patella. Leg in extension.



Measure thigh circumference 15 cm above mid-patella. Leg in extension.



Size	Calf circumference 15 cm below mid patella	Thigh circumference 15 cm above mid patella	Black
	cm	cm	Part#
XS	26 - 31	36 - 41	FP100001
S	31 - 36	41 - 46	FP100002
M	36 - 41	46 - 51	FP100003
L	41 - 46	51 - 56	FP100004
XL	46 - 50	56 - 60	FP100005