

FORMFIT® TENNIS ELBOW

Formfit Tennis Elbow is made with a patented, ergonomically-shaped plastic. There are two blue cushions with ribbed surface on each side that reduce traction on the tendons. The durable design and materials enable everyday use for both work and athletic activities.

INDICATIONS FOR USE

- Tendonitis of the elbow
- Prevention during work, daily life and sport
- Epicondylitis

FEATURES AND BENEFITS

- Adjustable strap
- Patented ergonomically-shaped plastic eliminates tourniquet effect associated with traditional designs.
- Twin blue cushions provide counterforce, reducing traction on the tendons and thus inflammation
- Ribbed surfaces reduce vibration to the tendon, increasing comfort
- Durable design and materials enable everyday use for both work and athletic activities

REQUIRED MEASUREMENTS

- Circumference of the elbow



ORDERING INFORMATION

Size	Circumference	Part#
Small	22-26cm	B-247000002
Medium	27-32cm	B-247000003
Large	33-38cm	B-247000004