

## FORMFIT® TENNIS ELBOW

### INDICATIONS FOR USE

- Epicondylitis

### FEATURES AND BENEFITS

- Patented ergonomically-shaped plastic eliminates tourniquet effect associated with traditional designs.
- Twin blue cushions provide counterforce, reducing traction on the tendons and thus inflammation.
- Ribbed surfaces reduce vibration to the tendon, increasing comfort.
- Durable design and materials enable everyday use for both work and athletic activities.



### REQUIRED MEASUREMENTS

- Circumference of the elbow

### ORDERING INFORMATION

Size	Circumference	Part#
Small	8.7-10.4" (22-26 cm)	B-247000002
Medium	10.4-12.8" (27-32 cm)	B-247000003
Large	12.8-15" (33-38 cm)	B-247000004

## AIRFORM® TENNIS ELBOW SUPPORT

### INDICATIONS FOR USE

- Epicondylitis

### FEATURES AND BENEFITS

- Focused pneumatic compression
- Conforms and comfortably cushions sensitive injured anatomy
- Easy to apply, simple design that fits most, reduces inventory needs



### REQUIRED MEASUREMENTS

- Universal: Fits elbow circumference 20 – 36 cm / 8-14"

### ORDERING INFORMATION

Description	Size	Part#
AirForm Tennis Elbow Support	Universal (fits elbow circumference 8-14" / 20 - 36 cm)	320000