

FORMFIT® THUMB SPICA

The Formfit Thumb Spica is designed to restrict wrist and thumb motion. It provides immobilisation and protection for patients suffering from tendonitis (e.g. DeQuervain's Tendonitis), select soft tissue injuries – including mild to moderate sprains and strains and Gamekeeper's Thumb –, as well as for post-cast healing. The lightweight, breathable material and contour fit make this comfortable product an excellent choice.



INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
 - e.g., DeQuervain's Tendonitis
- Post-cast healing
- Other selected soft tissue injuries
 - e.g., Gamekeeper's Thumb

FEATURES AND BENEFITS

- Breathable material
- Contoured fit
- Contact closure straps
- Customisable stays

REQUIRED MEASUREMENTS

- Wrist circumference

ORDERING INFORMATION

Size	Wrist Circumference	Left	Right
XSmall	13 - 16cm	3020	3120
Small	16 - 18.5cm	3030	3130
Medium	18.5 - 20.5cm	3050	3150
Large	20.5 - 23.5cm	3070	3170
XLarge	23.5 - 25.5cm	3080	3180

ACCESSORIES

GEL PAD FOR HOT/COLD THERAPY

Description	Left	Right
Form Fit Wrist (8" version)*	GB20L	GB20R
Form Fit Wrist & Forearm (10" version)	GB22L	GB22R



* Fits XSmall/Small. For Medium-Xlarge, please use Form Forearm 10" part numbers