

FORMFIT® WRIST & THUMB, UNIVERSAL

Formfit Universal Wrist and Thumb products are designed to restrict motion of the wrist and the thumb. Providing immobilization and protection from select soft tissue injuries, mild to moderate sprains and sprains, including gamekeeper's thumb and post cast healing. The lightweight, universal design also reduces inventory needs.

INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
- Post-cast healing
- Other selected soft tissue injuries

FEATURES AND BENEFITS

- Quick-Pull close strap
- Positional dorsal stay
- Customizable thumb stay
- Malleable palmar stay

REQUIRED MEASUREMENTS

- Wrist circumference



ORDERING INFORMATION

FORMFIT WRIST

Size	Wrist Circumference	Left	Right
Pediatric	3.5 – 5" (9 – 13 cm)	B-252603301	B-252603300
8" Universal	4.5 – 9.5" (11 – 24 cm)	B-252603109	B-252503109
8" X-Large	7.5 – 11" (19 – 28 cm)	B-252603106	B-252603105

FORMFIT WRIST & FOREARM UNIVERSAL

Size	Wrist Circumference	Left	Right
10" Universal	4.5 – 9.5" (11 – 24 cm)	B-252603210	B-252603209

FORMFIT THUMB SPICA UNIVERSAL

Size	Wrist Circumference	Left	Right
Pediatric	3.5 – 5" (9 – 13 cm)	B-253603301	B-253603300
8" Universal	4.5 – 9.5" (11 – 24 cm)	B-253603109	B-253503109
8" X-Large	7.5 – 11" (19 – 28 cm)	B-253603106	B-253603105