


## ICELOCK® 400 SERIES




### SMOOTH PIN SYSTEM

- No movement, therefore no noise and less wear
- Removable from inside of socket
- Very lightweight
- All impact levels



### DIAGRAM



| Part#     | Description                        | Ratings |        | Build Height |          | Weight |       |  |
|-----------|------------------------------------|---------|--------|--------------|----------|--------|-------|--|
| P20191*   | The Pin Standard Lamination        | N/A     |        | 35mm         | 1 3/8"   | 97g    | 3.4oz |  |
| Part#     | Description                        | Ratings |        | Build Height |          | Weight |       |  |
| P20191-X* | The Pin, Single/Double Lamination  | N/A     |        | 37mm         | 1 7/16"  | 95g    | 3.4oz |  |
| Part#     | Description                        | Ratings |        | Build Height |          | Weight |       |  |
| P20192*   | The Pin for BOA-System, Lamination | N/A     |        | 47mm         | 1 14/16" | 110g   | 3.9oz |  |
| Part#     | Description                        | Ratings |        | Build Height |          | Weight |       |  |
| P20194*   | The Pin, w/ 4-Hole Lamination      | 120kg   | 220lbs | 34mm         | 1 5/16"  | 180g   | 6.3oz |  |

\* Note: Reusable fixtures and lamination dummies not included in lock kits. Must be ordered separately.

## ICELOCK® 400 SERIES

### FIXTURES (REUSABLE)

| Part#  | Description               |
|--------|---------------------------|
| P20109 | Fixture for all lock kits |

### DUMMIES (REUSABLE)

| Part#  | Description                 |
|--------|-----------------------------|
| P20151 | Dummy for P20191 & P20191-X |
| P20152 | Dummy for P20190 & P20192   |
| P20154 | Dummy for P20194            |

### SPARE PARTS AND DRAPE FORMERS (REUSABLE)

| Part#    | Description                          |
|----------|--------------------------------------|
| P20141   | Coupling Plate for P20191            |
| P20141-X | Coupling Plate for P20191-X          |
| P20142   | Coupling Plate for P20192            |
| P20144   | 4-hole Adapter for P20194            |
| P20151   | Dummy for 191/191-X                  |
| P20161   | Funnel for P20191 & P20191-X         |
| P20162   | Funnel for P20192                    |
| L-693506 | Attachment Pin Smooth, Standard 43mm |
| L-693507 | Attachment Pin Smooth, 53mm          |
| L-693508 | Attachment Pin Smooth, 63mm          |
| L-693509 | Attachment Pin Smooth, 75mm          |
| P20177   | The Pin, Replacement Lock            |
| L-693504 | Metal Sleeve for Push Button         |
| L-693502 | Push Button, Adjustable              |
| P20251   | Insert Part for P20151               |