

## MIAMI LSO®

Miami LSO is a spinal orthosis designed for pain relief, protection of injured ligaments or muscles and post-surgical immobilisation. It provides low-profile, anatomical immobilisation, easy fit and wear, options for modular support and an advanced pulley system with a length-adjustable pulley cord.



### INDICATIONS FOR USE

Conditions requiring gross immobilisation of the trunk in the lumbar region such as:

- Stable non-displaced spinal fractures (L1-L5)
- Spinal Stenosis
- Herniated Disks
- Degenerative Spinal Pathologies
- Spondylolysis
- Spondylolisthesis

### FEATURES AND BENEFITS

- Slim, breathable material
- Proven Miami double-pulley system with length-adjustable pulley cord
- Intuitive fitting without trimming with fold-over belt arms & Precise-Size™ landmarks
- Optional fit tool for perfect fitting directly on the patient
- Universal belt sizing accommodates waists up to 178cm
- New Universal Posterior panel to accommodate more patients with less inventory
- Contoured posterior panels sized S-XL for a precise, anatomical fit
- Quick-Fit pocket enclosure and blue highlights for patient ease of use

### REQUIRED MEASUREMENTS

- Waist circumference

## MIAMI LSO®

## ORDERING INFORMATION

## MIAMI LSO BELT

Size	Waist Circumference	Part#
6" Universal (Slim)	66-127cm	MLSO-1E
8" Universal (Standard)	66-127cm	MLSO-2E



## OPTIONAL BELT EXTENDERS (2 PCS)

Size	Waist Circumference	Part#
6" (Slim) Belt Extension	127-177cm	MLSO-1EXT
8" (Standard) Belt Extension	127-177cm	MLSO-2EXT

## STANDARD POSTERIOR PANEL

Size	Part#
Standard Posterior Panel	MLSO-P56



## POSTERIOR PANEL WITH LATERAL SUPPORTS

Size	Waist Circumference	Part#
S	66-86cm	MLSO-P12
M	86-96cm	MLSO-P22
L	96-127cm	MLSO-P32
XL	127-178 cm	MLSO-P42



\*Non-lateral and non-lordotic panels available on request.

## ANTERIOR PANEL

Size	Part#
6" Panel (Slim)	MLSO-PA6
8" Panel (Standard)	MLSO-PA8



## QUICK FIT TOOL

Size	Part#
Quick fit tool	B-040100000

