




# Foot Comparison

For additional information, visit [go.ossur.com/feetanz](https://go.ossur.com/feetanz)



	BALANCE SOLUTIONS				DYNAMIC SOLUTIONS													JUNIOR SOLUTIONS				
	BALANCE™ FOOT S	BALANCE™ FOOT S TORSION	FLEX-FOOT ASSURE®	BALANCE™ FOOT J	PROPRIO FOOT®	PRO-FLEX® TERRA	PRO-FLEX® PIVOT	PRO-FLEX® XC	PRO-FLEX® XC TORSION	PRO-FLEX® LP	PRO-FLEX® LP TORSION	PRO-FLEX® LP ALIGN	PRO-FLEX® ST	PRO-FLEX® MODULAR	VARI-FLEX®	LP VARI-FLEX®	RE-FLEX SHOCK™	RE-FLEX ROTATE™	CHEETAH® XPLORE	FLEX-FOOT® JUNIOR	VARI-FLEX® JUNIOR	
PATIENT WEIGHT LIMIT	147kg				125kg	166kg	125kg	166kg	147kg	166kg	147kg	116kg	166kg					147kg	45kg	55kg		
HEEL HEIGHT	10mm				Adjustable up to 50mm	10mm						Adjustable up to 70mm	10mm		10mm or 19mm	10mm	10mm or 19mm		10mm	10mm		
BUILD HEIGHT†	126mm	204mm	125mm		180mm	179mm	155mm		222mm	73mm	147mm	121mm	160mm	168 - 327mm	170mm	68mm	210mm - 238mm T	216mm - 265mm T	348mm - 521mm	46mm	120mm	
WARRANTY	36 Months		36 Months	24 Months	36 Months																18 Months	
SIZE	22-30		19-30	22-30	22-30						22-28	22-30					One size per category		16-21	16-24		
IMPACT LEVELS	Low	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	
	Moderate																					
	High																					
HEEL HEIGHT ADJUSTABLE					•						•											
BUILT-IN ROTATION		•				•			•		•								•			
VERTICAL SHOCK PYLON		•							•		•							•	•			
UNITY COMPATIBLE*	•	•	•c	•	•c		•	•	•	•	•				•	•	•	•				
FOOT BLADE SHAPE																						
WEATHER RESISTANCE																						

SANDAL TOE: Narrow, anatomically correct foot cover for use with a sandal toe shoe.
 SPLIT TOE: A split toe feature allows the foot to adjust to the surface below, allowing amputees to walk naturally on uneven ground.
 WEATHERPROOF: A weatherproof device can be used in a wet and/or humid environment (but not under water).
 WATERPROOF: Based on the design and materials used, products identified as waterproof have been tested to confirm that water is not a danger to support, function, or durability. Please refer to the instructions for each product for waterproof specifications.
 \* Weatherproof only † Size 27 ‡ PDAC Approved ▼ Weatherproof when used with pyramid or T-connector. Waterproof when used with carbon fiber tube, glass tube or direct lamination. Please refer to product instructions for further details. T All options available, D/P Flexion added, C Compatible with select foot sizes only.
 \*\* The responsibility for accurate coding lies solely with the provider treating the patient. Ossur assumes no responsibility or liability for the provider's coding decisions. Ossur's coding suggestions rest on its best judgment and are subject to revision.