

CANVAS ROCKER BOTTOM CAST SHOE

FOR USE OVER AN EXISTING CAST

CANVAS ROCKER BOTTOM CAST SHOE



FEATURES & BENEFITS

Cushioned rocker bottom cast sandal

Contact closure, open toe

Broad design can accommodate bulky casts

ORDERING INFORMATION

PART #	SIZE	INSOLE LENGTH
308XS	XSmall	up to 17.8cm
308SM	Small	17.8 - 22.9cm
308MED	Medium	25.4 - 26.7cm
308LG	Large	27.9 - 30.5cm
308XLG	XLarge	30.5 - 34.9cm

MESH TOP POST-OP SHOE

MESH TOP POST-OP SHOE



FEATURES & BENEFITS

Breathable mesh upper with reinforced heel counter and padded heel collar, more stability overall to the rear of the foot while reducing slippage

Semi-rigid sole provides superior protection and control after osteotomies

Extended strap and dual side closure system accommodates extra bulky dressings

Compatible with DARCO Toe Guard

ORDERING INFORMATION - MENS

PART #	SIZE	SHOE SIZE
17013	Small	4 - 6
17015	Medium	6.5 - 8.5
17017	Large	9 - 11
17018	XLarge	11.5 - 14

ORDERING INFORMATION - WOMENS

PART #	SIZE	SHOE SIZE
17003	Small	5 - 6
17005	Medium	7 - 8
17007	Large	9+

INDICATIONS

Post-surgery use

Short term footwear option