

MIAMI JTO® THORACIC EXTENSION

The Miami JTO Thoracic Extension is designed in conjunction with leading spine surgeons, orthotists and nurses. The orthosis is used exclusively with all Miami J® collars when even greater stability at the cervical junctions is indicated.



INDICATIONS FOR USE

- | | |
|--|--|
| <p>Lower C-Spine disruptions</p> <ul style="list-style-type: none"> • Level one instabilities • Isolated compartment fracture • Lateral mass fracture • Spinous Process fracture | <p>Upper C-Spine disruptions</p> <ul style="list-style-type: none"> • For increased stability over the collar alone • Stable Jefferson or Hangman’s Fracture • Type I Odontoid Fracture |
|--|--|

FEATURES AND BENEFITS

- Provides more control than a collar alone can provide
- Innovative device designed to treat cervical and high-thoracic injuries
- Patented design with ‘floating back’ permits maximum immobilization without the risk of compensatory motion
- Snap on and off without tools
- Sorbatex II padding is antibacterial and clinically proven to inhibit microbial growth, which enhances comfort and helps assure patient compliance
- MRI, CT and X-ray lucent, MR compatible

REQUIRED MEASUREMENTS

- Universal, one size

ORDERING INFORMATION

MIAMI JTO THORACIC EXTENSION

Description	Part#
Miami JTO	MT-900
Miami JTO with replacement softgoods	MTR-900

Compatible with all sizes of the Miami J®, level 1 through 5 of the Miami J® Select, and the P3 Miami Jr®. Any model of Miami J collar is not included with Miami JTO Thoracic Extension and is purchased separately.

ACCESSORIES

Description	Part#
Replacement softgoods	MTP-950