

MIAMI LUMBAR® POSTEO

Designed with patient comfort and compliance in mind, the Miami Lumbar Posteo provides functional relief of vertebral fracture pain, and facilitates return to mobility. The lumbar compression system offers both pain relief via active vertebral offloading and better grip of patient's hips for greater ease of use and compliance. The Miami Lumbar Posteo includes an innovative approach to ease of use and consistent donning.



INDICATIONS FOR USE

- Back pain caused by stable vertebral compression fractures (VCF)

FEATURES AND BENEFITS

- Active offloading via thoracic extension and lumbar compression
- Integrated pulley system for customized fit
- Trimmable axillary straps with adjustable tension for optimal stabilization
- Cold-moldable aluminum PTE to fit different anatomies
- Extra cushioning around the aluminum strut for minimized skin pressure

REQUIRED MEASUREMENTS

- Back length measured from the top of the sacrum up to C7 vertebrae, following the spinal curvature

Please note that the Miami Lumbar Posteo includes a universal lumbar belt which can be trimmed down to the appropriate size. Calculation of the appropriate back length size, measured from the tip of the sacrum to the C7 vertebra, following the curvature of the spine.

ORDERING INFORMATION

MIAMI LUMBAR POSTEO

Size	Back Length	Part#
X-Small	14.5" – 17" (37 – 44 cm)	LMB-POSTEO-XS
Small	17" – 20" (44 – 51 cm)	LMB-POSTEO-SM
Medium	19" – 22" (49 – 56 cm)	LMB-POSTEO-MD
Large	21" – 24" (55 – 61 cm)	LMB-POSTEO-LG

ACCESSORIES

Description	Part#
Strap Pad	LMB-POSTEO ST PD

