

MIAMI TLSO 464

The Miami TLSO from Össur directly addresses the needs of patients and practitioners, simplifying preparation without compromising on fit or performance. This easy-to-use pre-assembled Miami TLSO brace offers three configurations in one.

Compliance: Easily applied in over the shoulder or underarm configuration, donning and doffing is further simplified with the conveniently located quick release buckles.

Adjustability: The universal belt designs allows for setting the belt length on the patient for a more accurate sizing experience.

Comfort: Truly shaped to fit the anatomy, the anatomically contoured design better immobilises the spine while maximizing patient comfort.

Ease of Use: The Miami TLSO brace requires minimal preparation and is designed to provide optimal patient comfort.



INDICATIONS FOR USE

Conditions requiring gross immobilisation of the thoracic and lumbar spine. Such conditions may include:

- Post-surgical immobilisation
- Spinal support TLSO attachment is up to T7
- Spinal stenosis
- Herniated disc
- Degenerative spinal fractures
- Spondylolithesis
- Spondylolysis

CONTRAINDICATIONS

- Unstable, displaced fractures

FEATURES AND BENEFITS

- Significantly superior comfort and ease of donning/doffing
- Fully assembled for fitting with ease
- Accommodates heights 155-193 cm (5'1" – 6'5")
- Posterior thoracic extension option with ergonomic axillary strapping
- Step-down to an LSO as one progresses through rehabilitation
- One-touch telescoping anterior-thoracic extension option with single or double-hinge
- The convenience of a one-size Universal or individualise with a uniquely-sized anatomical design
- Improve comfort and immobilisation with multiple strapping configurations and removable anchors

REQUIRED MEASUREMENTS

- Waist circumference

MIAMI TLSO 464

ORDERING INFORMATION

ATE/PTE, LORDOTIC

Size	Part#	Waist circumference cm	Waist circumference inch
Universal	MTLSO-2351	66 - 177	26 - 70
S	MTLSO-2111	66 - 86	26 - 34
M	MTLSO-2121	86 - 107	34 - 42
L	MTLSO-2131	107 - 127	42 - 50
XL	MTLSO-2141	127 - 177	50 - 70

*Non-lateral and non-lordotic available on request

*Belt extensions are required for waist circumferences larger than 127 cm / 50"

ATE/PTE, LORDOTIC, SLIM

Size	Part#	Waist circumference cm	Waist circumference inch
Universal	MTLSO-1351	66 - 177	26 - 70
S	MTLSO-1111	66 - 86	26 - 34
M	MTLSO-1121	86 - 107	34 - 42
L	MTLSO-1131	107 - 127	42 - 50
XL	MTLSO-1141	127 - 177	50 - 70

*Belt extensions are required for waist circumferences larger than 127 cm / 50"

OTHER ACCESSORIES

Description	Part#
5MM Hex Long Wrench	MTLSO-0004
Axillary strap accessory kit	MTLSO-0003
Quick fit tool	B-040100000
Belt extension, standard (Recommended for waist circumferences 127-167 cm (50"-66")). Fits size XL of Miami LSO.)	MLSO-2EXT
Belt extension, slim (Recommended for waist circumferences 127-167 cm (50"-66")). Fits size XL of Miami LSO.)	MLSO-1EXT