

# MIAMI TLSO™

MIAMI TLSO™





Watch the Miami TLSO Fitting Video



Miami TLSO™ 464 Universal



Miami TLSO™ 456 and Universal

The Miami TLSO from Össur directly addresses the needs of patients and practitioners, simplifying preparation without compromising on fit or performance. This easy-to-use pre-assembled Miami TLSO brace offers three configurations in one. Easily applied in over the shoulder or underarm configuration, donning and doffing is further simplified with the conveniently located quick release buckles.

## INDICATIONS

Indications that require gross thoracic and lumbar spine immobilisation in all three planes, in which the following conditions may be present:

- Herniated disc
- Degenerative Spinal Pathologies
- Post-Surgical Stabilisation
- Stable, non-displaced spinal fractures
- Spinal stenosis
- Spondylolisthesis
- Spondylolysis

## CONTRAINDICATIONS

Unstable, displaced fractures

## FEATURES AND BENEFITS

- Significantly superior comfort and ease of donning/doffing.
- Fully assembled for fitting with ease.
- Accommodates heights 155–193 cm
- Posterior thoracic extension option with ergonomic axillary strapping.
- Step-down to an LSO as one progresses through rehabilitation.
- One-touch telescoping anterior-thoracic extension option with single or double-hinge.
- The convenience of a one-size Universal or individualize with a uniquely-sized anatomical design.
- Improve comfort and immobilization with multiple strapping configurations and removable anchors.

### MIAMI TLSO 464 - PREASSEMBLED KIT

PART #	DESCRIPTION	SIZE	WAIST CIRC.
MTLSO-2111	Miami TLSO 464	Small	66-86 cm
MTLSO-2121	Miami TLSO 464	Medium	86-107 cm
MTLSO-2131	Miami TLSO 464	Large	107-127 cm
MTLSO-2141	Miami TLSO 464	Extra Large	127-177cm*

\*Belt extensions are required for waist circumferences larger than 127cm

### MIAMI TLSO 456 - PREASSEMBLED KIT

PART #	DESCRIPTION	SIZE	WAIST CIRC.
MTLSO-2112	Miami TLSO 456	Small	66-86 cm
MTLSO-2122	Miami TLSO 456	Medium	86-107 cm
MTLSO-2132	Miami TLSO 456	Large	107-127 cm
MTLSO-2142	Miami TLSO 456	Extra Large	127-177cm*

\*Belt extensions are required for waist circumferences larger than 127cm

### ACCESSORIES

PART #	DESCRIPTION
B-040100000	Ossur Fit Tool
MTLSO-0003	Axillary Strap Accessory Kit
MTLSO-0004	5mm Hex Long Wrench

### MIAMI TLSO 464 & 456 - UNIVERSAL

PART #	DESCRIPTION	SIZE
MTLSO-2351	Miami TLSO 464	Universal
MTLSO-2352	Miami MTLISO 456	Universal

PART #	DESCRIPTION
MLSO-2EXT	20 cm Standard Belt Extension*

\*Recommended for waist circumferences 127-177 cm