

PATRIOT® COLLAR

INDICATIONS FOR USE

- Post-trauma

FEATURES AND BENEFITS

- Multi-height adjustment – Adult (4 positions) and Pediatric (3 positions)
- Open rear ventilation panel for fluid drainage, airflow and cervical spine palpation
- Universal sizing optimises inventory requirements



REQUIRED MEASUREMENTS

- Neck circumference

ORDERING INFORMATION

Size	Neck circumference	Blue
Child	20 – 46 cm	PHP-PP
Adult	28 – 58 cm	PHP-PA