

# PHILADELPHIA® ADJUSTABLE TRACHEOTOMY COLLAR

## INDICATIONS FOR USE

- Post-trauma
- Herniated cervical disc
- Cervical spondylitis
- Post-operative
- Cervical spondylosis
- RA and OA cervical spine

## FEATURES AND BENEFITS

- Circumference and height adjustable
- Permanent lock on height adjustment
- Anatomically shaped Plastazote® foam (Not made with natural rubber latex, non-toxic and hypoallergenic to reduce skin irritation), leaving clearance for cervical veins
- Total cervical arch support maintains neutral alignment of the cervical spine
- Water-resistant for bathing and hydro-therapy.
- Well ventilated
- X-ray and CT lucent. MR conditional.



## REQUIRED MEASUREMENTS

- Neck height: 1 inch (or 25mm) of neck length adjustability accommodates varying patient neck length.

## ORDERING INFORMATION

### PHILADELPHIA ADJUSTABLE TRACHEOTOMY COLLAR

Size	Neck circumference
Small (2 tears)	10 – 13" (25.4 – 33 cm)
Medium (1 tear)	13 – 16" (33 – 40.6 cm)
Large (no tear)	16 – 19" (40.6 – 48.3 cm)

NOTE: The table above represents common sizing only. Sizing must be verified by the healthcare provider for each patient individually.

Part#	Size / Description
PHP-AA	Conforms to sizes 2, 3 and 4 of the original Philadelphia Tracheotomy Collar