

PAPOOSE®

FOR INFANTS

PAPOOSE - FOR INFANTS



Papoose is designed for the youngest and most challenging patients - newborn infants. It cradles the baby's body, securely and gently immobilizing the head and spine.

INDICATIONS

C-Spine precaution for infants

Immobilisation for pre and post c-spine surgery

Other conditions requiring gross immobilisation of the infant cervical spine

FEATURES & BENEFITS

X-ray and CT lucent, MR safe

ORDERING INFORMATION

PART #	DESCRIPTION	RANGE
PA-100	PAPOOSE W/ Front and Pads	Front and pads Infants up to 3 months
PAP-100	Replacement pad	

PHILADELPHIA® ADJUSTABLE COLLAR

PHILADELPHIA® ADJUSTABLE COLLAR



INDICATIONS

Situations requiring gross immobilisation of the cervical spine, these may include:

Pre- and post c-spine surgery

C-spine precaution for trauma patients

FEATURES & BENEFITS

Circumference adjustable and height adjustable with permanent lock.

Well ventilated and water-resistant for bathing and hydrotherapy.

X-ray and CT lucent. MR safe.

ORDERING INFORMATION

PART#	SIZE/ DESCRIPTION
PHP-AA	Conforms to sizes 2, 3 and 4 of the original Philadelphia Tracheotomy Collar

NOTE: The opposite table represents common sizing only. Sizing must be verified by the healthcare provider for each patient individually.

SIZING INFORMATION

CIRCUMFERENCE SIZE	RANGE
Small (remove 2 tabs)	10"-13" (25.4-33.0 cm)
Medium (remove 1 tab)	13"-16" (33.0-40.6 cm)
Large (no tear)	16"-19" (40.6-48.3 cm)

CLINICAL TIPS

The adjustable collar fits most patients but is not suited for long or thick neck shapes. For patients with long necks, the 5" height sized original model should be considered. For patients with neck circumference exceeding 48cm, the XLG sized original model should be considered. Refer to next page for sizing guidelines.