



Pro-Flex[®] Terra



“I’m rediscovering aspects of my mobility that I had forgotten about since my amputation.”

Albert Lin, Ph.D.

Explorer, Entrepreneur, Communicator



Tom Triol is a recently retired photographer and videographer for the Indiana State Police. He keeps busy walking, traveling, working around the house and yard, and birdwatching with his wife, Shelley. Tom needs a foot that allows him to enjoy his life pain and worry free.

Pro-Flex Terra's flexibility means that it remains flat on the ground at most angles. Tom receives more push-off power and reduced sound side impact due to the **1 Pre-Tensioned Carbon**. Furthermore, the **2 Single Anchor Point with Gliding Contact Pad** provides a "neutral range" rather than a "neutral point" on the foot, resulting in superior ground compliance. The **3 Pre-Compressed Foam** allows Tom to experience a smoother transition from heel to toe at any speed, leading to a more natural gait (which Shelley noticed right away).

"I can do any amount of yard work, we can go to concerts, we can enjoy functions with family and friends. I can help clean the house. I enjoy all of my hobbies without pain and without restriction. With Pro-Flex Terra, I can stay active and comfortable all day."

Tom Triol
Retired Police Photographer

Albert Lin is a National Geographics explorer and father of two. His show *Lost Cities with Albert Lin* can involve caving, rock scrambling, hiking, and fording rivers. Albert needs a foot that can keep up with him, a **Waterproof** foot that he can trust in any terrain.

Pro-Flex Terra's exceptional flexibility means Albert's foot is always in contact with the ground, no matter how steep the slope is. The **1 Pre-Tensioned Carbon** unit provides consistent push-off power at those extreme angles. During Albert's high impact activities, the friction generated by the **2 Single Anchor Point with Gliding Contact Pad** delivers increased energy return and the **3 Pre-Compressed Foam** bumper dampens the impact of every step. All of these innovations in Pro-Flex Terra make Albert feel like the foot is matching his energy and can accommodate his demanding lifestyle.

"There are moments when I should be looking down to see what's underneath me, but now I look up to see what is in front of me. I'm rediscovering aspects of my mobility that I had forgotten about since my amputation."

Albert Lin, Ph.D.
Explorer, Entrepreneur, Communicator



1 Pre-Tensioned Carbon

- Increased flexibility provides greater range of motion
- Increased comfort while standing
- Same force to achieve higher energy return

2 Single Anchor Point with Gliding Contact Pad

- Accommodates minor shifts in weight while standing
- Provides a neutral range during midstance by allowing the point of contact to move as weight shifts from heel to toe
- Provides multi-axial rotation for better ground adaptation



3 Pre-Compressed Foam

- Provides progressive stiffness across a range of impact levels
- Smoother transition from heel to toe
- Increased dampening at higher walking speeds



Learn more about our Pro-Flex Terra


We collaborate extensively with prosthetic users during product development. Building from that feedback, we created lifestyle features that tackle issues that patients face every day. From a snap-on footshell to optional accessories like a sandal fix or heel wedges, each combines with Pro-Flex Terra to provide great performance across a diverse range of daily requirements.



NEW! Footshell – easy-to-fasten buttons enable users to quickly remove the foot cover for cleaning or drying



Tread Design – Ensures traction and stability with or without the foot shell



Heel Wedges† – an optional heel wedge allows for extra stability while doing vigorous activities



SandalFix† – allows users to wear open heel footwear by securely attaching the heel of the foot cover to the sandal

Contact your Össur Clinical Specialist to learn more!

* Whether it's fresh water, salt water or chlorinated water, Waterproof Össur products will give you the maximum protection you need. The use of first-rate materials provides permanent protection against harmful ingress of water, even after submersion, whilst continuing to provide the functionality you expect.

† Optional accessories sold separately.


WWW.OSSUR.COM.AU
WWW.OSSUR.CO.NZ

Össur ANZ
TEL AU 1300 123 268
TEL NZ 0800 369 524
anzprosthetics@ossur.com

FAX AU +61 2 9475 1114
FAX NZ 0800 448 265

 **ÖSSUR**[®]
LIFE WITHOUT LIMITATIONS