

REBOUND® CARTILAGE

The success of an articular cartilage repair is dependent on patient compliance to the rehabilitation process. The protocol is based on restricted cartilage load combined with gradual increase in activity. Many rehabilitation protocols suggest crutches for partial weight-bearing and it is the standard of care for successful outcomes. However, patient compliance to partial weight-bearing instructions with crutches has been shown to be inconsistent. The Rebound Cartilage is a protective functional solution designed to support regeneration of the knee cartilage after repair procedures such as Microfracture, OATS and (M)ACI. Based on the clinically-proven 3-Point Leverage System and patented Cartilage Protection Straps, joint unloading is maintained in flexion for protection of the cartilage during the healing process.



INDICATIONS FOR USE

For knee conditions that may benefit from unicompartmental load reduction or range of movement restriction, such as:

- Articular cartilage defect repair
- Meniscal repair
- Avascular necrosis
- Condylar bone marrow lesions (bone bruises)
- Tibial plateau fractures
- Mild to severe unicompartmental knee osteoarthritis

FEATURES AND BENEFITS

- Easy-to-Stock - Just 4 configurations to fit most of your patients
- Comfortable - Skin-friendly liner feels great against the skin and reduces pressure points.
- Anti-Migration Technology - Strap design captures the thickest part of the calf to reduce brace migration during activity.
- Control - SmartDosing® dials put you in control of your pain by allowing you to fine-tune strap tension, while Cartilage Protection Straps maintain joint unloading in flexion.
- Tool-Free - Tool-free polycentric hinge with optional ROM stops for range-of-motion restriction.

REQUIRED MEASUREMENTS

- Circumference measurement: 6" (15 cm) below mid-patella

REBOUND® CARTILAGE

ORDERING INFORMATION

Rebound Cartilage	Left medial	Right medial	Left lateral	Right lateral
Universal 31 - 61.5 cm (12.2 - 24")	B-254610006	B-254510006	B-254620006	B-254520006

Rebound Cartilage - No Box	Left medial	Right medial	Left lateral	Right lateral
Universal 31 - 61.5 cm (12.2 - 24")	B-254610007	B-254510007	B-254620007	B-254520007

Rebound Cartilage - Custom	Left	Right
Calf circumference 25 - 61.5 cm (9.9 - 24")	B-154600001	B-154500001

SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



Description	Part#
SmartMeasure Measuring Kit	CMS0000
Constant Force Caliper	B-312000000

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

CUSTOM STANDARD COLORS

Red	Pink	Silver
Orange	Black	White
Yellow	Navy Blue	Champagne
Lime Green	Ocean Blue	Charcoal

CUSTOM PAINTED DESIGNS



Additional fees apply for custom colors, custom painted designs, and custom printed graphics.

ACCESSORIES

Accessory Parts	Part#
Suspension Strap Kit	UNL00211
OA Hinge Extension Stop Tree (5°, 10°, 15°, 30°)	OA40029
OA Hinge Flexion Stop Tree (0°, 30°, 60°, 90°)	OA40039
Wash Bag	UON60200

REBOUND® CARTILAGE

Sleeves	Black Ossur One*	Black - Sleeve with Mesh**	Black - Ossur Knee Undersleeve***
XXS	OSO60021	NA	NA
XS	OSO60022	SLV60011	B-250150101
S	OSO60023	SLV60012	B-250150102
M	OSO60024	SLV60013	B-250150103
L	OSO60025	SLV60014	B-250150104
XL	OSO60026	SLV60015	B-250150105
XXL	OSO60027	SLV60016	B-250150106
3XL	OSO60028	NA	NA

*Ossur One - 51% Tencel/30% Polyamide/19% Elastane

**Sleeve with Mesh - Main fabric: 92% Cotton/8% Spandex - Mesh: 100% Polyester - Band: 80% Polyester/20% Spandex

***Ossur Knee Undersleeve: 80% Nylon/20% Spandex

REPLACEMENT PARTS

Service Kit*	Right medial/Left lateral	Right lateral/Left medial
Black, Service Kit	CRT416901	CRT415901

*Includes Jersey knit thigh liner, Sensil calf liner, Popliteal pad standard, Hinge pad & Strap Kit Complete

Thigh Liners	Right medial/Left lateral	Right lateral/Left medial
Black, Jersey Knit	CRT526901	CRT525901

Calf Liners	Right medial/Left lateral	Right lateral/Left medial
Black, Sensil	CRT516952	CRT515952
Black, Jersey Knit	CRT516951	CRT515951

Strap Kit Complete*	Part #
Black, Strap Kit Complete	B-701290005

*Includes 4 straps (1 Thigh, 1 Calf, 2 CPS), 5 Alligator tabs, 2 CPS Strap pads

CPS Strap Kit*	Part #
Black, CPS Kit	B-700290005

*Includes 2 CPS Straps, 2 Alligator tabs, 2 CPS Strap pads

Buckle Clasp Kit	Right medial/Left lateral	Right lateral/Left medial
Black	B-708217001	B-708218001

Popliteal Pad	Part #
Black, Short	B-304010105
Black, Standard	B-304010106

Other Replacement Parts	Part #
Velcro Twist-Fit Attachment Kit	B-705290000
Hinge Pad	B-304120101
Alligator Kit-8 Pieces	B-704210001