

REBOUND® DUAL

REBOUND® DUAL - OFF THE SHELF



Watch the Rebound Dual Fitting Video

The Rebound DUAL offers simplicity without compromise. With its malleable yet supportive frame, Rebound DUAL is the go-to for knee ligament support.

INDICATIONS

The brace is indicated for combined and non-combined instabilities of the:

- Anterior cruciate ligament (ACL)
- Medial collateral ligament (MCL)
- Lateral collateral ligament (LCL)
- Posterior cruciate ligament (PCL)

Acculign version only: Conditions requiring unicompartamental load control.

FEATURES & BENEFITS

Cold-malleable aluminum frame, trimmable straps, and field-serviceable D-rings accommodate various patient anatomies

Acculign™ upper frame allows height and varus/valgus adjustment

Polycentric hinge with ROM stop options for range-of-motion restriction

Össur Hinge Alignment Guide, Blue is You™ touch points, and tactile numbered strap tabs guide the fitting process

Optional SmartDosing® Kit accessory offers fine-tuned unicompartamental offloading

Optional PCL Kit accessory for multi-ligament support

REBOUND® DUAL -CUSTOM



CUSTOM PAINT COLOURS:

- Black (default)
- Electric blue
- Metallic silver
- White
- Red

CUSTOM LINER COLOURS:

- Metallic Silver (default)
- Electric blue
- Yellow
- Purple
- Black

READY TO FIT MODEL AVAILABLE IN BLACK FRAME WITH BLACK LINER ONLY.

REBOUND® DUAL

ORDERING INFORMATION - OFF THE SHELF MODEL

PART #	LEG	SIZE	MEASURING	
			CALIPER*	CIRCUMFERENCE**
B-242510001	Right	XS	7.6 - 8.9cm	38.1 - 41.9cm
B-242510002	Right	Small	8.9 - 10.2cm	41.9 - 47cm
B-242510003	Right	Medium	10.2 - 11.4cm	47 - 52.1cm
B-242510004	Right	Large	11.4 - 12.7cm	52.1 - 57.2cm
B-242510005	Right	X-Large	12.7 - 14cm	57.2 - 62.2cm
B-242510006	Right	XX-Large	14 - 15.2cm	62.2 - 69.9cm
B-242510007	Right	XXX-Large	15.2 - 16.5 cm	69.9 - 76cm
B-242610001	Left	XS	7.6 - 8.9cm	38.1 - 41.9cm
B-242610002	Left	Small	8.9 - 10.2cm	41.9 - 47cm
B-242610003	Left	Medium	10.2 - 11.4cm	47 - 52.1cm
B-242610004	Left	Large	11.4 - 12.7cm	52.1 - 57.2cm
B-242610005	Left	X-Large	12.7 - 14cm	57.2 - 62.2cm
B-242610006	Left	XX-Large	14 - 15.2cm	62.2 - 69.9cm
B-242610007	Left	XXX-Large	15.2 - 16.5 cm	69.9 - 76cm

* Caliper measurement: M/L at joint width in a standard, weight-bearing position

** Thigh circumference 15cm above mid-patella

ORDERING INFORMATION - ACCESSORIES

PART #	MEASURING
B-155X10001	SmartMeasure App
B-701090060	SmartDosing Kit RM/LL
B-701090061	SmartDosing Kit RL/LM
B-701090088	Flexion/Extension Stop Set - Tool-Free Hinge
B-701090059	D-Ring Kit
B-701090001	Tibial Wedge Kit
B-701090009	Condyle Pad Kit

* Flexion Stops: 0, 10, 20, 30, 45, 60, 75, 90

** Extension Stops: 0, 5, 10, 20, 30, 45

REBOUND DUAL CUSTOM - MEASUREMENT CHART

Rebound DUAL Custom can be manufactured between one shape difference for thigh, calf and M/L.

SHAPE	M/L (CM)		* CALF (CM)		SHAPE	M/L (CM)		THIGH (CM)	
	MIX/	MAX	MIN	MAX		MIX/	MAX	MIN	MAX
A	7.5 - 9		27	29	A	7.5 - 9		38	42
B	9 - 10		29	34	B	9 - 10		42	47
C	10 - 11.5		34	39	C	10 - 11.5		47	52
D	11.5 - 13		39	45	D	11.5 - 13		52	57
E	13 - 14		45	50	E	13 - 14		57	62
F	14 - 15		50	54	F	14 - 15		62	70
G	15 - 17		54	60	G	15 - 17		70	76

MORE ACCESSORIES AVAILABLE AT: [OSSUR.COM/EN-AU/BRACING-AND-SUPPORTS/KNEE/REBOUND-DUAL](https://www.ossur.com/en-au/bracing-and-supports/knee/rebound-dual)