

# RHEO KNEE®

UNITY  
COMPATIBLE

## INTRODUCTION

The RHEO KNEE is a microprocessor-controlled knee designed to promote stability and safety. It enables natural and effortless gait – even on difficult terrain. The RHEO KNEE is easy to fit, and delivers reliable user outcomes.

## KEY FEATURES

- Auto adaptive real-time stance and swing control.
- Dependable and easy swing initiation.
- Automatic stumble recovery.
- Manual extension lock.
- Weatherproof design allows exposure to fresh water splashes.
- Accompanying Össur Logic app for iOS devices provides:
  - Simple practitioner setup and adjustment, and printable user activity reports.
  - Functional training exercise.
  - Selection of user preferences and monitoring of step count and battery charge level.



## RECOMMENDED WITH



Seal-In® X TF liner



Pro-Flex® Pivot®



Proprio Foot®



Pro-Flex® XC



Pro-Flex® XC  
Torsion

\*Note. Pro-Flex Pivot can not be used for high impact activities.




# RHEO KNEE®




User Information	
Amputation Level:	Transfemoral and Knee Disarticulation
Impact Level:	Low to Moderate
Maximum Patient Weight:	136kg (300lbs)
Knee Information	
Net weight:	Approx. 1.6 kg (3.5 lbs)
Build Height:	236mm (9 1/4")
Knee Flexion:	120°
Aluminum Frame Construction	
Weatherproof	

## CLEARANCE

RHEO KNEE



RHEO KNEE WITH THREADED TOP ADAPTER (A - 145310)



# RHEO KNEE®

## KITS

### RHEO KNEE KIT

Part#	Description
RKN130003	Rheo Knee incl. 3 year warranty

Includes:

- RHEO KNEE
- RHEO KNEE Power Supply w/Plug
- Terms of warranty
- Instructions for use

### RHEO KNEE KIT WITH THREADED TOP ADAPTER:

RHEO KNEE with threaded top adapter is the solution for knee disarticulation amputees and those with long residual limbs.

To order RHEO KNEE with threaded top adapter, order the RHEO KNEE Kit and add Össur's Female Pyramid Insert to the order.

Part#	Description	Height		Weight	
A-145310	Female Pyramid Insert / Threaded Top Adapter	9mm	3/8"	33g	1.16oz



### RHEO KNEE PROTECTOR

Part#	Description
BIC00110	Bionic Protector for RHEO KNEE



# RHEO KNEE®

## EXTENDED WARRANTY

### RHEO KNEE KIT EXTENDED WARRANTY\*


Part#	Description
RKN3GE3YR	Additional 3 Year Warranty (Incl. Service Check at 40 Months)

\* Warranty extensions must be purchased while the product is still under warranty. Maximum total warranty period is 6 years.

## UNITY

### UNITY KITS

Part#	Description	Applicable Foot Kits
UTF0001	Unity TF Valve Kit	All



## POWER SUPPLY

### DUAL CHARGER KIT

The Dual Charger Kit is an optional kit which allows for the simultaneous charging of Rheo Knee / Rheo Knee XC and Proprio Foot.

Part#	Description
CME10110	Power Supply and Cable Attachment



### POWER SUPPLY KIT

Part#	Description
CME50110	Power Supply w/Plug

