

## i-LIMB® ACCESS

A myoelectric, multi-articulating prosthetic hand that offers five individually powered digits and manually rotatable thumb.

- Individually motorised digits with compliant grip and stall-out ability
- Manually rotatable thumb that can be positioned between lateral and oppositional mode for creating different grip positions
- With our speed boost feature, digit speed can be increased by up to 30%
- App control provides instant access to Quick Grips™ at the touch of an icon.
- Muscle control uses specific muscle signals called triggers to instruct the hand to activate a specific grip
- 12 different automated grip options available.
- Auto-Grasp™ – anti drop safety feature for added reassurance
- Biosim™ and My i-Limb™ Apps
- Available in 4 sizes – extra small, small, medium and large
- Multiple covering options available



The i-Limb Access is also available with titanium strengthened fingers, which increase the carry load by 50%. Available in the sizes S / M and L.

### CONTROL METHODS



TRIGGERS  
MUSCLE CONTROL



QUICK GRIPS  
APP CONTROL

Device Weight			
	Extra Small	Small	Medium/Large
QWD	432g	468g	478g
WD	392g	428g	438g
Flexion Wrist	532g	568g	578g
Friction Wrist	427g	463g	473g

Technical Information	
Voltage	7.4 V (nominal)
Max. Current	5 A
Battery Capacity	Rechargeable lithium polymer; 7.4 V (nominal); 2000 mAh capacity; 1300 mAh capacity
Max. hand load limit (static limit)	40kg (Extra Small) 90kg (Small/Medium/Large)
Finger carry load (static limit)	20kg (Extra Small) 32kg (Small/Medium/Large)
Time from full open to full close	0.8 seconds

## i-LIMB® ACCESS

### i-LIMB ACCESS PACKAGES

i-LIMB ACCESS PACKAGES		
TBX50420	i-Limb Access	PK i-Limb NA Access QWD
TBX50421	i-Limb Access	PK i-Limb NA Access WD
TBX50424	i-Limb Access	PK i-Limb NA Access Flexion
TBX50425	i-Limb Access	PK i-Limb NA Access Friction
TBX50422	i-Limb Access - Titanium	PK i-Limb NA Access Titanium QWD
TBX50423	i-Limb Access - Titanium	PK i-Limb NA Access Titanium WD
TBX50426	i-Limb Access - Titanium	PK i-Limb NA Access Titanium Flexion
TBX50427	i-Limb Access - Titanium	PK i-Limb NA Access Titanium Friction

All packages include:

- Socket Kit
- Battery
- Mains Charger
- Car Charger
- 2 i-Limb Skin Active/TS
- 4 i-Limb Skin Contour/TS
- Years 1 & 2 Touch Care Warranty

TOUCH CARE WARRANTY EXTENSION*	
TBX50173	Year 3
TBX50174	Year 4
TBX50175	Year 5

\* Years 1 & 2 touch care warranty is included in the packages