

Össur Legs for balance users

Digital orderform

Liner

Socket

Adapter

Knee


Tube

Foot



ÖSSUR LEGS - OUR AIM IS TO MAKE DOING BUSINESS EASY

We make it easy to order a prosthetic leg, configured by you to fit your patients' needs. You are in control of every step of the way as Össur Legs can be tailored to fit your practice needs. Let us, help you improve:



Outcomes

In today's O&P environment, your greatest value is in strengthening the relationships essential to growing your business. Whether providing next-level care to your patients, networking with referral sources, or collecting outcomes data to support claims, there are more demands on your time than ever.

With Össur Legs, you'll spend less time purchasing and collecting parts, and lowers the administrative burden to comply to quality and regulatory matters. Instead, you can spend your time just where you want to, with your patients. Our high-quality products, are vigorously tested and verified



Efficiency

to work together. On the following pages you will find our product range to fulfill the specific needs of – balance users. Please contact us for our other product ranges.

7 reasons to use Össur Legs for Balance users:

- Components designed to be used together
- One place to select and order
- One place of contact
- All products comply to EU MDR
- Use of safe and compatible products
- Attractive package offering
- Clinical support



Patient Satisfaction

FOLLOW THE NEXT STEPS TO ORDER AN OSSUR LEG FOR BALANCE USERS

- 1** Select the configuration.
- 2** Fill out this form > select amount, size and/or other information.
- 3** Send this form by using the 'send button' at the bottom of the last page.

SELECT YOUR SOLUTION

Solution	Foot	Adapter(s)*	Knee	Liner	Socket	Discount
TF Solution + liner + Socket	✓	✓	✓	✓	✓	20%
TF Solution + liner	✓	✓	✓	✓		12%
TF Solution	✓	✓	✓			8%
TT Solution + liner + Socket	✓	✓	n.a.	✓	✓	20%
TT Solution + liner	✓	✓	n.a.	✓		11%
TT Solution	✓	✓	n.a.			7%

*A solution must at least contain 3 adapters.

SELECT THE LINER

Please select when not applicable

Sizing TT liners:

Standard: 16, 18, 20, 22, 23.5, 25, 26.5, 28, 30, 32, 34, 36, 40

Sizing Seal Rings:

20, 22, 24, 26, 28, 30, 32, 35, 38, 41, 44, 47, 51, 55, 60, 65

Sizing TF liners:

Standard: 25, 26.5, 28, 30, 32, 34, 36, 38, 40, 45, 50, 55
Conical: 23.5, 25, 26.5, 28, 30, 32, 34, 36, 38, 40

TRANSTIBIAL LOCKING LINERS

Iceross Dermo Locking	Size	Pieces
Dermo Locking 3mm TT Conical		
Dermo Locking 3mm TT Standard		
Dermo Locking 3mm TT Wave		
Dermo Locking 6mm TT Std		

Iceross Synergy Locking	Size	Pieces
Synergy Locking 3mm TT Standard		
Synergy Locking 3mm TT Wave		

Iceross Comfort Locking	Size	Pieces
Comfort Locking 3mm TT Standard		
Comfort Locking 3mm TT Wave		
Comfort Locking 6mm TT Standard		

TRANSTIBIAL CUSHION LINERS

Iceross Comfort Cushion	Size	Pieces
Comfort Cushion 3mm TT Standard		
Comfort Cushion 6mm TT Standard		

Iceross Dermo Cushion	Size	Pieces
Dermo Cushion 3mm TT Standard		
Dermo Cushion 3mm TT Wave		
Dermo Cushion 6mm TT Standard		

Iceross Synergy Cushion	Size	Pieces
Synergy Cushion 3mm TT Standard		
Synergy Cushion 3mm TT Wave		

TRANSTIBIAL SEAL-IN LINERS

Iceross Dermo Seal-In (HSM)	Size	Pieces
Dermo 3mm TT Standard		

Iceross Seal-In X5	Size	Pieces
Seal-In X5 3mm TT Standard		
Seal-In X5 3mm TT Wave		

Iceross Seal-In V	Size	Pieces
Seal-In V 3mm TT Standard		
Seal-In V 3mm TT HP		

Iceross Seal-In X*	Size	Pieces
Seal-In X 3mm TT		
Seal-In X 6mm TT		

*Order in combination with Seal-In Seal.

TRANSFEMORAL LINERS

Iceross Transfemoral	Size	Pieces
Transfemoral Locking Conical		
Transfemoral Locking Standard		

Iceross TF Seal-In (HSM)	Size	Pieces
Seal-In TF Conical		
Seal-In TF Standard		

Iceross Seal-In X5	Size	Pieces
Seal-In X5 TF Conical		
Seal-In X5 TF Standard		

Iceross Seal-In X*	Size	Pieces
Seal-In X TF Conical		
Seal-In X TF Standard		

*Order in combination with Seal-In Seal.

TRANSTIBIAL SEAL-IN X LOCKING

Iceross Seal-In X Locking*	Size	Pieces
Seal-In X Locking 3mm TT		
Seal-In X Locking 6mm TT		

*Order in combination with Seal-In Seal.

SEAL-IN SEALS

Iceross Seal-In Seals	Size	Pieces
Iceross Seal-In X-Classic S20-S65		
Iceross Seal-In X-Grip S20-S65		
Iceross Seal-In X-Volume S20-S65		

SELECT THE SOCKET SOLUTION

Please select when not applicable

DIRECT SOCKET TT & DIRECT SOCKET TF



Direct Socket Kits Basalt	Pieces
Direct Socket, 4" Basalt Fibers with standard distal connector, 200ml Black Resin	
Direct Socket, 5" Basalt Fibers with standard distal connector, 200 ml Black Resin	
Direct Socket, 7" Basalt Fibers with standard distal connector, 400ml Black Resin	
Direct Socket, 7" Basalt Fibers with large distal connector, 400ml Black Resin	
Direct Socket, 9" Basalt Fibers with large distal connector, 400ml Black Resin	

For Direct Socket TF we recommend the use of the Iceross Transfemoral Locking liner, TF Seal-In and Seal in X5. For Direct Socket TT we recommend the use of Seal-In X Locking.



Össur Brim (combine with DS TF)	Size	Pieces
Össur Brim, size 40.0	39 - 42cm	
Össur Brim, size 42.5	42 - 44cm	
Össur Brim, size 45.0	44 - 47cm	
Össur Brim, size 47.5	47 - 49cm	
Össur Brim, size 50.0	49 - 52cm	
Össur Brim, size 52.5	52 - 54cm	
Össur Brim, size 55.0	54 - 57cm	
Össur Brim, size 57.5	57 - 59cm	
Össur Brim, size 60.0	59 - 62cm	
Össur Brim, size 62.5	62 - 64cm	
Össur Brim, size 65.0	64 - 67cm	
Össur Brim, size 67.5	67 - 69cm	
Össur Brim, size 70.0	69 - 72cm	

Direct Socket can only be ordered in combination with a complete leg configuration as displayed in the table 'select your solution' at the previous bottom left page.

CONNECT TF



Connect TF Socket	Pieces
Medium Short, left side	
Medium Standard, left side	
Large Short, left side	
Large Standard, left side	
Medium Short, right side	
Medium Standard, right side	
Large Short, right side	
Large Standard, right side	

SELECT THE KNEE (WITH TF SOLUTIONS)

Locking Knee	
Amputation level	TF
Max. body weight	125kg
Build height	26mm
Pieces	

OFM1 Knee with pyramid	
Amputation level	TF
Max. body weight	136kg
Build height	158mm
Pieces	

KFM1 Knee with loop adapter	
Amputation level	KD
Max. body weight	136kg
Build height	150mm
Pieces	

Balance Knee OFM2	
Amputation level	TF
Max. body weight	125kg
Build height	115mm
Pieces	

OP4 Knee	
Amputation level	TF
Max. body weight	100kg
Build height	198mm
Pieces	

Total Knee 1900	
Amputation level	TF
Max. body weight	100kg
Build height	173mm
Pieces	

OHP3 with pyramid	
Amputation level	TF
Max. body weight	125kg
Build height	210mm
Pieces	

KHP3 Knee with loop adapter	
Amputation level	KD
Max. body weight	125kg
Build height	202mm
Pieces	

SELECT THE ADAPTERS (MIN. 3)

4-Hole socket adapter		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		

4-Prong socket adapter		
Weight limit	166kg	
Material	Stainless Steel	
Pieces		

3-Prong socket adapter		
Weight limit	166kg	
Material	Stainless Steel	
Pieces		

Male insert for prong		
Weight limit	100kg	166kg
Material	Stainless Steel	Titanium
Pieces		

Female insert for prong		
Weight limit	100kg	166kg
Material	Stainless Steel/Titanium	
Pieces		

Height adjustable pylon		
Weight limit	100kg	
Material	Aluminium/Stainless Steel	
Length	Standard 7,7-9,5cm	Long 9,4-12,8cm
Pieces		

4-Hole male pyramid		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		

14mm offset tube female adapter		
Weight limit	100kg	
Material	Aluminium	
Pieces		

4-Hole female pyramid		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		

4-Hole male pyramid with 5/16" / 8mm offset		
Weight limit	166kg	
Material	Aluminium	
Pieces		

Female pyramid tube clamp		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		

Knee rotation adapter		
Weight limit	125kg	
Material	Stainless Steel	
Pieces		

4-Hole pyramid 20/20 offset			
Weight limit	125kg		
Material	Aluminium		
Length	Standard	Left 10"	Right 10"
Pieces			

Female pylon		
Weight limit	100kg	
Material	Aluminium	
Length	Standard 250mm	Long 400mm
Pieces		

Female pylon		
Weight limit	166kg	
Material	Aluminium/ Titanium	
Length	Standard 250mm	Long 400mm
Pieces		

Icelock 562 Hybrid for Seal-In X Locking	Weight limit	Pieces
Icelock 562 Hybrid for unity	166kg	
Icelock 562 Hybrid for passive vacuum	166kg	

PLEASE SEE NEXT PAGE FOR MORE ADAPTER OPTIONS

SELECT THE ADAPTERS (MIN. 3)



Icelock 200	Weight limit	Pieces
Icelock 211 clutch	n.a.	
Icelock 214 clutch 4-Hole	120kg	



Icelock 500	Pieces
Icelock Expulsion Valve 551	
Icelock Expulsion Valve 552	
Icelock Expulsion Kit 544	



Icelock 600 without adapter	Pieces
Icelock Ratchet 621	
Icelock Lanyard 631	
Icelock Valve 641	
Icelock Smooth 651	

SELECT THE FOOT



Flex-Foot Balance	
Build height	62mm
Weight limit	136kg
Foot size 21-28	
Category 1-5	
Foot cover beige or brown	
Side: left or right	

Order with Unity for TT solution Order with Unity for TF solution

Category selection chart						
Weight kg	36-44	45-61	62-80	81-97	98-114	115-136
21	0	1	N.V.	N.V.	N.V.	N.V.
22	0	1	N.V.	N.V.	N.V.	N.V.
23	0	1	2	N.V.	N.V.	N.V.
24	0	1	2	N.V.	N.V.	N.V.
25	N.V.	1	2	3	4	5
26	N.V.	1	2	3	4	5
27	N.V.	N.V.	2	3	4	5
28	N.V.	N.V.	2	3	4	5



Flex-Foot Balance J	
Build height	120mm
Weight limit	136kg
Foot size 21-30	
Category 1-5	
Foot cover beige or brown	
Side: left or right	

Order with Unity for TT solution Order with Unity for TF solution

Category selection chart					
Weight kg	45-61	62-80	81-97	98-114	115-136
22	1	2	N.V.	N.V.	N.V.
23	1	2	3	4	N.V.
24	1	2	3	4	N.V.
25	1	2	3	4	5
26	1	2	3	4	5
27	1	2	3	4	5
28	N.V.	2	3	4	5
29	N.V.	N.V.	3	4	5
30	N.V.	N.V.	N.V.	4	5



Flex-Foot Assure	
Build height	125mm
Weight limit	136kg
Foot size 19-30	
Category 1-5	
Foot cover beige or brown	
Side: left or right	

Order with Unity for TT solution Order with Unity for TF solution

Category selection chart					
Weight kg	45-61	62-80	81-97	98-114	115-136
19-20*	1	2	N.V.	N.V.	N.V.
21-22*	1	2	N.V.	N.V.	N.V.
23-24	1	2	3	4	N.V.
25-26	1	2	3	4	5
27-28	1	2	3	4	5
29-30	N.V.	2	3	4	5

*Unity for Flex-Foot Assure is available from sizes 22 and in all categories.



Talux*	
Build height	182mm
Weight limit	147kg
Foot size 23-30	
Category 1-8	
Foot cover beige or brown	
Side: left or right	

*Talux is not available with Unity.

Category selection chart									
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
23-24	1	1	2	3	4	5	N.V.	N.V.	N.V.
25-26	N.V.	N.V.	2	3	4	5	6	7	N.V.
27-28	N.V.	N.V.	N.V.	3	4	5	6	7	8
29-30	N.V.	N.V.	N.V.	N.V.	4	5	6	7	8



Balance Foot S	
Build height	126mm
Weight limit	147kg
Foot size 22-28	
Category 1-8	
Foot cover beige or brown	
Side: left or right	

Order with Unity for TT solution Order with Unity for TF solution

Category selection chart									
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
22-24	1	1	2	3	4	5	N.V.	N.V.	N.V.
25-27	1	1	2	3	4	5	6	7	8
27-28	3	3	3	3	4	5	6	7	8



Balance Foot S Torsion	
Build height	204mm
Weight limit	147kg
Foot size 22-28	
Category 1-8	
Foot cover beige or brown	
Side: left or right	

Order with Unity for TT solution Order with Unity for TF solution

Category selection chart									
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
22-24	1	1	2	3	4	5	N.V.	N.V.	N.V.
25-27	1	1	2	3	4	5	6	7	8
27-28	3	3	3	3	4	5	6	7	8



YOU ARE **ALMOST READY**

Please follow the next steps to finalise the order:

- Check the specific conditions that apply to your Össur Leg configuration
- Please fill out the details down below and click the send button to submit the order or use the reference code when placing the order.

Conditions

The discount only applies to a configuration for one individual user/patient. However it is possible to order multiple products for trial purposes.

The invoice discount only relates to the purchase of the selected product configuration. In the situation where returning products results in incomplete configurations, the standard conditions agreed with your organisation apply. We will charge the difference of the discount.

Name: _____ Contact person: _____

Phone: _____ E-mail: _____

Customer number: _____ Reference: _____

Date: _____

Please first send me a quote:

Comments:

Clicking on the yellow button should open an external window to send the e-mail. If nothing happens, please attach this orderform to an e-mail and send it to info@ossur.com.

Thank you very much and please contact our customer care team if you have questions.