



## Össur Legs with dynamic TT solutions

Digital orderform





# ÖSSUR LEGS - OUR AIM IS TO MAKE DOING BUSINESS EASY

We make it easy to order a prosthetic leg, configured by you to fit your patients' needs. You are in control of every step of the way as Össur Legs can be tailored to fit your practice needs. Let us, help you improve:



**Outcomes**

In today's O&P environment, your greatest value is in strengthening the relationships essential to growing your business. Whether providing next-level care to your patients, networking with referral sources, or collecting outcomes data to support claims, there are more demands on your time than ever.

With Össur Legs, you'll spend less time purchasing and collecting parts, and lowers the administrative burden to comply to quality and regulatory matters. Instead, you can spend your time just where you want to, with your patients. Our high-quality products, are vigorously tested and verified



**Efficiency**

to work together. On the following pages you will find our product range to fulfill the specific needs of active, transtibial users. Please contact us for our other product ranges.

### 7 reasons to use Össur Legs with Proprio Foot & Pro-Flex Feet:

- Components designed to be used together
- One place to select and order
- One place of contact
- All products comply to EU MDR
- Use of safe and compatible products
- Attractive package offering
- Clinical support



**Patient Satisfaction**

## FOLLOW THE NEXT STEPS TO ORDER AN OSSUR LEG

- 1 Select the configuration.
- 2 Fill out this form > select amount, size and/or other information.
- 3 Send this form by using the 'send button' at the bottom of the last page.

## SELECT YOUR SOLUTION

Solution	Foot	Adapter(s)*	Liner	Socket	Discount
TT Solution + liner + Direct Socket	✓	✓	✓	✓	20%
TT Solution + liner	✓	✓	✓		11%
TT Solution	✓	✓			7%

\*A solution must at least contain 3 adapters.

## SELECT THE LINER

Please select when not applicable

### Sizing TT liners:

Standard: 16,18, 20, 22, 23.5, 25, 26.5, 28, 30, 32, 34, 36, 40

### Sizing Seal Rings:

20, 22, 24, 26, 28, 30, 32, 35, 38, 41, 44, 47, 51, 55, 60, 65

### TRANSTIBIAL LOCKING LINERS

Iceross Dermo Locking	Size	Pieces
Dermo Locking 3mm TT Conical		
Dermo Locking 3mm TT Standard		
Dermo Locking 3mm TT Wave		
Dermo Locking 6mm TT Std		

Iceross Synergy Locking	Size	Pieces
Synergy Locking 3mm TT Standard		
Synergy Locking 3mm TT Wave		

Iceross Comfort Locking	Size	Pieces
Comfort Locking 3mm TT Standard		
Comfort Locking 3mm TT Wave		
Comfort Locking 6mm TT Standard		

### TRANSTIBIAL CUSHION LINERS

Iceross Comfort Cushion	Size	Pieces
Comfort Cushion 3mm TT Standard		
Comfort Cushion 6mm TT Standard		

Iceross Dermo Cushion	Size	Pieces
Dermo Cushion 3mm TT Standard		
Dermo Cushion 3mm TT Wave		
Dermo Cushion 6mm TT Standard		

Iceross Synergy Cushion	Size	Pieces
Synergy Cushion 3mm TT Standard		
Synergy Cushion 3mm TT Wave		

### TRANSTIBIAL SEAL-IN LINERS

Iceross Dermo Seal-In (HSM)	Size	Pieces
Dermo 3mm TT Standard		

Iceross Seal-In X5	Size	Pieces
Seal-In X5 3mm TT Standard		
Seal-In X5 3mm TT Wave		

Iceross Seal-In V	Size	Pieces
Seal-In V 3mm TT Standard		
Seal-In V 3mm TT HP		

Iceross Seal-In X*	Size	Pieces
Seal-In X 3mm TT		
Seal-In X 6mm TT		

\*Order in combination with Seal-In Seal.

### TRANSTIBIAL SEAL-IN X LOCKING

Iceross Seal-In X Locking*	Size	Pieces
Seal-In X Locking 3mm TT		
Seal-In X Locking 6mm TT		

\*Order in combination with Seal-In Seal.

### SEAL-IN SEALS

Iceross Seal-In Seals	Size	Pieces
Iceross Seal-In X-Classic S20-S65		
Iceross Seal-In X-Grip S20-S65		
Iceross Seal-In X-Volume S20-S65		

## SELECT THE SOCKET SOLUTION

Please select when not applicable


### DIRECT SOCKET TT




Direct Socket Kits Basalt	Pieces
Direct Socket, 4" Basalt Fibers with standard distal connector, 200ml Black Resin	
Direct Socket, 5" Basalt Fibres with standard distal connector, 200 ml Black Resin	
Direct Socket, 7" Basalt Fibers with standard distal connector, 400ml Black Resin	

Direct Socket can only be ordered in combination with a complete leg configuration as displayed in the table 'select your solution' at the bottom left page.


## SELECT THE ADAPTERS (MIN. 3)




4-Hole socket adapter		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		




4-Prong socket adapter		
Weight limit	166kg	
Material	Stainless Steel	
Pieces		




3-Prong socket adapter		
Weight limit	166kg	
Material	Stainless Steel	
Pieces		




Male insert for prong		
Weight limit	100kg	166kg
Material	Stainless Steel	Titanium
Pieces		




Female insert for prong		
Weight limit	100kg	166kg
Material	Stainless Steel/Titanium	
Pieces		




Height adjustable pylon		
Weight limit	100kg	
Material	Aluminium/Stainless Steel	
Length	Standard 7,7-9,5cm	Long 9,4-12,8cm
Pieces		




4-Hole male pyramid		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		




14mm offset tube female adapter		
Weight limit	100kg	
Material	Aluminium	
Pieces		




4-Hole female pyramid		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		




4-Hole male pyramid with 5/16" / 8mm offset		
Weight limit	166kg	
Material	Aluminium	
Pieces		



Female pyramid tube clamp		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		



Female pylon		
Weight limit	100kg	
Material	Aluminium	
Length	Standard 250mm	Long 400mm
Pieces		



Female pylon		
Weight limit	166kg	
Material	Aluminium/ Titanium	
Length	Standard 250mm	Long 400mm
Pieces		




Icelock 562 Hybrid for Seal-In X Locking		
Weight limit	166kg	
Pieces		


Weight limit	166kg	
Pieces		


Weight limit	166kg	
Pieces		



Icelock 200	Weight limit	Pieces
Icelock 211 clutch	n.a.	
Icelock 214 clutch 4-Hole	120kg	



Icelock 500	Pieces
Icelock Expulsion Valve 551	
Icelock Expulsion Valve 552	
Icelock Expulsion Kit 544	
Icelock 544 Direct Socket Adapter Plate	
Icelock 544 Plate Expulsion	
Icelock 544 Plate Unity	



Icelock 600 without adapter	Pieces
Icelock Ratchet 621	
Icelock Lanyard 631	
Icelock Valve 641	
Icelock Smooth 651	

## SELECT THE FOOT



Proprio Foot	
Amputation level	Transtibial and Transfemoral
Impact level	Low to Moderate
Maximum patient Weight	125kg
Build height	180mm
Heel height Adjustability	Up to 50mm
Foot size 22-30	
Category 1-8	
Foot cover beige or brown	
Side: left or right	

Order with Unity

Category selection chart								
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-125
Low Impact Level	1	1	2	3	4	5	6	7
Moderate Impact Level	1	2	3	4	5	6	7	8

Unity can be ordered with Proprio Foot in sizes 25 and above.



Pro-Flex LP Align	
Build height	147mm
Weight limit	116kg Low impact 100kg Moderate impact
Foot size 22-30	
Category 1-9	
Foot cover beige or brown	
Side: left or right	

Order with Unity

Category selection chart							
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116
<b>Low Impact Category</b>							
Size 22-24	1	1	2	3	4	5	N/A
Sizes 25-27	1	1	2	3	4	5	6
Sizes 28	3	3	3	3	4	5	6
<b>Moderate Impact Category</b>							
Size 22-24	1	2	3	4	5	N/A	N/A
Sizes 25-27	1	2	3	4	5	6	N/A
Sizes 28	3	3	3	4	5	6	N/A

PLEASE SEE NEXT PAGE FOR MORE FOOT OPTIONS

# SELECT THE FOOT



Pro-Flex XC	
Build height	155mm
Weight limit	166kg
Foot size 22-30	
Category 1-9	
Foot cover beige or brown	
Side: left or right	

Order with Unity



Pro-Flex XC Torsion	
Build height	222mm
Weight limit	147kg
Foot size 22-30	
Category 1-9	
Foot cover beige or brown	
Side: left or right	

Order with Unity



Pro-Flex LP	
Build height	73mm
Weight limit	166kg
Foot size 22-30	
Category 1-9	
Foot cover beige or brown	
Side: left or right	

Order with Unity



Pro-Flex LP Torsion	
Build height	147mm
Weight limit	147kg
Foot size 22-30	
Category 1-9	
Foot cover beige or brown	
Side: left or right	

Order with Unity

Category selection chart											
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166	
<b>Low Impact Category</b>											
Sizes 22-24	1	1	2	3	4	5	6	N/A	N/A	N/A	
Sizes 25-27	1	1	2	3	4	5	6	7	8	9	
Sizes 28-30	3	3	3	3	4	5	6	7	8	9	
<b>Moderate Impact Category</b>											
Sizes 22-24	1	2	3	4	5	6	N/A	N/A	N/A	N/A	
Sizes 25-27	1	2	3	4	5	6	7	8	9	N/A	
Sizes 28-30	3	3	3	4	5	6	7	8	9	N/A	
<b>High Impact Category</b>											
Sizes 22-24	2	3	4	5	6	N/A	N/A	N/A	N/A	N/A	
Sizes 25-27	2	3	4	5	6	7	8	9	N/A	N/A	
Sizes 28-30	3	3	4	5	6	7	8	9	N/A	N/A	

Category selection chart										
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	
<b>Low Impact Category</b>										
Sizes 22-24	1	1	2	3	4	5	6	N/A	N/A	
Sizes 25-27	1	1	2	3	4	5	6	7	8	
Sizes 28-30	3	3	3	3	4	5	6	7	8	
<b>Moderate Impact Category</b>										
Sizes 22-24	1	2	3	4	5	6	N/A	N/A	N/A	
Sizes 25-27	1	2	3	4	5	6	7	8	N/A	
Sizes 28-30	3	3	3	4	5	6	7	8	N/A	
<b>High Impact Category</b>										
Sizes 22-24	2	3	4	5	6	N/A	N/A	N/A	N/A	
Sizes 25-27	2	3	4	5	6	7	8	N/A	N/A	
Sizes 28-30	3	3	4	5	6	7	8	N/A	N/A	

Category selection chart											
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166	
<b>Low Impact Category</b>											
Sizes 22-24	1	1	2	3	4	5	N/A	N/A	N/A	N/A	
Sizes 25-27	1	1	2	3	4	5	6	7	8	9	
Sizes 28-30	3	3	3	3	4	5	6	7	8	9	
<b>Moderate Impact Category</b>											
Sizes 22-24	1	2	3	4	5	N/A	N/A	N/A	N/A	N/A	
Sizes 25-27	1	2	3	4	5	6	7	8	9	N/A	
Sizes 28-30	3	3	3	4	5	6	7	8	9	N/A	
<b>High Impact Category</b>											
Sizes 22-24	2	3	4	5	N/A	N/A	N/A	N/A	N/A	N/A	
Sizes 25-27	2	3	4	5	6	7	8	9	N/A	N/A	
Sizes 28-30	3	3	4	5	6	7	8	9	N/A	N/A	

Category selection chart										
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	
<b>Low Impact Category</b>										
Sizes 22-24	1	1	2	3	4	5	N/A	N/A	N/A	
Sizes 25-27	1	1	2	3	4	5	6	7	8	
Sizes 28-30	3	3	3	3	4	5	6	7	8	
<b>Moderate Impact Category</b>										
Sizes 22-24	1	2	3	4	5	N/A	N/A	N/A	N/A	
Sizes 25-27	1	2	3	4	5	6	7	8	N/A	
Sizes 28-30	3	3	3	4	5	6	7	8	N/A	
<b>High Impact Category</b>										
Sizes 22-24	2	3	4	5	N/A	N/A	N/A	N/A	N/A	
Sizes 25-27	2	3	4	5	6	7	8	N/A	N/A	
Sizes 28-30	3	3	4	5	6	7	8	N/A	N/A	



Pro-Flex Pivot	
Build height	155mm
Weight limit	125kg
Foot size 22-30	
Category 1-7	
Foot cover beige or brown	
Side: left or right	

Order with Unity

Category selection chart								
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-125
<b>Low Impact Category</b>								
Sizes 22	1	1	2	3	4	N/A	N/A	N/A
Sizes 23-24	1	1	2	3	4	5	N/A	N/A
Sizes 25-26	2	2	2	3	4	5	6	7
Sizes 27-28	N/A	N/A	N/A	3	4	5	6	7
Sizes 29-30	N/A	N/A	N/A	N/A	4	5	6	7
<b>Moderate Impact Category</b>								
Sizes 22	1	2	3	4	N/A	N/A	N/A	N/A
Sizes 23-24	1	2	3	4	5	N/A	N/A	N/A
Sizes 25-26	3	3	3	4	5	6	7	N/A
Sizes 27-28	N/A	N/A	4	5	6	6	7	N/A
Sizes 29-30	N/A	N/A	N/A	5	6	7	8	N/A





## YOU ARE ALMOST READY

**Please follow the next steps to finalise the order:**

- Check the specific conditions that apply to your Össur Leg configuration
- Please fill out the details down below and click the send button to submit the order or use the reference code when placing the order.

**Conditions**

The discount only applies to a configuration for one individual user/patient. However it is possible to order multiple products for trial purposes.

The invoice discount only relates to the purchase of the selected product configuration. In the situation where returning products results in incomplete configurations, the standard conditions agreed with your organisation apply. We will charge the difference of the discount.

Name: \_\_\_\_\_ Contact person: \_\_\_\_\_

Phone: \_\_\_\_\_ E-mail: \_\_\_\_\_

Customer number: \_\_\_\_\_ Reference: \_\_\_\_\_

Date: \_\_\_\_\_

Please first send me a quote:

Comments:

\_\_\_\_\_

Clicking on the yellow button should open an external window to send the e-mail. If nothing happens, please attach this orderform to an e-mail and send it to [info@ossur.com](mailto:info@ossur.com).

Thank you very much and please contact our customer care team if you have questions.