

CUSTOM MADE SEAL-IN®

ORDER FORM

Name of workshop: _____

Name of CPO: _____

Product number: _____

Product name: _____

Amputation level TF TT KD Other


Liner profile Cylindrical Conical

Reason for order _____

Full description of liner _____

MEASUREMENT INSTRUCTIONS

- Distance from proximal edge of seal to lowest edge of socket brim (A) must be at least 2,5 cm (1") to ensure suspension.
- All measurements should be done on bear skin.



<p>Liner Size.</p> <p>Circumference at 4 cm distance from the distal end (cm).</p>	<p>B= ____ cm</p>	
<p>Distance from the distal end where the seal-in ring will be placed (cm).</p> <p>It is not possible to offset the seal-in ring more distally than standard.</p>	<p>C= ____ cm</p>	
<p>Circumference where the lower edge of the seal-in ring will be located.</p>	<p>D= ____ cm</p>	

This product is custom made and not to be CE market nor tested according to CE requirements. This product is covered with warranty until initial fit has been performed.

I agree with the above statement:

Signed _____ Date ____ / ____ / ____