



NAKED
PROSTHETICS™

PIPDriver™ | 2nd GENERATION

A robust, custom, functional solution.

The PIPDriver is a body-driven prosthesis designed for finger and partial-hand limb difference at the middle or distal phalanx of digits two through five.

The design of this device anatomically matches the PIP and DIP joints so it operates predictably and naturally.

IDEAL ANATOMICAL PRESENTATIONS FOR FIT:

This device is driven by an intact PIP joint with enough residuum to engage the cage.

For a person with finger or partial-hand limb difference who doesn't meet the typical requirements for the PIPDriver, it is still possible to obtain functional benefit from an added device. Our engineers customise the design to fit your patient, resulting in a functional device

Our intake form will walk you through the steps of ordering a device using a Naked Prosthetics sizing kit, submitting detailed pictures of both hands, and a video showing range of motion.



OTHER COMFORT ASPECTS INCLUDE:

- Design customised to fit comfortably and improve function for activities of daily living.
- Easy to clean, care for, don and doff.
- Made from strong and durable medical grade nylon 12.
- Cage-like structure protects residuum.
- Slender and sleek design allows for wearing of multiple devices.
- Conforming silicone rubber tip pads resemble the natural fingertip.
- Lightweight titanium fasteners.
- Conductive Tip option for touchscreen use.

STYLE + COATINGS

PIPDriver is available in two nail bed styles, Organic and Minimalistic.

To view our palette of colour coatings please visit npdevices.com.



ORGANIC



MINIMALIST