

Össur Balance Leg Trial Form

STEP 1 General Information

Prosthetist Name: _____ Acc Nr: _____

Patient Identifier: _____ Weight: _____ Kg Date of Birth: ____ / ____ / ____

PLEASE NOTE TRIAL WILL ONLY BE MADE AVAILABLE FOR A MAXIMUM TRIAL PERIOD OF 30 (THIRTY) DAYS.

STEP 2 Select the Icross TF Locking Liner specifications *The liner will be invoiced and can't be returned for credit

Standard Conical Size: _____


STEP 3 Select the Connect TF size and side

CONNECT[®] TF

Size Selection

1. Have the user sit down
2. Measure the circumference of the residual limb at 40mm from the distal end (A) **Note.** It is important to take this measurement with the user sitting
3. Measure the length of the residual limb from the distal end to perineum (B)
4. Select the correct size from the sizing guidelines (C) or size selection table below

Sizing Guidelines



Size Selection Table

Connect TF Size	Liner size (circumference 40mm from distal)	Residual limb length
Medium Short	26.5 – 32 (265-320mm)	200-240mm
Medium Standard	26.5 – 32 (265-320mm)	240-280mm
Large Short	34 – 40 (340-400mm)	200-240mm
Large Standard	34 – 40 (340-400mm)	240-280mm

Ordering Information

Part#	Description	Part#	Description
CTF0202L	Medium Short, left side	CTF0202R	Medium Short, right side
CTF0203L	Medium Standard, left side	CTF0203R	Medium Standard, right side
CTF0204L	Large Short, left side	CTF0204R	Large Short, right side
CTF0205L	Large Standard, left side	CTF0205R	Large Standard, right side

Note: A locking pin is included

STEP 4 The Balance Knee OFM2 is included

STEP 5 Complete Balance Foot S specifications **Foot will be credited after return

Category: _____ Foot Size: _____ Left/Right: _____
FOOT COVER Beige Brown

STEP 6 Select an Adapter. A female Pylon Long AL is included with all three options

OPTION 1

4-HOLE FEMALE PYRAMID WITH
20/20
OFFSET

OPTION 2

FEMALE PYRAMID 20/20
OFFSET 10° RIGHT

OPTION 3

4-HOLE FEMALE PYRAMID 20/20
OFFSET 10° LEFT

Please note delivery period is 14 days. Fitting support request from Össur specialist: Yes No

By completing this form it is understood that you accept that the trial period for all relevant products is valid for only 30 days. Upon the completion of the 30-day trial period, you have a further 7-day grace period within which to return all trial products. Should the product(s) not be returned to Össur within the grace period, penalties may be incurred.