

Össur Custom Solutions Ohio 6640 Riverside Drive, Suite 360 Dublin, OH 43017 Email: tb_orders@ossur.com

FACILITY INFORMATION (Billing)


Account #	Shipping same as billing address		
Clinician Name	Ship to Address		
Company Name			
Company Address			
City	State	Zip	
Telephone			
Purchase Order #			
Sales Quotation #	Ship Via	Next Day	2 Day
Please provide email address or fax number for order confirmation			

PATIENT INFORMATION

Patient Identifier	Insurance Company
--------------------	-------------------

AFFECTED SIDE(S)



Left	Right
-------------	--------------



Extra Small
Extra Small mounting plate with a 5, 5, 3, 2 digit configuration will be supplied as standard. If an alternative configuration is required, then please note digit strength may be compromised.



Medium
Medium mounting plate with a 6, 6, 5, 3 digit configuration will be supplied as standard. If an alternative configuration is required, then please note digit strength may be compromised.



INDICATE DIGITS AND SIZE

Left Prosthesis Configuration	Right Prosthesis Configuration
thumb	thumb
index	index
middle	middle
ring	ring
pinky	pinky

ACCESSORIES:

Covers:	Black	Clear	
Optional (Check all that apply):			
Touch Care Warranty Extension:	Year 3	Year 4	Year 5

NOTES:

FABRICATION PROCESS:

Please take impressions using the instructions below. Watch the video here: <https://www.youtube.com/watch?v=dr2AqyGaAQo>

Note: It is recommended that you refrigerate the silicone 24 hours ahead of taking the impressions.



Step 1:

Apply Vaseline to the affected hand and arm up to the elbow. If patient has hair, pat down with Vaseline and then wipe the excessive Vaseline off with a paper towel.

Step 2:

Take impression up to the elbow crease on the affected side.



Step 3:

When the silicone is completely cured, remove the mold. Have the patient wiggle their digit(s) or wrist gently while you pull the mold off. The silicone can be ripped slightly at the top for ease of removal.

SHIP IMPRESSIONS TO:

Össur Custom Solutions Ohio 6640 Riverside Drive, Suite 360 Dublin, OH 43017

Email: tb_orders@ossur.com



page 2 of 2

FOLLOW ÖSSUR ON



USA (800) 233-6263
WWW.OSSUR.COM

CANADA (800) 663-5982
WWW.OSSUR.CA

