

TT Socketless Össur Leg Order Form

tel: (888)839-6213, fax: (800)788-9878 | cs@ossur.com | www.ossur.com

Area Manager Information

Area Manager: _____ phone: _____ email: _____
 Notes: _____

BILLING & SHIPPING INFORMATION

Össur Account #: _____ PO# _____
 Ship To: _____
 Company: _____ Contact: _____
 Address: _____ Phone: _____ Fax: _____
 City/State/Zip: _____ Email: _____













Completion of this order form with the most accurate and up-to-date information, including all patient information, measurements, and construction information enables us to provide the highest quality product for you and your patient.

Required Patient Information

Last Name: _____ Height: ft. _____ in. Impact Level: High Med Low Amputation Level: TT Symes
 First Name: _____ Weight: lbs. K-Level: K1 K2 K3 K4 Left Right Bilat Order

1: Liners - Two Are ALWAYS Included with Each Leg (check one then fill in desired sizes and options)

Send liners first for measurements/scanning/casting

<p>Comfort Locking</p> 	<p>Dermo Locking</p> 	<p>Iceform Locking</p> 	<p>LINER SIZE</p> <p>Standard Wave 3 mm 6 mm</p>	<p>Seal-In V</p> 	<p>Seal-In V HP</p> 	<p>Seal-In X5</p> 	<p>LINER SIZE</p> <p>Standard Wave</p>
<p>Comfort w/ Sleeve</p> 	<p>Dermo w/ Sleeve</p> 	<p>Iceform w/ Sleeve</p> 	<p>LINER SIZE</p> <p>Standard Wave 3 mm 6 mm</p> <p>SLEEVE SIZE</p>	<p>Seal-In X with Rings</p> 	<p>Seal-In X Locking with Rings</p> 	<p>(Ring Options)</p>  <p>Classic Volume Grip</p>	<p>RING SIZE</p>
Iceform Sleeve	Iceross Sleeve	Genu Sleeve					

2. Suspension & Socket Adapter for Conventional Socket Fabrication (check all that apply)


Locking Options

Icelock 621 Ratchet & Extra Pin 	Icelock 651 Smooth & Extra Pin 	Icelock 562 Hybrid 	Icelock 600 XM 	Icelock 211 Clutch & Extra Pin Icelock 214 4-Hole Clutch & Extra Pin 	Icelock 234 Lanyard 
Works With					
Icelock 672 Male Pyramid Ti  Icelock 673 Male Pyramid Al	Icelock 674 4-Hole Adapter Al  Icelock 675 4-Hole Adapter Ti				

Expulsion Options






Icelock 544 Expulsion Kit 	Icelock 551 Expulsion Valve 	Unity Valve TT 	Icelock 544 Unity Plate 
Icelock 544 Expulsion Valve Kit for sleeves 	Works With		
	Icelock 544 Socket Adapter 	3-Prong Socket Adapter 4-Prong Socket Adapter	

Other (write in):

 Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

3. Direct Socket Suspension Options (Choose when fabricating your socket using a Direct Socket Kit)

Locking































621 Ratchet Pin 	651 Smooth Pin 	214 Clutch Lock 	562 Hybrid Lock with Unity (for TT only) 	562 Hybrid Lock with Passive Vacuum (for TT only) 
Requires: Modular Spacer Plate 				

Add Extra Pin: Long Standard Short

Expulsion

Icelock 551 Expulsion Valve 	Icelock 544 Expulsion Plate 	Unity TT Valve 	Icelock 544 Unity Plate 
Requires: Bottom Seal 	Requires: DS Adapter Plate 	Requires: Bottom Seal 	Requires: DS Adapter Plate 

TT Socketless Össur Leg Order Form




























Products		
4. Endoskeletal Components (check all that apply for a single prosthesis)		
4-Hole Adapters	Double Adapters	Single Adapters
 4-Hole Male <input type="checkbox"/> <input type="checkbox"/>	 Female Double Adapter 32mm <input type="checkbox"/> <input type="checkbox"/>	 Male Single Adapter Short <input type="checkbox"/> <input type="checkbox"/>
 4-Hole Female <input type="checkbox"/> <input type="checkbox"/>	 Female Double Adapter 45mm <input type="checkbox"/> <input type="checkbox"/>	 Male Single Adapter Long <input type="checkbox"/> <input type="checkbox"/>
 4-Hole Male w/ Rotation <input type="checkbox"/> <input type="checkbox"/>	 Female Double Adapter 60mm <input type="checkbox"/> <input type="checkbox"/>	 Female Single Adapter Short <input type="checkbox"/> <input type="checkbox"/>
 4-Hole Male w/ 8mm Ax Offset <input type="checkbox"/>	 Female Double Adapter 75mm <input type="checkbox"/> <input type="checkbox"/>	 Female Single Adapter Long <input type="checkbox"/>
 4-Hole Female w/ Rotation <input type="checkbox"/> <input type="checkbox"/>	Pylon Kits	Pylons
 4-Hole Female Pyramid 20/20 Offset 10° Right <input type="checkbox"/>	 Female Kit Short <input type="checkbox"/> <input type="checkbox"/>	 Female Short <input type="checkbox"/> <input type="checkbox"/>
 4-Hole Female Pyramid 20/20 Offset 10° Left <input type="checkbox"/>	 Female Kit Long <input type="checkbox"/> <input type="checkbox"/>	 Female Long <input type="checkbox"/> <input type="checkbox"/>
Prong Adapters	Height Adjustable Pylons	Tube Clamps
 4-Prong Socket Adapter <input type="checkbox"/>	 Standard <input type="checkbox"/>	 Male Tube Clamp <input type="checkbox"/> <input type="checkbox"/>
 3-Prong Socket Adapter <input type="checkbox"/>	 Long <input type="checkbox"/>	 Female Tube Clamp <input type="checkbox"/> <input type="checkbox"/>
 Male Pyramid Insert for Prong <input type="checkbox"/> <input type="checkbox"/>	Torsion Shock Adapter (Extra Cost)	 4-Hole Tube Clamp <input type="checkbox"/> <input type="checkbox"/>
 Female Pyramid Insert for Prong <input type="checkbox"/> <input type="checkbox"/>	 Torsion Shock Adapter <input type="checkbox"/>	<p>* Patients under 220lb will be supplied with Aluminium adaptors unless specified.</p> <p>* Patients between 220-365lb will be provided with Titanium adaptors.</p> <p>Please check this box if you would prefer titanium adaptors.</p>
 Female Pyramid Insert for Prong LP <input type="checkbox"/>		
Specify Other Components and/or Quantities Desired		


Össur may select different components to fit your dimensions if necessary. If you do not want Össur to make any changes, check this box:

TT Socketless Össur Leg Order Form

Products

5: Foot: With or Without Unity (check one and fill out the options below)

No Foot		Foot Options:		Foot Size:	cm.	Category (1-9)	Left	Right	Foot Shell:	Beige	Brown														
						Suggested L-Codes†							Suggested L-Codes†												
	Balance Foot S  with Unity	5972, 5986 5781*		Balance Foot S Torsion  with Unity	5972, 5984, 5988 5781*																				
	Balance J with Unity	5976 5781*		K2 Sensation with DP Flexion with Unity	5972, 5986 5972, 5968 5781*																				
	Pro-Flex ST 	5981, 5986		Assure with Unity	5981 5781*																				
	Pro-Flex LP Torsion with Unity	5981, 5984, 5988 5781*		Vari-Flex with Unity	5981‡ 5781*																				
	Pro-Flex Pivot with Unity	See Össur reimbursement guide		Pro-Flex LP  with Unity 	5981, 5986 5781*																				
	Pro-Flex LP Align	5981, 5986, 5990		Proprio Foot  with Unity 	5973‡ 5781*																				
	ReFlex Shock with Unity	5987‡ 5781*		ReFlex Rotate with Unity	5987‡, 5984‡ 5781*																				
	Pro-Flex XC  with Unity 	5981‡ 5781*		Pro-Flex XC Torsion with Unity	5987‡, 5984‡ 5781*																				
	Pro-Flex Modular 	5980‡, 5986‡	Desired build height (mm): <small>With foot cover on. No return or exchange after foot is cut.</small> <table border="1"> <thead> <tr> <th rowspan="2">Foot Size</th> <th colspan="2">Build Height</th> </tr> <tr> <th>Min</th> <th>Max</th> </tr> </thead> <tbody> <tr> <td>22-24cm</td> <td>157mm</td> <td>313mm</td> </tr> <tr> <td>25-27cm</td> <td>168mm</td> <td>327mm</td> </tr> <tr> <td>28-30cm</td> <td>176mm</td> <td>340mm</td> </tr> </tbody> </table>									Foot Size	Build Height		Min	Max	22-24cm	157mm	313mm	25-27cm	168mm	327mm	28-30cm	176mm	340mm
Foot Size	Build Height																								
	Min	Max																							
22-24cm	157mm	313mm																							
25-27cm	168mm	327mm																							
28-30cm	176mm	340mm																							

 Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

*Unity code added to other codes, e.g. Proprio Foot with Unity: 5973, 5781.

†The responsibility for accurate coding lies solely with the healthcare provider, Össur assumes no responsibility or liability for the provider's coding decisions.

‡PDAC approved.

Knit-Rite Soft-Socks (sold in packs of 3)

Knit-Rite socks sold in packs of three, not included in Össur Legs pricing.

Measurement Instructions

Socks are sized in flat measurement, half of the circumference measurement correspond to the size measurements listed in the table. If the sock is used with a liner please measure the circumference with the liner on the patient.



Össur Soft-Sock — Black Polyester/Lycra® (Sold in packs of 3)							
Lightweight	3-Ply	5-Ply	Size	Length	Top width	Toe width*	Fits to liner size
Qty.	Qty.	Qty.	Narrow Short	10"–14" (25–35 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Narrow Medium	16"–20" (40–50 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Narrow Long	22"–28" (55–70 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Regular X-Short	8"–12" (20–30 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Short	10"–14" (25–35 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Medium	16"–20" (40–50 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Long	22"–28" (55–70 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Wide Short	10"–14" (25–35 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34
Qty.	Qty.	Qty.	Wide Medium	16"–20" (40–50 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34
Qty.	Qty.	Qty.	X-Wide Short	10"–14" (25–35 cm)	10" (25 cm)	7" (17.5 cm)	34, 36, 38, 40
Qty.	Qty.	Qty.	X-Wide Medium	16"–20" (40–50 cm)	10.0–12.0" (25–30 cm)	7–8" (17.5–20 cm)	36, 38, 40, 42, 45

* measured 2 inch (5cm) above distal end

Össur Hole-In-Toe Soft-Sock - Black Polyester/Lycra® (Sold in packs of 3)							
Lightweight	3-Ply	5-Ply	Size	Length	Top width	Toe width*	Fits to liner size
Qty.	Qty.	Qty.	Narrow Short	10"–14" (25–35 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Narrow Medium	16"–20" (40–50 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Regular X-Short	8"–12" (20–30 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Short	10"–14" (25–35 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Medium	16"–20" (40–50 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Long	22"–28" (55–70 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Wide Short	10"–14" (25–35 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34
Qty.	Qty.	Qty.	Wide Medium	16"–20" (40–50 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34

* Measured 2 inch (5cm) above distal end

CLEAR FORM

5 of 5

FOLLOW ÖSSUR ON



ÖSSUR CUSTOM SOLUTIONS
 9545 S JOHN YOUNG PARKWAY, STE #100 | ORLANDO, FL 32819
 TEL (888) 839-6213, FAX (407) 367-3695
 CS@OSSUR.COM | OSSUR.COM

