



Unloader One® X

Improves mobility by
providing pain relief
from knee osteoarthritis





Designed with patient compliance in mind

The Unloader One X builds on the existing strengths and over 20 years experience of the Unloader One. This clinically proven OA Knee Brace* incorporates a host of features designed to improve patient comfort, ease of use and aesthetics while maintaining the same pain relief for a wider range of patients.

EASE OF USE



Buckle system

Improved usability, press anywhere to lock



Dual action buckles

For easier donning and doffing



'Blue is you™' contact points

Patient touch points throughout the brace are highlighted using 'Blue is you™' feature



Straps adjustable anteriorly

For ease of fitting




Strap Seam

Easily distinguish inside from outside




Injection Molded Hooks

Simple and safe liner attachment



Hinge alignment guide

To align & assist in correct placement of brace



Machine washable

Entire brace is machine washable (washbag included)

*data on file at Ossur

Unloader One® braces can...

- ... provide **pain relief** for patients with knee OA
- ... provide **functional improvements** for patients with knee OA
- ... provide health related **quality of life improvements** for patients with degenerative meniscal tears
- ... provide **pain relief** for patients with symptomatic degenerative meniscal tears
- ... provide a **reduction of unicompartmental load**
- ... **improve quality of life** for patients with knee OA
- ... **delay surgery** for those with knee OA

Scan the QR code for more information on Unloader One clinical studies and research



COMFORT

- 1 Skin-friendly liners feel great against the skin and reduce pressure points
- 2 3 Point Leverage System with Dynamic Force System® straps, rigid shells, and a single upright hinge biomechanically unloads the affected knee compartment
- 3 Tool-free SmartDosing® dials enable the patient to fine-tune the strap tension to control pain
- 4 Anti-migration strap technology design captures the widest part of the calf to reduce brace migration
- 5 Alignment guide on hinge facilitates correct positioning of the brace on the leg and color coded touch points aid in easy donning and doffing
- 6 Tool-free polycentric hinge with optional ROM stops for range-of-motion restriction
- 7 Slender frame with soft edges easily contours to individual anatomies and remains discreet whether worn under or over clothing
- 8 Weighing only 360g



AESTHETICS

- 1 Stylish, low-profile frame
- 2 Uniform black colour scheme for a holistic and sleek look
- 3 Blue is you™

Accessories

Sundial

For easy grip on boa dial



Finger Loops

Effortless opening of buckles



The Össur One Sleeve

Size	Calf circumference 15 cm (6") above mid patella	Black - Össur One*
XXS	31 - 36	OSO60021
XS	36 - 42	OSO60022
S	42 - 48	OSO60023
M	48 - 55	OSO60024
L	55 - 62	OSO60025
XL	62 - 70	OSO60026
XXL	70 - 78	OSO60027

* Össur One: 51% Tencel™/30% Polyamide/19% Elastane



Sensil Liners



The Sensil Liners provide a combination of comfort and grip that helps to make sure that the brace stays in position.

	Right Medial, Left Lateral	Right Lateral, Left Medial
Thigh Liners	UOX526801	UOX525801
Calf Liners, XS-M	UOX516952	UOX515952
Calf Liners, L-XL	UOX516962	UOX515962

Ordering Information

Hinge is placed on the side of the affected compartment. Short version recommended for patients with a leg length of ≤ 85 cm (33.5").

Standard	Size	Right Medial	Right Lateral	Left Medial	Left Lateral
XS - M	25 - 43cm	UOX0518713	UOX0528713	UOX0618713	UOX0628713
L - XL	43 - 62cm	UOX0518715	UOX0528715	UOX0618715	UOX0628715

Short	Size	Right Medial	Right Lateral	Left Medial	Left Lateral
XS - M	25 - 43cm	UOX1518713	UOX1528713	UOX1618713	UOX1628713
L - XL	43 - 62cm	UOX1518715	UOX1528715	UOX1618715	UOX1628715

* Circumference: 15cm (6") below mid-patella.

UOX SERVICE KIT*	Right Medial, Left Lateral	Right Lateral, Left Medial	Left Medial	Left Lateral
Black, XS-M	UOX416901	UOX415901	UOX0618713	UOX0628713
Black, L-XL	UOX416911	UOX415911	UOX0618715	UOX0628715

* Includes Doeskin Thigh Liner, Sensil Calf Liner, Thigh and Lower DFS Strap with Buckle, Calf and Upper DFS Strap with Buckle, Hinge Pad, Popliteal Pad, and Alligator Kit

Additional information

Scan the QR code

to watch the Unloader One X Classroom Catch-up Playlist on Youtube.



Scan the QR code

to watch the Unloader One X clinician fitting video.



WWW.OSSUR.COM

Össur Europe B.V.
De Schakel 70
5651 GH Eindhoven
P.O. Box 120
5690 AC Son en Breugel
The Netherlands

TEL NL 00800 022 3719
info-europe@ossur.com

