

# TT Solution Order Form Package

Thank you for ordering a TT Solution from Össur Solutions.

This package includes all the order forms needed to define the transtibial solution, order the check socket, and order the definitive socket.

## Step One: Order the TT Solution and Check Socket

- Complete the TT Solution Order Form, and send to Össur Solutions.
- Össur Solutions will then create and send you a quote detailing all items in the TT Solution.
- Approve the quote and provide a purchase order number to finalize the order.
- The new liners included with the Össur Solutions package can be shipped before the check or definitive socket for casting/scanning/measurements by checking the box next to the Liners in the Products Section AND a Purchase Order is provided on the form.
- If special fabrication will be needed for check socket and/or definitive socket, complete the Special Fabrication portion of the check and/or definitive order forms. This ensures those items are included in the quote.
- Measure the residual limb.

*Cast/measure/scan over the same type and size liners the patient will use. Complete the Measurements section of the attached TT Solution Order Form. Send cast and TT Solution Order Form to Össur Solutions (see form for address).*

OR

### Order Online Using Össur Portal

You can order your TT Solution using these order forms, or you can order using the Össur Portal. Scan here for more information and to request access.



## Step Two: Fit the Check Socket, Make any Necessary Modifications

## Step Three: Order the Definitive Socket


- Fill out the attached TT Definitive Socket Order Form.
- Send TT Definitive Socket Order Form and modified check socket to Össur Solutions (see form for address).

This form is available in electronic form you can fill out electronically and email or print.

Please contact Össur Customer Service for a copy.



Instagram



Facebook

ÖSSUR SOLUTIONS  
9545 S JOHN YOUNG PARKWAY, STE #100 | ORLANDO, FL 32819  
TEL (888) 839-6213, FAX (800) 831-3160  
CS@OSSUR.COM | OSSUR.COM

© ÖSSUR, 02.2025. CRE-20828

# TT Solution Order Form Package

Send all casts, sockets, files to: 9545 S John Young Parkway, Ste #100, Orlando, FL 32819

tel: (888)839-6213, fax: (800)831-3160 | [cs@ossur.com](mailto:cs@ossur.com) | [www.ossur.com](http://www.ossur.com)

## Area Manager Information

Area Manager:

phone:

email:

Notes:

## Billing & Shipping Information

Össur Account #:

Contact for ordering or delivery questions:

Contact:

Email:

Mobile Phone:

Ship To:

Company:

Address:

City/State/Zip:

Contact:

Phone:

Fax:

Email:

PO#:

Shipping Check Priority: Next Day Air **2nd Day**  
Include return shipping label (\$50 charge inside 48 states)

Completion of this order form with the most accurate and up-to-date information, including all patient information, measurements, and construction information enables us to provide the highest quality product for you and your patient.

## Required Patient Information

Last Name:

Height: ft. in.

Impact Level: High Med Low Amputation Level: TT Symes

First Name:

Weight: lbs.

K-Level: K1 K2 K3 K4 Left Right Bilat Order

## Check Socket Measurements

**Call to discuss before design \*may result in a delay in your order**

Measurements**		Shape	Socket Styles	Trim Lines	Check Socket Material
	MPT to distal length	Shape by cast Shape by app scan Fully Modified CAD file CAD file name:	Total Surface Bearing Patellar Tendon Bearing Modified PTB	Standard Supracondylar Supracondylar Suprapatella	Polyethylene terephthalate glycol (PETG) option
	Level Measurement	Liner Used	*Alignment Flexion: Adduction:		High Impact Polystyrene (Oritrans Stiff, Thermolyn) option
	10cm 5cm 0cm 5cm 10cm 15cm 20cm	*If a specific alignment isn't provided, ÖCS' standard alignment is 0° flexion and 0° adduction. Notes:			
	PML cm ML cm AP cm	Total Reduction* % Ply	*3% global reduction (if using Iceross liners) OR 4% for any other liner unless requested differently.  *The percent reduction needed for correct socket fit is highly dependent on practitioner's measurement technique as well as patient's tissue density. It's not uncommon to require a different percent reduction.		
		Notes:			

\*\*Length and circumferences are required if Össur Custom Solutions is designing the socket. If not provided and socket fit requires remake, customer is responsible for cost of remake. Cast sent to Össur Custom Solutions (optional) Check here to receive foam carving for this order.

# TT Solution Order Form Package

## Special Fabrications

Flexible Inner Socket Plastics <i>Suggested L-code: 5645</i>	Fab Custom Inner From Sheet Material <i>Suggested L-codes: 5655 or 5654</i>		Other <i>Suggested L-codes: 5655 or 5654</i>
Flexible EVA with Silicone ( <b>Proflex with Silicone</b> ) Flexible EVA without Silicone ( <b>Proflex without Silicone</b> ) Flexible EVA – Black ( <b>Optek Flex Comfort - Black</b> ) Flexible EVA – Optek Flex Comfort White + Black ( <b>BiLam</b> ) Modified Polyethylene ( <b>MPE</b> ) EVA ( <b>ÖssurFlex</b> )	Lightweight Polyethylene – Firm 25 Shore A ( <b>PeLite</b> ) 12mm (1/2") beige white 5mm (3/16") beige 3mm (1/8") beige	Stovepipe Slit in lateral side Slit in medial side	Lightweight Ethylene Vinyl Acetate (Keasy Cone) 6.5mm profile, 10mm distal

Fab Custom Pad	RevoFit
<b>Fab Custom Pad on Check Socket</b> <b>Fab Custom Pad on Laminated Socket</b> Lightweight Polyethylene – Medium Shore A ( <b>Plastizote</b> ) 12mm (1/2") pink white 6mm (1/4") pink 5mm (3/16") pink  Lightweight Polyethylene – Firm 25 Shore A ( <b>PeLite</b> ) 12mm (1/2") beige white 5mm (3/16") beige 3mm (1/8") beige  Location Distal end (1/2") Other (see notes below) <i>Suggested L-code: 5668</i>	<b>RevoFit on Check Socket</b> <i>Suggested L-code: 5783</i> <b>RevoFit on Laminated Socket</b> <i>Suggested L-code: 5783</i>  Dial Location Distal anterior medial Proximal lateral Other (see notes below)  Panel Locations Pretibial & posterior (three panels) with Pads Posterior only (one panel) with Pad through trim line  Pad on Panel Firm Ethylene Vinyl Acetate 55 Shore A ( <b>Puff</b> ) 12mm (1/2") black 10mm (3/8") black 6mm (1/4") black

Windows and Doors	Cosmetics
<b>Windows and Doors on Check Socket</b> <b>Windows and Doors on Laminated Socket</b> Fab Window(s) or cut out(s) Number: Fab panel Fab overlap door (draw window or door on cast/check socket, or use notes field below) Add strap Dacron Velcro	Decorative Lamination on Definitive Using Fred's Legs Sleeve Customer supplied Össur to purchase Fred's item:  Customer shaped foam cover (measurements in TT measurement chart) <i>Suggested L-code: 5704</i> Fab cosmetic distal end

## Notes

# TT Solution Order Form Package

## Products

### 1. Assembly (check one)

Partial Assembly: Customer cuts pylon to final height (default option)	Partial Assembly: Ossur cuts pylon to height using landmark to floor measurement

Full Assembly by Customer


















Include Ossur Leg Bag
-----------------------

### 2: Liners - Two Are ALWAYS Included with Each Leg (check one then fill in desired sizes and options)

Send liners first for measurements/scanning/casting
---

<i>Suggested L-code: 5673</i>				<i>Suggested L-code: 5679</i>			
Comfort Locking	Dermo Locking	Iceform Locking	LINER SIZE	Seal-In V	Seal-In V HP	Seal-In X5	LINER SIZE
			Standard Wave 3 mm 6 mm				Standard Wave
<i>Suggested L-codes: 5679 and 5685</i>				<i>Suggested L-code: 5679</i>	<i>Suggested L-code: 5673</i>	<i>Suggested L-code: 7700 (only billed with L-5679) (Ring Options)</i>	
Comfort w/ Sleeve	Dermo w/ Sleeve	Iceform w/ Sleeve	LINER SIZE	Seal-In X with Rings	Seal-In X Locking with Rings	Classic Volume Grip	
			Standard Wave 3 mm 6 mm				
Iceform Sleeve	Iceroos Sleeve	Genu Sleeve	SLEEVE SIZE	Ring Size:	20	22	24
				26	28	30	32
				35	38	41	44
				47	51	55	60
							65
							LINER SIZE

# TT Solution Order Form Package

Products					
3: Sockets (check options)					
Check socket ONLY package		Check socket + Definitive package		Definitive ONLY package	
 <p>Polyethylene terephthalate glycol <b>(PETG)</b> High Impact Polystyrene <b>(Orfitrans Stiff, Thermolyn)</b></p> <p>Flexible Inner Socket Flexible Inner Socket - Plastic Flexible Inner Socket - Sheet Material/Other (PeLite, Keasy Cone) Include second check socket</p>		 <p>Polyethylene terephthalate glycol <b>(PETG)</b> High Impact Polystyrene <b>(Orfitrans Stiff, Thermolyn)</b></p> <p>Flexible Inner Socket Flexible Inner Socket - Plastic Flexible Inner Socket - Sheet Material/Other (PeLite, Keasy Cone) Include second check socket</p>		 <p>Laminated no flexible inner Laminated, with flexible inner Flexible Inner Socket - Plastic Flexible Inner Socket - Sheet Material/Other (PeLite, Keasy Cone)</p>	
<p><b>NOTE:</b> Exact Material selection is on definitive order form and not needed until the definitive form is ready for submission</p>					
Direct Socket TT					
		Fiberglass black resin	Fiberglass neutral resin	Basalt	Carbon
4" Material Kit (small distal connector, 200 ml resin)					
5" Material Kit (standard distal connector, 200 ml resin)					
7" Material Kit (400 ml resin) Medium Distal      Large Distal (Most Common)					
4. Suspension & Socket Adapter for Conventional Socket Fabrication (check all that apply)					
Locking Options					
<p>Icelock 621 Ratchet &amp; Extra Pin </p> <p><i>Suggested L-code: 5671</i></p>		<p>Icelock 651 Smooth &amp; Extra Pin </p> <p><i>Suggested L-code: 5671</i></p>		<p>Icelock 562 Hybrid </p> <p><i>Suggested L-code: 5671</i></p>	<p>Icelock 600 XM </p> <p><i>Suggested L-code: 5671</i></p>
Works With					
<p>Icelock 672 Male Pyramid Ti </p> <p>Icelock 673 Male Pyramid Al</p>		<p>Icelock 674 4-Hole Adapter Al </p> <p>Icelock 675 4-Hole Adapter Ti </p>		<p><i>Suggested L-code: 5671</i> Icelock 211 Clutch &amp; Extra Pin Icelock 214 4-Hole Clutch &amp; Extra Pin </p> <p><i>Suggested L-code: 5671</i></p>	
Expulsion Options			Unity Elevated Vacuum Options		
<p>Icelock 544 Expulsion Kit </p> <p><i>Suggested L-code: 5647</i></p>			<p>Icelock 551 Expulsion Valve </p> <p><i>Suggested L-code: 5647</i></p>		<p>Unity Valve TT </p>
<p>Icelock 544 Expulsion Valve Kit for sleeves </p> <p><i>Suggested L-code: 5647</i></p>			<p>Works With</p> <p>Icelock 544 Socket Adapter </p>		<p>Icelock 544 Unity Plate </p>
			<p>3-Prong Socket Adapter 4-Prong Socket Adapter</p>		








#### 4. Suspension & Socket Adapter for Conventional Socket Fabrication Cont. (check all that apply)

Bulldog Locks	KISS Suspension
3GEN-A – standard push lock BDSG-63 – ball bearing safety lock 3GEN-AGG (push lock with grooves) HDT-AP-TL Heavy Duty Top Loading Attachment Plate only (Does not include top load lock, HD pyramid, HD pylon, or HD tube clamp. Those items can be ordered through Bulldog Tools directly or a distributor)	KS1-A (Select either a 544 socket adapter or a prong adapter above) KS2-A KS4-A
ESP	
Lyn Valve - AK RV Auto Expulsion PA0002/B	

#### Products

#### 5. Direct Socket Suspension Options (Choose when fabricating your socket using a Direct Socket Kit)


##### Locking

621 Ratchet Pin <i>Suggested L-code: 5671</i> 	651 Smooth Pin <i>Suggested L-code: 5671</i> 	214 Clutch Lock <i>Suggested L-code: 5671</i> 	562 Hybrid Lock with Unity (for TT only)  <i>Suggested L-code: 5671</i> 	562 Hybrid Lock with Passive Vacuum (for TT only)  <i>Suggested L-code: 5671</i> 
<b>Requires: Modular Spacer Plate</b> 				

**Add Extra Pin:**    Long    Standard    Short

##### Expulsion
























































































Icelock 551 Expulsion Valve  <i>Suggested L-code: 5647</i> Lyn Valve - AK RV Auto Expulsion PA0002/B  <b>Requires: Bottom Seal</b> 	Icelock 544 Expulsion Plate <i>Suggested L-code: 5647</i>  <b>Requires: DS Adapter Plate</b> 	Unity TT Valve   <b>Requires: Bottom Seal</b> 	Icelock 544 Unity Plate  <b>Requires: DS Adapter Plate</b> 
---	---	--	--

 Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

# TT Solution Order Form Package

## Products

### 6. Endoskeletal Components (check all that apply for a single prosthesis)

4-Hole Adapters			Double Adapters			Single Adapters		
	4-Hole Male	 		Female Double Adapter 32mm	 		Male Single Adapter Short	 
	4-Hole Female	 		Female Double Adapter 45mm	 		Male Single Adapter Long	 
	4-Hole Male w/ Rotation	 		Female Double Adapter 60mm	 		Female Single Adapter Short	 
	4-Hole Male w/ 8mm Ax Offset			Female Double Adapter 75mm	 		Female Single Adapter Long	
	4-Hole Female w/ Rotation	 	<b>Pylon Kits</b>			Other:		
	4-Hole Female Pyramid 20/20 Offset 10° Right			Female Kit Short	 	<b>Pylons</b>		
	4-Hole Female Pyramid 20/20 Offset 10° Left			Female Kit Long	 		Female Short	 
	Flexion Plate w/ 45mm Offset, 10° Flexion	 	<b>Height Adjustable Pylons</b>			<b>Tube Clamps</b>		
	Flexion Plate w/ 60mm Offset, 10° Flexion	 		Standard	 		Male Tube Clamp	 
	Flexion Plate w/ 75mm Offset, 15° Flexion	 		Long	 		Female Tube Clamp	 
<b>Prong Adapters</b>			<b>Torsion Shock Adapter (Extra Cost)</b>					
	4-Prong Socket Adapter			Torsion Shock Adapter	 	* Patients under 220lb will be supplied with Aluminium adaptors unless specified.		
	3-Prong Socket Adapter		Suggested L-codes: 5984 and 5988			* Patients between 220-365lb will be provided with Titanium adaptors.		
	Male Pyramid Insert for Prong	 				Please check this box if you would prefer titanium adaptors.		
	Female Pyramid Insert for Prong	 						
	Female Pyramid Insert for Prong LP							





























Össur may select different components to fit your dimensions if necessary.


If you do not want Össur to make any changes, check this box:

# TT Solution Order Form Package

## Products

**7: Foot: With or Without Unity** (check one and fill out the options below)

No Foot		Foot Options:		Foot Size:	cm.	Category (1-9)	Left	Right	Foot Shell:	Beige	Brown														
						Suggested L-Codes†																			
	Balance Foot S  with Unity			5972, 5986 5781*				Balance Foot S Torsion  with Unity	5972, 5984, 5988 5781*		Suggested L-Codes†														
	Balance J with Unity			5976 5781*				K2 Sensation with DP Flexion with Unity	5972, 5986 5972, 5968 5781*																
	Pro-Flex ST  with Unity 			5981, 5986 5781*				Assure with Unity	5981 5781*																
	Pro-Flex LP Torsion with Unity			5981, 5984, 5988 5781*				Vari-Flex with Unity	5981‡ 5781*																
	Pro-Flex Pivot with Unity			5999, See Össur reimbursement guide				Pro-Flex LP  with Unity 	5981, 5986 5781*																
	Pro-Flex LP Align			5981, 5986, 5990				Proprio Foot  with Unity 	5973‡ 5781*																
	ReFlex Shock with Unity			5987‡ 5781*				ReFlex Rotate with Unity	5987‡, 5984‡ 5781*																
	Pro-Flex XC  with Unity 			5981‡ 5781*				Pro-Flex XC Torsion with Unity	5987‡, 5984‡ 5781*																
	Pro-Flex Modular 			5980‡, 5986‡			Desired build height (mm): With foot cover on. No return or exchange after foot is cut. <table border="1"> <thead> <tr> <th rowspan="2">Foot Size</th> <th colspan="2">Build Height</th> </tr> <tr> <th>Min</th> <th>Max</th> </tr> </thead> <tbody> <tr> <td>22-24cm</td> <td>157mm</td> <td>313mm</td> </tr> <tr> <td>25-27cm</td> <td>168mm</td> <td>327mm</td> </tr> <tr> <td>28-30cm</td> <td>176mm</td> <td>340mm</td> </tr> </tbody> </table>					Foot Size	Build Height		Min	Max	22-24cm	157mm	313mm	25-27cm	168mm	327mm	28-30cm	176mm	340mm
Foot Size	Build Height																								
	Min	Max																							
22-24cm	157mm	313mm																							
25-27cm	168mm	327mm																							
28-30cm	176mm	340mm																							

 Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

\*Unity code added to other codes, e.g. Proprio Foot with Unity: 5973, 5781.

†The responsibility for accurate coding lies solely with the healthcare provider, Össur assumes no responsibility or liability for the provider's coding decisions.

‡PDAC approved.



## Knit-Rite Soft-Socks (sold in packs of 3)

Knit-Rite socks sold in packs of three, not included in Össur Legs pricing.

### Measurement Instructions

Socks are sized in flat measurement, half of the circumference measurement correspond to the size measurements listed in the table. If the sock is used with a liner please measure the circumference with the liner on the patient.



Össur Soft-Sock — Black Polyester/Lycra® (Sold in Packs of 3)							
Lightweight	3-Ply	5-Ply	Size	Length	Top width	Toe width*	Fits to liner size
Qty.	Qty.	Qty.	Narrow Short	10"–14" (25–35 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Narrow Medium	16"–20" (40–50 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Narrow Long	22"–28" (55–70 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Regular X-Short	8"–12" (20–30 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Short	10"–14" (25–35 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Medium	16"–20" (40–50 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Long	22"–28" (55–70 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Wide Short	10"–14" (25–35 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34
Qty.	Qty.	Qty.	Wide Medium	16"–20" (40–50 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34
Qty.	Qty.	Qty.	X-Wide Short	10"–14" (25–35 cm)	10" (25 cm)	7" (17.5 cm)	34, 36, 38, 40
Qty.	Qty.	Qty.	X-Wide Medium	16"–20" (40–50 cm)	10.0–12.0" (25–30 cm)	7–8" (17.5–20 cm)	36, 38, 40, 42, 45

\* measured 2 inch (5cm) above distal end

Össur Hole-In-Toe Soft-Sock - Black Polyester/Lycra® (Sold in Packs of 3)							
Lightweight	3-Ply	5-Ply	Size	Length	Top width	Toe width*	Fits to liner size
Qty.	Qty.	Qty.	Narrow Short	10"–14" (25–35 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Narrow Medium	16"–20" (40–50 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Regular X-Short	8"–12" (20–30 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Short	10"–14" (25–35 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Medium	16"–20" (40–50 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Long	22"–28" (55–70 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Wide Short	10"–14" (25–35 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34
Qty.	Qty.	Qty.	Wide Medium	16"–20" (40–50 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34

\* Measured 2 inch (5cm) above distal end

## Billing & Shipping Information

Össur Account #:

CONTACT FOR ORDERING OR DELIVERY QUESTIONS:

Contact:

Email:

Mobile Phone:

Ship To:

Company:

Address:

City/State/Zip:

Contact:

Phone:

Fax:

Email:

PO#:

Shipping Check Priority: Next Day Air **2nd Day**

Completion of this order form with the most accurate and up-to-date information, including all patient information, measurements, and construction information enables us to provide the highest quality product for you and your patient.

## Patient Information

Last Name:

First Name:

Weight:

lbs.

Impact Level: High Medium Low

Amputation Level: BK SYMES

K-Level: K1 K2 K3 K4

Left: Right:

Bilat Order:

## Definitive Socket

■ Call to discuss before fabrication


From Socket

From Cast

From CAD file dated:

Fully Modified CAD file

### TT Measurements for Cover



## Lamination

Carbon

Basalt

Heavy Duty Lamination

Decorative Lamination from

Color (PRS Pigment 1-16)

Fred's Legs Sleeve

Customer supplied

Össur to purchase Fred's item:

## Alignment

Transfer

Transfer and zero out

Use Lines

Other

Use Angles

Flexion

Adduction

## Modifications

### Socket Sizing

### Prosthetic Height (optional)

Increase

Decrease

Ply:

%:

No Changes

Transfer

Extend:

Shorten:

## Special Fabrications

### Flexible Inner Socket Plastics

*Suggested L-code: 5645*

### Fab Custom Inner From Sheet Material

*Suggested L-codes: 5655 or 5654*

### Other

*Suggested L-codes: 5655 or 5654*

Flexible EVA with Silicone ( <b>Proflex with Silicone</b> ) Flexible EVA without Silicone ( <b>Proflex without Silicone</b> ) Flexible EVA – Black ( <b>Optek Flex Comfort - Black</b> ) Flexible EVA – Optek Flex Comfort White + Black ( <b>BiLam</b> ) Modified Polyethylene ( <b>MPE</b> ) EVA ( <b>ÖssurFlex</b> )	Lightweight Polyethylene – Firm 25 Shore A ( <b>PeLite</b> ) 12mm (1/2") beige white 5mm (3/16") beige 3mm (1/8") beige	Stovepipe Slit in lateral side Slit in medial side	Lightweight Ethylene Vinyl Acetate (Keasy Cone) 6.5mm profile, 10mm distal
--	--	--	--

### Fab Custom Pad

### RevoFit

*Suggested L-code: 5783*

Lightweight Polyethylene – Medium Shore A ( <b>Plastizote</b> ) 12mm (1/2") pink white 6mm (1/4") pink 5mm (3/16") pink  Lightweight Polyethylene – Firm 25 Shore A ( <b>PeLite</b> ) 12mm (1/2") beige white 5mm (3/16") beige 3mm (1/8") beige  Location <i>Suggested L-code: 5668</i> Distal end (1/2") Other (see notes below)	Dial Location Distal anterior medial Proximal lateral Other (see notes below)  Panel Locations Pretibial & posterior (three panels) with Pads Posterior only (one panel) with Pad through trim line  Pad on Panel Firm Ethylene Vinyl Acetate 55 Shore A ( <b>Puff</b> ) 12mm (1/2") black 10mm (3/8") black 6mm (1/4") black  See markings on cast/socket Other (see notes below)
---	---

### Windows and Doors

Fab Window(s) or cut out(s) Number:  
 Fab panel  
 Fab overlap door (draw window or door on cast/check socket, or use notes field below)  
 Add strap  
 Dacron Velcro

### Cosmetics

Decorative Lamination on Definitive Using Fred's Legs Sleeve  
 Customer supplied  
 Össur to purchase Fred's item:  
 Customer shaped foam cover (measurements in TF measurement chart) *Suggested L-code: 5704*  
 Fab cosmetic distal end

### Notes