



TT Össur Leg Order Form Package

Thank you for ordering a TT Össur Leg from Össur Custom Solutions. This package includes all the order forms needed to define the transtibial leg, order the check socket, and order the definitive socket.

Step One: Order the Complete Leg and Check Socket

- Complete the TT Össur Leg Order Form, and send to Össur Custom Solutions.
- Össur Custom Solutions will then create and send you a quote detailing all items in the Össur Leg.
- Approve the quote and provide a purchase order number to finalize the order.
- The new liners included with the Össur Leg package can be shipped before the check or definitive socket for casting/scanning/measurements by checking the box next to the Liners in the Products Section AND a Purchase Order is provided on the form.
- If special fabrication will be needed for check socket and/or definitive socket, complete the Special Fabrication portion of the check and/or definitive order forms. This ensures those items are included in the quote.
- Measure the residual limb.

Cast/measure/scan over the same type and size liners the patient will use. Complete the Measurements section of the attached TT Össur Leg Order Form. Send cast and TT Össur Leg order form to Össur Custom Solutions (see form for address).

OR

Order Online Using Össur Portal

You can order your Össur Leg using these order forms, or you can order using the Össur Portal. Scan here for more information and to request access.



Step Two: Fit the Check Socket, Make any Necessary Modifications

Step Three: Order the Definitive Socket

- Fill out the attached TT Össur Leg Definitive Socket Order Form.
- Send TT Össur Leg Definitive Socket Order Form and modified check socket to Össur Custom Solutions (see form for address).

This form is available in electronic form you can fill out electronically and email or print. Please contact Össur Customer Service for a copy.

FOLLOW ÖSSUR ON



ÖSSUR CUSTOM SOLUTIONS
9545 S JOHN YOUNG PARKWAY, STE #100 | ORLANDO, FL 32819
TEL (888) 839-6213, FAX (407) 367-3695
CS@OSSUR.COM | OSSUR.COM

TT Össur Leg Order Form Package

Send all casts, sockets, files to: 9545 S John Young Parkway, Ste #100, Orlando, FL 32819
 tel: (888)839-6213, fax: (800)788-9878 | cs@ossur.com | www.ossur.com

Area Manager Information		
Area Manager:	phone:	email:
Notes:		

Billing & Shipping Information		
Össur Account #:	Ship To:	
Contact for ordering or delivery questions:	Company:	
Contact:	Address:	
	City/State/Zip:	
	Contact:	
Email:	Phone:	Fax:
Mobile Phone:	Email:	

PO#:	Shipping Check Priority: Next Day Air 2nd Day Include return shipping label (\$50 charge inside 48 states)
------	----------------------------------------------------------------------------------------------------------------------

Completion of this order form with the most accurate and up-to-date information, including all patient information, measurements, and construction information enables us to provide the highest quality product for you and your patient.

Required Patient Information		
Last Name:	Height: ft. in.	Impact Level: High Med Low Amputation Level: TT Symes
First Name:	Weight: lbs.	K-Level: K1 K2 K3 K4 Left Right Bilat Order

Check Socket Measurements Call to discuss before design

Measurements**	Shape	Socket Styles	Trim Lines
	Shape by cast Shape by app scan Fully Modified CAD file CAD file name:	Total Surface Bearing Patellar Tendon Bearing Modified PTB	Standard Supracondylar Supracondylar Suprapatellar
	Liner Used	*Alignment Flexion: Adduction:	
*If a specific alignment isn't provided, ÖCS' standard alignment is 0° flexion and 0° adduction.			
Notes:			
Total Reduction* %		*3% global reduction (if using Iceross liners) OR 4% for any other liner unless requested differently.	
Ply		*The percent reduction needed for correct socket fit is highly dependent on practitioner's measurement technique as well as patient's tissue density. It's not uncommon to require a different percent reduction.	
Notes:			

**Length and circumferences are required if Össur Custom Solutions is designing the socket. If not provided and socket fit requires remake, customer is responsible for cost of remake. Cast sent to Össur Custom Solutions (optional) Check here to receive foam carving for this order.

TT Össur Leg Order Form Package

Special Fabrications

Flexible Inner Socket Plastics <i>Suggested L-code: 5645</i>	Fab Custom Inner From Sheet Material <i>Suggested L-codes: 5655 or 5654</i>	Other <i>Suggested L-codes: 5655 or 5654</i>
Flexible EVA with Silicone (Proflex with Silicone) Flexible EVA without Silicone (Proflex without Silicone) Flexible EVA – Black (Optek Flex Comfort - Black) Flexible EVA – Optek Flex Comfort White + Black (BiLam) Modified Polyethylene (MPE) EVA (ÖssurFlex)	Lightweight Polyethylene – Firm 25 Shore A (PeLite) 12mm (1/2") beige white 5mm (3/16") beige 3mm (1/8") beige	Stovepipe Slit in lateral side Slit in medial side Lightweight Ethylene Vinyl Acetate (Keasy Cone) 6.5mm profile, 10mm distal

Fab Custom Pad

RevoFit

<p>Fab Custom Pad on Check Socket</p> <p>Fab Custom Pad on Laminated Socket</p> <p>Lightweight Polyethylene – Medium Shore A (Plastizote) 12mm (1/2") pink white 6mm (1/4") pink 5mm (3/16") pink</p> <p>Lightweight Polyethylene – Firm 25 Shore A (PeLite) 12mm (1/2") beige white 5mm (3/16") beige 3mm (1/8") beige</p> <p>Location Distal end (1/2") Other (see notes below) <i>Suggested L-code: 5668</i></p>	<p>RevoFit on Check Socket <i>Suggested L-code: 5783</i></p> <p>RevoFit on Laminated Socket <i>Suggested L-code: 5783</i></p> <p>Dial Location Distal anterior medial Proximal lateral Other (see notes below)</p> <p>Panel Locations Pretibial & posterior (three panels) with Pads Posterior only (one panel) with Pad through trim line</p> <p>Pad on Panel Firm Ethylene Vinyl Acetate 55 Shore A (Puff) 12mm (1/2") black 10mm (3/8") black 6mm (1/4") black</p>
<p>Windows and Doors</p>	<p>Cosmetics</p> <p>Decorative Lamination on Definitive Using Fred's Legs Sleeve Customer supplied Össur to purchase Fred's item: Customer shaped foam cover (measurements in TT measurement chart) <i>Suggested L-code: 5704</i> Fab cosmetic distal end</p>

<p>Windows and Doors on Check Socket</p> <p>Windows and Doors on Laminated Socket</p> <p>Fab Window(s) or cut out(s) Number: Fab panel Fab overlap door (draw window or door on cast/check socket, or use notes field below) Add strap Dacron Velcro</p>	<p>Cosmetics</p> <p>Decorative Lamination on Definitive Using Fred's Legs Sleeve Customer supplied Össur to purchase Fred's item: Customer shaped foam cover (measurements in TT measurement chart) <i>Suggested L-code: 5704</i> Fab cosmetic distal end</p>
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Notes

TT Össur Leg Order Form Package

Products

1. Assembly (check one)

Partial Assembly: Customer cuts pylon to final height (default option)

Partial Assembly: Össur cuts pylon to height using landmark to floor measurement



Full Assembly by Customer



Include Össur Leg Bag

2: Liners - Two Are ALWAYS Included with Each Leg (check one then fill in desired sizes and options)

Send liners first for measurements/scanning/casting

Suggested L-code: 5673

Comfort Locking



Dermo Locking



Iceform Locking



LINER SIZE

Standard
Wave
3 mm
6 mm

Suggested L-code: 5679

Seal-In V



Seal-In V HP



Seal-In X5



LINER SIZE

Standard
Wave

Suggested L-codes: 5679 and 5685

Comfort w/
Sleeve



Dermo w/
Sleeve



Iceform w/
Sleeve



LINER SIZE

Standard
Wave
3 mm
6 mm

SLEEVE SIZE

Iceform Sleeve

Iceross Sleeve

Genu Sleeve

Suggested L-code: 5679

Seal-In X
with
Rings



Suggested L-code: 5673

Seal-In X
Locking
with Rings



*Suggested L-code: 7700
(only billed with L-5679)*

(Ring Options)



Classic
Volume
Grip

Ring Size:	20	22	24	26	28
	30	32	35	38	41
	47	51	55	60	65

TT Össur Leg Order Form Package








Products









3: Sockets (check options)

Check socket ONLY package	Check socket + Definitive package	Definitive ONLY package
 <p>Polyethylene terephthalate glycol (PETG) High Impact Polystyrene (Orfitrans Stiff, Thermolyn)</p> <p>Flexible Inner Socket Flexible Inner Socket - Plastic Flexible Inner Socket - Sheet Material/Other (PeLite, Keasy Cone) Include second check socket</p>	 <p>Polyethylene terephthalate glycol (PETG) High Impact Polystyrene (Orfitrans Stiff, Thermolyn)</p> <p>Flexible Inner Socket Flexible Inner Socket - Plastic Flexible Inner Socket - Sheet Material/Other (PeLite, Keasy Cone) Include second check socket</p>	 <p>Laminated no flexible inner Laminated, with flexible inner Flexible Inner Socket - Plastic Flexible Inner Socket - Sheet Material/Other (PeLite, Keasy Cone)</p>
<p>NOTE: Exact Material selection is on definitive order form and not needed until the definitive form is ready for submission</p>		

Direct Socket TT	Fiberglass black resin	Fiberglass neutral resin	Basalt	Carbon
4" Material Kit (small distal connector, 200 ml resin)				
5" Material Kit (standard distal connector, 200 ml resin)				
7" Material Kit (400 ml resin) Medium Distal Large Distal (Most Common)				

4. Suspension & Socket Adapter for Conventional Socket Fabrication (check all that apply)

Locking Options					
<p>Icelock 621 Ratchet & Extra Pin </p> <p><i>Suggested L-code: 5671</i></p> 	<p>Icelock 651 Smooth & Extra Pin </p> <p><i>Suggested L-code: 5671</i></p> 	<p>Icelock 562 Hybrid </p> <p><i>Suggested L-code: 5671</i></p> 	<p>Icelock 600 XM</p> <p><i>Suggested L-code: 5671</i></p> 	<p><i>Suggested L-code: 5671</i></p> <p>Icelock 211 Clutch & Extra Pin Icelock 214 4-Hole Clutch & Extra Pin</p> 	<p>Icelock 234 Lanyard</p> <p><i>Suggested L-code: 5671</i></p> 
Works With					
<p>Icelock 672 Male Pyramid Ti </p> <p>Icelock 673 Male Pyramid Al</p> 	<p>Icelock 674 4-Hole Adapter Al Icelock 675 4-Hole Adapter Ti </p> 				

Expulsion Options	Unity Elevated Vacuum Options		
<p>Icelock 544 Expulsion Kit</p> <p><i>Suggested L-code: 5647</i></p> 	<p>Icelock 551 Expulsion Valve </p> <p><i>Suggested L-code: 5647</i></p> 	<p>Unity Valve TT </p> 	<p>Icelock 544 Unity Plate</p> 
Works With			
<p>Icelock 544 Expulsion Valve Kit for sleeves</p> <p><i>Suggested L-code: 5647</i></p> 	<p>Icelock 544 Socket Adapter</p> 	<p>3-Prong Socket Adapter 4-Prong Socket Adapter</p>	









4. Suspension & Socket Adapter for Conventional Socket Fabrication Cont. (check all that apply)

Bulldog Locks	KISS Suspension
3GEN-A – standard push lock BDSG-63 – ball bearing safety lock 3GEN-AGG (push lock with grooves) HDT-AP-TL Heavy Duty Top Loading Attachment Plate only (Does not include top load lock, HD pyramid, HD pylon, or HD tube clamp. Those items can be ordered through Bulldog Tools directly or a distributor)	KS1-A (Select either a 544 socket adapter or a prong adapter above) KS2-A KS4-A
ESP	
Lyn Valve - AK RV Auto Expulsion PA0002/B	

Products

5. Direct Socket Suspension Options (Choose when fabricating your socket using a Direct Socket Kit)

Locking


621 Ratchet Pin <i>Suggested L-code: 5671</i> 	651 Smooth Pin <i>Suggested L-code: 5671</i> 	214 Clutch Lock <i>Suggested L-code: 5671</i> 	562 Hybrid Lock with Unity (for TT only)  <i>Suggested L-code: 5671</i> 	562 Hybrid Lock with Passive Vacuum (for TT only)  <i>Suggested L-code: 5671</i> 
Requires: Modular Spacer Plate				
				

Add Extra Pin: Long Standard Short































Expulsion

Unity

Icelock 551 Expulsion Valve  <i>Suggested L-code: 5647</i> Lyn Valve - AK RV Auto Expulsion PA0002/B  Requires: Bottom Seal 	Icelock 544 Expulsion Plate <i>Suggested L-code: 5647</i>  Requires: DS Adapter Plate 	Unity TT Valve   Requires: Bottom Seal 	Icelock 544 Unity Plate  Requires: DS Adapter Plate 
-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

 Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

TT Össur Leg Order Form Package





























Products		
6. Endoskeletal Components (check all that apply for a single prosthesis)		
4-Hole Adapters	Double Adapters	Single Adapters
 4-Hole Male <input type="checkbox"/> <input type="checkbox"/>	 Female Double Adapter 32mm <input type="checkbox"/> <input type="checkbox"/>	 Male Single Adapter Short <input type="checkbox"/> <input type="checkbox"/>
 4-Hole Female <input type="checkbox"/> <input type="checkbox"/>	 Female Double Adapter 45mm <input type="checkbox"/> <input type="checkbox"/>	 Male Single Adapter Long <input type="checkbox"/> <input type="checkbox"/>
 4-Hole Male w/ Rotation <input type="checkbox"/> <input type="checkbox"/>	 Female Double Adapter 60mm <input type="checkbox"/> <input type="checkbox"/>	 Female Single Adapter Short <input type="checkbox"/> <input type="checkbox"/>
 4-Hole Male w/ 8mm Ax Offset <input type="checkbox"/>	 Female Double Adapter 75mm <input type="checkbox"/> <input type="checkbox"/>	 Female Single Adapter Long <input type="checkbox"/>
 4-Hole Female w/ Rotation <input type="checkbox"/> <input type="checkbox"/>	Pylon Kits	Pylons
 4-Hole Female Pyramid 20/20 Offset 10° Right <input type="checkbox"/>	 Female Kit Short <input type="checkbox"/> <input type="checkbox"/>	 Female Short <input type="checkbox"/> <input type="checkbox"/>
 4-Hole Female Pyramid 20/20 Offset 10° Left <input type="checkbox"/>	 Female Kit Long <input type="checkbox"/> <input type="checkbox"/>	 Female Long <input type="checkbox"/> <input type="checkbox"/>
Prong Adapters	Height Adjustable Pylons	Tube Clamps
 4-Prong Socket Adapter <input type="checkbox"/>	 Standard <input type="checkbox"/>	 Male Tube Clamp <input type="checkbox"/> <input type="checkbox"/>
 3-Prong Socket Adapter <input type="checkbox"/>	 Long <input type="checkbox"/>	 Female Tube Clamp <input type="checkbox"/> <input type="checkbox"/>
 Male Pyramid Insert for Prong <input type="checkbox"/> <input type="checkbox"/>	Torsion Shock Adapter	 4-Hole Tube Clamp <input type="checkbox"/> <input type="checkbox"/>
 Female Pyramid Insert for Prong <input type="checkbox"/> <input type="checkbox"/>	 <i>Suggested L-codes: 5984 and 5988</i> Torsion Shock Adapter (not included in package price) <input type="checkbox"/>	<p>* Patients under 220lb will be supplied with Aluminium adaptors unless specified.</p> <p>* Patients between 220-365lb will be provided with Titanium adaptors.</p> <p>Please check this box if you would prefer titanium adaptors.</p>
 Female Pyramid Insert for Prong LP <input type="checkbox"/>		
Specify Other Components and/or Quantities Desired		


Össur may select different components to fit your dimensions if necessary. If you do not want Össur to make any changes, check this box:

TT Össur Leg Order Form Package

Products

7: Foot: With or Without Unity (check one and fill out the options below)

No Foot		Foot Options:		Foot Size:	cm.	Category (1-9)	Left	Right	Foot Shell:	Beige	Brown														
						Suggested L-Codes†							Suggested L-Codes†												
	Balance Foot S  with Unity	5972, 5986 5781*		Balance Foot S Torsion  with Unity	5972, 5984, 5988 5781*																				
	Balance J with Unity	5976 5781*		K2 Sensation with DP Flexion with Unity	5972, 5986 5972, 5968 5781*																				
	Pro-Flex ST  with Unity 	5981, 5986 5781*		Assure with Unity	5981 5781*																				
	Pro-Flex LP Torsion with Unity	5981, 5984, 5988 5781*		Vari-Flex with Unity	5981‡ 5781*																				
	Pro-Flex Pivot with Unity	5999, See Össur reimbursement guide		Pro-Flex LP  with Unity 	5981, 5986 5781*																				
	Pro-Flex LP Align	5981, 5986, 5990		Proprio Foot  with Unity 	5973‡ 5781*																				
	ReFlex Shock with Unity	5987‡ 5781*		ReFlex Rotate with Unity	5987‡, 5984‡ 5781*																				
	Pro-Flex XC  with Unity 	5981‡ 5781*		Pro-Flex XC Torsion with Unity	5987‡, 5984‡ 5781*																				
	Pro-Flex Modular 	5980‡, 5986‡	Desired build height (mm): <small>With foot cover on. No return or exchange after foot is cut.</small> <table border="1"> <thead> <tr> <th rowspan="2">Foot Size</th> <th colspan="2">Build Height</th> </tr> <tr> <th>Min</th> <th>Max</th> </tr> </thead> <tbody> <tr> <td>22-24cm</td> <td>157mm</td> <td>313mm</td> </tr> <tr> <td>25-27cm</td> <td>168mm</td> <td>327mm</td> </tr> <tr> <td>28-30cm</td> <td>176mm</td> <td>340mm</td> </tr> </tbody> </table>									Foot Size	Build Height		Min	Max	22-24cm	157mm	313mm	25-27cm	168mm	327mm	28-30cm	176mm	340mm
Foot Size	Build Height																								
	Min	Max																							
22-24cm	157mm	313mm																							
25-27cm	168mm	327mm																							
28-30cm	176mm	340mm																							

 Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

*Unity code added to other codes, e.g. Proprio Foot with Unity: 5973, 5781.

†The responsibility for accurate coding lies solely with the healthcare provider, Össur assumes no responsibility or liability for the provider's coding decisions.

‡PDAC approved.

Billing & Shipping Information

Össur Account #:

CONTACT FOR ORDERING OR DELIVERY QUESTIONS:

Contact:

Email:

Mobile Phone:

Ship To:

Company:

Address:

City/State/Zip:

Contact:

Phone:

Email:

Fax:

PO#: Shipping Check Priority: Next Day Air **2nd Day**

Completion of this order form with the most accurate and up-to-date information, including all patient information, measurements, and construction information enables us to provide the highest quality product for you and your patient.


Patient Information

Last Name:	First Name:	Weight:	lbs.	Impact Level: High	Medium	Low
Amputation Level: BK SYMES	K-Level: K1 K2 K3 K4	Left:	Right:	Bilat Order:		

Definitive Socket ■ Call to discuss before fabrication

From Socket From Cast From CAD file dated:
Fully Modified CAD file

TT Measurements for Cover



Lamination

Carbon	Basalt	Heavy Duty Lamination
Decorative Lamination from		Color (PRS Pigment 1-16)
Fred's Legs Sleeve		
Customer supplied		
Össur to purchase Fred's item:		

Alignment

Transfer	Use Angles
Transfer and zero out	Flexion
Use Lines	Adduction
Other	

Modifications

Socket Sizing				Prosthetic Height (optional)		
Increase	Decrease	Ply:	%:	Transfer	Extend:	Shorten:
No Changes						

Special Fabrications

Flexible Inner Socket Plastics <i>Suggested L-code: 5645</i>	Fab Custom Inner From Sheet Material <i>Suggested L-codes: 5655 or 5654</i>		Other <i>Suggested L-codes: 5655 or 5654</i>
Flexible EVA with Silicone (Proflex with Silicone) Flexible EVA without Silicone (Proflex without Silicone) Flexible EVA – Black (Optek Flex Comfort - Black) Flexible EVA – Optek Flex Comfort White + Black (BiLam) Modified Polyethylene (MPE) EVA (ÖssurFlex)	Lightweight Polyethylene – Firm 25 Shore A (PeLite) 12mm (1/2") beige white 5mm (3/16") beige 3mm (1/8") beige	Stovepipe Slit in lateral side Slit in medial side	Lightweight Ethylene Vinyl Acetate (Keasy Cone) 6.5mm profile, 10mm distal

Fab Custom Pad	RevoFit <i>Suggested L-code: 5783</i>
Lightweight Polyethylene – Medium Shore A (Plastizote) 12mm (1/2") pink white 6mm (1/4") pink 5mm (3/16") pink Lightweight Polyethylene – Firm 25 Shore A (PeLite) 12mm (1/2") beige white 5mm (3/16") beige 3mm (1/8") beige Location <i>Suggested L-code: 5668</i> Distal end (1/2") Other (see notes below)	Dial Location Distal anterior medial Proximal lateral Other (see notes below) Panel Locations Pretibial & posterior (three panels) with Pads Posterior only (one panel) with Pad through trim line Pad on Panel Firm Ethylene Vinyl Acetate 55 Shore A (Puff) 12mm (1/2") black 10mm (3/8") black 6mm (1/4") black See markings on cast/socket Other (see notes below)
Windows and Doors	
Fab Window(s) or cut out(s) Number: Fab panel Fab overlap door (draw window or door on cast/check socket, or use notes field below) Add strap Dacron Velcro	

Cosmetics	Notes
Decorative Lamination on Definitive Using Fred's Legs Sleeve Customer supplied Össur to purchase Fred's item: Customer shaped foam cover (measurements in TF measurement chart) <i>Suggested L-code: 5704</i> Fab cosmetic distal end	