

DARCO APB™ ALL PURPOSE BOOT

DARCO APB™ ALL PURPOSE BOOT



The DARCO APB™ embodies nearly perfect four season postoperative protection with a closed-toe design that keeps dressings and toes warm and dry. The boot is designed to provide unparalleled versatility by accommodating both bulky dressings and most fiberglass casts.

FEATURES & BENEFITS

Post surgical healing.

Fits over casts and other bulky dressings.

Foot/ankle trauma.

Wounds or ulcerations present under metatarsal heads and toes.

ORDERING INFORMATION

PART #	SIZE	MENS SHOE	WOMENS SHOE	LENGTH	WIDTH
APQ0B	X-Small	-	4.0 - 6.5	22cm	9cm
APQ1B	Small	6.0 - 8	7.0 - 8.0	24cm	9.5cm
APQ2B	Medium	8.5 - 10	8.5 - 10	26cm	10.5cm
APQ3B	Large	10.5 - 12	10.5 -13	28cm	11cm
APQ4B	X-Large	12.5 - 14	13+	30cm	12cm