

# DARCO HEELWEDGE™ SHOE

FOR HEEL WEIGHT RELIEF

## DARCO HEELWEDGE™ SHOE



The DARCO HeelWedge™ is clinically proven to off-load pressure from the heel by shifting weight to the mid and forefoot to promote faster healing after surgery, trauma or when wounds or ulcerations are present on the heel.

### INDICATIONS

Rearfoot trauma.

Wounds or ulcerations present on the heel area.

Post Surgical healing for either soft tissue or bony structure of the heel.

### FEATURES & BENEFITS

Reduces Heel pressure by over 26%!

Square Toe Design acts as a bumper to protect the toes and provides better universal Left/Right fit.

Ankle Strap seats foot firmly in shoe to reduce heel slippage and eliminate friction.

Removable Forefoot Closure provides security while eliminating buckle pressure.

Zoned Outsole has an aggressive tread under the mid-foot where traction is needed most.

Removable Insole features twice the padding of standard insoles and can be modified as needed.

## ORDERING INFORMATION

PART #	SIZE	MENS SHOE	WOMENS SHOE	LENGTH	WIDTH
HQ0B	XSmall	-	4.0 - 6.5	22cm	9cm
HQ1B	Small	6.0 - 8	7.0 - 8.0	24cm	9.5cm
HQ2B	Medium	8.5 - 10	8.5 - 10	26cm	10.5cm
HQ3B	Large	10.5 - 12	10.5 -13	28cm	11cm
HQ4B	XLarge	12.5 - 14	13+	30cm	12cm