

DARCO MEDSURG™ SHOE

DARCO MEDSURG™ SHOE



The DARCO MedSurg Shoe is designed to provide the foot with solid protection and accommodate bulky bandages with ease. It features a square-toe design to provide added protection for the toes. The sole has a metashank design that provides added support to the metatarsal region when walking. The DARCO MedSurg Shoe has a velcro closure across the front which can fit up to a 17 inch circumference from large bandages and dressing. It also eliminates any added pressure caused by more traditional D ring closures that can irritate the forefoot area.

Very aggressive tread at the midfoot provides superior traction. The reduced friction tread at the toe and heel to allow toe drag or heel scuff without tripping up the patient. The strap around the front of the ankle prevents the product from slipping off the heel. The DARCO MedSurg has a rocker bottom sole which may help reduce some pressure on the bottom surface of your foot.

INDICATIONS

After foot surgery or when K-wires are present

When large bandages or dressing are required

FEATURES & BENEFITS

Square toe design acts as a bumper to provide additional protection when K-wires are present.

MetaShank™ Protection provides more rigid control under the metatarsal heads.

Strapless Closure eliminates buckle pressure on the forefoot and can expand to accommodate dressings up to 17 inches in circumference.

Rocker sole may help reduce plantar pressure on the forefoot and heel.

Compatible with DARCO Toe Guard

CLINICAL TIPS

For female patients with bulky dressings, the male shoe can be used. It provides more width e.g. use male medium instead of female large.

Where toe protection is required, add the DARCO toe guard, see page 16.

ORDERING INFORMATION - MENS

PART #	SIZE	SHOE SIZE	LENGTH	WIDTH
MQM1B	Small	6 - 8	26cm	10cm
MQM2B	Medium	8.5 - 10	27.5cm	10.5cm
MQM3B	Large	10.5 - 12	29cm	11cm
MQM4B	X-Large	12.5 - 14	30.5cm	12cm

ORDERING INFORMATION - WOMENS

PART #	SIZE	SHOE SIZE	LENGTH	WIDTH
MQW1B	Small	4 - 6	23cm	9cm
MQW2B	Medium	6.5 - 8	24.5cm	10cm
MQW3B	Large	8.5 - 10	26cm	10cm

ORDERING INFORMATION - PAEDIATRIC

PART #	SIZE	SHOE SIZE
MQPB	Paediatric	Children 12 - 1