

DARCO ORTHOWEDGE™ SHOE

FOR FOREFOOT WEIGHT RELIEF

DARCO ORTHOWEDGE™ SHOE



The DARCO OrthoWedge™ is clinically proven to reduce weight bearing pressure on the forefoot that may help with faster healing after surgery, trauma or when forefoot wounds or ulcerations are present.

INDICATIONS

Post Surgical healing

Forefoot trauma

Wounds or ulcerations present under metatarsal heads and toes.

Any situation that requires a pressure off-load of the forefoot.

FEATURES & BENEFITS

15° Wedge Sole effectively shifts body weight to the midfoot and heel. Reduces forefoot pressure by as much as 57%!

Square Toe Design acts as a bumper to protect the toes when k-wires are present and provides better Left/Right fit.

Ankle Strap seats foot firmly in shoe to reduce heel slippage and eliminate friction.

Removable Forefoot Closure provides security while eliminating buckle pressure.

Zoned Outsole has an aggressive tread under the mid-foot where traction is needed most.

Removable Insole features twice the padding of standard insoles and can be modified as needed.

Compatible with DARCO Toe Guard

POST-OPERATIVE
SHOE OPTIONS

ORDERING INFORMATION

PART #	SIZE	MENS SHOE	WOMENS SHOE	LENGTH	WIDTH
OQ0B	X-Small	-	4.0 - 6.5	22cm	9cm
OQ1B	Small	6.0 - 8	7.0 - 8.0	24cm	9.5cm
OQ2B	Medium	8.5 - 10	8.5 - 10	26cm	10.5cm
OQ3B	Large	10.5 - 12	10.5 - 13	28cm	11cm
OQ4B	X-Large	12.5 - 14	13+	30cm	12cm