

Waterproof Solutions

Transfemoral Order Form

Step 1: General Information

Prosthetist Name: _____ Acc Nr: _____
 Patient Identifier: _____ Weight: _____ kg Date of Birth: ____/____/____

Step 2: Select Suspension System

Seal-In®

Iceross Seal-In® X TF

Standard Conical

Size: _____

Iceross Seal-In® Ring

Ring size: _____

Classic Volume Grip

Iceross® TF Seal-In®

Standard Conical

Size: _____

Plus

Expulsion Valve System

Icelock® Expulsion Valve 552

Or

Unity System

Unity TF Valve Kit

Locking

Iceross® TF Locking Liner

Standard Conical

Size: _____

Plus Lock System

Icelock® 672 Titanium (L-672000)

Icelock® 621 (L-621000)

Icelock® 272 Male Pyramid TI (L-272000)

Icelock® 631 Lanyard (L-631000)

Icelock® 651 Smooth (L-651000)

Step 3: Select Knee

Capital™ Knee



Waterproof

Plus

Adapters

Option 1

3-Prong Socket Adapter (A-114030)

Female Insert for prong TI (A-145300)

Female Pylon TI Kit (Pyramid, Tube, Tube Clamp) WP-Long (A-764121)

Option 2

4-Hole Socket Adapter TI (A-125100)

Female Pylon TI Kit (Pyramid, Tube, Tube Clamp), WP-Long (A-764121)

4-Hole Female Pyramid Adapter TI (A-245100)

or

4-Hole Female Pyramid Adapter Rotation TI (A-245300)



Step 4: Select Waterproof Foot + Foot Cover

Balance™ Foot S
Balance™ Foot S Torsion
Flex-Foot Assure®*
Proprio Foot®**
Proprio Foot® with Unity®
Pro-Flex® LP
Pro-Flex® LP with Unity®
Pro-Flex® XC
Pro-Flex® XC with Unity®
Vari-Flex®*

Pro-Flex® ST
Pro-Flex® Modular
Vari-Flex® Modular (Only available direct lamination for higher cat)
Cheetah® Xplore
Flex-Run™
Cheetah® Xceed

FOOT COVER Beige Brown
Foot Size: _____
Left/Right: _____

- Waterproof only when used with tube kit.
- Proprio Foot only allows for full immersion into 1m deep water for 30 minutes and is not recommended to be used in chlorinated water.

Additional Waterproof Adapters

Waterproof Socket Adapters

| | |
|-------------------------------|------------|
| Male Socket Adapter | (A-135100) |
| 3-Prong Socket Adapter | (A-114030) |
| 4-Prong Socket Adapter | (A-135100) |
| 4-Hole Socket Adapter TI | (A-125100) |
| Male Insert For Prong TI | (A-135300) |
| Female Insert For Prong TI | (A-145300) |
| Female Insert For Prong LP TI | (A-145310) |

Waterproof 4-Hole Adapters

| | |
|---|------------|
| 4-Hole Male Pyramid Adapter TI | (A-235100) |
| 4-Hole Male Pyramid Adapter Rotation TI | (A-235300) |
| 4-Hole Male Pyramid Adapter 8mm Offset TI | (A-235400) |
| 4-Hole Female Pyramid Adapter TI | (A-245100) |
| 4-Hole Female Pyramid Adapter Rotation TI | (A-245300) |
| 4-Hole Tube Clamp TI | (A-325100) |

Waterproof Tube Clamps And Single Adapters

| | |
|------------------------------|------------|
| Male Tube Clamp TI | (A-335100) |
| Female Tube Clamp TI | (A-345100) |
| Female Single Adapter TI | (A-445100) |
| Male Single Adapter TI-Short | (A-435110) |
| Male Single Adapter TI-Long | (A-435120) |
| Male Pyramid Adapter | (A-835300) |
| Female Pyramid Adapter | (A-845300) |

Waterproof Double Adapters

| | |
|----------------------------------|------------|
| Female Double Adapter 32mm | (A-545110) |
| Female Double Adapter 45mm | (A-545120) |
| Female Double Adapter 60mm | (A-545130) |
| Female Double Adapter 75mm | (A-545140) |
| Male Double Adapter TI | (A-535100) |
| Female Double Adapter with Shift | (A-645210) |

Waterproof Pylons And Pylon Kits

| | |
|--|------------|
| Female Pylon TI AL-Short | (A-748010) |
| Female Pylon TI AL-Long | (A-748020) |
| Female Pylon Titanium, WP-Short | (A-748030) |
| Female Pylon Titanium, WP-Long | (A-748040) |
| Female Pylon Titanium Kit (Pyramid, Tube, Tube Clamp), WP-Short | (A-764111) |
| Female Pylon Titanium Kit (Pyramid, Tube, Tube Clamp), WP-Long | (A-764121) |

Disclaimer: Please refer to the instructions (IFU) for each product exact waterproof specifications.

Warning: For liners and sleeves it's specifically important to check the liner suspension during use in winter and dry the socket and the liner after leaving the water.



WWW.OSSUR.COM/EN-ZA

Össur South Africa
Unit 1, Stonewood Business Park
9 Tanzanite Crescent
Kraaifontein
Cape Town
7550

T 0800 996 220
F 0800 981 744
E ossursa@ossur.com

