

# Bracing & Supports

Catalogue

---





# ÖSSUR BRACING & SUPPORTS CATALOGUE

## TABLE OF CONTENTS

|                                                   |    |
|---------------------------------------------------|----|
| ORDERING INFORMATION .....                        | 10 |
| INTRODUCTION TO SOLUTIONS.....                    | 12 |
| SPINE .....                                       | 15 |
| MIAMI J® .....                                    | 17 |
| MIAMI J® SELECT .....                             | 19 |
| MIAMI JTO® THORACIC EXTENSION .....               | 20 |
| OCCIAN® BACK .....                                | 21 |
| MIAMI Jr® .....                                   | 22 |
| PHILADELPHIA® ADJUSTABLE TRACHEOTOMY COLLAR ..... | 23 |
| PHILADELPHIA® TRACHEOTOMY COLLAR .....            | 24 |
| NECLOC® .....                                     | 25 |
| NECLOC® KIDS .....                                | 26 |
| MIAMI LSO® .....                                  | 27 |
| MIAMI TLSO 464 .....                              | 29 |
| MIAMI TLSO 456 .....                              | 31 |
| MIAMI LUMBAR® POSTEO .....                        | 33 |
| PAPOOSE® .....                                    | 34 |
| FORMFIT® PRO BACK .....                           | 35 |
| FORMFIT® BACK VS.....                             | 36 |
| FORMFIT® BACK SUPPORT AIR .....                   | 37 |
| FORMFIT® BACK SUPPORT .....                       | 38 |
| HIP .....                                         | 39 |
| REBOUND® HIP .....                                | 41 |
| UNLOADER® HIP .....                               | 43 |
| UPPER EXTREMITY .....                             | 45 |
| FORMFIT® SHOULDER BRACE .....                     | 47 |
| PREMIUM SHOULDER IMMOBILISER .....                | 48 |
| HIGH ARM SLING .....                              | 48 |
| PAEDIATRIC UNIVERSAL ARM SLING .....              | 49 |
| FORMFIT® PRO ELBOW .....                          | 50 |
| REBOUND® POST-OP ELBOW .....                      | 51 |
| FORMFIT® TENNIS ELBOW .....                       | 53 |
| FORMFIT® PRO WRIST .....                          | 54 |
| FORMFIT® WRIST / FORMFIT® WRIST & FOREARM .....   | 55 |
| SOFTPRO® FUNCTIONAL RESTING WHFO .....            | 56 |
| SOFTPRO® GRIP WHFO .....                          | 57 |
| FORMFIT® UNIVERSAL WRIST & THUMB .....            | 58 |
| FORMFIT® THUMB SPICA .....                        | 59 |
| WRIST / THUMB BRACE .....                         | 60 |
| ELASTIC THUMB SPICA .....                         | 60 |
| SPECTRA® WRIST BRACE .....                        | 61 |

# ÖSSUR BRACING & SUPPORTS CATALOGUE

|                                                                  |            |
|------------------------------------------------------------------|------------|
| <b>KNEE</b> .....                                                | <b>63</b>  |
| UNLOADER ONE® X .....                                            | 65         |
| UNLOADER ONE® .....                                              | 71         |
| UNLOADER ONE® CUSTOM .....                                       | 74         |
| UNLOADER ONE® PLUS .....                                         | 75         |
| FORMFIT® PRO KNEE OA .....                                       | 77         |
| REBOUND® DUAL .....                                              | 78         |
| REBOUND® DUAL ST .....                                           | 82         |
| REBOUND® DUAL CUSTOM .....                                       | 85         |
| REBOUND® CARTILAGE .....                                         | 87         |
| CTI® OA .....                                                    | 90         |
| REBOUND® PCL .....                                               | 92         |
| REBOUND® ACL .....                                               | 96         |
| REBOUND® DUAL BASIC AND REBOUND® DUAL BASIC ST .....             | 100        |
| OA EASE .....                                                    | 102        |
| CTi®3 .....                                                      | 104        |
| CTI® CUSTOM .....                                                | 106        |
| FLEX™ .....                                                      | 109        |
| C180™ ROCKET .....                                               | 111        |
| REBOUND® POST-OP KNEE .....                                      | 112        |
| FORMFIT POST-OP KNEE .....                                       | 113        |
| REBOUND® KNEE BRACE .....                                        | 114        |
| REBOUND® KNEE IMMOBILIZER UNIVERSAL .....                        | 115        |
| FORMFIT® PRO KNEE .....                                          | 116        |
| FORMFIT® TRACKER .....                                           | 118        |
| FORMFIT® PRO KNEE FLITE .....                                    | 119        |
| FORMFIT® HINGED LATERAL J .....                                  | 120        |
| FORMFIT® KNEE HINGED .....                                       | 121        |
| FORMFIT® KNEE MCL .....                                          | 122        |
| FORMFIT® KNEE ROM BRACE .....                                    | 123        |
| NEOPRENE HINGED KNEE SUPPORT .....                               | 124        |
| NEOPRENE KNEE SUPPORT WITH STABILISED PATELLA .....              | 124        |
| NEOPRENE WRAPAROUND HINGED KNEE SUPPORT .....                    | 125        |
| NEOPRENE KNEE SLEEVE .....                                       | 126        |
| MEASUREMENT FORM CTi® CUSTOM BRACE .....                         | 127        |
| MEASUREMENT FORM UNLOADER ONE® CUSTOM / REBOUND® ACL & PCL ..... | 130        |
| <b>FOOT AND ANKLE</b> .....                                      | <b>131</b> |
| REBOUND® DIABETIC WALKER .....                                   | 133        |
| DH OFFLOADING WALKER® .....                                      | 135        |
| REBOUND® AIR WALKER .....                                        | 136        |
| FORMFIT® WALKER AIR .....                                        | 138        |
| FORMFIT® WALKER .....                                            | 139        |
| PAEDIATRIC WALKER .....                                          | 140        |
| AFO DYNAMIC .....                                                | 141        |
| AFO LIGHT .....                                                  | 143        |

# ÖSSUR BRACING & SUPPORTS CATALOGUE

---

|                                        |            |
|----------------------------------------|------------|
| PROLITE AFO .....                      | 144        |
| EXTRA STRONG PROLITE .....             | 144        |
| REBOUND FOOT-UP® .....                 | 145        |
| FOOT-UP® .....                         | 146        |
| FORMFIT® PRO ANKLE .....               | 147        |
| FORMFIT® ANKLE BRACE .....             | 149        |
| FORMFIT® ANKLE WITH SPEEDLACE .....    | 150        |
| FORMFIT® ANKLE STIRRUPS .....          | 151        |
| EXOFORM® DORSAL NIGHT SPLINT .....     | 152        |
| FORMFIT® NIGHT SPLINT .....            | 153        |
| DH OFFLOADING POST-OP STYLE SHOE ..... | 154        |
| MESH TOP POST-OP SHOE .....            | 155        |
| SOFT TOP POST-OP SHOE .....            | 155        |
| CANVAS ROCKER BOTTOM CAST SHOE .....   | 156        |
| <b>CASTING &amp; TRACTION .....</b>    | <b>157</b> |
| TECHFORM® PREMIUM .....                | 159        |
| NOTES .....                            | 160        |



## WELCOME TO THE ÖSSUR® UK BRACING AND SUPPORTS PRODUCT CATALOGUE 2024

---

Össur's diverse range of bracing and supports products offer expert treatment for a range of conditions. Combining scientifically advanced products with high standards of quality assurance and customer care, our destination is Life Without Limitations. Össur's extensive orthotics range offers some of the most effective orthopaedic solutions available today.

Your business is important to us. We are committed to offering products and services that will enhance your operation, from highly innovative products backed by exceptional service and educational support, to more common items that represent great value for money.

We hope this catalogue will guide you through our comprehensive range of products.

Should you need any assistance, please do not hesitate to contact our Customer Care Team on 03450 065 065, email us at [ukinfo@ossur.com](mailto:ukinfo@ossur.com) or visit [www.ossur.com](http://www.ossur.com).

Yours sincerely



Peter Gill, UK Sales Director, Össur North Europe

## CUSTOMER CARE TEAM

Our Customer Care team really understand the UK healthcare industry and the NHS in particular – we've been working and changing with it for years.

Tel: [03450 065 065](tel:03450065065)

Enquiries: [ukinfo@ossur.com](mailto:ukinfo@ossur.com)

Orders: [ossuruk@ossur.com](mailto:ossuruk@ossur.com)



If you would like to know more about consolidated ordering you can contact our friendly and experienced Customer Care team or your Area Sales Manager who will help set you up if you like the concept.

Our sales team are well trained and are there to help you in any way they can – training staff, fitting products, solving logistical problems and demonstrating new technology.

As a global company we are able to innovate and advance our industry, investing heavily in research and development, we bring a continuous stream of new products, new ideas and materials to our customers in the UK.

These innovations can save cost, some as alternative treatments to surgery and others by reducing inpatient stay.

We're a pretty friendly bunch who work to a set of values that guide our actions, rather than a rule book.

### HONESTY

Honest with our customers  
and honest with each other.

### FRUGALITY

Use resources responsibly,  
ours and the environment's.


### COURAGE

Do the right thing and  
not the easy thing.

Össur UK are working towards implementation of GS1 as part of the NHS eProcurement Strategy.



## SALES AREA CONTACTS



**1 Scotland, Northern Ireland & North East**  
Tel: +44 (0) 7551 170897

**2 Pennines & Yorkshire**  
Tel: +44(0) 7551 170897

**3 North West & Republic of Ireland**  
Tel: +44 (0) 7968 950286

**4 South Wales & South West**  
Tel: +44 (0) 7827 950195

**5 South Central & East Mids**  
Tel: +44(0) 7408 802712

**6 North London & Anglia**  
Tel: +44 (0) 7408 809622

**7 South East, The Channel Islands & Private Health**  
Tel: +44(0) 7795 105748

## ORDERING INFORMATION

---

- Please quote your account number when placing an order
- Please quote model numbers and any other specification required when ordering
- Please quote your purchase order number
- Orders placed by telephone then faxed, posted or emailed should be clearly marked with 'CONFIRMATION ORDER' to avoid risk of duplication
- Össur® participate in the NHS PURCHASING CARD SCHEME and can accept orders / payment by this method

### RESPONSE TO ORDERS

We are proud of our service to our customers and maintain large stocks of products to expedite delivery. Orders for stock items are usually processed and despatched the same day.

### METHOD OF DISPATCH

Your order will be dispatched by the most economical means available, with regard to weight / volume and the degree of urgency. Unless otherwise instructed, orders will be sent by standard 2 - 3 day delivery. Timed deliveries are available on request. Please contact Customer Care for pricing.

- Next working day before 2pm only\*
- Next working day before 5pm\*

*\* where available at an additional cost*

### CUSTOM PRODUCTS

Össur® UK Ltd., through our own manufacturing facility, have the expertise and resources to create one-off custom requirements for products such as the Unloader One Custom, Rebound PCL Custom, Rebound ACL Custom, CTi Custom. Where this is the case, please allow 10 working days for delivery.

### INCORRECT DELIVERY

On receipt of the goods the customer is required to check the consignment. In the unlikely event that the type or quantity of goods is incorrect then you are requested to notify our Customer Care Department within 7 working days of the date of delivery. After 7 working days Össur® UK Ltd cannot accept responsibility for an alleged incorrect delivery.

### RETURN OF GOODS POLICY

You can return your unwanted products to us free of charge for up to 365 days after you purchased them. Authorisation must be issued by Össur® UK Ltd prior to goods being returned, please contact Customer Care for a return reference number and ensure sections 1-5 of Return Advice Note (on the back of the dispatch note) are fully completed. Soiled, custom manufactured, modified or obsolete goods may not be returned. Goods accepted for return must be in a resalable condition, guaranteed infection free and in the original packaging. Össur® UK Ltd reserves the right to charge delivery costs and / or re-stocking fee of up to 20% on returned goods.

### 30-DAY MONEY BACK GUARANTEE ON THE UNLOADER ONE AND UNLOADER ONE X

A product that is risk-free; if the product doesn't work for a patient, you get a full refund.

### ÖSSUR NO QUIBBLES GUARANTEE.

At Össur we know that accidents can happen. To put your mind at ease and to avoid unexpected costs, we are pleased to announce that The Unloader One X is covered by our Össur Care, No Quibbles Guarantee.

### TERMS AND CONDITIONS OF SALE

Our full Terms and Conditions of Sale can be provided upon request.

The product range and availability may change during the period that the current catalogue runs.

The returns policy does not apply to custom made or modified goods. Statutory rights are not affected. Please check with Customer Care at time of sale if in doubt.

Telephone Orders: 03450 065 065

Facsimile Orders: 0330 808 3926

E-Mail: [ossuruk@ossur.com](mailto:ossuruk@ossur.com) for orders, [ukinfo@ossur.com](mailto:ukinfo@ossur.com) for enquiries

### SINGLE PATIENT USE

It is a condition of Sale that any orthopaedic product purchased from Össur® UK Ltd is restricted to 'Single use Only' by the originally-fitted patient in order to protect the care provider and the patient against potentially adverse consequences of infectious disease transmission or material instability resulting from previous use by the original wearer. Any expressed or implied warranties are voided if the orthosis is re-used for another patient. Any orthosis purchased from Össur® UK Ltd must be prescribed and applied by a qualified clinician to ensure that it meets the needs of the particular patient and accomplishes the desired results.





## ÖSSUR OA SOLUTIONS

### Enhance quality of life and mobility for people living with osteoarthritis.

These product solutions have been built on decades of experience and are designed to enhance the quality of life of people living with osteoarthritis by reducing pain and improving mobility. Through innovation leadership, OA expertise and partnerships with medical professionals, we offer comprehensive, proven and effective treatment options.

## UNLOADER®

The Össur Unloader® brand name has become synonymous with osteoarthritis bracing, representing OA Solutions' technology as the industry standard, with more than 15 clinical studies demonstrating its efficacy.

Össur strives to maintain its position as the innovation leader for OA Solutions by continuing to develop simple, yet effective clinically proven treatment solutions that increase people's quality of life.



## ÖSSUR INJURY SOLUTIONS

For people recovering from fractures, ligament injuries or in need of post-operative treatment solutions.

These product solutions have been built on decades of experience and are designed to enhance the healing process of bone and soft tissue injuries. Through innovation leadership, proven biomechanical approaches, and partnerships with leading medical professionals, we strive to offer comprehensive and effective treatment options.

Össur Injury Solutions are globally recognised, primarily under the brand names CTi®, Rebound® and Miami®, all of which are developed to assist healing and increase people's quality of life.

Our expertise in Injury Solutions has driven us to take the next step forward by developing indication-specific products which are designed to optimise healing while maintaining function and mobility. These Functional Healing products incorporate Össur proprietary technology and are designed to enhance clinically relevant outcomes.



FUNCTIONAL HEALING –  
WHERE MOBILITY MEANS RECOVERY.









## TABLE OF CONTENTS

---

|                                                   |    |
|---------------------------------------------------|----|
| SPINE .....                                       | 15 |
| MIAMI J® .....                                    | 17 |
| MIAMI J® SELECT .....                             | 19 |
| MIAMI JTO® THORACIC EXTENSION .....               | 20 |
| OCCIAN® BACK .....                                | 21 |
| MIAMI Jr® .....                                   | 22 |
| PHILADELPHIA® ADJUSTABLE TRACHEOTOMY COLLAR ..... | 23 |
| PHILADELPHIA® TRACHEOTOMY COLLAR .....            | 24 |
| NECLOC® .....                                     | 25 |
| NECLOC® KIDS .....                                | 26 |
| MIAMI LSO® .....                                  | 27 |
| MIAMI TLSO 464 .....                              | 29 |
| MIAMI TLSO 456 .....                              | 31 |
| MIAMI LUMBAR® POSTEO .....                        | 33 |
| PAPOOSE® .....                                    | 34 |
| FORMFIT® PRO BACK .....                           | 35 |
| FORMFIT® BACK VS .....                            | 36 |
| FORMFIT® BACK SUPPORT AIR .....                   | 37 |
| FORMFIT® BACK SUPPORT .....                       | 38 |



## MIAMI J®

Known as the superior c-spine immobiliser through multiple independent studies, the Miami J's patented design meets the specialised needs of individuals with cervical injuries. Bioengineered to minimise pressure points in key known areas: Chin, occiput, trapezius, and clavicle. Miami J offers features that enhance compliance, while providing immobilisation, inhibition of skin deterioration, and patient comfort.



### INDICATIONS FOR USE

Situations requiring gross immobilisation of the cervical spine, these may include:

- Pre and post c-spine surgery
- C-spine precaution for trauma patients


### FEATURES AND BENEFITS

- Sorbatex padding is antibacterial and clinically proven to inhibit microbial growth, which enhances comfort and helps assure patient compliance
- Comfort adjustment buttons for occiput and sternal pad is removable when supine, lying and swallowing
- MRI, CT and X-ray lucent, MR compatible

### REQUIRED MEASUREMENTS

- Sizing based on 6 phenotypes to ensure correct fit and optimum immobilisation

Phenotype-driven sizing system makes size selection simple while drastically reducing margin of error in fitting a patient. Each collar has a bioengineered front with a specific back designed to work together for superior immobilisation. For paediatric sizing, please see Miami Jr in this catalogue.



# MIAMI J®

## ORDERING INFORMATION

### MIAMI J

| Size        | Collar set* | Collar  | Collar Back | Collar Front |
|-------------|-------------|---------|-------------|--------------|
| Super Short | MJR-200S    | MJ-200S | MJ-400BK    | MJ-200SFRT   |
| Stout       | MJR-200L    | MJ-200L | MJ-200LBK   | MJ-200LFRT   |
| XS          | MJR-250     | MJ-250  | MJ-250BK    | MJ-250FRT    |
| Short       | MJR-300     | MJ-300  | MJ-300BK    | MJ-300FRT    |
| Regular     | MJR-400     | MJ-400  | MJ-400BK    | MJ-400FRT    |
| Tall        | MJR-500     | MJ-500  | MJ-500BK    | MJ-500FRT    |
| Bariatric   | MJR-200B    |         |             |              |

\*Collar Sets come with one additional set of replacement pads.

Miami J Front and back panels may be interchanged to accommodate patient anatomy, e.g. a Short front may be used with a Tall back. Choose collar size based on patient neck height, not patient body size. Super Short collars are packaged with regular back panels

## ACCESSORIES

### REPLACEMENT PADS

| Size        | Part#    |
|-------------|----------|
| Super Short | MJP-100  |
| Stout       | MJP-200L |
| XS          | MJP-250  |
| Short       | MJP-100  |
| Regular     | MJP-100  |
| Tall        | MJP-100  |
| Bariatric   | MJP-200B |

## MIAMI J® SELECT

The Miami J Select cervical collar by Össur offers the range of height adjustability in the market with a single, universal design that replicates the comfort, immobilisation and quality you expect from a Miami J product. Featuring a patented, first-to-market locking mechanism to encourage patient compliance and anti-microbial Sorbatex padding, Miami J Select is designed to support successful patient outcomes.



### INDICATIONS FOR USE

- Situations requiring gross immobilisation of the cervical spine, these may include:
- Pre and post c-spine surgery
- C-spine precaution for trauma patients

### FEATURES AND BENEFITS

- Easy height adjustability with a tamper-proof locking mechanism – ‘Left to Lock’
- Anti-microbial & biocompatible Sorbatex padding prevents bacterial build-up
- Universal design that doesn’t compromise immobilisation or comfort
- 7 precise fitting sizes with opportunity to fine-tune
- Miami J® inspired chin tray
- One-piece, anatomically inspired collar back with Flex Edge® designed to reduce occipital pressure
- Easy to don and doff with patented Reproducible Fit option
- Intuitive blue Patient Touch-Points to send patients home with confidence – Blue is You™
- Simplified comfort adjustments with integrated and patented Sternal-Relief Dial
- Large tracheal opening – improved function for ICU patients
- Easy to clean and replace multi-coloured padding
- X-ray and CT lucent, MR-safe
- Compatible with Occian® Back & Miami JTO® accessories

### ORDERING INFORMATION

#### COLLAR SET

| Description | Part#    |
|-------------|----------|
| Collar set  | MJSR-101 |

#### MIAMI J SELECT COLLAR

| Description           | Part#   |
|-----------------------|---------|
| Miami J Select Collar | MJS-101 |

### OTHER ACCESSORIES

| Description                                         | Part#    |
|-----------------------------------------------------|----------|
| Miami J Select Replacement Pads                     | MJSP-101 |
| Miami J Select Tabs Reproducible Fit Tabs (10 pack) | MJST-101 |

## MIAMI JTO® THORACIC EXTENSION

The Miami JTO Thoracic Extension is designed in conjunction with leading spine surgeons, orthotists and nurses. The orthosis is used exclusively with all Miami J® collars when even greater stability at the cervical junctions is indicated.

Miami JTO Thoracic Extension is a versatile, tool-free device with a patented chin strut, providing more control than a collar alone can provide, designed to treat cervical and high-thoracic injuries.



### INDICATIONS FOR USE

- Low c-spine and high t-spine precaution for trauma patients
- Immobilisation for pre and post low c-spine and high t-spine surgery
- Other conditions requiring gross immobilisation of the low-cervical and/or high-thoracic-spine
- Conditions requiring gross immobilisation of the low cervical and/or high-thoracic-spine

### FEATURES AND BENEFITS

- Patented design with 'floating back' permits maximum immobilisation without the risk of compensatory motion
- Snap on and off without tools
- Sorbatex II padding is antibacterial and clinically proven to inhibit microbial growth, which enhances comfort and helps assure patient compliance
- MRI, CT and X-ray lucent

### REQUIRED MEASUREMENTS

- Universal, one size

### ORDERING INFORMATION

#### MIAMI JTO THORACIC EXTENSION

| Description                          | Part#   |
|--------------------------------------|---------|
| Miami JTO                            | MT-900  |
| Miami JTO with replacement softgoods | MTR-900 |

Compatible with all sizes of the Miami J®, Miami J® Select and the P3 Miami Jr®. Miami J collar is not included with Miami JTO Thoracic Extension. Miami J collar needs to be purchased also

### ACCESSORIES

| Description            | Part#   |
|------------------------|---------|
| Replacement softgoods* | MTP-950 |

\*Contains 2 x 4" Shoulder Strap Pads / 4 x 8" Waist Strap Pads / 4 x Snap Rivet Sets

## OCCIAN® BACK

The Occian Back is designed exclusively for the specialised needs of the multi-system trauma patient in the ICU. Universally sized, specifically engineered and designed with advanced, Intuitech™ pressure-relieving memory foam, the Occian Back helps caregivers manage c-spine patients who are at high risk of occipital breakdown.



### INDICATIONS FOR USE

- Patients not awake and alert within 24 - 48 hours after admission
- Patients with pre-existing skin injuries and/or lacerations to the head/neck
- Burn patients who must wear a collar
- Elective surgery patients desiring night time post-operative comfort

### FEATURES AND BENEFITS

- Intuitech pressure relieving memory foam
- Maintains c-spine stability
- Provides unique pressure-relieving comfort by replacing the back section of a Miami J Collar.
- Easy to clean
- MRI, CT and X-ray lucent

### REQUIRED MEASUREMENTS

- Universal, one size

### ORDERING INFORMATION

| Description        | Part#   |
|--------------------|---------|
| Occian Back        | ACB-100 |
| Occian Back Junior | ACB-200 |

Compatible with all sizes of the Miami J®, Miami J® Advanced, Miami J® Select and the Miami Jr®. Collar front panel not included.



## MIAMI Jr®

Miami Jr is engineered to fit the unique anatomies of patients aged 12 and under. The Miami Jr line, which spans the growth cycle, is sized according to age. Össur is the only company that has a paediatric line of collars based on the anatomical differences between adults and children. We realise that children are not just “little adults” - they require innovative design.



## INDICATIONS FOR USE

Situations requiring gross immobilisation of the cervical spine, these may include:

- pre and post c-spine surgery
- c-spine precaution for trauma patients

## FEATURES AND BENEFITS

- Designed exclusively for children 12 years and under
- Sorbatex II pads make the collar “skin friendly” and protect skin integrity during extended wear
- MRI, CT and X-ray lucent, MR compatible

## REQUIRED MEASUREMENTS

- Age

## ORDERING INFORMATION

### MIAMI JR

| Size          | Collar set* | Collar |
|---------------|-------------|--------|
| 0 - 6 months  | MJR-P0      | MJ-P0  |
| 6 - 24 months | MJR-P1      | MJ-P1  |
| 2 - 6 years   | MJR-P2      | MJ-P2  |
| 6 - 12 years  | MJR-P3      | MJ-P3  |

\*Collar sets come with one additional set of replacement pads.

## ACCESSORIES

### REPLACEMENT PADS

| Size          | Part#  |
|---------------|--------|
| 0 - 6 months  | MJP-P0 |
| 6 - 24 months | MJP-P1 |
| 2 - 6 years   | MJP-P2 |
| 6 - 12 years  | MJP-P3 |

## PHILADELPHIA® ADJUSTABLE TRACHEOTOMY COLLAR

Philadelphia tracheotomy collars have been an industry standard for decades. Philadelphia Adjustable is constructed from the same high-quality materials as the original, and combines convenient size adjustability with immobilisation and comfort.

### INDICATIONS FOR USE

Situations requiring gross immobilisation of the cervical spine, these may include:

- Pre and post c-spine surgery
- C-spine precaution for trauma patients

### FEATURES AND BENEFITS

- Circumference adjustable and height adjustable with permanent lock
- Anatomically shaped Plastazote® foam leaving clearance for cervical veins
- Well ventilated and water-resistant for bathing and hydro-therapy
- MRI, CT and X-ray lucent, MR compatible

### REQUIRED MEASUREMENTS

- Universal, one size



## ORDERING INFORMATION

| Description                                | Part#  |
|--------------------------------------------|--------|
| Philadelphia Adjustable Tracheotomy Collar | PHP-AA |

## PHILADELPHIA® TRACHEOTOMY COLLAR


The Philadelphia Tracheotomy Collar is a hospital and rehabilitation collar that is designed to maintain neutral alignment of the cervical spine after injury to the neck. The collar contains a large tracheal opening, allows carotid pulse monitoring and emergency tracheotomies.

### INDICATIONS FOR USE

- C-spine precaution for trauma patients
- Immobilisation for pre c-spine surgery
- Other conditions requiring gross immobilisation of the mid-cervical spine

### FEATURES AND BENEFITS

- Anterior and posterior rigid plastic reinforcement supports limit movement
- Total cervical arch support maintains neutral alignment of the cervical spine
- Water-resistant for bathing and hydro-therapy
- Easily customised for individual requirements
- Does not contain natural rubber latex, non-toxic and hypoallergenic Plastazote® foam reduces skin irritation



### REQUIRED MEASUREMENTS

- Neck circumference and height from tip of chin to sternal notch

## ORDERING INFORMATION

| Size | Neck circumference | 6cm height | 8cm height | 11cm height | 13cm height |
|------|--------------------|------------|------------|-------------|-------------|
| XS   | < 25cm             | PHP-T2XS   |            |             |             |
| S    | 25 - 33cm          | PHP-T2S    | PHP-T3S    | PHP-T4S     | PHP-T5S     |
| M    | 33 - 41cm          | PHP-T2M    | PHP-T3M    | PHP-T4M     | PHP-T5M     |
| L    | 41 - 48cm          | PHP-T2L    | PHP-T3L    | PHP-T4L     | PHP-T5L     |
| XL   | > 48cm             | PHP-T2XL   | PHP-T3XL   | PHP-T4XL    | PHP-T5XL    |

| Size       | Neck circumference | Part#  |
|------------|--------------------|--------|
| Infant     | 15 - 20cm          | PHP-TI |
| Paediatric | 20 - 28cm          | PHP-TP |

## ACCESSORIES

### PHILADELPHIA COOLMAX COTTON LINER

| Size  | Neck circumference | Part#   |
|-------|--------------------|---------|
| Child | 20 - 28cm          | PHP-LCP |
| S     | 25 - 33cm          | PHP-LCS |
| M     | 33 - 41cm          | PHP-LCM |
| L     | 41 - 48cm          | PHP-LCL |

### PHILADELPHIA TRACHEOTOMY ATTACHMENTS

| Size  | Stabiliser |
|-------|------------|
| Adult | PHP-OS     |

# NECLOC®

## INDICATIONS FOR USE

Situations requiring gross immobilisation of the cervical spine, these may include:

- Pre and post c-spine surgery
- C-spine precaution for trauma patients

## FEATURES AND BENEFITS

- Ergonomic two-piece clamshell design and retaining strap secure flexion during application
- Five unique phenotype-driven sizes make it easy to select the correct size
- Minimal manipulation of patient during application
- X-ray and CT lucent. MR safe




## REQUIRED MEASUREMENTS

- Sizing based on 5 phenotypes to ensure correct fit and optimum immobilisation

## SIZING

- Phenotype-driven sizing system makes size selection simple while drastically reducing margin of error in fitting a patient. Each collar has bioengineered size specific front and back sections designed to work together for superior immobilisation. For paediatric sizing, please see NecLoc Kids.



## ORDERING INFORMATION

| Size  | Part#    |
|-------|----------|
| Stout | NL-200EL |
| XS    | NL-250E  |
| S     | NL-300E  |
| M     | NL-400E  |
| L     | NL-500E  |

## NECLOC® KIDS

### INDICATIONS FOR USE

Situations requiring gross immobilisation of the cervical spine, these may include:

- Pre and post c-spine surgery
- C-spine precaution for trauma patients

### FEATURES AND BENEFITS

- Ergonomic two-piece clamshell design and retaining strap secure flexion during application
- Three paediatric sizes ensure optimal airway and cervical spine immobilisation
- Minimal manipulation of patient during application to children
- Meets the unique developing anatomical requirements of children
- X-ray and CT lucent. MR Safe



### REQUIRED MEASUREMENTS

- Age

### ORDERING INFORMATION

| Age          | Part# |
|--------------|-------|
| 0 - 2 years  | NL-P1 |
| 2 - 6 years  | NL-P2 |
| 6 - 12 years | NL-P3 |

## MIAMI LSO®

Miami LSO is a spinal orthosis designed for pain relief, protection of injured ligaments or muscles and post-surgical immobilisation. It provides low-profile, anatomical immobilisation, easy fit and wear, options for modular support and an advanced pulley system with a length-adjustable pulley cord.



### INDICATIONS FOR USE

Conditions requiring gross immobilisation of the trunk in the lumbar region such as:

- Stable non-displaced spinal fractures (L1-L5)
- Spinal Stenosis
- Herniated Disks
- Degenerative Spinal Pathologies
- Spondylolysis
- Spondylolisthesis

### FEATURES AND BENEFITS

- Slim, breathable material
- Proven Miami double-pulley system with length-adjustable pulley cord
- Intuitive fitting without trimming with fold-over belt arms & Precise-Size™ landmarks
- Optional fit tool for perfect fitting directly on the patient
- Universal belt sizing accommodates waists up to 178cm
- New Universal Posterior panel to accommodate more patients with less inventory
- Contoured posterior panels sized S-XL for a precise, anatomical fit
- Quick-Fit pocket enclosure and blue highlights for patient ease of use

### REQUIRED MEASUREMENTS

- Waist circumference



# MIAMI LSO®

## ORDERING INFORMATION

### MIAMI LSO BELT

| Size                    | Waist Circumference | Part#   |
|-------------------------|---------------------|---------|
| 6" Universal (Slim)     | 66-127cm            | MLSO-1E |
| 8" Universal (Standard) | 66-127cm            | MLSO-2E |



### OPTIONAL BELT EXTENDERS (2 PCS)

| Size                         | Waist Circumference | Part#     |
|------------------------------|---------------------|-----------|
| 6" (Slim) Belt Extension     | 127-177cm           | MLSO-1EXT |
| 8" (Standard) Belt Extension | 127-177cm           | MLSO-2EXT |

### STANDARD POSTERIOR PANEL

| Size                     | Part#    |
|--------------------------|----------|
| Standard Posterior Panel | MLSO-P56 |



### POSTERIOR PANEL WITH LATERAL SUPPORTS

| Size | Waist Circumference | Part#    |
|------|---------------------|----------|
| S    | 66-86cm             | MLSO-P12 |
| M    | 86-96cm             | MLSO-P22 |
| L    | 96-127cm            | MLSO-P32 |
| XL   | 127-178 cm          | MLSO-P42 |



\*Non-lateral and non-lordotic panels available on request.

### ANTERIOR PANEL

| Size                | Part#    |
|---------------------|----------|
| 6" Panel (Slim)     | MLSO-PA6 |
| 8" Panel (Standard) | MLSO-PA8 |



### QUICK FIT TOOL

| Size           | Part#       |
|----------------|-------------|
| Quick fit tool | B-040100000 |



## MIAMI TLSO 464

The Miami TLSO from Össur directly addresses the needs of patients and practitioners, simplifying preparation without compromising on fit or performance. This easy-to-use pre-assembled Miami TLSO brace offers three configurations in one.

**Compliance:** Easily applied in over the shoulder or underarm configuration, donning and doffing is further simplified with the conveniently located quick release buckles.

**Adjustability:** The universal belt designs allows for setting the belt length on the patient for a more accurate sizing experience.

**Comfort:** Truly shaped to fit the anatomy, the anatomically contoured design better immobilises the spine while maximising patient comfort.

**Ease of Use:** The Miami TLSO brace requires minimal preparation and is designed to provide optimal patient comfort.



## INDICATIONS FOR USE

Conditions requiring gross immobilisation of the thoracic and lumbar spine. Such conditions may include:

- Post-surgical immobilisation
- Spinal support TLSO attachment is up to T7
- Spinal stenosis
- Herniated disc
- Degenerative spinal fractures
- Spondylolithesis
- Spondylolysis

## CONTRAINDICATIONS

- Unstable, displaced fractures

## FEATURES AND BENEFITS

- Significantly superior comfort and ease of donning/doffing
- Fully assembled for fitting with ease
- Accommodates heights 155-193 cm
- Posterior thoracic extension option with ergonomic axillary strapping
- Step-down to an LSO as one progresses through rehabilitation
- One-touch telescoping anterior-thoracic extension option with single or double-hinge
- The convenience of a one-size Universal or individualise with a uniquely-sized anatomical design
- Improve comfort and immobilisation with multiple strapping configurations and removable anchors

## REQUIRED MEASUREMENTS

- Waist circumference

## MIAMI TLSO 464

### ORDERING INFORMATION

#### ATE/PTE, LORDOTIC

| Size      | Part#      | Waist circumference |
|-----------|------------|---------------------|
| Universal | MTLSO-2351 | 66 - 177cm          |
| S         | MTLSO-2111 | 66 - 86cm           |
| M         | MTLSO-2121 | 86 - 107cm          |
| L         | MTLSO-2131 | 107 - 127cm         |
| XL        | MTLSO-2141 | 127 - 177cm         |

\*Non-lateral and non-lordotic available on request

\*Belt extensions are required for waist circumferences larger than 127cm

#### ATE/PTE, LORDOTIC, SLIM

| Size      | Part#      | Waist circumference |
|-----------|------------|---------------------|
| Universal | MTLSO-1351 | 66 - 177            |
| S         | MTLSO-1111 | 66 - 86             |
| M         | MTLSO-1121 | 86 - 107            |
| L         | MTLSO-1131 | 107 - 127           |
| XL        | MTLSO-1141 | 127 - 177           |

\*Belt extensions are required for waist circumferences larger than 127cm

### OTHER ACCESSORIES

| Description                                                                                              | Part#       |
|----------------------------------------------------------------------------------------------------------|-------------|
| 5MM Hex Long Wrench                                                                                      | MTLSO-0004  |
| Axillary strap accessory kit                                                                             | MTLSO-0003  |
| Quick fit tool                                                                                           | B-040100000 |
| Belt extension, standard<br>(Recommended for waist circumferences 127-167cm. Fits size XL of Miami LSO.) | MLSO-2EXT   |
| Belt extension, slim<br>(Recommended for waist circumferences 127-167cm. Fits size XL of Miami LSO.)     | MLSO-1EXT   |

## MIAMI TLSO 456

The Miami TLSO 456 from Össur directly addresses the needs of patients and practitioners, simplifying preparation without compromising on fit or performance. This easy-to-use pre-assembled Miami TLSO brace offers the convenience of a one-size Universal or a uniquely-sized anatomical design.

**Compliance:** Blue is You™ patient touch-points help send patients home with confidence.

**Adjustability:** The universal belt design allows for setting the belt length on the patient for a more accurate sizing experience. Heat moldable posterior panels allow a customised fit, while the modular design allows step-down support to an LSO.

**Comfort:** Featuring the Miami LSO™ Advanced Compression pulley system with lumbar compression you can feel, and an anatomically contoured design that better immobilises the spine while maximising patient comfort.

**Ease of Use:** Easily applied with over-the-shoulder straps, the Miami TLSO brace requires minimal preparation and is designed to provide optimal patient comfort.



## INDICATIONS FOR USE

Conditions requiring gross immobilisation of the thoracic and lumbar spine. Such conditions may include:

- Post-surgical immobilisation
- Spinal support TLSO attachment is up to T7
- Spinal stenosis
- Herniated disc
- Degenerative spinal fractures
- Spondylolithesis
- Spondylolysis

## CONTRAINDICATIONS

- Unstable, displaced fractures

## FEATURES AND BENEFITS

- Significantly superior comfort and ease of donning/doffing (when compared to Miami Lumbar TLSO)
- Fully assembled for fitting with ease
- Accommodates heights 155-193cm
- Step-down to an LSO as one progresses through rehabilitation
- The convenience of a one-size Universal or individualise with a uniquely-sized anatomical design

## REQUIRED MEASUREMENTS

- Waist circumference

## MIAMI TLSO 456

### ORDERING INFORMATION

#### PTE LORDOTIC

| Size      | Part#      | Waist circumference |
|-----------|------------|---------------------|
| Universal | MTLSO-2352 | 66 - 177cm          |
| S         | MTLSO-2112 | 66 - 86cm           |
| M         | MTLSO-2122 | 86 - 107cm          |
| L         | MTLSO-2132 | 107 - 127cm         |
| XL        | MTLSO-2142 | 127 - 177cm         |

\*Belt extensions are required for waist circumferences larger than 127cm

#### PTE LORDOTIC, SLIM

| Size      | Part#      | Waist circumference |
|-----------|------------|---------------------|
| Universal | MTLSO-1352 | 66 - 177cm          |
| S         | MTLSO-1112 | 66 - 86cm           |
| M         | MTLSO-1122 | 86 - 107cm          |
| L         | MTLSO-1132 | 107 - 127cm         |
| XL        | MTLSO-1142 | 127 - 177cm         |

\*Non-lateral and non-lordotic available on request.

\*Belt extensions are required for waist circumferences larger than 127cm

### OTHER ACCESSORIES

| Description                                                                                             | Part#       |
|---------------------------------------------------------------------------------------------------------|-------------|
| Quick fit tool                                                                                          | B-040100000 |
| Belt extension, standard<br>(Recommended for waist circumferences 127-167cm. Fits size XL of Miami LSO) | MLSO-2EXT   |
| Belt extension, slim<br>(Recommended for waist circumferences 127-167cm. Fits size XL of Miami LSO)     | MLSO-1EXT   |

## MIAMI LUMBAR® POSTEO

Designed with patient comfort and compliance in mind, the Miami Lumbar Posteo provides functional relief of vertebral fracture pain, and facilitates return to mobility. The lumbar compression system offers both pain relief via active vertebral offloading and better grip of patient's hips for greater ease of use and compliance. The Miami Lumbar Posteo includes an innovative approach to ease of use and consistent donning.



### INDICATIONS FOR USE

- Back pain caused by stable vertebral compression fractures (VCF)
- Post-operative immobilisation

### FEATURES AND BENEFITS

- Active offloading via thoracic extension and lumbar compression
- Integrated pulley system for customised fit
- Trimmable axillary straps with adjustable tension for optimal stabilisation
- Cold-moldable aluminum PTE to fit different anatomies
- Extra cushioning around the aluminum strut for minimised skin pressure

### REQUIRED MEASUREMENTS

- Back length measured from the top of the sacrum up to C7 vertebrae, following the spinal curvature

### SIZING

Please note that the Miami Lumbar Posteo includes a universal lumbar belt which can be trimmed down to the appropriate size.

### ORDERING INFORMATION

#### MIAMI LUMBAR POSTEO

| Size | Back Length | Part#         |
|------|-------------|---------------|
| XS   | 37 - 44cm   | LMB-POSTEO-XS |
| S    | 44 - 51cm   | LMB-POSTEO-SM |
| M    | 49 - 56cm   | LMB-POSTEO-MD |
| L    | 55 - 61cm   | LMB-POSTEO-LG |

## PAPOOSE®

Papoose is an innovative design for newborns up to 3 months (Broselow Grey). This true cervical spine immobiliser cradles the baby's body, securely while gently immobilising the head and spine in the proper anatomical airway and spinal alignment to help achieve optimal patient outcome while allowing for continued intimacy between caregiver and child.



### INDICATIONS FOR USE

- C-Spine precaution for infants
- Immobilisation for pre and post c-spine surgery
- Other conditions requiring gross immobilisation of the infant cervical spine

### FEATURES AND BENEFITS

- X-ray & CT Lucent. MRI Safe

### REQUIRED MEASUREMENTS

- Universal, one size

### ORDERING INFORMATION

#### PAPOOSE

| Description                                 | Part#   |
|---------------------------------------------|---------|
| Papoose without front collar                | PO-100  |
| Papoose without front with replacement pads | POR-100 |
| Papoose with front collar                   | PA-100  |
| Papoose with front with replacement pads    | PAR-100 |

### ACCESSORIES

#### REPLACEMENT PADS

| Description                              | Part#   |
|------------------------------------------|---------|
| Papoose replacement pads (with front)    | PAP-100 |
| Papoose replacement pads (without front) | POP-100 |



## FORMFIT® PRO BACK

Össur Formfit® Pro Back provides relief for patients suffering from pains and strains in the lower back. Manufactured with our dynamic, 3D MotionTech™ knitting technology, our belts feature superior moisture wicking capabilities and are highly breathable. The patented back panel design utilises our StayBilizer™ technology to provide optimum shape and fit, while conveniently placed finger loops make donning easy. The SmartPad™ provides a targeted support and cushioning and can be removed if desired.



MotionTech™



StayBilizer™



SmartPad™



CoolVent™



## INDICATIONS FOR USE

- For lower back conditions that may benefit from compression and gross motion restriction, such as:
- Lower back pain (Lumbago)
- Lumbosacral sprain/strain
- Degenerative joint/disc disease


## CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

## FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit
- Removable lumbar SmartPad™ for targeted support and cushioning
- Lightweight, malleable aluminum stays are shaped to the anatomy of the lower back

## REQUIRED MEASUREMENTS



Measure circumference of the trunk at a level midway between the lowest ribs and the upper iliac crest. Standing.



## ORDERING INFORMATION

| Size | Waist circumference | Part#    |
|------|---------------------|----------|
| XS   | 70 - 80cm           | FP500001 |
| S    | 80 - 90cm           | FP500002 |
| M    | 90 - 100cm          | FP500003 |
| L    | 100 - 115cm         | FP500004 |
| XL   | 115 - 130cm         | FP500005 |

## FORMFIT® BACK VS

Formfit® Back VS is a lumbar support belt with versatile double strap. V-shaped metal stays in the back of the belt support the lumbosacral joint. The height of the of the second strap can be adjusted to apply the desired level of compression on the lower back.



### INDICATIONS FOR USE

For lower back conditions that may benefit from gross motion restriction, such as:

- Muscular low back pain
- Lumbosacral strains
- Chronic low back pain
- Sciatica

### CONTRAINDICATIONS

- Conditions requiring more rigid or higher immobilization (> D12) or a surgical procedure, unless medically indicated.
- Pelvic floor dysfunction

### FEATURES AND BENEFITS

- Removable, versatile double strap, adjustable in height
- Patented V-shaped stays
- Breathable materials
- Wide finger loops and hook & loop areas

### REQUIRED MEASUREMENTS

- Waist circumference

### ORDERING INFORMATION

#### FORMFIT BACK VS - 21 CM

| Size | Waist circumference | Part#  |
|------|---------------------|--------|
| XS   | 70 - 80 cm          | 209112 |
| S    | 80 - 90 cm          | 209113 |
| M    | 90 - 100 cm         | 209115 |
| L    | 100 - 110 cm        | 209117 |
| XL   | 110 - 125 cm        | 209118 |
| XXL  | 125 - 140 cm        | 209119 |

#### FORMFIT BACK VS - 26 CM

| Size | Waist circumference | Part#  |
|------|---------------------|--------|
| XS   | 70 - 80 cm          | 209122 |
| S    | 80 - 90 cm          | 209123 |
| M    | 90 - 100 cm         | 209125 |
| L    | 100 - 110 cm        | 209127 |
| XL   | 110 - 125 cm        | 209128 |
| XXL  | 125 - 140 cm        | 209129 |

### ACCESSORIES

| Description              | Part# |
|--------------------------|-------|
| Hot/Cold Therapy Gel Pad | GB28  |

## FORMFIT® BACK SUPPORT AIR

The Form Fit Back Support with Air provides compression and support for those suffering from lower back pain. The integrated pneumatic system is intuitive and conveniently located so that adjustments can be made as demands on the body change throughout the day. The double wrap design improves compression, while durable finger loops make applying the brace easy. Made from breathable materials, the Form Fit Back Support provides compression and support for those suffering from lower back pain. The double wrap design improves compression, while the durable finger loops make applying the brace easy.

### INDICATIONS FOR USE

Lower back conditions that may benefit from gross motion restriction, such as;

- Muscular low back pain
- Lumbosacral strains
- Chronic low back pain
- Sciatica

### FEATURES AND BENEFITS

- Self-contained pneumatic system is low-profile and easy to use
- Mesh panels and ventilation holes increase breathability and enhance comfort
- Overlapping double pull strap design provides additional support to the lumbar and abdominal areas
- Flexible stays provide optimum support in the lumbar and abdominal areas
- Gripper strips prevent migration
- Adjustable air bladder provides a custom fit and offers focused compression in the lumbar area

### REQUIRED MEASUREMENTS

- Waist circumference



## ORDERING INFORMATION

| Size | Waist circumference | Standard |
|------|---------------------|----------|
| S    | 63 - 76cm           | 209133   |
| M    | 76 - 86cm           | 209135   |
| L    | 86 - 97cm           | 209137   |
| XL   | 97 - 107cm          | 209138   |
| XXL  | 107 - 142cm         | 209139   |

## ACCESSORIES

| Description | Part# |
|-------------|-------|
| Gel pad     | GB28  |

## FORMFIT® BACK SUPPORT

Made from breathable materials, the Formfit Back Support provides compression and support for those suffering from lower back pain. The double wrap design improves compression, while the durable finger loops make applying the brace easy

### INDICATIONS FOR USE

- Muscular low back pain
- Lumbosacral strains
- Chronic low back pain
- Sciatica

### FEATURES AND BENEFITS

- Anti-migration system reduces slippage
- Double wrap design improves compression
- Finger loops for easy application
- Hot/Cold Gel pad available separately

### REQUIRED MEASUREMENTS

- Waist circumference



## ORDERING INFORMATION

| Size | Waist Circumference | Standard |
|------|---------------------|----------|
| S    | 25-30"              | 209143   |
| M    | 30-34"              | 209145   |
| L    | 34-38"              | 209147   |
| XL   | 38-42"              | 209148   |
| XXL  | 42-56"              | 209149   |
| XXXL | 56-60"              | 209151   |

## ACCESSORIES

| Description | Part# |
|-------------|-------|
| Gel pad     | GB28  |



HIP

---



# TABLE OF CONTENTS

---

|                     |    |
|---------------------|----|
| HIP .....           | 39 |
| REBOUND® HIP .....  | 41 |
| UNLOADER® HIP ..... | 43 |

## REBOUND® HIP

Rebound Hip delivers gentle controlling forces and effective functional ROM restriction in extension / flexion and abduction / adduction to promote optimal post-surgical outcomes. The low profile, lightweight design supports patient comfort and compliance. The universal frame and softgood are easy to stock, assemble and quickly fit.



### INDICATIONS FOR USE

Hip conditions that could benefit from motion restriction such as:

- Postoperative after (total) hip arthroplasty
- Postoperative after revision of (total) hip arthroplasty
- Postoperative hip arthroscopy
- Post-op following FAI arthroscopic surgery
- Post-op following gluteus medius repair surgery
- Post-op following hamstring repair surgery
- Inoperable hip abnormalities
- Preventively when there is risk of hip dislocation
- Conventional treatment after dislocation/reduction

### FEATURES AND BENEFITS

- Low-profile, lightweight design that supports patient comfort and compliance
- Quick-release pegs easily adjust the desired range-of-motion needed for the patient's optimal rehab
- Quick-fit with sliding mechanism on frame to adjust for patient's height
- Optional, adjustable ROM hinge enables adduction/abduction from 0° to 10°, 20° or 30°
- OR kit offers patient protection of ankle and perineum skin during intraoperative traction on a surgical table

### REQUIRED MEASUREMENTS

- Circumference at iliac crest height
- Length from greater trochanter to mid patella

### ORDERING INFORMATION

#### REBOUND HIP

| Description                       | Part#    |
|-----------------------------------|----------|
| Rebound Hip                       | 98013XRB |
| Rebound Hip with adjustable hinge | 98012XRB |

X = 0 (Preset to left) or 5 (Preset to right)

\*Frame: Universal; Softgood: Universal fitting to 122cm hip circumference



# REBOUND® HIP

## ACCESSORIES

| Description                                           | Part#       |
|-------------------------------------------------------|-------------|
| Rebound Hip alligator kit                             | 980501      |
| Rebound Hip v2 softgood kit                           | 980601      |
| Rebound Hip v2 belt only kit                          | 980701      |
| OR Kit                                                | 980301      |
| Belt Extender (Fitting 122 - 170cm hip circumference) | 980121RB    |
| Rotation Strap                                        | RHK-100024  |
| Fitting Tool                                          | B-040100000 |

## UNLOADER<sup>®</sup> HIP

Unloader Hip is designed to improve activities of daily living by optimising load dispersion for patients suffering from mild and moderate osteoarthritis of the hip, thus contributing to hip stability and improved mobility.

Unloader Hip delivers proprioceptive support during the gait cycle by providing compression with a proven pulley system design. If required it can decrease the internal rotation and adduction moment with its unique and clinically proven Dynamic Rotation Straps™.



### INDICATIONS FOR USE

Conditions benefitting from hip repositioning and/or compression around the greater trochanters, this may include:

- Hip osteoarthritis (OA)
- Femoroacetabular impingement (FAI)
- Labrum tears

### FEATURES AND BENEFITS

- The Dynamic Rotation Straps with SmartDosing dial allows individual adjustment of the external rotation force, to achieve optimum pain relief for each patient
- Easy donning and doffing thanks to pocket wrap closure system and zipper
- Pulley system for improved stability
- Premium quality Lycra<sup>®</sup> for breathability and a snug fit


### REQUIRED MEASUREMENTS

- Measure the hip circumference at the anterior superior iliac-spine.

### ORDERING INFORMATION

#### UNLOADER HIP

| Size | Hip Circumference | Woman       | Men         |
|------|-------------------|-------------|-------------|
| S    | 80 – 88cm         | B-465350030 | B-465350031 |
| M    | 88 – 96cm         | B-465350032 | B-465350033 |
| L    | 96 – 104cm        | B-465350034 | B-465350035 |
| XL   | 104 – 114cm       | B-465350036 | B-465350037 |
| XXL  | 114 – 124cm       | B-465350038 | B-465350039 |



# UNLOADER® HIP

## ACCESSORIES/ OPTIONS

### REPLACEMENT TROUSERS

| Size | Woman       | Men         |
|------|-------------|-------------|
| S    | B-465350077 | B-465350078 |
| M    | B-465350079 | B-465350080 |
| L    | B-465350081 | B-465350082 |
| XL   | B-465350083 | B-465350084 |
| XXL  | B-465350085 | B-465350086 |


### OTHER ACCESSORIES

| Description        | Part#       |
|--------------------|-------------|
| Rotation strap     | B-465350064 |
| Washbag Kit        | WB000001    |
| Trochanter pad     | B-465350066 |
| Thigh strap S      | B-465350067 |
| Thigh strap M      | B-465350068 |
| Thigh strap L      | B-465350069 |
| Thigh strap XL     | B-465350070 |
| Thigh strap XXL    | B-465350071 |
| Bilateral kit* S   | B-465350072 |
| Bilateral kit* M   | B-465350073 |
| Bilateral kit* L   | B-465350074 |
| Bilateral kit* XL  | B-465350075 |
| Bilateral kit* XXL | B-465350076 |

\*The bilateral kit consists of a rotation control strap, a thigh strap and a trochanter pad

## UPPER EXTREMITY

---



## TABLE OF CONTENTS

---

|                                                 |    |
|-------------------------------------------------|----|
| UPPER EXTREMITY .....                           | 45 |
| FORMFIT® SHOULDER BRACE .....                   | 47 |
| PREMIUM SHOULDER IMMOBILISER .....              | 48 |
| HIGH ARM SLING .....                            | 48 |
| PAEDIATRIC UNIVERSAL ARM SLING .....            | 49 |
| FORMFIT® PRO ELBOW .....                        | 50 |
| REBOUND® POST-OP ELBOW .....                    | 51 |
| FORMFIT® TENNIS ELBOW .....                     | 53 |
| FORMFIT® PRO WRIST .....                        | 54 |
| FORMFIT® WRIST / FORMFIT® WRIST & FOREARM ..... | 55 |
| SOFTPRO® FUNCTIONAL RESTING WHFO .....          | 56 |
| SOFTPRO® GRIP WHFO .....                        | 57 |
| FORMFIT® UNIVERSAL WRIST & THUMB .....          | 58 |
| FORMFIT® THUMB SPICA .....                      | 59 |
| WRIST / THUMB BRACE .....                       | 60 |
| ELASTIC THUMB SPICA .....                       | 60 |
| SPECTRA® WRIST BRACE .....                      | 61 |

## FORMFIT® SHOULDER BRACE

The Formfit Shoulder Brace is designed to provide immobilisation following shoulder injury/surgery. With its streamlined fit, sporty look and breathability it offers excellent comfort and ease of use. Same as its predecessor, the SmartSling®, the Formfit Shoulder Brace also provides abduction of the shoulder (when the abduction pillow is attached), or external rotation of the shoulder (when the external rotation pillow is attached).

Optional axilla pillow helps distribute weight over a larger surface area.

### INDICATIONS FOR USE

- Immobilisation for rotator cuff repairs
- Bankart repairs
- Capsular shifts
- Glenohumeral dislocations/subluxations
- Soft tissue repairs/strains
- Shoulder instabilities/injuries
- Post-surgical rehabilitation

### FEATURES AND BENEFITS

- Breathable material
- Quick-release buckles
- Offloading strap
- Extra padding
- Universal thumb loop

### REQUIRED MEASUREMENTS

- Olecranon to end of closed fist



Formfit Shoulder Brace

Axilla Pillow

## ORDERING INFORMATION

### FORMFIT SHOULDER BRACE

| Size | Length      | Part   |
|------|-------------|--------|
| S    | up to 29cm  | SSS-SO |
| M    | 29 - 40cm   | SSR-SO |
| L    | 40 - 44.5cm | SSL-SO |

### FORMFIT SHOULDER BRACE WITH ABDUCTION

| Size | Length      | Part |
|------|-------------|------|
| S    | up to 29cm  | SSS  |
| M    | 29 - 40cm   | SSR  |
| L    | 40 - 44.5cm | SSL  |

### FORMFIT SHOULDER BRACE WITH EXTERNAL ROTATION

| Size | Length      | Part   |
|------|-------------|--------|
| S    | up to 29cm  | SSS-ER |
| M    | 29 - 40cm   | SSR-ER |
| L    | 40 - 44.5cm | SSL-ER |

| Description   | Part  |
|---------------|-------|
| Axilla Pillow | SSAXP |



## PREMIUM SHOULDER IMMOBILISER

### INDICATIONS FOR USE

For shoulder conditions that may benefit from immobilising the shoulder, such as:

- Pre or post-operative
- Soft tissue repairs/strains
- Shoulder instabilities/injuries

### FEATURES AND BENEFITS

- Effectively controls abduction and rotation
- Humerus cuff, wrist cuff and shoulder strap may be placed in desired position for optimum fit/comfort
- One size for reduced stockholding
- Fits sizes S-XL

### REQUIRED MEASUREMENTS

- Universal, one size (for upper waist diameter up to 106cm)



### ORDERING INFORMATION

| Description                     | Part#  |
|---------------------------------|--------|
| Premium Shoulder Immobiliser    | 0380   |
| Premium Shoulder Immobiliser XL | 0380XL |

## HIGH ARM SLING

### INDICATIONS FOR USE

- Post upper limb trauma
- Post upper limb surgery
- Hand surgery

### FEATURES AND BENEFITS

- Maintains upper arm in raised position
- Breathable nylon/foam laminate construction
- Easy to apply with touch and close fastening
- Waist strap included



### ORDERING INFORMATION

| Description    | Part# |
|----------------|-------|
| High Arm Sling | 0860  |



# PAEDIATRIC UNIVERSAL ARM SLING

## INDICATIONS FOR USE

- Post upper limb trauma
- Post upper limb surgery

## FEATURES AND BENEFITS

- Prong buckle closure correct fit without slippage

## REQUIRED MEASUREMENTS

- Universal, one size



## ORDERING INFORMATION

| Description                                         | Part#  |
|-----------------------------------------------------|--------|
| Paediatric Universal Arm Sling Packed in boxes of 6 | 204020 |

## FORMFIT® PRO ELBOW

Formfit® Pro Elbow provides relief from elbow pain. This 3D compressive elbow sleeve is manufactured using our MotionTech™ technology and features two integrated epicondyle SmartPads™ for more targeted support and cushioning. The RadialStrap™ provides additional compression and support.



MotionTech™



CoolVent™



SmartPad™



FoldLock™

### INDICATIONS FOR USE

- For elbow conditions that may benefit from compression, such as:
- Medial and lateral epicondylitis
- Elbow pain and swelling


### CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

### FEATURES AND BENEFITS

- 3D MotionTech™ for compression and dynamic fit keep
- Epicondyle SmartPads™ for targeted support and cushioning
- Patented FoldLock™ system for easier and more secure application
- RadialStrap™ for increased compression and support

### REQUIRED MEASUREMENTS



Measure circumference  
12cm below elbow.

### ORDERING INFORMATION

| Size | Elbow circumference | Part#    |
|------|---------------------|----------|
| XS   | 17 - 20cm           | FP300001 |
| S    | 20 - 23cm           | FP300002 |
| M    | 23 - 26cm           | FP300003 |
| L    | 26 - 29cm           | FP300004 |
| XL   | 29 - 32cm           | FP300005 |



## REBOUND® POST-OP ELBOW

The Rebound® Post-Op Elbow brace by Össur combines intuitive features for fitters and users with a design that aims to improve the post-op experience. The industry's only elbow immobiliser that combines a drop-lock hinge and four telescoping paddles, individual post-op protocols are easily catered for. The Rebound Post-Op Elbow, lightest of its kind, promotes user comfort with straps that stay in place, anatomical left and right fit with mouldable stability cuffs and blue patient touch-points for easy education.



### INDICATIONS FOR USE

Indications requiring elbow immobilisation and/or range of movement restriction. This may include:

- Stable fractures of the elbow, distal humerus, proximal radius or ulna
- Tendon and ligament injuries or repairs
- Collateral ligament reconstructions
- Elbow hyperextension
- Triceps tendon repair
- Chronic elbow injuries
- Tennis elbow release

### FEATURES AND BENEFITS

- Designed to be the most lightweight brace on the market
- Telescoping paddles relocate to individualise length and keep incision sites clear
- Easy education for users with intuitive patient touchpoints, including Blue is You™ buckles and numbered strapping system
- Drop-Lock hinge enables secure range of movement with bump-resistant pull adjustment stops and optional ROM lock-outs
- Mouldable stability cuffs designed to be individually customised for improved fit and support
- Strap guides maintain strap location for added security and comfort
- Universal fit for both left and right elbows

### REQUIRED MEASUREMENTS

- Universal, one size

# REBOUND® POST-OP ELBOW

## ORDERING INFORMATION

| Side  | Part#     |
|-------|-----------|
| Left  | B-8391680 |
| Right | B-8391580 |

## ACCESSORIES

| Description             | Part#     |
|-------------------------|-----------|
| Elbow Wrist Control Kit | B-8390003 |
| Elbow ROM Locks Kit     | B-8390005 |

## REPLACEMENT PARTS

| Description               | Part#     |
|---------------------------|-----------|
| Elbow Liner Kit           | B-8390002 |
| Elbow Sling Kit           | B-8390004 |
| Strap Tab Kit             | B-8390006 |
| Adhesive Hook Sticker Kit | B-8390007 |
| Anti-Migration Pad Kit    | B-8390008 |
| Elbow Strap Kit           | B-8390009 |

## FORMFIT® TENNIS ELBOW

Formfit Tennis Elbow is made with a patented, ergonomically-shaped plastic. There are two blue cushions with ribbed surface on each side that reduce traction on the tendons. The durable design and materials enable everyday use for both work and athletic activities.

### INDICATIONS FOR USE

- Tendonitis of the elbow
- Prevention during work, daily life and sport
- Epicondylitis

### FEATURES AND BENEFITS

- Adjustable strap
- Patented ergonomically-shaped plastic eliminates tourniquet effect associated with traditional designs.
- Twin blue cushions provide counterforce, reducing traction on the tendons and thus inflammation
- Ribbed surfaces reduce vibration to the tendon, increasing comfort
- Durable design and materials enable everyday use for both work and athletic activities

### REQUIRED MEASUREMENTS

- Circumference of the elbow



## ORDERING INFORMATION

| Size | Circumference | Part#       |
|------|---------------|-------------|
| S    | 22-26cm       | B-247000002 |
| M    | 27-32cm       | B-247000003 |
| L    | 33-38cm       | B-247000004 |

## FORMFIT® PRO WRIST

Formfit Pro Wrist is constructed using our superior signature 3D MotionTech knit, and ideal for injuries that may benefit from compression. Malleable aluminum stays can be adjusted to conform to an individual's anatomy. CoolVent™ technology in the basal thumb area enhances breathability and comfort.



MotionTech™



CoolVent™



StayBilizer™



FoldLock™



### INDICATIONS FOR USE

- For wrist conditions that may benefit from compression, such as:
- Rheumatoid Arthritis (RA) and other chronic pain conditions
- Wrist sprain and strains
- Overuse and pain conditions
- Wrist tendonitis


### CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

### FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit has been proven to be superior to the leading brands at moisture wicking and breathability, keeping patients cool during wear
- Malleable aluminum stay for a customised fit
- CoolVent™ technology in thumb base area for extra breathability and comfort

### REQUIRED MEASUREMENTS



Measure wrist circumference at the level of the radial styloid with the hand outstretched. Wrist in 0 degrees of flexion.

### ORDERING INFORMATION

| Size | Wrist circumference | Left     | Right    |
|------|---------------------|----------|----------|
| XS   | 14 - 15cm           | FP400201 | FP400101 |
| S    | 15 - 16cm           | FP400202 | FP400102 |
| M    | 16 - 17cm           | FP400203 | FP400103 |
| L    | 17 - 18cm           | FP400204 | FP400104 |
| XL   | 18 - 19cm           | FP400205 | FP400105 |
| XXL  | 19 - 20cm           | FP400206 | FP400106 |



## FORMFIT® WRIST / FORMFIT® WRIST & FOREARM

Formfit Wrist, and Formfit Wrist & Forearm, provides immobilisation and protection for patients suffering from tendonitis, select soft tissue injuries – including mild to moderate sprains and strains, as well as for post-cast healing. The lightweight, breathable material and contour fit make this comfortable product an excellent choice.

### INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
- Post-cast healing
- Other selected soft tissue injuries

### FEATURES AND BENEFITS

- Breathable material
- Contoured fit
- Contact closure straps
- Customisable stays

### REQUIRED MEASUREMENTS

- Wrist circumference



## ORDERING INFORMATION

### FORMFIT WRIST & FOREARM

|      |                     | 8"     |        | 10"    |        |
|------|---------------------|--------|--------|--------|--------|
| Size | Wrist Circumference | Left   | Right  | Left   | Right  |
| XS   | 13 - 16cm           | 317082 | 317072 | 617082 | 617072 |
| S    | 16 - 18.5cm         | 317083 | 317073 | 617083 | 617073 |
| M    | 18.5 - 20.5cm       | 317085 | 317075 | 617085 | 617075 |
| L    | 20.5 - 23.5cm       | 317087 | 317077 | 617087 | 617077 |
| XL   | 23.5 - 25.5cm       | 317088 | 317078 | 617088 | 617078 |

## ACCESSORIES

### GEL PAD FOR HOT/COLD THERAPY

| Description                            | Left  | Right |
|----------------------------------------|-------|-------|
| Form Fit Wrist (8" version)*           | GB20L | GB20R |
| Form Fit Wrist & Forearm (10" version) | GB22L | GB22R |



\* Fits XS/S. For M-XL, please use Form Forearm 10" part numbers





## SOFTPRO® FUNCTIONAL RESTING WHFO

The SoftPro® Functional Wrist/Hand/Finger Orthosis (WHFO) is clinically indicated for providing support and protection of a weak or injured wrist, hand, fingers, including thumb. The “bend to fit” endoskeleton frame can be hand shaped to accommodate the desired angle of support for the wrist, hand, fingers, and thumb. The endoskeleton frame can be modified in the sagittal plane as well as medially, laterally, and rotationally to accommodate most angles of the wrist/hand/and fingers without the use of tools. Each endoskeleton frame is closed cell foam. The endoskeleton frame is a closed cell foam, covered by a laminated foam/cloth cover that can be laundered as necessary.



### INDICATIONS FOR USE

- Support and protection of a weak or injured wrist, hand, fingers, including thumb

### FEATURES AND BENEFITS

- Extra long wrist strap holds the wrist in the proper position and applies a gentle stretch
- Maintains/supports the wrist/hand/fingers in a comfortable functional position
- “Bend to Fit” – Easily adjust to correct alignment & support joints without the use of a heat gun

### REQUIRED MEASUREMENTS

- Forearm-fingertip length
- Hand width at metacarpals

## ORDERING INFORMATION

| Part#       | Description               | Width of MP Joints | Length |
|-------------|---------------------------|--------------------|--------|
| SP31FFR-C-L | SoftPro Rest Child Left   | 4.5 - 5cm          | 18cm   |
| SP31FFR-C-R | SoftPro Rest Child Right  | 4.5 - 5cm          | 18cm   |
| SP31FFR-Y-L | SoftPro Rest Youth Left   | 5 - 6.5cm          | 20cm   |
| SP31FFR-Y-R | SoftPro Rest Youth Right  | 5 - 6.5cm          | 20cm   |
| SP31FFR-S-L | SoftPro Rest Small Left   | 6.5 - 7.5cm        | 23cm   |
| SP31FFR-S-R | SoftPro Rest Small Right  | 6.5 - 7.5cm        | 23cm   |
| SP31FFR-M-L | SoftPro Rest Medium Left  | 7.5 - 9cm          | 25cm   |
| SP31FFR-M-R | SoftPro Rest Medium Right | 7.5 - 9cm          | 25cm   |
| SP31FFR-L-L | SoftPro Rest Large Left   | 9 - 10cm           | 28cm   |
| SP31FFR-L-R | SoftPro Rest Large Right  | 9 - 10cm           | 28cm   |

## SOFTPRO® GRIP WHFO

The SoftPro Grip WHFO “bend to fit” endoskeleton frame can be hand shaped to accommodate the desired angle of support for the wrist, hand, fingers, and thumb. The endoskeleton frame can be modified in the sagittal plane as well as medially, laterally, and rotationally to accommodate most angles of the wrist/hand/and fingers without the use of tools. Each endoskeleton frame is closed cell foam. The endoskeleton frame is a closed cell foam, covered by a laminated foam/ cloth cover that can be laundered as necessary.

### INDICATIONS FOR USE

- Provides necessary support and protection of a weak or injured wrist and/or hand and/or fingers, including the thumb

### FEATURES AND BENEFITS

- Two incrementally larger palmar rolls included to gradually increase extension of the MCP and PIP joints
- Maintains/supports the wrist/hand/fingers in a comfortable position
- “Bend to Fit” – Easily adjust to correct alignment & support joints without the use of a heat gun

### REQUIRED MEASUREMENTS

- Forearm-fingertip length
- Hand width at metacarpals



## ORDERING INFORMATION

| Part#      | Description               | Width of MP Joints | Length |
|------------|---------------------------|--------------------|--------|
| SP31FG-C-L | SoftPro Grip Child Left   | 4.5 - 5cm          | 18cm   |
| SP31FG-C-R | SoftPro Grip Child Right  | 4.5 - 5cm          | 18cm   |
| SP31FG-Y-L | SoftPro Grip Youth Left   | 5 - 6.5cm          | 20cm   |
| SP31FG-Y-R | SoftPro Grip Youth Right  | 5 - 6.5cm          | 20cm   |
| SP31FG-S-L | SoftPro Grip Small Left   | 6.5 - 7.5cm        | 23cm   |
| SP31FG-S-R | SoftPro Grip Small Right  | 6.5 - 7.5cm        | 23cm   |
| SP31FG-M-L | SoftPro Grip Medium Left  | 7.5 - 9cm          | 25cm   |
| SP31FG-M-R | SoftPro Grip Medium Right | 7.5 - 9cm          | 25cm   |
| SP31FG-L-L | SoftPro Grip Large Left   | 9 - 10cm           | 28cm   |
| SP31FG-L-R | SoftPro Grip Large Right  | 9 - 10cm           | 28cm   |

## FORMFIT® UNIVERSAL WRIST & THUMB

Formfit Universal Wrist and Thumb products are designed to restrict motion of the wrist and the thumb. Providing immobilisation and protection from select soft tissue injuries, mild to moderate strains and sprains, including gamekeeper's thumb and post cast healing. The lightweight, universal design also reduces inventory needs.

### INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
- Post-cast healing
- Other selected soft tissue injuries

### FEATURES AND BENEFITS

- Breathable material
- Quick-Pull close strap
- Positional dorsal stay
- Customisable thumb stay
- Malleable palmar stay

### REQUIRED MEASUREMENTS

- Wrist circumference



## ORDERING INFORMATION

### FORMFIT WRIST

| Size         | Wrist Circumference | Left        | Right       |
|--------------|---------------------|-------------|-------------|
| Paediatric   | 3.5 - 5"            | B-252603301 | B-252603300 |
| 8" Universal | 4.5 - 9.5"          | B-252603109 | B-252503109 |
| 8" XL        | 7.5 - 11"           | B-252603106 | B-252603105 |

### FORMFIT WRIST & FOREARM UNIVERSAL

| Size          | Wrist Circumference | Left        | Right       |
|---------------|---------------------|-------------|-------------|
| 10" Universal | 4.5 - 9.5"          | B-252603210 | B-252603209 |

### FORMFIT THUMB SPICA UNIVERSAL

| Size         | Wrist Circumference | Left        | Right       |
|--------------|---------------------|-------------|-------------|
| Paediatric   | 3.5 - 5"            | B-253603301 | B-253603300 |
| 8" Universal | 4.5 - 9.5"          | B-253603109 | B-253503109 |
| 8" XL        | 7.5 - 11"           | B-253603106 | B-253603105 |

## FORMFIT® THUMB SPICA

The Formfit Thumb Spica is designed to restrict wrist and thumb motion. It provides immobilisation and protection for patients suffering from tendonitis (e.g. DeQuervain's Tendonitis), select soft tissue injuries – including mild to moderate sprains and strains and Gamekeeper's Thumb –, as well as for post-cast healing. The lightweight, breathable material and contour fit make this comfortable product an excellent choice.

### INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
  - e.g., DeQuervain's Tendonitis
- Post-cast healing
- Other selected soft tissue injuries
  - e.g., Gamekeeper's Thumb

### FEATURES AND BENEFITS

- Breathable material
- Contoured fit
- Contact closure straps
- Customisable stays

### REQUIRED MEASUREMENTS

- Wrist circumference



## ORDERING INFORMATION

| Size | Wrist Circumference | Left | Right |
|------|---------------------|------|-------|
| XS   | 13 - 16cm           | 3020 | 3120  |
| S    | 16 - 18.5cm         | 3030 | 3130  |
| M    | 18.5 - 20.5cm       | 3050 | 3150  |
| L    | 20.5 - 23.5cm       | 3070 | 3170  |
| XL   | 23.5 - 25.5cm       | 3080 | 3180  |

## ACCESSORIES

### GEL PAD FOR HOT/COLD THERAPY

| Description                            | Left  | Right |
|----------------------------------------|-------|-------|
| Form Fit Wrist (8" version)*           | GB20L | GB20R |
| Form Fit Wrist & Forearm (10" version) | GB22L | GB22R |



\* Fits XS/S. For M-XL, please use Form Forearm 10" part numbers

## WRIST / THUMB BRACE

A simple yet effective wrist brace offering superior comfort.

### INDICATIONS FOR USE

- Stabilising hyper-or hypotonic wrist/thumb flexors
- Wrist / Thumb over use injury
- Wrist / Thumb trauma

### FEATURES AND BENEFITS

- Simple touch-and-close fastenings
- Full length malleable aluminium thumb stay

### REQUIRED MEASUREMENTS

- Wrist circumference



## ORDERING INFORMATION

|      |                     | Elastic |        | Neoprene |        |
|------|---------------------|---------|--------|----------|--------|
| Size | Wrist circumference | Left    | Right  | Left     | Right  |
| XS   | 11 - 13cm           | 150/ET  | 155/ET | 160/NT   | 165/NT |
| S    | 13 - 15cm           | 151/ET  | 156/ET | 161/NT   | 166/NT |
| M    | 15 - 17cm           | 152/ET  | 157/ET | 162/NT   | 167/NT |
| L    | 17 - 19cm           | 153/ET  | 158/ET | 163/NT   | 168/NT |
| XL   | 19 - 21cm           | 154/ET  | 159/ET | 164/NT   | 169/NT |

## ELASTIC THUMB SPICA

The Elastic Thumb Spica is designed to protect the thumb from unwanted motion whilst allowing full range of motion of the fingers. Made from elastic, Thumb Spica is simple and very easy to fit. The thumb spica is designed to restrict movement of the thumb post injury, such as a broken thumb, which aids in the recuperation process. It also offers moderate and flexible control following CMC ligament sprains or strains.

### KEY BENEFITS:

Constructed from elastic, simple design with spiral thumb stay, restricts movement of first MCP and CMC joints, touch-and-close fastenings

### INDICATIONS FOR USE

- CMC ligament sprains or strains requiring moderate, flexible control



## ORDERING INFORMATION

| Size | Wrist circumference | Left   | Right  |
|------|---------------------|--------|--------|
| S    | 13 - 17cm           | 191/TS | 196/TS |
| M    | 17 - 21cm           | 192/TS | 197/TS |
| L    | 21 - 25cm           | 193/TS | 198/TS |
| XL   | 25 - 29cm           | 194/TS | 199/TS |

## SPECTRA® WRIST BRACE

The Spectra wrist brace is designed to immobilise the wrist at a neutral angle. Its single-pull lacing closure makes it easier to apply than a conventional wrist brace, meeting the needs of patients with arthritis and other hand dexterity issues.

### INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
- Carpal tunnel and other repetitive stress injuries
- Post cast healing

### FEATURES AND BENEFITS

- Wicking liner moves moisture away from the skin for increased comfort
- Single-pull lace closure enables easy application
- Adjustable webspace strap allows size adjustment for increased comfort
- Removable/adjustable palmar stay provides adjustability of immobilisation angle
- Suede soft good material is durable and breathable

### REQUIRED MEASUREMENTS

- Wrist circumference



## ORDERING INFORMATION

| Size | Wrist circumference | Left    | Right   |
|------|---------------------|---------|---------|
| XS   | 12 - 14cm           | 6602210 | 6602315 |
| S    | 14 - 16cm           | 6602422 | 6602535 |
| M    | 16 - 18cm           | 6602624 | 6602734 |
| L    | 18 - 20cm           | 6602755 | 6602793 |
| XL   | 20 - 22cm           | 6602826 | 6602932 |





KNEE





## TABLE OF CONTENTS

|                                                                  |     |
|------------------------------------------------------------------|-----|
| KNEE .....                                                       | 63  |
| UNLOADER ONE® X .....                                            | 65  |
| UNLOADER ONE® .....                                              | 71  |
| UNLOADER ONE® CUSTOM .....                                       | 74  |
| UNLOADER ONE® PLUS .....                                         | 75  |
| FORMFIT® PRO KNEE OA .....                                       | 77  |
| REBOUND® DUAL .....                                              | 78  |
| REBOUND® DUAL ST .....                                           | 82  |
| REBOUND® DUAL CUSTOM .....                                       | 85  |
| REBOUND® CARTILAGE .....                                         | 87  |
| CTI® OA .....                                                    | 90  |
| REBOUND® PCL .....                                               | 92  |
| REBOUND® ACL .....                                               | 96  |
| REBOUND® DUAL BASIC AND REBOUND® DUAL BASIC ST .....             | 100 |
| OA EASE .....                                                    | 102 |
| CTi®3 .....                                                      | 104 |
| CTI® CUSTOM .....                                                | 106 |
| FLEX™ .....                                                      | 109 |
| C180™ ROCKET .....                                               | 111 |
| REBOUND® POST-OP KNEE .....                                      | 112 |
| FORMFIT POST-OP KNEE .....                                       | 113 |
| REBOUND® KNEE BRACE .....                                        | 114 |
| REBOUND® KNEE IMMOBILIZER UNIVERSAL .....                        | 115 |
| FORMFIT® PRO KNEE .....                                          | 116 |
| FORMFIT® TRACKER .....                                           | 118 |
| FORMFIT® PRO KNEE FLITE .....                                    | 119 |
| FORMFIT® HINGED LATERAL J .....                                  | 120 |
| FORMFIT® KNEE HINGED .....                                       | 121 |
| FORMFIT® KNEE MCL .....                                          | 122 |
| FORMFIT® KNEE ROM BRACE .....                                    | 123 |
| NEOPRENE HINGED KNEE SUPPORT .....                               | 124 |
| NEOPRENE KNEE SUPPORT WITH STABILISED PATELLA .....              | 124 |
| NEOPRENE WRAPAROUND HINGED KNEE SUPPORT .....                    | 125 |
| NEOPRENE KNEE SLEEVE .....                                       | 126 |
| MEASUREMENT FORM CTi® CUSTOM BRACE .....                         | 127 |
| MEASUREMENT FORM UNLOADER ONE® CUSTOM / REBOUND® ACL & PCL ..... | 130 |

## UNLOADER ONE® X

Unloader One X improves mobility by providing pain relief from knee osteoarthritis, degenerative meniscal tears, and other unicompartmental knee conditions that may benefit from unloading.

It features Össur's unique, proprietary, and clinically proven 3-point Leverage System that generates unicompartmental unloading of the knee joint.

Available in custom and OTS versions to help you fit your patient's different leg anatomies.



### INDICATIONS FOR USE


- Moderate to severe unicompartmental knee osteoarthritis
- Degenerative meniscal tears
- Other knee conditions such as:
  - Articular cartilage defect repair
  - Meniscal repair
  - Avascular necrosis
  - Condylar bone marrow lesions
  - Tibial plateau fractures

### FEATURES AND BENEFITS

- Relieves Pain - 3 Point Leverage System with Dynamic Force System® straps, rigid shells, and a single upright hinge biomechanically unloads the affected knee compartment.
- Control - Tool-free SmartDosing® dials enable the patient to fine-tune the strap tension to control pain.
- Tool-Free - Tool-free polycentric hinge with optional ROM stops for range-of-motion restriction.
- Comfortable - Skin-friendly liners feel great against the skin and reduce pressure points. Anti-migration strap technology design captures the widest part of the calf to reduce brace migration.
- Easy-to-Use - Alignment guide on hinge facilitates correct positioning of the brace on the leg and color coded touch points aid in easy donning and doffing.
- Easy-to-Wear - Slender frame with soft edges easily contours to individual anatomies and remains discreet whether worn under or over clothing.
- Easy Care - Entire brace is machine washable (washbag included).
- Field serviceable.


# UNLOADER ONE® X

## REQUIRED MEASUREMENTS: UNLOADER ONE X OTS



Measure calf circumference 15 cm below mid-patella. Leg in extension.

Measure leg length from floor to groin.



| Size   | Calf circumference 15 cm below mid patella |
|--------|--------------------------------------------|
| XS - M | 25 - 43                                    |
| L - XL | 43 - 62                                    |

| Model    | Leg Length      |
|----------|-----------------|
| Standard | > 85 cm (33.5") |
| Short    | ≤ 85 cm (33.5") |

## COLOURS

### OTS COLOUR



Black Brace /  
Black Components

Actual colours may vary.

## ORDERING INFORMATION

Hinge is placed on the side of the affected compartment. Short version recommended for patients with a leg length of ≤ 85 cm.

### OTS


| Standard | Right Medial | Right Lateral | Left Medial | Left Lateral |
|----------|--------------|---------------|-------------|--------------|
| XS - M   | UOX0518713   | UOX0528713    | UOX0618713  | UOX0628713   |
| L - XL   | UOX0518715   | UOX0528715    | UOX0618715  | UOX0628715   |

| Short  | Right Medial | Right Lateral | Left Medial | Left Lateral |
|--------|--------------|---------------|-------------|--------------|
| XS - M | UOX1518713   | UOX1528713    | UOX1618713  | UOX1628713   |
| L - XL | UOX1518715   | UOX1528715    | UOX1618715  | UOX1628715   |


# UNLOADER ONE® X

## REPLACEMENT PARTS



# UNLOADER ONE® X

## ACCESSORIES



## UNLOADER ONE® X

| THIGH LINERS   | THIGH, Right Medial, Left Lateral | THIGH, Right Lateral, Left Medial | THIGH EXTENSION, Right Medial, Left Lateral | THIGH EXTENSION, Right Lateral, Left Medial |
|----------------|-----------------------------------|-----------------------------------|---------------------------------------------|---------------------------------------------|
| Black, Sensil  | UOX526801                         | UOX525801                         | NA                                          | NA                                          |
| Black, Doeskin | UOX526901                         | UOX525901                         | UOX546901                                   | UOX545901                                   |

| CALF LINERS         | Right Medial, Left Lateral | Right Lateral, Left Medial |
|---------------------|----------------------------|----------------------------|
| Black, Sensil XS-M  | UOX516952                  | UOX515952                  |
| Black, Sensil L-XL  | UOX516962                  | UOX515962                  |
| Black, Doeskin XS-M | UOX516951                  | UOX515951                  |
| Black, Doeskin L-XL | UOX516961                  | UOX515961                  |

| THIGH AND LOWER DFS STRAP* | Right Medial, Left Lateral | Right Lateral, Left Medial |
|----------------------------|----------------------------|----------------------------|
| Black                      | UOX406902                  | UOX405902                  |

\*Includes Thigh and DFS Strap with buckle

| CALF AND UPPER DFS STRAP* | Right Medial, Left Lateral | Right Lateral, Left Medial |
|---------------------------|----------------------------|----------------------------|
| Black                     | UOX406901                  | UOX405901                  |

\*Includes Calf and DFS Strap with buckle

| UOX SERVICE KIT* | Right Medial, Left Lateral | Right Lateral, Left Medial |
|------------------|----------------------------|----------------------------|
| Black, XS-M      | UOX416901                  | UOX415901                  |
| Black, L-XL      | UOX416911                  | UOX415911                  |

\* Includes Doeskin Thigh Liner, Sensil Calf Liner, Thigh and Lower DFS Strap with Buckle, Calf and Upper DFS Strap with Buckle, Hinge Pad, Popliteal Pad, and Alligator Kit

| ANTI-MIGRATION PILLAR | Right Medial, Left Lateral | Right Lateral, Left Medial |
|-----------------------|----------------------------|----------------------------|
| Black                 | UOX700015                  | UOX700035                  |

## OTHER REPLACEMENT PARTS AND ACCESSORIES

| Description                  | Black      |
|------------------------------|------------|
| Alligator Kit*               | UOX900001  |
| Hinge Pad                    | UOX700091  |
| Popliteal Pad                | UOX700012  |
| Popliteal Pad - Short        | UOX700022  |
| Sensil Strap Pad             | UON0181    |
| Sensil Strap Pad - 10 Pieces | UOX700113  |
| Smartdosing Dial             | UOX1000001 |
| Strap Wrap Set               | UOX700014  |

\* Includes 4 Alligator Clips: 2 DFS, 1 Thigh and 1 Calf

| Description                                      | Part #     |
|--------------------------------------------------|------------|
| Ease of Use Kit*                                 | UOX800001  |
| OA Hinge Extension Stop Tree (5°, 10°, 15°, 30°) | OA40029    |
| OA Hinge Flexion Stop Tree (0°, 30°, 60°, 90°)   | OA40039    |
| Sundial                                          | UOX1200001 |
| Wash Bag                                         | UON60200   |

\* Includes Sundial and 2v buckle finger-loops



# UNLOADER ONE® X

## OTHER REPLACEMENT PARTS AND ACCESSORIES

### SLEEVE WITH MESH\*

| Size | Circumference 15 cm above mid patella | Part #   |
|------|---------------------------------------|----------|
| XS   | < 43 cm                               | SLV60011 |
| S    | 43 – 46 cm                            | SLV60012 |
| M    | 46 – 51 cm                            | SLV60013 |
| L    | 51 – 53 cm                            | SLV60014 |
| XL   | 53 – 56 cm                            | SLV60015 |
| XXL  | 56 – 61 cm                            | SLV60016 |

\* Sleeve with Mesh: Main fabric: 92% Cotton/8% Spandex - Mesh: 100% Polyester - Band: 80% Polyester/20% Spandex



## REQUIRED MEASUREMENTS: UNLOADER ONE X CUSTOM

Custom made to measure: Smartmeasure, 3D Scan (.AOP, .STL, .OBJ, .CPV) or Patient Cast

### SMARTMEASURE®


A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



| Description            | Part#       |
|------------------------|-------------|
| Constant Force Caliper | B-312000000 |

## COLOURS

### CUSTOM STANDARD COLOURS



Actual colours may vary.

## ORDERING INFORMATION

Hinge is placed on the side of the affected compartment. Short version recommended for patients with a leg length of ≤ 85 cm.

### CUSTOM

| Description           | Right | Left  |
|-----------------------|-------|-------|
| Unloader One X Custom | UOX0R | UOX0L |

## UNLOADER ONE®

Unloader One® is designed to relieve osteoarthritis knee pain and wearing it may help improve mobility. It features Össur's unique, patented 3-point leverage system, which combines a single hinge with dual Dynamic Force System (DFS) straps to provide unloading of the knee joint. Clinical studies demonstrate that this system may relieve osteoarthritis knee pain and improve function and quality of life while being worn. Colour-coded Quick Fit buckles make donning and doffing simple, while flexible shells and breathable liners are designed to improve comfort and reduce migration, thus helping to ensure compliance.



### INDICATIONS FOR USE

- Moderate to severe unicompartmental knee osteoarthritis
- Degenerative meniscal tears
- Other knee conditions that may benefit from unicompartmental unloading, such as Articular cartilage defect repair, Avascular necrosis, Tibial plateau fracture, Bone marrow lesions (bone bruises)


### FEATURES AND BENEFITS

- Dual Dynamic Force Straps® designed to reduce pressure on sensitive bony surfaces of the knee, resulting in reduction of pain
- Lightweight frame, perforated for flexibility and ventilation, with flexible, low profile upright that contours to patient anatomy and allows use underneath clothing
- Tool-free SmartDosing® dial enables fine-tuned pain management and can be adjusted even when wearing the brace under clothing
- Colour-coded quick fit buckles enable easy application and removal
- Proprietary Sensil® silicone liners designed to reduce migration
- Available in standard and short version, for a perfect fit on shorter legs, too

### REQUIRED MEASUREMENTS

- Circumference measurement: 15cm below mid-patella

| Size | Measurements  |
|------|---------------|
| S    | 31 - 35.5cm   |
| M    | 35.5 - 42.5cm |
| L    | 42.5 - 51cm   |
| XL   | 51 - 61.5cm   |



# UNLOADER ONE®

## ORDERING INFORMATION

### UNLOADER ONE WITH SMARTDOSING

| Size | Calf circumference | Left medial | Right medial | Left lateral | Right lateral |
|------|--------------------|-------------|--------------|--------------|---------------|
| S    | 31 - 35.5cm        | B-240618712 | B-240518712  | B-240628712  | B-240528712   |
| M    | 35.5 - 42.5cm      | B-240618713 | B-240518713  | B-240628713  | B-240528713   |
| L    | 42.5 - 51cm        | B-240618714 | B-240518714  | B-240628714  | B-240528714   |
| XL   | 51 - 61.5cm        | B-240618715 | B-240518715  | B-240628715  | B-240528715   |

### UNLOADER ONE SHORT WITH SMARTDOSING

| Size | Calf circumference | Left medial | Right medial | Left lateral | Right lateral |
|------|--------------------|-------------|--------------|--------------|---------------|
| XS   | 25 - 31cm          | B-241618711 | B-241518711  | B-241628711  | B-241528711   |
| S    | 31 - 35.5cm        | B-241618712 | B-241518712  | B-241628712  | B-241528712   |
| M    | 35.5 - 42.5cm      | B-241618713 | B-241518713  | B-241628713  | B-241528713   |
| L    | 42.5 - 51cm        | B-241618714 | B-241518714  | B-241628714  | B-241528714   |
| XL   | 51 - 61.5cm        | B-241618715 | B-241518715  | B-241628715  | B-241528715   |

### UNLOADER ONE WITH RATCHET

| Size | Calf circumference | Left medial | Right medial | Left lateral | Right lateral |
|------|--------------------|-------------|--------------|--------------|---------------|
| S    | 31 - 35.5cm        | B-240619712 | B-240519712  | B-240629712  | B-240529712   |
| M    | 35.5 - 42.5cm      | B-240619713 | B-240519713  | B-240629713  | B-240529713   |
| L    | 42.5 - 51cm        | B-240619714 | B-240519714  | B-240629714  | B-240529714   |
| XL   | 51 - 61.5cm        | B-240619715 | B-240519715  | B-240629715  | B-240529715   |

### UNLOADER ONE SHORT WITH RATCHET

| Size | Calf circumference | Left medial | Right medial | Left lateral | Right lateral |
|------|--------------------|-------------|--------------|--------------|---------------|
| XS   | 25 - 31cm          | B-241619711 | B-241519711  | B-241629711  | B-241529711   |
| S    | 31 - 35.5cm        | B-241619712 | B-241519712  | B-241629712  | B-241529712   |
| M    | 35.5 - 42.5cm      | B-241619713 | B-241519713  | B-241629713  | B-241529713   |
| L    | 42.5 - 51cm        | B-241619714 | B-241519714  | B-241629714  | B-241529714   |
| XL   | 51 - 61.5cm        | B-241619715 | B-241519715  | B-241629715  | B-241529715   |

Standard version is recommended for patients above 1.65m  
Short version is recommended for patients under 1.65m tall.  
Hinge is placed on the side of the affected compartment

## ACCESSORIES

### ACCESSORY KITS

|                                      | SMARTDOSING |             | RATCHET     |             |
|--------------------------------------|-------------|-------------|-------------|-------------|
| Description                          | RM/LL       | RL/LM       | RM/LL       | RL/LM       |
| Strap & Strap Pad & Buckle Clasp Kit | B-705277003 | B-705278003 | B-701217008 | B-701218008 |
| Strap & Strap Pad Kit                | B-705270003 |             | B-701210006 |             |
| Strap Kit                            | B-701270013 |             | B-701210004 |             |
| Strap Kit Short                      | B-701240012 |             |             |             |
| Pad Kit                              | B-702210003 |             | B-702210003 |             |

### REPLACEMENT PARTS AND ACCESSORIES

| Description                            | Part#    |
|----------------------------------------|----------|
| Extension Stop Kit (5°, 10°, 15°, 30°) | OAH40029 |
| Flexion Stop Kit (0°, 30°, 60°, 90°)   | OAH40039 |

# UNLOADER ONE®

## ACCESSORIES

### ACCESSORY KITS

| Description                                     | Right medial/Left lateral | Right lateral/Left medial |
|-------------------------------------------------|---------------------------|---------------------------|
| Buckle Clasp Kit                                | B-708217003               | B-708218003               |
| Alligator Clip Kit (included in all strap kits) | B-704210003               |                           |
| Extension Stop Kit                              | B-301040046               |                           |
| Twist-Fit Hook & Loop Attachment Kit            | B-721240003               |                           |

### UNDERSLEEVE

| Size | Thigh Circumference | Part#       |
|------|---------------------|-------------|
| XS   | 25.4 - 30.9cm       | B-250150201 |
| S    | 30.9 - 35.6cm       | B-250150202 |
| M    | 35.6 - 42.5cm       | B-250150203 |
| L    | 42.5 - 50.8cm       | B-250150204 |
| XL   | 50.8 - 61cm         | B-250150205 |
| XXL  | 61 - 81.3cm         | B-250150206 |

When using an undersleeve we advice to use the suspension strap kit to minimise migration

### SENSIL LINERS

|      | Tibia       |             | Femur       |             |
|------|-------------|-------------|-------------|-------------|
| Size | RM/LL*      | RL/LM**     | RM/LL*      | RL/LM**     |
| XS   | B-419700311 | B-419800311 | B-409700321 | B-409800321 |
| S    | B-419700312 | B-419800312 | B-409700322 | B-409800322 |
| M    | B-419700313 | B-419800313 | B-409700323 | B-409800323 |
| L    | B-419700314 | B-419800314 | B-409700324 | B-409800324 |
| XL   | B-419700315 | B-419800315 | B-409700325 | B-409800325 |

\* Right medial/Left lateral / \*\* Right lateral/Left medial

### DOESKIN LINERS

|      | Tibia       |             | Femur       |             |
|------|-------------|-------------|-------------|-------------|
| Size | RM/LL*      | RL/LM**     | RM/LL*      | RL/LM**     |
| XS   | B-424703721 | B-424803721 | B-412701721 | B-412801721 |
| S    | B-424703722 | B-424803722 | B-412701722 | B-412801722 |
| M    | B-424703723 | B-424803723 | B-412701723 | B-412801723 |
| L    | B-424703724 | B-424803724 | B-412701724 | B-412801724 |
| XL   | B-424703725 | B-424803725 | B-412701725 | B-412801725 |

When using a doeskin tibia liner we advice to use a suspension strap kit to minimise migration / \* Right medial/Left lateral / \*\* Right lateral/Left medial

### JERSEY LINERS

|      |       | CALF                      |                           | THIGH                     |                           |
|------|-------|---------------------------|---------------------------|---------------------------|---------------------------|
| Size | Color | Right medial/Left lateral | Right lateral/Left medial | Right medial/Left lateral | Right lateral/Left medial |
| XS   | Black | UNL015111                 | UNL016111                 | UNL015121                 | UNL016121                 |
| S    | Black | UNL015112                 | UNL016112                 | UNL015122                 | UNL016122                 |
| M    | Black | UNL015113                 | UNL016113                 | UNL015123                 | UNL016123                 |
| L    | Black | UNL015114                 | UNL016114                 | UNL015124                 | UNL016124                 |
| XL   | Black | UNL015115                 | UNL016115                 | UNL015125                 | UNL016125                 |

### SENSIL AMS (ANTI-MIGRATION SYSTEM) WRAPS

|      | Black                     |                           | Grey                      |                           |
|------|---------------------------|---------------------------|---------------------------|---------------------------|
| Size | Right medial/Left lateral | Right lateral/Left medial | Right medial/Left lateral | Right lateral/Left medial |
| XS   | UNL035111                 | UNL036111                 | B-715217031               | B-715218031               |
| S    | UNL035112                 | UNL036112                 | B-715217032               | B-715218032               |
| M    | UNL035113                 | UNL036113                 | B-715217033               | B-715218033               |
| L    | UNL035114                 | UNL036114                 | B-715217034               | B-715218034               |
| XL   | UNL035115                 | UNL036115                 | B-715217035               | B-715218035               |

### SUSPENSION STRAP KIT

| Description          | Part#    |
|----------------------|----------|
| Suspension Strap Kit | UNL00212 |

# UNLOADER ONE® CUSTOM

## ORDERING INFORMATION

### UNLOADER ONE CUSTOM WITH SMARTDOSING

| Colour | Thigh circumference | Calf circumference | Left        | Right       |
|--------|---------------------|--------------------|-------------|-------------|
| Black  | 34.3 - 73.5cm       | 25 - 61.5cm        | B-115600112 | B-115500112 |
| Grey   | 34.3 - 73.5cm       | 25 - 61.5cm        | B-115601712 | B-115501712 |

### UNLOADER ONE CUSTOM

| Colour | Thigh circumference | Calf circumference | Left        | Right       |
|--------|---------------------|--------------------|-------------|-------------|
| Black  | 34.3 - 73.5cm       | 25 - 61.5cm        | B-115600111 | B-115500111 |
| Grey   | 34.3 - 73.5cm       | 25 - 61.5cm        | B-115601711 | B-115501711 |

\*Please allow 10 working days for delivery



### SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace, including the Unloader One. Utilises a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



### DIGITAL MEASUREMENT SYSTEM (DMS)

This system requires minimal physical measurements and two digital photos, which can then be sent to customer service. This system is a quick, easy and reliable method of measuring for an Unloader One® Custom without the need for training.

\*For custom accessories, please contact customer care on 03450 065 065.

### CUSTOM STANDARD COLOURS

|            |            |           |
|------------|------------|-----------|
| Red        | Pink       | Silver    |
| Orange     | Black      | White     |
| Yellow     | Navy Blue  | Champagne |
| Lime Green | Ocean Blue | Charcoal  |

## UNLOADER ONE® PLUS

Many sufferers of OA knee pain are overweight. For these patients, containing excess soft tissue is essential for comfortable and effective bracing. That's why our new Unloader One Plus model features a wider thigh strap than the standard version, plus a gastroc strap that helps secure the brace in place. And because bigger legs require more support, we included a sturdier upright with an Adjustable Dynamic Joint, enabling additional leverage for optimal pain management. Multiple, peer-reviewed clinical and biomechanical studies have proven the efficacy of this technology.



### INDICATIONS FOR USE

- Moderate to severe unicompartmental knee osteoarthritis
- Degenerative meniscal tears
- Other knee conditions that may benefit from unicompartmental unloading such as: Articular cartilage defect repair, Avascular necrosis, Tibial plateau fracture, Bone marrow lesions (bone bruises)

### FEATURES AND BENEFITS

- 3-Point Leverage System provides clinically relevant pain relief of affected compartment
- Single SmartDosing dial allows fine-tuned pain management
- Helps keep patients mobile during daily activities
- Discreet hinge and sleeve design ensure a snug fit underneath clothing
- Lightweight and easy to use

### REQUIRED MEASUREMENTS

- Thigh circumferences of 53.3cm to 73.7cm and calf circumferences of 38.1cm to 50.8cm

# UNLOADER ONE® PLUS

## ORDERING INFORMATION

### UNLOADER ONE PLUS WITH RATCHET

| Size | Thigh Circumference | Calf circumference | Left medial | Right medial | Left lateral | Right lateral |
|------|---------------------|--------------------|-------------|--------------|--------------|---------------|
| L    | 53 - 64cm           | 36 - 43cm          | B-245610004 | B-245510004  | B-245620004  | B-245520004   |
| XL   | 64 - 81cm           | 43 - 51cm          | B-245610005 | B-245510005  | B-245620005  | B-245520005   |

Standard version is recommended for patients between 1,65 and 1,92 metres tall.

## ACCESSORIES

### SENSIL AMS (ANTI-MIGRATION SYSTEM) WRAPS

| Size | Right medial/Left lateral | Right lateral/Left medial |
|------|---------------------------|---------------------------|
| M    | B-715237033               | B-715238033               |
| L    | B-715237034               | B-715238034               |
| XL   | B-715237035               | B-715238035               |

### LINERS

|      | Tibia, silicone           |                           | Femur, doeskin            |                           |
|------|---------------------------|---------------------------|---------------------------|---------------------------|
| Size | Right medial/Left lateral | Right lateral/Left medial | Right medial/Left lateral | Right lateral/Left medial |
| M    | B-413700313               | B-413800313               | B-414700323               | B-414800323               |
| L    | B-413700314               | B-413800314               | B-414700324               | B-414800324               |
| XL   | B-413700314               | B-413800315               | B-414700325               | B-414800325               |
| XXL  |                           |                           | B-414700326               | B-414800326               |

### ACCESSORY KITS

| Description                                     | Right medial/Left lateral | Right lateral/Left medial |
|-------------------------------------------------|---------------------------|---------------------------|
| Strap & Strap Pad & Buckle Clasp Kit            | B-701237008               | B-701238008               |
| Strap & Strap Pad Kit                           | B-701230006               |                           |
| Strap Kit                                       | B-701230004               |                           |
| Pad Kit                                         | B-702230003               |                           |
| Extra Pads                                      | B-702230013               |                           |
| Buckle Clasp Kit                                | B-708237003               | B-708238003               |
| Alligator Clip Kit (included in all strap kits) | B-704210003               |                           |
| Gastroc Kit                                     | B-705230001               |                           |
| Elastic Thigh Strap 8"                          | B-701237083               | B-701238083               |
| Elastic Thigh Strap 10"                         | B-701237103               | B-701238103               |



## FORMFIT® PRO KNEE OA

Össur Formfit® Pro Knee OA is a unique, compressive knee sleeve with a Dynamic Force Strap system, designed to provide pain relief from the early stages of osteoarthritis (OA). The dynamic, 3D MotionTech™ knit is supportive and meets certified medical grade compression requirements, and has been proven to be superior to the leading brands at moisture wicking and breathability, keeping your knee cool during wear. CoolVent™ technology in the popliteal area enhances breathability and comfort, while the distinct patella graphics help with intuitive placing of the sleeve. The integrated StayBilizer™ hinge, constructed from lightweight, flexible EVA material, is suitable for sports, such as football.



MotionTech™



CoolVent™



StayBilizer™



FoldLock™



### INDICATIONS FOR USE

For osteoarthritis and other knee conditions which may benefit from unloading and/or compression, such as:

- Medial or lateral knee pain
- Early or mild unicompartmental osteoarthritis of the knee
- Symptomatic, unicompartmental meniscus degeneration

### CONTRAINDICATIONS


- Must not be used by individuals for whom compression is contraindicated.

### FEATURES AND BENEFITS

- Patented DFS system.
- Certified medical grade compression (20-30 mmHg).
- Lightweight CoolVent™ knit in popliteal area for added breathability and comfort.
- Patented FoldLock™ system for easier and more secure application.
- Sports-compatible StayBilizer™ hinge made from light, flexible EVA material for optimum shape and fit.
- Patella graphics for intuitive placement.
- Available in 5 standard sizes and 3 comfort sizes.

### REQUIRED MEASUREMENTS

| Size       | Calf circumference<br>15 cm below mid patella | Thigh circumference<br>15 cm above mid patella | Right Medial,<br>Left Lateral | Right Lateral,<br>Left Medial |
|------------|-----------------------------------------------|------------------------------------------------|-------------------------------|-------------------------------|
| XS         | 26 - 31 cm                                    | 36 - 41 cm                                     | FP110401                      | FP110301                      |
| S          | 31 - 36 cm                                    | 41 - 46 cm                                     | FP110402                      | FP110302                      |
| M          | 36 - 41 cm                                    | 46 - 51 cm                                     | FP110403                      | FP110303                      |
| M Comfort  | 36 - 41 cm                                    | 51 - 56 cm                                     | FP110413                      | FP110313                      |
| L          | 41 - 46 cm                                    | 51 - 56 cm                                     | FP110404                      | FP110304                      |
| L Comfort  | 41 - 46 cm                                    | 56 - 60 cm                                     | FP110414                      | FP110314                      |
| XL         | 46 - 50 cm                                    | 56 - 60 cm                                     | FP110405                      | FP110305                      |
| XL Comfort | 46 - 50 cm                                    | 60 - 64 cm                                     | FP110415                      | FP110315                      |



Measure calf circumference  
15 cm below mid-patella.  
Leg in extension.



Measure thigh circumference  
15 cm above mid-patella.  
Leg in extension.

## REBOUND® DUAL

The Rebound® DUAL offers simplicity without compromise. With its malleable yet supportive frame, Rebound DUAL is the go-to for knee ligament support.



### INDICATIONS FOR USE

The brace is indicated for combined and non-combined instabilities of the:

- Anterior cruciate ligament (ACL)
- Medial collateral ligament (MCL)
- Lateral collateral ligament (LCL)
- Posterior cruciate ligament (PCL)
- Acculign version only: Conditions requiring unicompartmental load control.

### FEATURES AND BENEFITS

- Cold-malleable aluminum frame, trimmable straps, and field-serviceable D-rings accommodate various patient anatomies
- Acculign™ upper frame allows height and varus/valgus adjustment
- Tool-free polycentric hinge with ROM stop options for range-of-motion restriction
- Össur Hinge Alignment Guide, Blue is You™ touch points, and tactile numbered strap tabs guide the fitting process
- Optional SmartDosing® Kit accessory offers fine-tuned unicompartmental offloading
- Optional PCL kit accessory for multi-ligament support

### REQUIRED MEASUREMENTS

- Caliper measurement: Medio-laterally at joint width in a standard, weight-bearing position
- Circumference: 15 cm above mid-patella

# REBOUND® DUAL

## ORDERING INFORMATION

### REBOUND DUAL

| Size | Caliper        | Circumference  | Right       | Left        |
|------|----------------|----------------|-------------|-------------|
| XS   | 7.6 - 8.9 cm   | 38.1 - 41.9 cm | B-242510001 | B-242610001 |
| S    | 8.9 - 10.2 cm  | 41.9 - 47 cm   | B-242510002 | B-242610002 |
| M    | 10.2 - 11.4 cm | 47 - 52.1 cm   | B-242510003 | B-242610003 |
| L    | 11.4 - 12.7cm  | 52.1 - 57.2 cm | B-242510004 | B-242610004 |
| XL   | 12.7 - 14cm    | 57.2 - 62.2 cm | B-242510005 | B-242610005 |
| XXL  | 14 - 15.2cm    | 62.2 - 69.9 cm | B-242510006 | B-242610006 |
| 3XL  | 15.2 - 16.5 cm | 69.9 - 76 cm   | B-242510007 | B-242610007 |

M/L at joint width in a standard, weight-bearing position 2 Thigh circumference: 6" above mid-patella

### REBOUND DUAL CUSTOM

| Measuring System                                                 | Right       | Left        |
|------------------------------------------------------------------|-------------|-------------|
| SmartMeasure, DMS, Coordinate Cast System, Computer Aided Design | B-155510001 | B-155610001 |

## SMARTMEASURE®


A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



| Description            | Part#       |
|------------------------|-------------|
| Constant Force Caliper | B-312000000 |

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

## CUSTOM BRACE COLOURS



# REBOUND® DUAL

## ACCESSORIES

### SERVICE KITS & LINER KITS

| Size | Softgoods Kit, Black* | Liner Kit, Black** |
|------|-----------------------|--------------------|
| XS   | B-701090091           | B-701090110        |
| S    | B-701090092           | B-701090111        |
| M    | B-701090093           | B-701090112        |
| L    | B-701090094           | B-701090113        |
| XL   | B-701090095           | B-701090114        |
| XXL  | B-701090096           | B-701090115        |
| 3XL  | B-701090097           | B-701090116        |

\*Includes Doeskin thigh liner, Active Grip calf liner, straps, gastroc pad, tibial pad, numbered alligator tabs, hook strip

\*\*Includes Doeskin thigh liner and Sensil calf liner

### CALF LINERS

| Size | Active Grip, Black |
|------|--------------------|
| XS   | B-701090160        |
| S    | B-701090161        |
| M    | B-701090162        |
| L    | B-701090163        |
| XL   | B-701090164        |
| XXL  | B-701090165        |
| 3XL  | B-701090166        |

### THIGH LINERS

| Size | Doeskin, Black |
|------|----------------|
| XS   | B-701090150    |
| S    | B-701090151    |
| M    | B-701090152    |
| L    | B-701090153    |
| XL   | B-701090154    |
| XXL  | B-701090155    |
| 3XL  | B-701090156    |

### STRAP KITS & PAD KITS

| Size   | Replacement Strap Kit | PCL Strap Kit | Active Grip Strap Pad Kit |
|--------|-----------------------|---------------|---------------------------|
| XS-L   | B-701090105           | B-701090057   | B-701090039               |
| XL-3XL | B-701090106           | B-701090058   | B-701090040               |

# REBOUND® DUAL

## ACCESSORIES

### OTHER

| Description                                  | Part#       |
|----------------------------------------------|-------------|
| SmartDosing Kit RM/LL                        | B-701090060 |
| SmartDosing Kit RL/LM                        | B-701090061 |
| Flexion/Extension Stop Set - Tool-Free Hinge | B-701090088 |
| *Flexion Stop Set - Custom & Legacy          | B-701090004 |
| **Extension Stop Set - Custom & Legacy       | B-701090008 |
| D-Ring Kit                                   | B-701090059 |
| Tibial Wedge Kit                             | B-701090001 |
| Condyle Pad Kit                              | B-701090009 |

\*Flexion Stops: 0,10,20,30,45,60,75,90

\*\*Extension Stops: 0,5,10,20,30,45

### SLEEVES

| Size  | Circumference<br>15 cm above<br>mid patella | Neoprene<br>undersleeve | Neoprene<br>oversleeve |
|-------|---------------------------------------------|-------------------------|------------------------|
| S     | < 38 cm                                     | 20251                   | 20051                  |
| M     | 38 – 45 cm                                  | 20252                   | 20052                  |
| L     | 45 – 49.5 cm                                | 20253                   | 20053                  |
| XL    | 49.5 – 54.5 cm                              | 20254                   | 20054                  |
| XXL   | 54.5 – 62 cm                                | 20255                   | 20055                  |
| XXXL  | 62 – 71 cm                                  | 20256                   | 20056                  |
| XXXXL | 71 – 86 cm                                  | 20257                   | 20057                  |

### SLEEVE WITH MESH\*

| Size | Circumference<br>15 cm above<br>mid patella | Sleeve with<br>Mesh |
|------|---------------------------------------------|---------------------|
| XS   | < 43 cm                                     | SLV60011            |
| S    | 43 – 46 cm                                  | SLV60012            |
| M    | 46 – 51 cm                                  | SLV60013            |
| L    | 51 – 53 cm                                  | SLV60014            |
| XL   | 53 – 56 cm                                  | SLV60015            |
| XXL  | 56 – 61 cm                                  | SLV60016            |

\* Sleeve with Mesh: Main fabric: 92% Cotton/8% Spandex - Mesh: 100% Polyester - Band: 80% Polyester/20% Spandex



## REBOUND® DUAL ST

The Rebound® DUAL ST offers simplicity without compromise. With its malleable yet supportive frame, Rebound DUAL ST is the go-to for knee ligament support.



### INDICATIONS FOR USE

The brace is indicated for combined and non-combined instabilities of the:

- Anterior cruciate ligament (ACL)
- Medial collateral ligament (MCL)
- Lateral collateral ligament (LCL)
- Posterior cruciate ligament (PCL)

### FEATURES AND BENEFITS

- Cold-malleable aluminum frame, trimmable straps, and field-serviceable D-rings accommodate various patient anatomies
- Step-through design moves the rigid frame cuff behind the calf
- Tool-free polycentric hinge with ROM stop options for range-of-motion restriction
- Össur Hinge Alignment Guide, Blue is You™ touch points, and tactile numbered strap tabs guide the fitting process
- Optional PCL Kit accessory for multi-ligament support

### REQUIRED MEASUREMENTS



Measure knee caliper medio-laterally at joint width in standing, weight-bearing position.



Measure thigh circumference 15 cm above mid-patella. Leg in extension.

| Size | Caliper        | Circumference  | Right       | Left        |
|------|----------------|----------------|-------------|-------------|
| XS   | 7.6 - 8.9 cm   | 38.1 - 41.9 cm | B-242520001 | B-242620001 |
| S    | 8.9 - 10.2 cm  | 41.9 - 47 cm   | B-242520002 | B-242620002 |
| M    | 10.2 - 11.4 cm | 47 - 52.1 cm   | B-242520003 | B-242620003 |
| L    | 11.4 - 12.7 cm | 52.1 - 57.2 cm | B-242520004 | B-242620004 |
| XL   | 12.7 - 14 cm   | 57.2 - 62.2 cm | B-242520005 | B-242620005 |
| XXL  | 14 - 15.2 cm   | 62.2 - 69.9 cm | B-242520006 | B-242620006 |
| 3XL  | 15.2 - 16.5 cm | 69.9 - 76 cm   | B-242520007 | B-242620007 |

## REBOUND® DUAL ST

## ORDERING INFORMATION

## REBOUND DUAL ST CUSTOM

| Measuring    | Right       | Left        |
|--------------|-------------|-------------|
| SmartMeasure | B-155520001 | B-155620002 |

## SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilises a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



| Description            | Part#       |
|------------------------|-------------|
| Constant Force Caliper | B-312000000 |

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

## ACCESSORIES

## LINER KITS

| Size | Liner Kit, Black* |
|------|-------------------|
| XS   | B-701090120       |
| S    | B-701090121       |
| M    | B-701090122       |
| L    | B-701090123       |
| XL   | B-701090124       |
| XXL  | B-701090125       |
| 3XL  | B-701090126       |

\*Includes Doeskin thigh liner and Sensil calf liner

## THIGH LINERS

| Size | Doeskin, Black | Doeskin, Silver |
|------|----------------|-----------------|
| XS   | B-701090170    | B-701090067     |
| S    | B-701090171    | B-701090068     |
| M    | B-701090172    | B-701090069     |
| L    | B-701090173    | B-701090070     |
| XL   | B-701090174    | B-701090071     |
| XXL  | B-701090175    | B-701090072     |
| 3XL  | B-701090176    | B-701090073     |

## CALF LINERS

| Size | Active Grip, Black | Active Grip, Silver |
|------|--------------------|---------------------|
| XS   | B-701090180        | B-701090074         |
| S    | B-701090181        | B-701090075         |
| M    | B-701090182        | B-701090076         |
| L    | B-701090183        | B-701090077         |
| XL   | B-701090184        | B-701090078         |
| XXL  | B-701090185        | B-701090079         |
| 3XL  | B-701090186        | B-701090080         |



# REBOUND® DUAL ST

## ACCESSORIES

### STRAP KITS & PAD KITS

| Size   | Replacement Strap Kit | PCL Strap Kit | Active Grip Strap Pad Kit |
|--------|-----------------------|---------------|---------------------------|
| XS-L   | B-701090107           | B-701090065   | B-701090039               |
| XL-3XL | B-701090108           | B-701090066   | B-701090040               |

### OTHER

| Description                                  | Part#       |
|----------------------------------------------|-------------|
| Flexion/Extension Stop Set - Tool-Free Hinge | B-701090088 |
| *Flexion Stop Set - Custom & Legacy          | B-701090004 |
| **Extension Stop Set - Custom & Legacy       | B-701090008 |
| D-Ring Kit                                   | B-701090059 |
| Condyle Pad Kit                              | B-701090009 |

\*Flexion Stops: 0,10,20,30,45,60,75,90

\*\*Extension Stops: 0,5,10,20,30,45

### SLEEVES

| Size  | Circumference<br>15 cm above<br>mid patella | Neoprene<br>undersleeve | Neoprene<br>oversleeve |
|-------|---------------------------------------------|-------------------------|------------------------|
| S     | < 38 cm                                     | 20251                   | 20051                  |
| M     | 38 – 45 cm                                  | 20252                   | 20052                  |
| L     | 45 – 49.5 cm                                | 20253                   | 20053                  |
| XL    | 49.5 – 54.5 cm                              | 20254                   | 20054                  |
| XXL   | 54.5 – 62 cm                                | 20255                   | 20055                  |
| XXXL  | 62 – 71 cm                                  | 20256                   | 20056                  |
| XXXXL | 71 – 86 cm                                  | 20257                   | 20057                  |

### SLEEVE WITH MESH\*

| Size | Circumference<br>15 cm above<br>mid patella | Sleeve with<br>Mesh |
|------|---------------------------------------------|---------------------|
| XS   | < 43 cm                                     | SLV60011            |
| S    | 43 – 46 cm                                  | SLV60012            |
| M    | 46 – 51 cm                                  | SLV60013            |
| L    | 51 – 53 cm                                  | SLV60014            |
| XL   | 53 – 56 cm                                  | SLV60015            |
| XXL  | 56 – 61 cm                                  | SLV60016            |

\* Sleeve with Mesh: Main fabric: 92% Cotton/8%  
Spandex - Mesh: 100% Polyester - Band: 80%  
Polyester/20% Spandex



## REBOUND® DUAL CUSTOM

The Rebound® DUAL family is a series of double upright aluminum frame knee ligament braces. Available in step-through, anterior calf frame, or a custom upper/lower frame combination, Rebound DUAL offers a solution for you.



### INDICATIONS FOR USE

The brace is indicated for combined and non-combined instabilities of the:

- Anterior cruciate ligament (ACL)
- Medial collateral ligament (MCL)
- Lateral collateral ligament (LCL)
- Posterior cruciate ligament (PCL)
- Acculign version only: Conditions requiring unicompartmental load control

### FEATURES AND BENEFITS

- Cold-malleable aluminum frame, trimmable straps, and field-serviceable D-rings accommodate various patient anatomies
- Acculign™ upper frame allows height and varus/valgus adjustment
- Polycentric hinge with ROM stop options for range-of-motion restriction
- Össur Hinge Alignment Guide, Blue is You™ touch points, and tactile numbered strap tabs guide the fitting process
- Custom-made for your patient
- Optional SmartDosing® Kit accessory offers fine-tuned unicompartmental offloading
- Optional PCL Kit accessory for multi-ligament support

### REQUIRED MEASUREMENTS

- Thigh circumference 20 cm above mid patella
- Thigh circumference 15 cm above mid patella
- Thigh circumference 7.5 cm above mid patella
- Knee circumference
- Knee circumference (at) mid-patella
- Calf circumference 7.5 cm below mid patella
- Calf circumference 15 cm below mid patella
- Calf circumference 20 cm below mid patella

### ORDERING INFORMATION

#### STANDARD

| Type  | Part#       |
|-------|-------------|
| Left  | B-155610001 |
| Right | B-155510001 |

M/L at joint width in a standard, weight-bearing position 2 Thigh Circumference: 6" above mid-patella

# REBOUND® DUAL CUSTOM

## ACCESSORIES

### SLEEVES

| Size  | Circumference<br>15 cm above mid<br>patella | Neoprene<br>undersleeve | Neoprene<br>oversleeve | Sport<br>undersleeve | Padded sports<br>oversleeve |
|-------|---------------------------------------------|-------------------------|------------------------|----------------------|-----------------------------|
| S     | < 38 cm                                     | 20251                   | 20051                  | 20682                | 20448                       |
| M     | 38 – 45 cm                                  | 20252                   | 20052                  | 20683                | 20448                       |
| L     | 45 – 49.5 cm                                | 20253                   | 20053                  | 20684                | 20449                       |
| XL    | 49.5 – 54.5 cm                              | 20254                   | 20054                  | 20685                | 20449                       |
| XXL   | 54.5 – 62 cm                                | 20255                   | 20055                  |                      |                             |
| XXXL  | 62 – 71 cm                                  | 20256                   | 20056                  |                      |                             |
| XXXXL | 71 – 86 cm                                  | 20257                   | 20057                  |                      |                             |

### EXTENSION STOP KIT

| Colour | Part#       |
|--------|-------------|
| Black  | B-701090008 |

Includes: 5°, 10°, 20°, 30°, 40° extension stops. (0° and 10° are included with brace).

### FLEXION STOP KIT

| Size      | Colour | Part#       |
|-----------|--------|-------------|
| Universal | Black  | B-701090004 |

Includes: 45°, 60°, 75°, 90° flexion stops

### PCL STRAP KIT

| Size   | Part#       |
|--------|-------------|
| XS-XL  | B-701090057 |
| XL-3XL | B-701090058 |

### SMARTDOSING® KITS

| Description                | Part#       |
|----------------------------|-------------|
| Right Medial, Left Lateral | B-701090060 |
| Right Lateral, Left Medial | B-701090061 |

### TIBIAL WEDGE KIT

| Size      | Part#       |
|-----------|-------------|
| Universal | B-701090001 |

### LINER KIT

| Part#       |
|-------------|
| B-701090050 |

### CONDYLE PAD KIT

| Part#       |
|-------------|
| B-701090009 |

### FIELD SERVICABLE D-RING KIT

| Part#       |
|-------------|
| B-701090059 |

### STRAP KIT

| Size   | Part#       |
|--------|-------------|
| XS-L   | B-701090105 |
| XL-3XL | B-701090106 |

### LINERS

| Size | Doeskin calf liner | ActiveGrip thigh liner |
|------|--------------------|------------------------|
| XS   | B-701090028        | B-701090041            |
| S    | B-701090029        | B-701090042            |
| M    | B-701090030        | B-701090043            |
| L    | B-701090031        | B-701090044            |
| XL   | B-701090032        | B-701090045            |
| XXL  | B-701090033        | B-701090046            |
| 3XL  | B-701090034        | B-701090047            |

### THERMOFORMED LINER KITS (STANDARD)

| Size | Part#       |
|------|-------------|
| XS   | B-701090110 |
| S    | B-701090111 |
| M    | B-701090112 |
| L    | B-701090113 |
| XL   | B-701090114 |
| XXL  | B-701090115 |
| 3XL  | B-701090116 |

## REBOUND® CARTILAGE

The success of an articular cartilage repair is dependent on patient compliance to the rehabilitation process. The protocol is based on restricted cartilage load combined with gradual increase in activity. Many rehabilitation protocols suggest crutches for partial weight-bearing and it is the standard of care for successful outcomes.

However, patient compliance to partial weight-bearing instructions with crutches has been shown to be inconsistent. The Rebound Cartilage is a protective functional solution designed to support regeneration of the knee cartilage after repair procedures such as Microfracture, OATS and (M)ACI. Based on the clinically-proven 3-Point Leverage System and proprietary Cartilage Protection Straps, joint unloading is maintained in flexion for protection of the cartilage during the healing process.



### INDICATIONS FOR USE

For knee conditions that may benefit from unicompartmental load reduction or range of movement restriction, such as:

- Articular cartilage defect repair
- Meniscal repair
- Avascular necrosis
- Condylar bone marrow lesions (bone bruises)
- Tibial plateau fractures
- Mild to severe unicompartmental knee osteoarthritis

### FEATURES AND BENEFITS

- Easy-to-Stock - Just 4 configurations to fit most of your patients
- Comfortable - Skin-friendly liner feels great against the skin and reduces pressure points.
- Anti-Migration Technology - Strap design captures the thickest part of the calf to reduce brace migration during activity.
- Control - SmartDosing® dials put you in control of your pain by allowing you to fine-tune strap tension, while Cartilage Protection Straps maintain joint unloading in flexion.
- Tool-Free - Tool-free polycentric hinge with optional ROM stops for range-of-motion restriction.

### REQUIRED MEASUREMENTS

- Circumference measurement: 15 cm below mid-patella

# REBOUND® CARTILAGE

## ORDERING INFORMATION

| Rebound Cartilage      | Left medial | Right medial | Left lateral | Right lateral |
|------------------------|-------------|--------------|--------------|---------------|
| Universal 31 - 61.5 cm | B-254610006 | B-254510006  | B-254620006  | B-254520006   |

| Rebound Cartilage - No Box | Left medial | Right medial | Left lateral | Right lateral |
|----------------------------|-------------|--------------|--------------|---------------|
| Universal 31 - 61.5 cm     | B-254610007 | B-254510007  | B-254620007  | B-254520007   |

| Rebound Cartilage - Custom      | Left        | Right       |
|---------------------------------|-------------|-------------|
| Calf circumference 25 - 61.5 cm | B-154600001 | B-154500001 |

## SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilises a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



| Description                | Part#       |
|----------------------------|-------------|
| SmartMeasure Measuring Kit | CMS0000     |
| Constant Force Caliper     | B-312000000 |

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

### CUSTOM STANDARD COLORS



### CUSTOM PAINTED DESIGNS



## ACCESSORIES

| Accessory Parts                                  | Part#    |
|--------------------------------------------------|----------|
| Suspension Strap Kit                             | UNL00211 |
| OA Hinge Extension Stop Tree (5°, 10°, 15°, 30°) | OAH40029 |
| OA Hinge Flexion Stop Tree (0°, 30°, 60°, 90°)   | OAH40039 |
| Wash Bag                                         | UON60200 |

### SLEEVE WITH MESH\*

| Size | Circumference 15 cm above mid patella | Part #   |
|------|---------------------------------------|----------|
| XS   | < 43 cm                               | SLV60011 |
| S    | 43 – 46 cm                            | SLV60012 |
| M    | 46 – 51 cm                            | SLV60013 |
| L    | 51 – 53 cm                            | SLV60014 |
| XL   | 53 – 56 cm                            | SLV60015 |
| XXL  | 56 – 61 cm                            | SLV60016 |

\* Sleeve with Mesh: Main fabric: 92% Cotton/8% Spandex - Mesh: 100% Polyester - Band: 80% Polyester/20% Spandex



# REBOUND® CARTILAGE

## REPLACEMENT PARTS

| Service Kit*       | Right medial/Left lateral | Right lateral/Left medial |
|--------------------|---------------------------|---------------------------|
| Black, Service Kit | CRT416901                 | CRT415901                 |

\*Includes Jersey knit thigh liner, Sensil calf liner, Popliteal pad standard, Hinge pad & Strap Kit Complete

| Thigh Liners       | Right medial/Left lateral | Right lateral/Left medial |
|--------------------|---------------------------|---------------------------|
| Black, Jersey Knit | CRT526901                 | CRT525901                 |

| Calf Liners        | Right medial/Left lateral | Right lateral/Left medial |
|--------------------|---------------------------|---------------------------|
| Black, Sensil      | CRT516952                 | CRT515952                 |
| Black, Jersey Knit | CRT516951                 | CRT515951                 |

| Strap Kit Complete*       | Part #      |
|---------------------------|-------------|
| Black, Strap Kit Complete | B-701290005 |

\*Includes 4 straps (1 Thigh, 1 Calf, 2 CPS), 5 Alligator tabs, 2 CPS Strap pads

| CPS Strap Kit* | Part #      |
|----------------|-------------|
| Black, CPS Kit | B-700290005 |

\*Includes 2 CPS Straps, 2 Alligator tabs, 2 CPS Strap pads

| Buckle Clasp Kit | Right medial/Left lateral | Right lateral/Left medial |
|------------------|---------------------------|---------------------------|
| Black            | B-708217001               | B-708218001               |

| Popliteal Pad   | Part #      |
|-----------------|-------------|
| Black, Short    | B-304010105 |
| Black, Standard | B-304010106 |

| Other Replacement Parts         | Part #      |
|---------------------------------|-------------|
| Velcro Twist-Fit Attachment Kit | B-705290000 |
| Hinge Pad                       | B-304120101 |
| Alligator Kit-8 Pieces          | B-704210001 |

## CTi® OA

CTi braces provide the ultimate combination of stabilisation and protection of the knee joint. CTi uses Accutrak® hinges with extension stops to glide with the knee and breathable liners coated with bio-inert Sensil® silicone to ensure the CTi remains properly positioned on the leg, for optimal functional benefit. Allowing for varus or valgus correction, the custom made CTi OA is an excellent, dual hinged hybrid OA and ligament knee brace for your more active patients.



### INDICATIONS FOR USE

- For knee conditions that may benefit from increased A-P and M-L stability around the knee, such as: ACL, MCL, LCL, PCL, rotary and combined instabilities
- Also treats mild to severe unicompartmental osteoarthritis

### FEATURES AND BENEFITS

- Reversible field-serviceable d-rings and trimmable straps facilitate fast adjustments and customised fit
- Corrosion resistant materials, ideal for water sports
- Anatomically-correct Accutrak® hinges with extension stops
- Proprietary Sensil® silicone liners designed to reduce migration

### REQUIRED MEASUREMENTS

- Custom-to-measurement
- 5 point knee measurement (Smartmeasure)





## CTi® OA

### ORDERING INFORMATION

#### CTi OA

| Description | Left        | Right       |
|-------------|-------------|-------------|
| CTi OA      | B-116630010 | B-116530010 |

For accessories see CTi Custom.



#### SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace, including the CTi OA. Utilises a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



#### DIGITAL MEASUREMENT SYSTEM (DMS)

This system simply requires minimal physical measurements and two digital photos, which can then be sent to customer service. This system is a quick, easy and reliable method of measuring for a CTi OA without the need for training.

Please contact customer care for a free DMS Kit. Item number - 89100

## REBOUND® PCL

Unlike the anterior cruciate ligament (ACL), the PCL has the intrinsic ability to heal and regain functionality after injury. Physiologically correct positioning of the PCL during conservative or post-surgical rehabilitation is important, as the PCL experiences variable tensile forces during knee flexion, and healing of the PCL in an elongated position can lead to chronic instability and disability. Insufficient healing after a PCL rupture leads to significant pathologies, such as chronic anterior knee pain or knee osteoarthritis.

The Rebound PCL applies a physiologically correct dynamic force to the PCL during rehabilitation. The dynamic load on the tibia is generated by applying an anterior directed dynamic force on the calf area, and opposing counterforces on the anterior aspect of the leg.

**FUNCTIONAL HEALING - WHERE MOBILITY MEANS RECOVERY**




### INDICATIONS FOR USE

- Non-surgical treatment of PCL ruptures

Post-surgical rehabilitation for:

- Partial PCL tears
- PCL reconstruction/augmentation
- PCL revision



LaPrade RF, Smith SD, Wilson KJ, Wijdicks CA. Quantification of functional brace forces for posterior cruciate ligament injuries on the knee joint: an in vivo investigation. Knee Surg Sports Traumatol Arthrosc [Internet]. 2014 Aug 22 [cited 2014 Dec 11]; Available from: <http://www.ncbi.nlm.nih.gov/pubmed/25145947>

### FEATURES AND BENEFITS

- The Dynamic Tension System™ (DTS) is designed to facilitate healing by increasing the load on the tibia as the knee goes through flexion, and providing an anterior drawer force necessary to decrease load on the healing posterior cruciate ligament (PCL).
- Three shear knob options permit customisation of the biochemical correct posterior tibial load for each patient, and allow for adjustment, e.g. in case of swelling.
- A polycentric hinge with extension and flexion stops allows for range of motion adjustment. Lightweight brace with proprietary Sensil® silicone calf liner and doeskin thigh liner to reduce migration, while providing comfort to the user.
- Posterior frame for optimised support, with open anterior for ease of donning/doffing.

### REQUIRED MEASUREMENTS

- Caliper measurement: Medio-laterally at joint width in a standard, weight-bearing position
- Calf circumference 15 cm below mid patella
- Thigh circumference 15 cm above mid patella

# REBOUND® PCL

## ORDERING INFORMATION

### REBOUND PCL

| Size | Caliper M/L*   | Calf Circ.** | Thigh Circ.*** | Left          | Right         |
|------|----------------|--------------|----------------|---------------|---------------|
| S    | 8.9 - 10.2cm   | 29 - 34.5cm  | 42 - 45cm      | B-124501161   | B-124501158   |
| M    | 10.2 - 11.4cm  | 34.5 - 38cm  | 45 - 51cm      | B-124501162   | B-124501159   |
| L    | 11.4 - 12.7cm  | 38 - 42cm    | 51 - 58.5cm    | B-124501163   | B-124501160   |
| XL   | 12.7 - 14.0 cm | 42 - 45.7 cm | 58.5 - 66.7cm  | B-124501155XL | B-124501154XL |

\* M/L measurement: In a standard, weight bearing position - Using a Caliper: M/L at joint width, Using a Constant Force Caliper: M/L at knee axis line

\*\* Calf circumference measurement: 15cm below knee centre

\*\*\* Thigh circumference measurement: 15cm above knee centre

### DOESKIN THIGH LINER

| Size | Left     | Right    |
|------|----------|----------|
| XS   | PCL60031 | PCL60011 |
| S    | PCL60032 | PCL60012 |
| M    | PCL60033 | PCL60013 |
| L    | PCL60034 | PCL60014 |
| XL   | PCL60035 | PCL60015 |
| XXL  | PCL60036 | PCL60016 |

### DOESKIN CALF LINER

| Size | Left     | Right    |
|------|----------|----------|
| XS   | PCL60171 | PCL60071 |
| S    | PCL60172 | PCL60072 |
| M    | PCL60173 | PCL60073 |
| L    | PCL60174 | PCL60074 |
| XL   | PCL60175 | PCL60075 |
| XXL  | PCL60176 | PCL60077 |

### SILICONE CALF LINER

| Size | Left     | Right    |
|------|----------|----------|
| XS   | PCL60121 | PCL60021 |
| S    | PCL60122 | PCL60022 |
| M    | PCL60123 | PCL60023 |
| L    | PCL60124 | PCL60024 |
| XL   | PCL60125 | PCL60025 |
| XXL  | PCL60126 | PCL60026 |

### STANDARD THIGH LINER KIT

| Size | Left        | Right       |
|------|-------------|-------------|
| XS   | B-124501170 | B-124501164 |
| S    | B-124501171 | B-124501165 |
| M    | B-124501172 | B-124501166 |
| L    | B-124501173 | B-124501167 |
| XL   | B-124501174 | B-124501168 |
| XXL  | B-124501175 | B-124501169 |

Includes: Doeskin thigh liner, doeskin femoral cuff liner (softgood on lower thigh strap), doeskin pre-femoral liner (softgood on upper thigh strap), hook and loop (4x round coins, 8x rec. tabs)

# REBOUND® PCL

## ORDERING INFORMATION

### STANDARD CALF LINER KIT

| Size | Left        | Right       |
|------|-------------|-------------|
| XS   | B-124601176 | B-124501176 |
| S    | B-124601177 | B-124501177 |
| M    | B-124601178 | B-124501178 |
| L    | B-124601179 | B-124501179 |
| XL   | B-124601180 | B-124501180 |
| XXL  | B-124601181 | B-124501181 |

Includes: Silicone calf liner/AMS, doeskin tibia cuff liner (softgood on lower calf strap), hook and loop (2x round coins, 4x rec. tabs)

### DOESKIN CALF LINER KIT

| Size | Left        | Right       |
|------|-------------|-------------|
| XS   | B-124601182 | B-124501182 |
| S    | B-124601183 | B-124501183 |
| M    | B-124601184 | B-124501184 |
| L    | B-124601185 | B-124501185 |
| XL   | B-124601186 | B-124501186 |
| XXL  | B-124601187 | B-124501187 |

Includes: Doeskin calf liner/AMS, doeskin tibia cuff liner (softgood on lower calf strap), tibia strap pad, hook and loop (2x round coins, 4x rec. tabs)



### SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilises a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



### DIGITAL MEASUREMENT SYSTEM (DMS)

This system simply requires minimal physical measurements and two digital photos, which can then be sent to customer service. This system is a quick, easy and reliable method of measuring for a custom knee brace without the need for training.

Please contact customer care for a free DMS Kit. Item number - 89100

### REBOUND PCL CUSTOM

| Description              | Part#       |
|--------------------------|-------------|
| Rebound PCL Custom Left  | B-124501155 |
| Rebound PCL Custom Right | B-124501154 |

\*Please allow 10 working days for delivery

### CUSTOM STANDARD COLOURS

|            |            |           |
|------------|------------|-----------|
| Red        | Pink       | Silver    |
| Orange     | Black      | White     |
| Yellow     | Navy Blue  | Champagne |
| Lime Green | Ocean Blue | Charcoal  |

# REBOUND® PCL

## ACCESSORIES

### STRAPS

| Size     | Part#    |
|----------|----------|
| S - L    | PCL00061 |
| XL - XXL | PCL00062 |

### OTHER ACCESSORIES

| Description        | Part#       |
|--------------------|-------------|
| Extension Stop Kit | B-705130042 |
| Flexion Stop Kit   | B-705130043 |
| Shear Knob Kit     | B-465350087 |
| Condyle Pad Kit    | B-705130046 |



Shear Knob Kit

## REBOUND® ACL

The Rebound ACL brace is an evolution in bracing, designed to comfortably apply an anatomically correct, dynamic load to the anterior cruciate ligament (ACL). As an element in either conservative or surgical management, the Rebound ACL brace complements clinical protocols and is conducive to an ideal biomechanical joint environment as a foundation for long-term successful patient outcomes. The innovative technology in the Rebound® Dynamic Tension System (DTS) offers individualised adjustment options to reduce strain on the healing ACL.




FUNCTIONAL HEALING – WHERE MOBILITY MEANS RECOVERY

### INDICATIONS FOR USE

Conditions benefiting from an agonistic force that varies throughout flexion in proportion to ACL tension applied to the thigh and calf in the sagittal plane. These may include:

- Anterior cruciate ligament (ACL) injury
- Reconstruction
- Non-surgical treatment



### FEATURES AND BENEFITS


- Dynamic technology for knee stability, with a Dynamic Tension System™ that applies an agonistic, dynamic force to the ACL that is flexion-angle dependent (LaPrade et al., 2017) and designed to reduce ACL deficiency (positive anterior-drawer test), promoting knee stability.
- Specific load adaptation according to the patient's individual anatomy and rehabilitation requirements. Clinical and biomechanical studies indicate that the Rebound ACL brace reduces strain on the ACL during activity and significantly reduces strain on the meniscus in ACL-deficient knees (Tomescu et al., 2017). Biocompatible, hi-tech materials, including Blue is You™ patient touchpoints and ActiveGrip™ liners, promote suspension and comfort.
- Lightweight, aircraft-grade aluminium frame
- Available in sizes S-XL
- Flexion & extension range-of-motion lock-outs.
- Posterior frame optimised for fit and support, with open anterior for ease of donning/doffing Clinically-tested comfort, validated by wearers and with load mapping data collected during rehabilitation activities, the Rebound® Dynamic Tension System (DTS) provides a comfortable load to the posterior thigh and anterior tibia.

# REBOUND® ACL


## REQUIRED MEASUREMENTS



Measure knee caliper medio-laterally at joint width in standing, weight-bearing position.



Measure calf circumference 15 cm below mid-patella. Leg in extension.



Measure thigh circumference 15 cm above mid-patella. Leg in extension.

## ORDERING INFORMATION

### STANDARD

| Size | Left        | Right       | Caliper       | Calf circumference 15 cm below mid patella | Thigh circumference 15 cm above mid patella |
|------|-------------|-------------|---------------|--------------------------------------------|---------------------------------------------|
| SM   | B-124601191 | B-124501191 | 8.9 - 10.2cm  | 29.2 - 34.3cm                              | 41.9 - 47.0cm                               |
| MD   | B-124601192 | B-124501192 | 10.2 - 11.4cm | 34.3 - 39.4cm                              | 47.0 - 52.1cm                               |
| LG   | B-124601193 | B-124501193 | 11.4 - 12.7cm | 39.4 - 44.5cm                              | 52.1 - 57.2cm                               |
| XL   | B-124601194 | B-124501194 | 12.7 - 14.0cm | 44.5 - 49.5cm                              | 57.2 - 62.2cm                               |

Recommended for patients 1.56m and taller

### REPLACEMENT PARTS

| Description                                                         | Part#       |
|---------------------------------------------------------------------|-------------|
| Extension stop kit (0°, 5°, 10°, 20°, 30°, 40°)                     | B-124702035 |
| Flexion stop kit (0°, 15°, 25°, 35°, 45°, 55°, 65°, 75°, 90°, 105°) | B-124702034 |
| Shear Knob Kit                                                      | B-705130051 |
| SM/MD Strap Kit                                                     | B-124702030 |
| LG/XL Strap Kit                                                     | B-124702031 |
| Condyle Pad Kit                                                     | B-701090009 |
| Tibial Pad                                                          | B-124702032 |
| Tibial Hex Pad                                                      | B-124702033 |

# REBOUND® ACL

## ACCESSORIES

### ANTI-MIGRATION WRAP

| SIZE     | RIGHT       | LEFT        |
|----------|-------------|-------------|
| S/M      | B-124501195 | B-124601195 |
| S/M HEX  | B-124501196 | B-124601196 |
| L/XL     | B-124501197 | B-124601197 |
| L/XL HEX | B-124501198 | B-124601198 |

### UPPER DOESKIN LINER KIT

| SIZE | RIGHT       | LEFT        |
|------|-------------|-------------|
| S    | B-124702004 | B-124702000 |
| M    | B-124702005 | B-124702001 |
| L    | B-124702006 | B-124702002 |
| XL   | B-124702007 | B-124702003 |

### LOWER ACTIVEGRIP LINER KIT

| SIZE | RIGHT       | LEFT        |
|------|-------------|-------------|
| S    | B-124702014 | B-124702010 |
| M    | B-124702015 | B-124702011 |
| L    | B-124702016 | B-124702012 |
| XL   | B-124702017 | B-124702013 |

### LOWER DOESKIN LINER KIT

| SIZE | RIGHT       | LEFT        |
|------|-------------|-------------|
| S    | B-124702024 | B-124702020 |
| M    | B-124702025 | B-124702021 |
| L    | B-124702026 | B-124702022 |
| XL   | B-124702027 | B-124702023 |



## REBOUND® ACL

### ORDERING INFORMATION FOR REBOUND ACL CUSTOM

#### REBOUND ACL CUSTOM

| Part#       | Type  | Min / Max Caliper | Min / Max Circ. above | Min / Max Circ. below | Measuring          |
|-------------|-------|-------------------|-----------------------|-----------------------|--------------------|
| B-124701000 | Right | 8.9 - 14cm        | 41.9 - 62.2cm         | 29.2 - 49.5cm         | SmartMeasure, Cast |
| B-124701001 | Left  |                   |                       |                       |                    |

Caliper = M/L width at knee center  
Circumference = 15cm above / below mid-patella  
Ideal for patients 1.56m and taller

#### SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilises a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



| Description            | Part#       |
|------------------------|-------------|
| Constant Force Caliper | B-312000000 |

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required.

#### REQUIRED MEASUREMENTS

- Custom made to measure ( Smartmeasure, 3D-Scan or patient cast)

## REBOUND® DUAL BASIC AND REBOUND® DUAL BASIC ST

The Rebound DUAL Basic and Rebound DUAL Basic ST offer the value and simplicity you need to navigate the complex coverage landscape. Within the versatility of the Rebound DUAL collection, Össur offers high quality solutions from a single source, to help users transition to a successful recovery without compromise.



Rebound DUAL Basic

Rebound DUAL Basic ST  
(Step-Through)

### INDICATIONS FOR USE

The Rebound DUAL Basic is indicated for combined and non-combined instabilities of the:

- Anterior cruciate ligament (ACL)
- Medial collateral ligament (MCL)
- Lateral collateral ligament (LCL)
- Posterior cruciate ligament (PCL)
- (Only Rebound DUAL Basic) Acculign version only: Conditions requiring unicompartmental load control

### FEATURES AND BENEFITS

- (Only Rebound DUAL Basic ST) Step-Through Design: advantage in comfort & fitting with ease
- (Only Rebound DUAL Basic) Multi-purpose length-adjustable upper frame accommodates the user's height, helping to ensure the perfect fit and enabling adjustment for varus/valgus alignment
- Össur Hinge alignment guide provides easy & accurate re-fitting for users, ensuring correct placement of the brace relative to the knee
- Lightweight aircraft aluminium frame keeps a low profile
- Doeskin liners for comfort & no-slip Activegrip™ strap pad
- Trimmable straps with tactile strap tabs
- Low-profile polycentric hinge with stop options for stability & range-of-motion restriction
- Step-through or Anterior frame design
- Suitable for low-high impact and sliding sports
- Optional PCL kit accessory for multi-ligament support
- Available in pre-fabricated sizes XS-3XL

### REQUIRED MEASUREMENTS

- Knee caliper
- Thigh circumference 15 cm above mid patella

# REBOUND® DUAL BASIC AND REBOUND® DUAL BASIC ST

## ORDERING INFORMATION

**X** = 5 (right)  
 6 (left)  
**Y** = 1 (Rebound DUAL Basic)  
 2 (Rebound DUAL Basic ST)

| Description                                  | Size | Caliper       | Thigh Circumference | Part #                |
|----------------------------------------------|------|---------------|---------------------|-----------------------|
| Rebound DUAL Basic or Rebound DUAL Basic ST* | XS   | 7.6 - 8.9cm   | 38.1 - 41.9cm       | B-242 <b>XY</b> 0001B |
|                                              | S    | 8.9 - 10.2cm  | 41.9 - 47cm         | B-242 <b>XY</b> 0002B |
|                                              | M    | 10.2 - 11.4cm | 47 - 52.1cm         | B-242 <b>XY</b> 0003B |
|                                              | L    | 11.4 - 12.7cm | 52.1 - 57.2cm       | B-242 <b>XY</b> 0004B |
|                                              | XL   | 12.7 - 14cm   | 57.2 - 62.2cm       | B-242 <b>XY</b> 0005B |
|                                              | XXL  | 14 - 15.2cm   | 62.2 - 69.9cm       | B-242 <b>XY</b> 0006B |
|                                              | 3XL  | 15.2 - 16.5cm | 69.9 - 76cm         | B-242 <b>XY</b> 0007B |

## REPLACEMENT LINER KITS

| Measuring | Basic Liner Part # | Basic ST Liner Part # |
|-----------|--------------------|-----------------------|
| XS        | B-701090130        | B-701090140           |
| S         | B-701090131        | B-701090141           |
| M         | B-701090132        | B-701090142           |
| L         | B-701090133        | B-701090143           |
| XL        | B-701090134        | B-701090144           |
| XXL       | B-701090135        | B-701090145           |
| 3XL       | B-701090136        | B-701090146           |

## REPLACEMENT PARTS

| Size                            | Part #      |
|---------------------------------|-------------|
| RDUAL Basic Strap Kit XS-L      | B-701090137 |
| RDUAL Basic Strap Kit XL-3XL    | B-701090138 |
| RDUAL Basic ST Strap Kit XS-L   | B-701090147 |
| RDUAL Basic ST Strap Kit XL-3XL | B-701090148 |
| Condyle Pad Kit                 | B-701090009 |

## ACCESSORIES

| Size                                | Part #      |
|-------------------------------------|-------------|
| RDUAL Basic ROM Stop Kit            | B-701090109 |
| RDUAL Basic PCL Strap Kit XS-L      | B-701090057 |
| RDUAL Basic PCL Strap Kit XL-3XL    | B-701090058 |
| RDUAL Basic ST PCL Strap Kit XS-L   | B-701090065 |
| RDUAL Basic ST PCL Strap Kit XL-3XL | B-701090066 |

\* Rebound DUAL Basic ST is 35.5cm

## OA EASE

The Formfit OA Ease is a single upright brace with a Dynamic Force System™ strap designed to relieve pain from mild-to-moderate unicompartmental knee osteoarthritis (OA). Its wraparound design is easy to apply and remove using simple rubber pull-tabs, colour-coded guide dots, and strap tension indicator lines. This low-profile knee brace is comfortable, easy to use, and adapts to a variety of leg anatomies. The Formfit OA Ease can also be machine washed for effortless care.



### INDICATIONS FOR USE


Take control of mild-to-moderate OA pain with the Formfit OA Ease, an easy-to-wear knee brace with a simple wraparound design that can be applied and removed with ease.

- Mild-to-moderate unicompartmental knee osteoarthritis and/or degenerative meniscal tears
- Other conditions that may benefit from unicompartmental unloading

### FEATURES AND BENEFITS

- Easy Wraparound Design - Simple rubber pull-tabs, colour-coded guide dots and a convenient hinge alignment guide make application and removal simple
- Relieves OA Pain - 3-Point Leverage System with Dynamic Force System™ strap and shells unload the affected knee compartment
- Machine Washable - Lightweight and breathable fabric feels great against the skin and can be machine washed for effortless care
- Tool-Free - Tool-free polycentric hinge with available ROM stop options for range-of-motion restriction
- Control – Fine-tune strap tension for individual pain relief, with marked dosing indicator lines
- Extension kits and liners are available to further enhance fit and comfort

### REQUIRED MEASUREMENTS



Measure calf circumference 15cm below mid-patella. Leg in extension.



Measure thigh circumference 15cm above mid-patella. Leg in extension.

# OA EASE

## ORDERING INFORMATION

| Size | Calf circ. 15cm<br>below mid patella | Thigh circ. 15cm<br>above mid patella | Right Medial | Left Medial | Right Lateral | Left Lateral |
|------|--------------------------------------|---------------------------------------|--------------|-------------|---------------|--------------|
| XS   | 28 - 32cm                            | 36 - 41cm                             | FOA0518711   | FOA0618711  | FOA0528711    | FOA0628711   |
| S*   | 32 - 36cm                            | 41 - 46cm                             | FOA0518712   | FOA0618712  | FOA0528712    | FOA0628712   |
| M    | 36 - 41cm                            | 46 - 53cm                             | FOA0518713   | FOA0618713  | FOA0528713    | FOA0628713   |
| L*   | 41 - 46cm                            | 53 - 60cm                             | FOA0518714   | FOA0618714  | FOA0528714    | FOA0628714   |
| XL   | 46 - 51cm                            | 60 - 67cm                             | FOA0518715   | FOA0618715  | FOA0528715    | FOA0628715   |
| XXL* | 51 - 56cm                            | 67 - 75cm                             | FOA0518716   | FOA0618716  | FOA0528716    | FOA0628716   |

\*All sizes are available but for consolidated stock option this size can be achieved by combining the smaller size with accessories below.

## ACCESSORIES

| Description                                             | Part#     |
|---------------------------------------------------------|-----------|
| ActiveGrip Calf Liner Kit 2MM                           | FOA000101 |
| ActiveGrip Calf Liner Kit 6MM                           | FOA000102 |
| Doeskin Thigh Liner Kit 2MM                             | FOA000201 |
| Doeskin Thigh Liner Kit 6MM                             | FOA000202 |
| Small Extension Kit 2 Pieces                            | FOA000301 |
| Small Extension Kit 10 Pieces                           | FOA000302 |
| Large Extension Kit XS-SM (Right Medial/Left Lateral)   | FOA000401 |
| Large Extension Kit MD-LG (Right Medial/Left Lateral)   | FOA000402 |
| Large Extension Kit XL-XXXL (Right Medial/Left Lateral) | FOA000403 |
| Large Extension Kit XS-SM (Right Lateral/Left Medial)   | FOA000404 |
| Large Extension Kit MD-LG (Right Lateral/Left Medial)   | FOA000405 |
| Large Extension Kit XL-XXXL (Right Lateral/Left Medial) | FOA000406 |
| Wash Bag                                                | WB000001  |
| Suspension Strap                                        | OAW00007  |
| OA Hinge Extension Stop Kit (5°, 10°, 15°, 30°)         | OA40029   |
| OA Hinge Flexion Stop Kit (0°, 30°, 60°, 90°)           | OA40039   |

## REPLACEMENT PARTS

| Description                                | Part#       |
|--------------------------------------------|-------------|
| Alligator Clip (Right Medial/Left Lateral) | FOW00501    |
| Alligator Clip (Right Lateral/Left Medial) | FOW00502    |
| Hinge Pad                                  | B-304120101 |

## CTi®3

CTi®3 brings a new vision to double upright knee braces. Innovative engineering makes CTi3 an excellent option for most of your patient population. Suitable for both everyday activities and sport, CTi3 offers great quality with its enhanced construction and anatomical fit for elevated comfort.



### INDICATIONS

- For knee conditions that may benefit from increased AP and ML stability and protection, such as injuries to the ACL, MCL, LCL, PCL, rotary and combined instabilities.

### FEATURES AND BENEFITS

- Enhanced construction – Lightweight, high strength carbon composite frames provide the optimal stiffness to weight ratio
- Anatomical fit – Designed to fit a wider range of leg shapes and includes an Accutrak® hinge to closely follow the knee's anatomical movement
- Elevated comfort – Unique liners enhance comfort with their special tread block pattern, which provides a balance of cushion and ventilation
- Easy-to-use accessories and replacement parts, including MX kits, replacement strap sets, liner sets and flexion stops

### REQUIRED MEASUREMENTS

- Knee caliper
- Thigh circumference 20 cm (8") above mid patella

### ORDERING INFORMATION

#### CTi3

| Size | Caliper  | Circumference | Left      | Right     |
|------|----------|---------------|-----------|-----------|
| S    | 8-10 cm  | 40-47 cm      | CTI-23002 | CTI-22002 |
| M    | 10-11 cm | 45-55 cm      | CTI-23003 | CTI-22003 |
| L    | 11-13 cm | 52-61 cm      | CTI-23004 | CTI-22004 |
| XL   | 13-14 cm | 58-67 cm      | CTI-23005 | CTI-22005 |

Includes: Condyle Pad Set and Extension Stop Set, 0° Extension Stop and Standard Condyle Pad are pre-installed.

Note: M/L measurements are the preferred and most accurate form of measurement for CTi. Circumferential measurements are approximate.

# CTi®3

## ACCESSORIES & REPLACEMENT PARTS

### CTi3 MX KIT

| Size | Left      | Right     |
|------|-----------|-----------|
| S    | CTI-23132 | CTI-22132 |
| M    | CTI-23133 | CTI-22133 |
| L    | CTI-23134 | CTI-22134 |
| XL   | CTI-23135 | CTI-22135 |

Includes: Impact Guard, Sleeve with Mesh, CTi Decal Sheet and CTi Product Bag.

### IMPACT GUARD

| Size | Left      | Right     |
|------|-----------|-----------|
| S-M  | CTI-23112 | CTI-22112 |
| L-XL | CTI-23115 | CTI-22115 |

### LINER SET

| Size | Left      | Right     |
|------|-----------|-----------|
| S    | CTI-23322 | CTI-22322 |
| M    | CTI-23323 | CTI-22323 |
| L    | CTI-23324 | CTI-22324 |
| XL   | CTI-23325 | CTI-22325 |

### HOOK SET

| Size | Left      | Right     |
|------|-----------|-----------|
| S    | CTI-23332 | CTI-22332 |
| M    | CTI-23333 | CTI-22333 |
| L    | CTI-23334 | CTI-22334 |
| XL   | CTI-23335 | CTI-22335 |

### STRAP SET

| Size | Part #    |
|------|-----------|
| S    | CTI-24222 |
| M    | CTI-24223 |
| L    | CTI-24224 |
| XL   | CTI-24225 |

### MISCELLANEOUS

| Description                                | Part #    |
|--------------------------------------------|-----------|
| Guard Elastic Cord                         | CTI-24130 |
| Guard Screws                               | CTI-24140 |
| Tibial Wedge Pad                           | 22077     |
| Condyle Pad Set (Thin, Standard, Thick)    | CTI-24320 |
| Extension Stop Set (0°, 5°, 10°, 20°, 30°) | CTI-24430 |
| Flexion Stop Set (0°, 60°, 75°, 90°)       | CTI-24440 |
| Hinge Sticker Set                          | CTI-24450 |
| D-Ring Set                                 | CTI-24510 |
| CTi3 PCL Strap Kit                         | CTI-24520 |

### SLEEVE WITH MESH

| Size | Part #   |
|------|----------|
| S    | SLV60012 |
| M    | SLV60013 |
| L    | SLV60014 |
| XL   | SLV60015 |

## CTi® CUSTOM

CTi braces provide the ultimate combination of stabilisation and protection of the knee joint. This truly custom made brace, with rock-solid CTi carbon frame, uses Accutrac® hinges with extension stops to glide with the knee and breathable liners coated with bio-inert Sensil® silicone to ensure the CTi remains properly positioned on the leg, for optimal functional benefit. CTi custom braces are available in different models and can be further customised with a variety of options, colours, designs and accessories.



### INDICATIONS FOR USE

For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, rotary and combined instabilities

### FEATURES AND BENEFITS

- Reversible field-serviceable d-rings and trimmable straps facilitate fast adjustments and customised fit (optional)
- Corrosion resistant materials, ideal for water sports
- Anatomically-correct Accutrac® hinges with extension stops
- Proprietary Sensil® silicone liners designed to reduce migration
- Available in 3 different models: Standard model for moderate to high impact levels, Pro Sport model for the impact levels and Vapor model for low to moderate impact levels

### REQUIRED MEASUREMENTS

- Custom made to measure (CMS, Smartmeasure, 3D-Scan or patient cast)


#### CUSTOM STANDARD COLOURS



#### CUSTOM PAINTED DESIGNS



No charge for custom standard colours. Additional fees required for other custom colours, custom painted designs and custom printed graphics.





# CTI® CUSTOM

## ORDERING INFORMATION

### CTI CUSTOM

| Description | Left        | Right       |
|-------------|-------------|-------------|
| Standard    | B-116600010 | B-116500010 |
| Pro Sport   | B-117600010 | B-117500010 |
| Vapor       | B-118600010 | B-118500010 |

\*Please allow 10 working days for delivery



### SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace, including the CTi. Utilises a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



### DIGITAL MEASUREMENT SYSTEM (DMS)

This system simply requires minimal physical measurements and two digital photos, which can then be sent to customer service. This system is a quick, easy and reliable method of measuring for a CTi without the need for training.

Please contact customer care for a free DMS Kit. Item number - 89100

## ACCESSORIES

### CALF LINERS

| Description      | Part#       |
|------------------|-------------|
| CTi Custom liner | B-411500115 |

### STRAPS KITS

| Size       | Part#       |
|------------|-------------|
| S, M       | B-701130104 |
| L, XL, XXL | B-701130108 |

### PCL KITS

| Size       | Part#       |
|------------|-------------|
| S, M       | B-705000132 |
| L, XL, XXL | B-705000235 |



# CTI® CUSTOM

## ACCESSORIES

### STRAP PADS

| Size    | Part#       |
|---------|-------------|
| S       | B-304010057 |
| M, L    | B-304010058 |
| XL, XXL | B-304010059 |



Strap pads

### AMS (ANTI-MIGRATION SYSTEM) WRAPS

| Description     | Part# |
|-----------------|-------|
| AMS Wrap EVA    | 22095 |
| AMS Wrap Neo    | 22094 |
| AMS Wrap Neo XL | 22152 |



AMS Wrap

### MX GEAR

| Description                                                       | Part#       |
|-------------------------------------------------------------------|-------------|
| Patella Cup                                                       | B-711130011 |
| Gear Guards                                                       | B-705000003 |
| MX Kit S (patella cup, gear guards, undersleeve, stickers)        | B-711140002 |
| MX Kit M (patella cup, gear guards, undersleeve, stickers)        | B-711140003 |
| MX Kit L (patella cup, gear guards, undersleeve, stickers)        | B-711140004 |
| MX Kit XL (patella cup, gear guards, undersleeve, stickers)       | B-711140005 |
| MX Kit XXL (patella cup, gear guards, Sleeve with Mesh, stickers) | B-711140006 |

MX Kit



### SLEEVES

| Size  | Circumference<br>15 cm above<br>mid patella | Softsleeve | Neoprene<br>undersleeve | Neoprene<br>oversleeve |
|-------|---------------------------------------------|------------|-------------------------|------------------------|
| S     | < 38 cm                                     | 20643      | 20251                   | 20051                  |
| M     | 38 – 45 cm                                  | 20644      | 20252                   | 20052                  |
| L     | 45 – 49.5 cm                                | 20644      | 20253                   | 20053                  |
| XL    | 49.5 – 54.5 cm                              | 20645      | 20254                   | 20054                  |
| XXL   | 54.5 – 62 cm                                | 20645      | 20255                   | 20055                  |
| XXXL  | 62 – 71 cm                                  | 20646      | 20256                   | 20056                  |
| XXXXL | 71 – 86 cm                                  | 20647      | 20257                   | 20057                  |

### SLEEVE WITH MESH\*

| Size | Circumference<br>15 cm above<br>mid patella | Sleeve with<br>Mesh |
|------|---------------------------------------------|---------------------|
| XS   | < 43 cm                                     | SLV60011            |
| S    | 43 – 46 cm                                  | SLV60012            |
| M    | 46 – 51 cm                                  | SLV60013            |
| L    | 51 – 53 cm                                  | SLV60014            |
| XL   | 53 – 56 cm                                  | SLV60015            |
| XXL  | 56 – 61 cm                                  | SLV60016            |



\* Sleeve with Mesh: Main fabric: 92% Cotton/8% Spandex - Mesh: 100% Polyester - Band: 80% Polyester/20% Spandex

### OTHER ACCESSORIES

| Description                                       | Part#       |
|---------------------------------------------------|-------------|
| CTi D-Ring Kit Standard 2x 1" and 2x 1.5" d-rings | B-704132600 |
| CTi D-Ring Kit 2" 2x 2" d-rings                   | B-704132602 |
| CTi Condyle Pad Kit                               | B-704130150 |
| Extension Stop Kit                                | 22130       |
| Flexion Stop Kit                                  | 19254       |

## FLEX™

Flex is a post-op/functional rehab solution offering more versatility in both shape and function than most other options. It's heat-moldable, carbon composite frame accommodates oedema and/or changes in leg shape during recovery from surgery, while its Accutrac® range-of-motion hinges provide flexion and extension control. Gel Fit™ condylar pads improve suspension and are comfortable for tender knees recovering from trauma.

### INDICATIONS FOR USE

- ACL, MCL, LCL, PCL, rotary and functional support
- Long-term use for post-op, rehab and functional support
- Recommended for low to moderate activity

### FEATURES AND BENEFITS

- Accutrac® hinges with easy to use flexion/extension stops
- Heat-mouldable, carbon composite frame accommodates oedema and/or changes in leg shape during rehabilitation
- Gel Fit™ Condylar pad for improved suspension and comfort
- Total Support System™ for ACL, MCL, LCL, PCL, rotary and combined instabilities
- Adjustable Tibial Strap captured the tibia and allows adjustments in hinge depth for a customised fit
- Large, "open window" frame design accommodates post-op dressing and oedema
- Corrosion resistant materials, ideal for use in water



### REQUIRED MEASUREMENTS

- Caliper measurement: Medio-laterally at joint width in a standard, weight-bearing position
- Circumference: Just below the distal border of the patella

## ORDERING INFORMATION

|      |               |               | Standard |          | PCL          |              |
|------|---------------|---------------|----------|----------|--------------|--------------|
| Size | Caliper       | Circumference | Left     | Right    | Left         | Right        |
| S    | 8.9 - 10.2cm  | 27.9 - 32.1cm | FLEX-SL  | FLEX-SR  | FLEXPC-L-SL  | FLEXPC-L-SR  |
| M    | 10.2 - 11.4cm | 32.1 - 36.2cm | FLEX-ML  | FLEX-MR  | FLEXPC-L-ML  | FLEXPC-L-MR  |
| L    | 11.4 - 12.4cm | 36.2 - 39.1cm | FLEX-LL  | FLEX-LR  | FLEXPC-L-LL  | FLEXPC-L-LR  |
| XL   | 12.4 - 13.3cm | 39.1 - 42.2cm | FLEX-XLL | FLEX-XLR | FLEXPC-L-XLL | FLEXPC-L-XLR |

## ACCESSORIES

### REFURB KIT\*

| Size | Gel Left   | Gel Right  | Neoprene Left | Neoprene Right |
|------|------------|------------|---------------|----------------|
| S    | F/O-SETLS  | F/O-SETRS  | F/O-SETLSNEO  | F/O-SETRSNEO   |
| M    | F/O-SETLM  | F/O-SETRM  | F/O-SETLMNEO  | F/O-SETRMNEO   |
| L    | F/O-SETLL  | F/O-SETRL  | F/O-SETLLNEO  | F/O-SETRLNEO   |
| XL   | F/O-SETLXL | F/O-SETRXL | F/O-SETLXLNEO | F/O-SETRXLNEO  |

\*incl. liners, straps, strap pads, condyle pads

## FLEX™

### LINERS

| Size | Tibia Left | Tibia Right | Femur Left | Femur Right |
|------|------------|-------------|------------|-------------|
| S    | 44048      | 44049       | 44046      | 44047       |
| M    | 44040      | 44041       | 44042      | 44043       |
| L    | 44044      | 44045       | 44038      | 44039       |
| XL   | 44052      | 44053       | 44050      | 44051       |

### ACCESSORIES

#### STRAPS AND STRAP PAD KITS

| Size | Part# |
|------|-------|
| S    | 44071 |
| M    | 44069 |
| L    | 44070 |
| XL   | 44068 |

#### OTHER ACCESSORIES

| Description                                       | Part#       |
|---------------------------------------------------|-------------|
| Flexion / Extension Stop Kit                      | 44035       |
| Flexion Stop Kit (10,20,30°)                      | B-305070010 |
| Hinge Cover Kit                                   | 44083       |
| Condyle Pad Gel (1/4in)                           | 22183       |
| Condyle Pad Kit Neoprene (1/8in + 1/16in + cover) | 44066       |
| Rigid Condyle Pad Left                            | 44027       |
| Rigid Condyle Pad Right                           | 44028       |

## C180™ ROCKET

With a condensed size and Accutrac™ hinges that are modified to more accurately track young knees, the lightweight C180 Rocket provides support for children with ligament instabilities.

### INDICATIONS FOR USE

- For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, rotary and combined instabilities

### FEATURES AND BENEFITS

- Anatomically-correct Accutrac® hinges with easy-to-use extension stops
- Semi-rigid frame consisting of woven carbon fibre infused with acrylic resin, is heat-mouldable for customisation

### REQUIRED MEASUREMENTS

- Caliper measurement: Medio-laterally at joint width in a standard, weight-bearing position
- Circumference: Just below the distal border of the patella



## ORDERING INFORMATION

| Size | Caliper      | Circumference | Left     | Right    |
|------|--------------|---------------|----------|----------|
| S    | 7.6 - 8.9cm  | 24.1 - 27.9cm | C180KSL  | C180KSR  |
| L    | 8.3 - 9.5cm  | 26 - 29.8cm   | C180KLL  | C180KLR  |
| XL   | 8.9 - 10.2cm | 27.9 - 32.1cm | C1802XSL | C1802XSR |

## ACCESSORIES

| Description         | Part#       |
|---------------------|-------------|
| Motocross kit S / L | 10833       |
| Motocross kit XL    | B-705130072 |
| Patella Cup S / L   | 25029       |
| Patella Cup XL      | B-711130000 |

## REBOUND® POST-OP KNEE

The simple, versatile post-operative knee bracing solution that is easy for everyone. Whether adjusting ROM, setting the drop lock, or altering the brace length, you'll appreciate the simplicity of the Rebound Post-Op knee brace.



### INDICATIONS FOR USE

- ACL, PCL, MCL, LCL repairs/injuries
- Meniscal repairs
- Tibial plateau fractures
- Patellar tendon repairs
- Osteochondral repairs
- Condylar fractures
- Sprains/strains of the knee
- HTOs (High Tibial Osteotomy)

### FEATURES AND BENEFITS

- Tool-free design and multi-purpose Smart-Fit™ packaging
- Available in universal, left or right contour and paediatric
- Telescoping frame from 47cm-68cm (Adult) and 39.5cm - 54cm (paediatric)
- OneTouch™ easy to adjust ROM hinge with lock-out capability
- Colour-coding and numbered frame for patient education
- Quick-fit, gripping buckles for easy application

### REQUIRED MEASUREMENTS

- Circumference: 15cm above mid-patella

### ORDERING INFORMATION

#### REBOUND POST-OP

| Description                | Thigh circumference | Cool      | Full foam |
|----------------------------|---------------------|-----------|-----------|
| Paediatric                 | 30.5 - 63.5cm       | B-8381470 |           |
| Universal, Left contoured  | 41 - 81cm           | B-8381680 | B-8382680 |
| Universal, Right contoured | 41 - 81cm           | B-8381580 | B-8382580 |
| Universal                  | 41 - 81cm           | B-8381480 | B-8382480 |
| Universal XL               | 41 - 96,5cm         | B-8381490 | B-8382490 |

#### ACCESSORIES

| Description                                       | Part#     | Description               | Part#     |
|---------------------------------------------------|-----------|---------------------------|-----------|
| ROM stop locks                                    | B-8380005 | Standard buckle & pad kit | B-8381096 |
| Alligator hook & loop kit                         | B-8380006 | AMS Kit Paediatric        | B-8381093 |
| Adhesive Hook & Loop kit Paediatric               | B-8380008 | AMS Kit Regular           | B-8381094 |
| Adhesive Hook & Loop kit universal / universal XL | B-8380007 | AMS Kit XL                | B-8381095 |

#### REPLACEMENT SOFTGOODS

| Description  | Cool      | Full foam |
|--------------|-----------|-----------|
| Paediatric   | B-8381072 |           |
| Universal*   | B-8381082 | B-8382082 |
| Universal XL | B-8381092 | B-8382092 |

\*For universal, left contoured, right contoured

## FORMFIT POST-OP KNEE

### INDICATIONS FOR USE

- Post-op/post-injury ROM control or immobilisation
- Indicated for ligament and fracture-related surgeries and when controlled range of motion is necessary

### FEATURES AND BENEFITS

- Breakaway struts for easy, tool free length adjustment
- OneTouch™ easy to adjust ROM hinge with lock-out capability from -10° to 30° (10° incremental)
- Extension control from -10° to 90° (10° incremental)
- Flexion control from 0° to 120° (10° incremental)
- Available in cool or foam models, 3 lengths, and two sizes (Regular and XL)



### REQUIRED MEASUREMENTS



Measure thigh circumference  
15 cm above mid-patella.  
Leg in extension.

| Size    | Thigh circumference 15 cm above mid patella |
|---------|---------------------------------------------|
| XL      | 69 - 89cm                                   |
| Regular | 0 - 69cm                                    |

### ORDERING INFORMATION

| Size            | Thigh circumference | Cool       | Full Foam  |
|-----------------|---------------------|------------|------------|
| Knee PO 61cm    | UP TO 69cm          | FPOK403100 | FPOK403200 |
| Knee PO XL 61cm | 69 - 89cm           | FPOK403102 | FPOK403202 |
| Knee PO 71cm    | UP TO 69cm          | FPOK403101 | FPOK403201 |
| Knee PO XL 71cm | 69 - 89cm           | FPOK403103 | FPOK403203 |
| Knee PO 41cm    | 69cm                | -          | FPOK403300 |
| Knee PO XL 41cm | 69 - 89cm           | -          | FPOK403301 |

### SOFT GOODS

| Description                        | Cool   | Full Foam |
|------------------------------------|--------|-----------|
| Replacement Softgood, Regular      | 403175 | 403275    |
| Replacement Softgood, XL           | 403176 | 403277    |
| Replacement Softgood, Regular Long | 403276 | 403278    |

### REPLACEMENT SOFTGOODS

| Size                                  |           |
|---------------------------------------|-----------|
| ROM stop locks (set of 20 lock clips) | B-8380005 |
| Strap Reduction Kit Cool              | 403199    |
| Strap Reduction Kit Full Foam         | 403299    |

## REBOUND® KNEE BRACE

With six variations to choose from the Rebound hinged knee brace product line is a comprehensive portfolio of high-quality braces designed to serve the needs of a broad range of patients.

### INDICATIONS FOR USE

- ROM Hinge Models: Mild to moderate ACL and / or PCL, MCL and LCL instabilities, sprains and strains
- Non-ROM Hinge Models: Mild medial / lateral knee instabilities

### FEATURES AND BENEFITS

- Highly breathable aerospacer material reduces heat retention and improves comfort.
- Low profile hinges with extension and flexion stop kit options (ROM version) and malleable arms to provide a custom fit along with knee movement control.
- Contoured pre-tensed Cooltech™ fabric designed to reduce irritation and material bulging in the popliteal area.
- Removable straps allow application from right or left side, and strap pull-tab allows for consistent donning by patient.
- Universal buttresses are included with all sleeve versions for customizable patella support.

### REQUIRED MEASUREMENTS

- Circumference: 15cm above mid-patella



## ORDERING INFORMATION

### REBOUND ROM KNEE

| Size  | Circumference | Sleeve Long | Sleeve Short | Wrap Long | Wrap Short |
|-------|---------------|-------------|--------------|-----------|------------|
| XS    | 29 - 35cm     | 708152      | 708052       | 707152    | 707052     |
| S     | 35 - 41cm     | 708153      | 708053       | 707153    | 707053     |
| M     | 41 - 46cm     | 708155      | 708055       | 707155    | 707055     |
| L     | 46 - 52cm     | 708157      | 708057       | 707157    | 707057     |
| XL    | 52 - 57cm     | 708158      | 708058       | 707158    | 707058     |
| XXL   | 57 - 63cm     | 708159      | 708059       | 707159    | 707059     |
| XXXL  | 63 - 72cm     | 708160      | 708060       | 707160    | 707060     |
| XXXXL | 72 - 81cm     |             |              | 707161    | 707061     |

### REBOUND KNEE

| Size  | Circumference | Sleeve Short | Wrap Short |
|-------|---------------|--------------|------------|
| XS    | 29 - 35cm     | 704052       | 703052     |
| S     | 35 - 41cm     | 704053       | 703053     |
| M     | 41 - 46cm     | 704055       | 703055     |
| L     | 46 - 52cm     | 704057       | 703057     |
| XL    | 52 - 57cm     | 704058       | 703058     |
| XXL   | 57 - 63cm     | 704059       | 703059     |
| XXXL  | 63 - 72cm     | 704060       | 703060     |
| XXXXL | 72 - 81cm     |              | 703060     |



## REBOUND® KNEE IMMOBILIZER UNIVERSAL

Rebound® Knee Immobilizer Universal features a universal 3-panel design that can be adjusted to accommodate most leg circumferences. With an MR\* Conditional rating, patients can be safely scanned in an MR system, while maintaining immobilisation.

\*See IFU for details.



### INDICATIONS FOR USE

- Emergency immobilisation
- Pre- and post-surgical knee immobilisation
- Traumatic knee injuries

### FEATURES AND BENEFITS

- Adjustable medial/lateral stays and rigid posterior stays provide motion restriction of the knee
- Universal 3-piece design simplifies inventory needs
- MR\* Conditional: safely scan patients in an MR system while maintaining immobilisation
- Min thigh circumference: 26cm
- Max thigh circumference: 67.31cm

### ORDERING INFORMATION

| Length | Part#  |
|--------|--------|
| 30 cm  | 210012 |
| 35 cm  | 210014 |
| 40 cm  | 210016 |
| 45 cm  | 210018 |
| 50 cm  | 210020 |
| 55 cm  | 210022 |
| 60 cm  | 210024 |

### ACCESSORIES

| Description | Part#       |
|-------------|-------------|
| PCL Pad     | B-705000036 |



PCL Pad



## FORMFIT® PRO KNEE

Össur Formfit® Pro Knee is a compressive knee sleeve that features a SmartPad™ buttress for targeted support and cushioning. Our patented design realigns the patella, suitable for conditions such as Patella Tendonitis (Jumper's Knee) and Chondromalacia Patella. The dynamic, 3D MotionTech™ knit is supportive, meets certified medical grade compression requirements and has been proven to be superior to the leading brands at moisture wicking and breathability, keeping your knee cool during wear. Lightweight, removable spiral stays ensure optimum shape and fit, while distinct patella graphics make donning intuitive.



MotionTech™



CoolVent™



StayBilizer™



SmartPad™

### INDICATIONS FOR USE

For knee conditions that may benefit from compression, such as;

- Patella tendonitis
- Chondromalacia Patella
- Mild sprains and strains
- Knee pain
- Feeling of instability


### CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.


### FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit with superior moisture wicking and breathability, keeping your knee cool during wear
- Certified medical grade compression (20-30 mmHg)
- Lightweight CoolVent™ knit in popliteal area for added breathability and comfort
- Patented SmartPad™ for targeted support and cushioning
- Patella graphics for intuitive placement
- Patented FoldLock™ system for easier and more secure application
- Lightweight, removable spiral steel stays

### REQUIRED MEASUREMENTS



Measure calf circumference  
15 cm below mid-patella.  
Leg in extension.



Measure thigh circumference  
15 cm above mid-patella.  
Leg in extension.

| Size | Calf circumference 15 cm below mid patella | Thigh circumference 15 cm above mid patella |
|------|--------------------------------------------|---------------------------------------------|
| XS   | 26 - 31cm                                  | 36 - 41cm                                   |
| S    | 31 - 36cm                                  | 41 - 46cm                                   |
| M    | 36 - 41cm                                  | 46 - 51cm                                   |
| L    | 41 - 46cm                                  | 51 - 56cm                                   |
| XL   | 46 - 50cm                                  | 56 - 60cm                                   |

# FORMFIT® PRO KNEE

## ORDERING INFORMATION

### BLACK

| Size | Part#    | Size | Part#    |
|------|----------|------|----------|
| XS   | FP100001 | XS   | FP101001 |
| S    | FP100002 | S    | FP101002 |
| M    | FP100003 | M    | FP101003 |
| L    | FP100004 | L    | FP101004 |
| XL   | FP100005 | XL   | FP101005 |

### BLUE



## FORMFIT® TRACKER

The Formfit® Tracker with micro-adjustable ReflexWing™ buttress realigns the knee cap, providing dynamic lateral support for common knee cap conditions. Its FlyWeight™ fabric delivers lightweight breathability and a low-profile design that feels great against the skin and minimises migration. PowerLock™ straps work in tandem with heat mouldable CustomFit™ hinges to provide a customisable fit for your patient.

### INDICATIONS FOR USE

For knee conditions that may benefit from realignment of the knee cap and/or patients that would benefit from an added lateral support of the knee cap, such as:


- lateral knee cap subluxation or dislocation
- lateral maltracking

### FEATURES AND BENEFITS

- CustomFit™ hinges for customised fit and support
- Anti-Migration Tech™ minimises migration
- PowerLock™ straps for secure lateral support
- ReflexWing™ buttress for comfortable knee cap realignment
- FlyWeight™ fabric for soft feel and breathability
- Durable and machine washable



## REQUIRED MEASUREMENTS & ORDERING INFORMATION



Calf circumference 15 cm  
below mid patella



Thigh circumference 15 cm  
above mid patella

| Size | Calf circumference 15 cm below mid patella | Thigh circumference 15 cm above mid patella | Left   | Right  |
|------|--------------------------------------------|---------------------------------------------|--------|--------|
| XS   | 25 - 30cm                                  | 33 - 39cm                                   | 719162 | 719152 |
| S    | 30 - 35cm                                  | 39 - 47cm                                   | 719163 | 719153 |
| M    | 35 - 40cm                                  | 47 - 53cm                                   | 719164 | 719154 |
| L    | 40 - 46cm                                  | 53 - 60cm                                   | 719165 | 719155 |
| XL   | 46 - 51cm                                  | 60 - 67cm                                   | 719166 | 719156 |
| XXL  | 51 - 56cm                                  | 67 - 75cm                                   | 719167 | 719157 |
| 3XL  | 56 - 61cm                                  | 75 - 81cm                                   | 719168 | 719158 |

## ACCESSORIES

| Description          | Part#  |
|----------------------|--------|
| ReflexWing™ Buttress | 719100 |

## FORMFIT® PRO KNEE FLITE

Össur Formfit® Pro Knee Flite is a performance 3D knitted sleeve that is micro-engineered from our custom MotionTech™ yarn to provide breathability, moisture wicking, and a skintight fit. Our unique FoldLock™ and StayBilizer™ technologies combine to ensure the brace stays in place for maximum benefit.



MotionTech™



CoolVent™



StayBilizer™



FoldLock™



### INDICATIONS FOR USE

For knee conditions that may benefit from compression, such as:

- Mild sprains and strains
- Knee pain
- Feeling of instability


### FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit
- CoolVent™ knit in popliteal area for added breathability
- Patented FoldLock™ system for easy and secure application
- Lightweight, removeable spiral steel stays

### CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

### REQUIRED MEASUREMENTS



Measure calf circumference  
15 cm below mid-patella.  
Leg in extension.



Measure thigh circumference  
15 cm above mid-patella.  
Leg in extension.



| Size | Calf circumference<br>15 cm below mid patella | Thigh circumference<br>15 cm above mid patella | Blue     | Black    |
|------|-----------------------------------------------|------------------------------------------------|----------|----------|
| XS   | 26 - 31 cm                                    | 36 - 41 cm                                     | FP141001 | FP140001 |
| S    | 31 - 36 cm                                    | 41 - 46 cm                                     | FP141002 | FP140002 |
| M    | 36 - 41 cm                                    | 46 - 51 cm                                     | FP141003 | FP140003 |
| L    | 41 - 46 cm                                    | 51 - 56 cm                                     | FP141004 | FP140004 |
| XL   | 46 - 50 cm                                    | 56 - 60 cm                                     | FP141005 | FP140005 |

## FORMFIT® HINGED LATERAL J

Formfit Knee Hinged Lateral J is a sleeve design knee brace with a breathable Lycra-material and open popliteal area with two removeable hinges and J-shaped buttress providing support to the patella. The knee brace is designed for patients that would benefit from realignment of the patella, and/or patients that would benefit from an added lateral support of the patella.

### INDICATIONS FOR USE

Patients that would benefit from realignment of the patella, and/or patients that would benefit from an added lateral support of the patella.

- Lateral patellar subluxations and dislocations.
- Lateral patellofemoral malalignment.

### FEATURES AND BENEFITS

- Breathable material
- Removable hinges
- Easy pull-on
- Lateral J buttress
- Open popliteal area

### REQUIRED MEASUREMENTS

- Thigh Circumference: 15cm above mid-patella
- Calf Circumference: 15cm below mid-patella



## ORDERING INFORMATION

| Size | Thigh circumference | Knee circumference | Calf circumference | Right leg | Left leg |
|------|---------------------|--------------------|--------------------|-----------|----------|
| XS   | 33 - 39.4cm         | 30.5 - 33cm        | 25.4 - 30.5cm      | 709152    | 709162   |
| S    | 39.4 - 47cm         | 33 - 35.6cm        | 30.5 - 35.6cm      | 709154    | 709164   |
| M    | 47 - 53.3cm         | 35.6 - 38.1cm      | 35.6 - 40.6cm      | 709155    | 709165   |
| L    | 53.3 - 59.7cm       | 38.1 - 43.2cm      | 40.6 - 45.7cm      | 709157    | 709167   |
| XL   | 59.7 - 67.3cm       | 43.2 - 48.3cm      | 45.7 - 50.8cm      | 709158    | 709168   |
| XXL  | 67.3 - 74.9cm       | 48.3 - 53.3cm      | 50.8 - 55.9cm      | 709159    | 709169   |

## FORMFIT® KNEE HINGED

The Formfit Knee brace is for patients suffering from mild to moderate strains, sprains or ligament tears (ACL/PCL/MCL/LCL). It provides excellent stability, protection and compression during injury recovery. Made from strong and breathable materials it offers comfort while remaining highly durable. The open popliteal area reduces material bulging and potential irritation, while the malleable rigid hinge arms provide an intimate fit.

The Formfit Knee brace comes in either a wrap or a sleeve version.

### INDICATIONS FOR USE

- Mild to moderate ACL and PCL tears
- MCL and LCL tears
- Combined mild knee sprains and strains
- Instabilities of the knee

### FEATURES AND BENEFITS

- Breathable material
- Easy pull-on
- Malleable hinges
- Universal buttress
- Open popliteal area

### REQUIRED MEASUREMENTS

- Knee circumference at mid-patella



## ORDERING INFORMATION

| Size | Knee Circumference | Sleeve | Wrap   |
|------|--------------------|--------|--------|
| XS   | 33cm to 39cm       | 504252 |        |
| S    | 39cm to 47cm       | 504253 | 503253 |
| M    | 47cm to 53cm       | 504255 | 503255 |
| L    | 53cm to 60cm       | 504257 | 503257 |
| XL   | 60cm to 67cm       | 504258 | 503258 |
| XXL  | 67cm to 75cm       | 504259 | 503259 |
| XXXL | 75cm to 81cm       |        | 503260 |

## FORMFIT® KNEE MCL

The Formfit Knee MCL provides external stability and protection for low to medium contact/impact activities and while recovering from injury. A unique, single hinged, hybrid wrap/sleeve design sets this brace apart from the competition. The lateral, range of motion hinge keeps the medial side free of obstruction. The Dynamic Force Strap aids in the resistance of valgus forces during knee extension, while the adjustable tibial strap improves suspension by minimising slippage.

### INDICATIONS FOR USE

- Mild to Moderate Medial Collateral Ligament (MCL) sprains (Grade I and Grade II)

### FEATURES AND BENEFITS

- Single, steel ROM lateral hinge (0° or 20° extension stop, 70° flexion stop)
- Dynamic Force Strap resist valgus forces
- Hybrid wrap/sleeve for easy application
- Adjustable tibial strap for improved suspension



### REQUIRED MEASUREMENTS

- Thigh Circumference: 15cm above mid-patella

### ORDERING INFORMATION

| Size | Thigh Circumference | Formfit Knee MCL Left | Formfit Knee MCL Right |
|------|---------------------|-----------------------|------------------------|
| XS   | 33 - 39cm           | B-237600001           | B-237500001            |
| S    | 39 - 47cm           | B-237600002           | B-237500002            |
| M    | 47 - 53cm           | B-237600003           | B-237500003            |
| L    | 53 - 60cm           | B-237600004           | B-237500004            |
| XL   | 60 - 67cm           | B-237600005           | B-237500005            |



# FORMFIT® KNEE ROM BRACE

## INDICATIONS FOR USE

- Mild to moderate ACL and PCL tears
- MCL and LCL tears
- Combined mild knee sprains and strains
- Instabilities of the knee

## FEATURES AND BENEFITS

- With 6 options to choose from, you can confidently meet the unique needs of any patient
- When you need to limit range-of-motion, our ROM hinge provides the control you need. When you don't, you can opt for a more streamlined, polycentric hinge
- ROM braces available in 12" or 16" Length
- A removable, universal buttress is included with every sleeve version of the Formfit Knee Brace for additional, customisable patella support, if necessary



## REQUIRED MEASUREMENTS

- Knee circumference at mid-patella

## ORDERING INFORMATION

| Size | Knee Circumference | Sleeve Short | Wrap Short | Sleeve Long | Wrap Long |
|------|--------------------|--------------|------------|-------------|-----------|
| XS   | 33cm to 39cm       | 508252       | 507252     | 508352      | 507352    |
| S    | 39cm to 47cm       | 508253       | 507253     | 508353      | 507353    |
| M    | 47cm to 53cm       | 508255       | 507255     | 508355      | 507355    |
| L    | 53cm to 60cm       | 508257       | 507257     | 508357      | 507357    |
| XL   | 60cm to 67cm       | 508258       | 507258     | 508358      | 507358    |
| XXL  | 67cm to 75cm       | 508259       | 507259     | 508359      | 507359    |
| XXXL | 75cm to 81cm       |              | 507260     |             | 507360    |

## NEOPRENE HINGED KNEE SUPPORT

### INDICATIONS FOR USE

For patients suffering from knee conditions that may benefit from additional knee support, such as:

- Patellar tendonitis
- Subluxation dislocation
- Post-operation rehabilitation
- Post arthroscopic surgery
- Chondromalacia
- Bursitis
- Arthritis
- Tendonitis

### FEATURES AND BENEFITS

- Dual action hinges are moulded into the brace to maintain a low profile
- Accordion molding in the popliteal area reduces bunching behind the knee
- Lycra liner offers a breathable contact area for comfortable and extensive use
- Varied areas of thickness give additional protection

### REQUIRED MEASUREMENTS

- Knee circumference at mid-patella



### ORDERING INFORMATION

| Size | Knee circumference | Part#     |
|------|--------------------|-----------|
| S    | 29 - 33cm          | 302153BLK |
| M    | 33 - 37cm          | 302155BLK |
| L    | 37 - 41cm          | 302157BLK |
| XL   | 41 - 44cm          | 302158BLK |
| XXL  | 44 - 48cm          | 302159BLK |

## NEOPRENE KNEE SUPPORT WITH STABILISED PATELLA

### INDICATIONS FOR USE

For patients suffering from knee conditions that may benefit from additional knee and patella support, such as: patellar tendonitis, subluxation dislocation, post-operation rehabilitation, post arthroscopic surgery, chondromalacia, bursitis, arthritis and tendonitis.

### FEATURES AND BENEFITS

- Moulded material to maintain a low profile
- Accordion moulding in the popliteal area eliminates bunching

### REQUIRED MEASUREMENTS

- Knee circumference at mid-patella



### ORDERING INFORMATION

| Size | Knee circumference | Part#     |
|------|--------------------|-----------|
| S    | 11.5 - 13"         | 302563BLK |
| M    | 13 - 14.5"         | 302565BLK |
| L    | 14.5 - 16"         | 302567BLK |
| XL   | 16 - 17.5"         | 302568BLK |
| XXL  | 17.5 - 19"         | 302569BLK |

# NEOPRENE WRAPAROUND HINGED KNEE SUPPORT

## INDICATIONS FOR USE:

For patients with unstable or injured knees that may benefit from additional knee support, such as: post-operation rehabilitation, tendonitis, bursitis, arthritis, sprains and strains.

## FEATURES AND BENEFITS

- Dual action hinges are removable for easy cleaning
- Open popliteal area eliminates bunching and increases comfort
- Front contact closures allow for hassle-free application
- Wraparound style accommodates all leg shape variations

## REQUIRED MEASUREMENTS

- Knee circumference at mid-patella



## ORDERING INFORMATION

| Size  | Knee circumference | Part#  |
|-------|--------------------|--------|
| S     | 33 - 35.5cm        | 212003 |
| M     | 35.5 - 38cm        | 212005 |
| L     | 38 - 40.5cm        | 212007 |
| XL    | 40.5 - 43cm        | 212008 |
| XXL   | 43 - 45.5cm        | 212009 |
| XXXL  | 45.5 - 58.5cm      | 212010 |
| XXXXL | 58.5 - 63.5cm      | 212011 |

# NEOPRENE KNEE SLEEVE

A neoprene sleeve-design knee support providing compression and warmth

## INDICATIONS FOR USE

- For patients suffering from weak and/or injured knees that may benefit from additional knee support.

## FEATURES AND BENEFITS

- Unique accordion design enhances comfort by eliminating bunching in the popliteal area
- Open or closed patella

## REQUIRED MEASUREMENTS

- Knee circumference at mid-patella



## ORDERING INFORMATION

|      |               | 0.3cm Neoprene |                | 0.6cm Neoprene |                |
|------|---------------|----------------|----------------|----------------|----------------|
| Size | Circumference | Open Patella   | Closed Patella | Open Patella   | Closed Patella |
| S    | 29 - 33cm     | 402013BLK      | 412013BLK      | 302013BLK      | 312013BLK      |
| M    | 33 - 37cm     | 402015BLK      | 412015BLK      | 302015BLK      | 312015BLK      |
| L    | 37 - 41cm     | 402017BLK      | 412017BLK      | 302017BLK      | 312017BLK      |
| XL   | 41 - 44cm     | 402018BLK      | 412018BLK      | 302018BLK      | 312018BLK      |
| XXL  | 44 - 48cm     | 402019BLK      | 412019BLK      | 302019BLK      | 312019BLK      |

# MEASUREMENT FORM CTi® CUSTOM BRACE

## CUSTOM MEASUREMENT SYSTEM

Date: \_\_\_\_\_ Customer #: \_\_\_\_\_ Purchase order #: \_\_\_\_\_

### Required patient information

Patient ID: \_\_\_\_\_ Diagnosis: \_\_\_\_\_

Height: \_\_\_\_\_ cm. Weight: \_\_\_\_\_ kg. Date of Birth: \_\_\_\_\_ Sex(m/f): \_\_\_\_\_ ☐ Left leg ☐ Right leg

☐ Bilats - l/r

### Customer:

Contact: \_\_\_\_\_ Prescribing consultant: \_\_\_\_\_

Phone: \_\_\_\_\_ Clinic: \_\_\_\_\_

Company: \_\_\_\_\_

Please fill out all mandatory fields  
(marked in bold)

#### ☐ VAPOR

(For low to medium impact levels)

#### ☐ STANDARD\*

(For medium to high impact levels)

#### ☐ PRO SPORT

(For highest impact levels)

### OPTIONS

#### D-ring attachment

☐ Screwed ☐ Riveted\*

☐ Internal ☐ External\*

#### Thigh Liner

☐ None\* ☐ Sensil

☐ Evazote/Foam

#### Calf Liner

☐ Sensil\* ☐ Evazote/Foam

#### Extension Stops

☐ 0° ☐ 10°\* ☐ 20° ☐ 30° ☐ 40°

### SPECIAL MODIFICATIONS

(Additional charges may apply)

#### ☐ Flexion Stop Kit

\_\_\_\_\_° (0-90° installed)

#### ☐ PCL System

(Helps reduce posterior drawer)

☐ ACL Cable Kit (Additional stability in ext)

☐ Hyperextension Straps (5th strap added)

☐ Ski Boot Attachment (Brace-to-boot)

☐ Other: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

### SPECIAL NOTES:

☐ Please contact me about this order

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\*This configuration is delivered by default, if nothing is filled in.

### SPECIAL CONDITIONS:

☐ OA (Osteoarthritis unloading)

☐ Medial

☐ Lateral \_\_\_\_\_°

(Choose 1° to 7° of unloading, 4°\*)

☐ Prominent VMO

☐ Prominent VLO

☐ Hockey Quads

☐ Prominent Tibialis Anterior

☐ Prominent Fibular Head

☐ RTM

(Rounded Tibial Member)

☐ Osgood-Schlatters Disease

☐ Super Short

(For 1.60m and under)

☐ Shorten Frame Component

Tibial ☐ 1.3cm ☐ 2.6cm ☐ 3.9cm

Femoral ☐ 1.3cm ☐ 2.6cm ☐ 3.9cm

### COLOURS:

Select finish:

☐ MATTE\* ☐ GLOSS

Select colour:

☐ Black\* ☐ Yellow

☐ Ocean Blue ☐ Orange

☐ White ☐ Lime Green

☐ Silver ☐ Pink

☐ Red ☐ Charcoal

☐ Champagne ☐ Navy Blue

### ACCESSORIES:

(Bolded are free with initial order only)

#### ☐ AMS WRAP

☐ Neoprene\* ☐ Evazote

#### ☐ Sofsleeve

☐ Neoprene Undersleeve

☐ Neoprene Oversleeve

☐ Padded Sports Oversleeve

☐ Sleeve with Mesh

☐ MX Patella Cup

☐ MX Gear Guards

☐ MX Complete Kit  
(Patella Cup, Gear Guards,  
Sleeve with Mesh)

☐ Other: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

### CUSTOM PAINT:

(Additional fee plus extra delivery time)

Select custom finish:

☐ MATTE\* ☐ GLOSS

Select custom colour:

☐ Hibiscus ☐ Sports Equipment

☐ Flag ☐ Camouflage

☐ Flames ☐ Single Colour

☐ Metal ☐ 2-colour Fade

☐ Sponge ☐ 2-colour Half & Half

☐ Celestial ☐ Sports Fan

☐ Special

(Provide design file with order)

# MEASUREMENT FORM CTi® CUSTOM BRACE

## CUSTOM MEASUREMENT SYSTEM

### CUSTOM MEASURING SYSTEM

☐ Left leg ☐ Right leg

1: \_\_\_\_\_ Hinge Width

2: \_\_\_\_\_ Hinge Depth

**LATERAL MEASUREMENTS**

3: \_\_\_\_\_ Contour

4: \_\_\_\_\_ A-P

5: \_\_\_\_\_ Contour

6: \_\_\_\_\_ A-P

7: \_\_\_\_\_ Contour

8: \_\_\_\_\_ A-P

9: \_\_\_\_\_ Tibia

10: \_\_\_\_\_ Tibia

11: \_\_\_\_\_ Contour

12: \_\_\_\_\_ A-P

13: \_\_\_\_\_ Tibia

**MEDIAL CONTOURS**

14: \_\_\_\_\_ Contour

15: \_\_\_\_\_ Contour

**JOINT SPACE**

16: \_\_\_\_\_ Contour

17: \_\_\_\_\_ Contour

**CIRCUMFERENCES (in centimeters)**

18. Level A \_\_\_\_\_

19. Level B \_\_\_\_\_

20. Level C \_\_\_\_\_

21. Level E \_\_\_\_\_

22. Distal Border of Patella \_\_\_\_\_

23.

Caliper measurement at joint space (1.6 - 2.0 cm less than #1)  
Take measurement in a weight bearing position at full extension

**SHORTEN TIBIA COMPONENT/ SKI BOOT MEASUREMENT**

\_\_\_\_\_ cm from the distal border of the patella to the top of the boot. Ski boot attachment is not available if the measurement is less than 18cm.

### CMS INSTRUCTIONS

Patient should not be wearing shoes for the measurement procedure.

There should be no space between the block and the back of the knee. The system will not be accurate if the knee is measured in a flexed position.

1. Slide the lateral pad to the knee. Tighten to patient's individual comfort level. This measurement will be the actual width of the patient's knee brace.

2. Pre-set hinge depth at approximately 3.5cm, adjust if needed.

3. – 13. Patient should remain still with foot dorsiflexed throughout the lateral measurements and recording process.

9., 10. and 13. Measure using the horizontal bar of the tool measure the position of the highest point of the tibialis anterior.

• 9. As reference take the #1 hinge width / 2 - 1 = ...  
(for example 13cm / 2 - 1 = 5.5cm)

• 10. should be 0.5cm less than 9.

• 13. should be 1cm less than 9.


14. – 17. The patient may now relax throughout measuring the medial contours from the contracted position. Please note that with larger legs this measurement may be a negative value. Please read the numbers carefully and record negative values if necessary, such as -0.5 or -1.5

18. – 22. Take the measurements with a soft measuring tape ensuring that the tape is snug but not too loose/tight.


# MEASUREMENT FORM CTi® CUSTOM BRACE

CUSTOM MEASUREMENT SYSTEM


## SPECIAL CONDITIONS REFERENCES




Prominent VLO




Prominent VMO




Hockey Quads




Prominent Tibialis Anterior




Prominent Fibular Head



RTM  
(Rounded Tibial Member)




Osgood-Schlatters Disease



Shorten Femoral Component  
Shorten Tibial Component



Super Short



Osteoarthritis Unloading

To find out more about CTi custom, please scan the QR code



# MEASUREMENT FORM UNLOADER ONE® CUSTOM / REBOUND® ACL & PCL

FOR USE WITH DMS KIT

## PATIENT INFORMATION:


|                |                          |          |  |
|----------------|--------------------------|----------|--|
| Name:          |                          |          |  |
| Date of Birth: | Contact No.              |          |  |
| Height:        |                          | Weight : |  |
| Diagnosis:     | Prescribing consultant : |          |  |

## MEASUREMENTS:

| M-L Caliper               | Left | Right |
|---------------------------|------|-------|
| Circumferences            |      |       |
| 20cm Above mid patella:   |      |       |
| 15cm Above mid patella:   |      |       |
| 7.5cm Above mid patella : |      |       |
| 7.5cm Below mid patella:  |      |       |
| 15cm Below mid patella:   |      |       |
| 20cm Below mid patella:   |      |       |

## BRACE DETAILS:

|                           |        |                          |         |                          |
|---------------------------|--------|--------------------------|---------|--------------------------|
| Left / Right / Bilateral: |        |                          |         |                          |
| Unloader® One:            | Medial | <input type="checkbox"/> | Lateral | <input type="checkbox"/> |
| Unloader® Plus:           | Medial | <input type="checkbox"/> | Lateral | <input type="checkbox"/> |
| Rebound®:                 | ACL    | <input type="checkbox"/> | PCL     | <input type="checkbox"/> |
| Colour Options:           |        |                          |         |                          |
| Notes:                    |        |                          |         |                          |



## FILE NAME:

|       | Anterior | Lateral |
|-------|----------|---------|
| Left  |          |         |
| Right |          |         |

**Product Warranty** - Custom Made Knee Braces – Össur warrants the hinge components to be free of defects in material and workmanship for the lifetime of the patient. With regard to components of the brace that wear with use, which include straps, liners, covers, sleeves, return cords and hinge coatings, these components shall be warranted free from defect in material and workmanship for 6 months. Prefabricated Knee Braces – Össur warrants the frame and hinge components to be free of defects in material and workmanship for two years. With regard to components of the brace that wear with use, including straps, liners, covers, sleeves, return cords and hinge coatings, these components shall be warranted free from defect in material and workmanship for 6 months. Return Policy - Upon receipt of the Össur knee brace, the provider shall inspect it to ensure that it meets the requirements set forth in the order. If the brace does not meet the order requirements, the provider must call the Össur Customer Care Department at 03450 065 065, Monday to Friday, 8:30am to 5:00pm. Prior to returning a brace for any reason, including product warranty or patient satisfaction issues, the provider is required to obtain a Return Authorisation (RA) number from Customer Care. Final acceptance of all returned items remains the discretion of the Össur Quality Control Department. Satisfaction guarantee and product warranty time frames begin on the date the brace is delivered to the patient. Proof of brace delivery, signed and dated by the patient, is required for processing of guarantee and warranty returns. Satisfaction guarantee does not apply to the Össur post-operative braces. The stabilising effect of a knee brace will depend on the type of brace, patient age, height, weight, muscle tone, severity of the injury/OA, the type of injury and the medical treatment involved. A knee brace is designed to reduce the potential and severity of a knee injury and is not a guarantee against injury. There are no warranties expressed or implied made by Össur on their braces except the warranties against defects in materials or workmanship. Liability is limited to replacement or repair of defective parts or workmanship. Liability for incidental or consequential damage, personal injury or property damage as a result of defect is specifically denied. Össur neither assumes nor authorises any other person to assume for it any other liability for its braces.



## FOOT AND ANKLE



## TABLE OF CONTENTS

---

|                                        |            |
|----------------------------------------|------------|
| <b>FOOT AND ANKLE.....</b>             | <b>131</b> |
| REBOUND® DIABETIC WALKER .....         | 133        |
| DH OFFLOADING WALKER® .....            | 135        |
| REBOUND® AIR WALKER .....              | 136        |
| FORMFIT® WALKER AIR .....              | 138        |
| FORMFIT® WALKER .....                  | 139        |
| PAEDIATRIC WALKER .....                | 140        |
| AFO DYNAMIC .....                      | 141        |
| AFO LIGHT .....                        | 143        |
| PROLITE AFO .....                      | 144        |
| EXTRA STRONG PROLITE .....             | 144        |
| REBOUND FOOT-UP® .....                 | 145        |
| FOOT-UP® .....                         | 146        |
| FORMFIT® PRO ANKLE .....               | 147        |
| FORMFIT® ANKLE BRACE .....             | 149        |
| FORMFIT® ANKLE WITH SPEEDLACE .....    | 150        |
| FORMFIT® ANKLE STIRRUPS .....          | 151        |
| EXOFORM® DORSAL NIGHT SPLINT .....     | 152        |
| FORMFIT® NIGHT SPLINT .....            | 153        |
| DH OFFLOADING POST-OP STYLE SHOE ..... | 154        |
| MESH TOP POST-OP SHOE .....            | 155        |
| SOFT TOP POST-OP SHOE .....            | 155        |
| CANVAS ROCKER BOTTOM CAST SHOE .....   | 156        |

## REBOUND® DIABETIC WALKER

Offloading has proven to be a key factor in foot ulcer treatment and entails reducing pressure of the plantar area in general and the ulcer site in particular.

The Rebound Diabetic Walker is uniquely designed for effective treatment of plantar foot ulcers, with easy access and re-fitting during wound care, while improving the patient's mobility during the healing process.

The Rebound Diabetic Walker has major advantages over diabetic shoes and shoe inserts in its healing effectiveness and over total contact casting due to its ease of use, speed of application and low cost.



### INDICATIONS FOR USE

- Plantar foot ulcers and wound prevention

### FEATURES AND BENEFITS

- Effective tri-laminate insole with removable shock-absorbing hexagons that promote ulcer offloading and a Plastazote, heat-moldable top layer that evenly distributes the pressure and minimises friction
- Unique pressure relief valve for preventing over-inflation, while balancing compression with customised fit
- Positioning Portals allow the clinician to check the fit of the walker at the heel bottom, in order to ensure appropriate positioning of the foot in the walker
- Anatomically designed shell that minimises shear forces by immobilising the foot and ankle at 90 degrees while providing protection
- Flex Edge™ technology for relief of pressure points
- Rocker Bottom outsole for a smooth roll over while allowing for individual natural toe out
- Soft and breathable, low friction, machine washable liner for hygiene promotion and prevention of infection
- Tamper evident strapping for compliance monitoring

### REQUIRED MEASUREMENTS

- Shoe size



# REBOUND® DIABETIC WALKER

## ORDERING INFORMATION

### REBOUND DIABETIC WALKER

| Size | Shoe size UK | Shoe size EU | Part#       |
|------|--------------|--------------|-------------|
| S    | 4 - 6.5      | 36 - 39      | B-242500002 |
| M    | 7 - 10       | 39.5 - 44    | B-242500003 |
| L    | 10.5 - 12    | 44.5 - 46.5  | B-242500004 |
| XL   | 12.5 - 14.5  | 47 - 50      | B-242500005 |

## ACCESSORIES

### REPLACEMENT SOFTGOODS

| Size | Liners      | Insoles     |
|------|-------------|-------------|
| S    | B-418500102 | B-418500202 |
| M    | B-418500103 | B-418500203 |
| L    | B-418500104 | B-418500204 |
| XL   | B-418500105 | B-418500205 |



### WALKER COVERS

| Size | 1. Hygiene covers (2 pcs) | 2. Durable weather cover | 3. Disposable weather cover (20 pcs) |
|------|---------------------------|--------------------------|--------------------------------------|
| S    | B-418500602               | B-418500502              | B-418500512                          |
| M    | B-418500603               | B-418500503              |                                      |
| L    | B-418500604               | B-418500504              | B-418500514                          |
| XL   | B-418500605               | B-418500505              |                                      |



### DIABETIC SOCKS (2 PCS)

| Size  | Part#       |
|-------|-------------|
| S, M  | B-418500703 |
| L, XL | B-418500704 |

### EVENUP

| Size | Part# |
|------|-------|
| S    | 10350 |
| M    | 10351 |
| L    | 10352 |



## DH OFFLOADING WALKER®

Assists in the healing of plantar ulcers with a patented pressure relief insole

### INDICATIONS FOR USE

- Treatment of plantar ulcers

### FEATURES AND BENEFITS

- Shear force reduction by stabilising foot and ankle at 90°
- Wide low profile rocker bottom provides a smooth roll-through and increased stability for altered gaits
- Designed around Össur's unique Active Offloading Insole
- The softgood's patented inner bootie prevents the foot from migrating on the insole
- Provides a safe environment for the simultaneous use of wound care material
- Create an Instant TCC using a single roll of cast tape for patient compliance



### REQUIRED MEASUREMENTS

- Shoe size

## ORDERING INFORMATION

| Size | Shoe Size UK | Shoe Size EU | Part#     |
|------|--------------|--------------|-----------|
| XS   | 1.5 - 3.5    | 33 - 35.5    | DH0200BLK |
| S    | 4 - 6.5      | 36 - 39      | DH0400BLK |
| M    | 7 - 10       | 39.5 - 44    | DH0600BLK |
| L    | 10.5 - 12    | 44.5 - 46.5  | DH0800BLK |
| XL   | 12.5 - 14.5  | 47 - 50      | DH1000BLK |

## ACCESSORIES

### REPLACEMENT PARTS

| Size | Softgoods | Insoles |
|------|-----------|---------|
| XS   | 10958     | 11742C  |
| S    | 10958     | PP04C   |
| M    | 10959     | PP06C   |
| L    | 10960     | 11740C  |
| XL   | 10960     | 11741C  |



## REBOUND® AIR WALKER

Rebound Air Walker is a circumferential walker with fully integrated air pump and release valve that provides easy, tool-free inflation and deflation for fine-tuned, consistent compression. It features a host of details designed to provide added protection, simplified application and usage.



### INDICATIONS FOR USE

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- Tendon or ligamentous sprain
- During postoperative or rehab phase. If used with the Achilles Wedge Kit, Achilles ruptures can be treated as well.

- Rebound Air Walker Low Top
- Foot and ankle conditions that may benefit from motion restriction, such as:
- Stable foot fractures
- Tendon or ligamentous sprain
- During postoperative or rehab phase.

### FEATURES AND BENEFITS

- Linked medial and lateral air bladders with a single inflation point provide balanced compression and are shaped to cradle the Achilles, holding the heel in place
- Rocker bottom promotes a natural, stable gait, while a trimmable, semi-soft toe bumper provides added protection
- Moisture-wicking softgood helps keep skin dry and comfortable
- Straps fasten in front for easy application, adjustment and removal
- EVA material provides greater shock absorption in the heel
- Fully-integrated air pump and release valve provides easy, tool-free inflation and deflation for fine-tuned, consistent compression
- Flex Edge™ overmolding conforms to each patient's unique anatomy to help limit the formation of pressure points
- Ventilated panels wick moisture to prevent microbial buildup and improve patient comfort

### REQUIRED MEASUREMENTS

- Shoe size



# REBOUND® AIR WALKER

## ORDERING INFORMATION

### REBOUND AIR WALKER

| Size | Shoe Size UK | Shoe Size EU | Standard Liner / Standard | Standard Liner / Low Top | Slim Liner / Standard | Slim Liner / Low Top |
|------|--------------|--------------|---------------------------|--------------------------|-----------------------|----------------------|
| XS   | 1.5 - 3.5    | 33 - 35.5    | B-242900001E              | B-242900061E             | B-242900001SE         | B-242900061SE        |
| S    | 4 - 6.5      | 36 - 39      | B-242900002E              | B-242900062E             | B-242900002SE         | B-242900062SE        |
| M    | 7 - 10       | 39.5 - 44    | B-242900003E              | B-242900063E             | B-242900003SE         | B-242900063SE        |
| L    | 10.5 - 12    | 44.5 - 46.5  | B-242900004E              | B-242900064E             | B-242900004SE         | B-242900064SE        |
| XL   | 12.5 - 14.5  | 47 - 50      | B-242900005E              | B-242900065E             | B-242900005SE         | B-242900065SE        |

## ACCESSORIES

### REPLACEMENT SOFTGOODS

| Size | Standard Liner / Standard | Standard Liner / Low Top | Slim Liner / Standard | Slim Liner / Low Top |
|------|---------------------------|--------------------------|-----------------------|----------------------|
| XS   | B-418000101               | B-418000161              | B-418030101           | B-418030161          |
| S    | B-418000102               | B-418000162              | B-418030102           | B-418030162          |
| M    | B-418000103               | B-418000163              | B-418030103           | B-418030163          |
| L    | B-418000104               | B-418000164              | B-418030104           | B-418030164          |
| XL   | B-418000105               | B-418000165              | B-418030105           | B-418030165          |

### ACHILLES HEEL WEDGE KIT

| Size | Description                | Part#    |
|------|----------------------------|----------|
| S    | Wedge kit for walkers XS-S | AWK0001  |
| M    | Wedge Kit for walkers M    | AWK0005  |
| L    | Wedge kit for walkers L-XL | AWK0009  |
| N/A  | 5 layer Wedge kit          | AWK00015 |

### EVENUP

| Size | Part# |
|------|-------|
| S    | 10350 |
| M    | 10351 |
| L    | 10352 |

### WALKER COVERS

| Size | Hygiene covers (2 pcs) | Durable weather cover | Disposable weather cover (20 pcs) |
|------|------------------------|-----------------------|-----------------------------------|
| S    | B-418500602            | B-418500502           | B-418500512                       |
| M    | B-418500603            | B-418500503           |                                   |
| L    | B-418500604            | B-418500504           | B-418500514                       |
| XL   | B-418500605            | B-418500505           |                                   |

### OTHER ACCESSORIES

| Description                      | Part#       |
|----------------------------------|-------------|
| Adhesive Insole for Wedge        | WHW-W0107   |
| Strap Extension                  | B-502100101 |
| Rebound Hightop Sock SM/MD SM/MD | B-418500703 |
| Rebound Hightop Sock LG/XL LG/XL | B-418500704 |
| Rebound Lowtop Sock SM/MD SM/MD  | B-418500903 |
| Rebound Lowtop Sock LG/XL LG/XL  | B-418500904 |

## FORMFIT® WALKER AIR

The Formfit Walker Air provides immobilisation for patients suffering from soft tissue injury; grades 2 & 3 strains and sprains, and stable fractures. The patented pneumatic system is designed to increase stability while decreasing pain and swelling. The intuitive pump and clearly marked decompression valve make air/pressure adjustments easy.

### INDICATIONS FOR USE

For foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- Soft tissue injuries
- During postoperative or rehab phase

### FEATURES AND BENEFITS

- Available in low top or standard height
- Integrated pneumatic system
- Low, wide rocker bottom promotes a stable, natural gait
- Shock absorbing sole for reduced impact at heel strike
- Breathable, seamless liner

### REQUIRED MEASUREMENTS

- Shoe size



## ORDERING INFORMATION

| Size | Shoe Size UK | Shoe Size EU | Standard   | Low Top    |
|------|--------------|--------------|------------|------------|
| XS   | 1.5 - 3.5    | 33 - 35.5    | A-W0200BLK | A-W0300BLK |
| S    | 4 - 6.5      | 36 - 39      | A-W0400BLK | A-W0500BLK |
| M    | 7 - 10       | 39.5 - 44    | A-W0600BLK | A-W0700BLK |
| L    | 10.5 - 12    | 44.5 - 46.5  | A-W0800BLK | A-W0900BLK |
| XL   | 12.5 - 14.5  | 47 - 50      | A-W1000BLK | A-W1100BLK |

## ACCESSORIES

### REPLACEMENT SOFTGOODS

| Size | Standard | Low Top |
|------|----------|---------|
| XS   | 12349    | 12356   |
| S    | 10675    | 12357   |
| M    | 10678    | 12358   |
| L    | 10681    | 12359   |
| XL   | 12347    | 12360   |

### OTHER ACCESSORIES

| Description                                                         | Part#     |
|---------------------------------------------------------------------|-----------|
| 10cm softgood width extension                                       | WKSG-1BLK |
| Tibia and calf plates for use when additional splinting is required | WKSP100   |
| Hot/Cold therapy insert for safe and reusable cryotherapy           | WKGEL100  |
| 5 layer Wedge kit                                                   | AWK0005   |
| Adhesive insole for heel wedge                                      | WHW-W0107 |



## FORMFIT® WALKER

The Formfit Walker provides immobilisation for patients suffering from soft tissue injury; grades 2 & 3 strains and sprains; and stable fractures.

### INDICATIONS FOR USE

For foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- Soft tissue injuries
- During postoperative or rehab phase

### FEATURES AND BENEFITS

- Available in low top or standard height
- Low, wide rocker bottom promotes a stable, natural gait
- Shock absorbing sole for reduced impact at heel strike
- Breathable, seamless liner

### REQUIRED MEASUREMENTS

- Shoe size



## ORDERING INFORMATION

| Size | Shoe Size UK | Shoe Size EU | Standard | Low Top  |
|------|--------------|--------------|----------|----------|
| XS   | 1.5 - 3.5    | 33 - 35.5    | W0200BLK | W0300BLK |
| S    | 4 - 6.5      | 36 - 39      | W0400BLK | W0500BLK |
| M    | 7 - 10       | 40 - 44      | W0600BLK | W0700BLK |
| L    | 10.5 - 12    | 44.5 - 46.5  | W0800BLK | W0900BLK |
| XL   | 12.5 - 14.5  | 47 - 50      | W1000BLK | W1100BLK |

## ACCESSORIES

### REPLACEMENT SOFTGOODS

| Size | Standard    | Low Top     |
|------|-------------|-------------|
| XS   | WKSG200BLK  | WKSG0300BLK |
| S    | WKSG400BLK  | WKSG500BLK  |
| M    | WKSG600BLK  | WKSG700BLK  |
| L    | WKSG800BLK  | WKSG900BLK  |
| XL   | WKSG1000BLK | WKSG1100BLK |

### OTHER ACCESSORIES

| Description                                                         | Part#     |
|---------------------------------------------------------------------|-----------|
| 10cm softgood width extension                                       | WKSG-1BLK |
| Tibia and calf plates for use when additional splinting is required | WKSP100   |
| Hot/Cold therapy insert for safe and reusable cryotherapy           | WKGEL100  |
| 5 layer Wedge kit                                                   | AWK0005   |
| Adhesive insole for heel wedge                                      | WHW-W0107 |

## PAEDIATRIC WALKER

Paediatric Walkers are designed to stabilise and protect the lower leg, ankle and foot. They are easy to operate and comfortable to wear, making them ideal for an active child who requires stability following a sprain, stress fracture, or surgery.

### INDICATIONS FOR USE

- Soft tissue injuries
- Grade 2 and 3 sprains
- Stable fractures
- Post-operative ankle stabilisation
- Trauma and rehab use

### FEATURES AND BENEFITS

- High quality rigid outer skeleton for protection
- Dense foam lining within the sole to provide comfort and support
- Shock absorbing sole to improve comfort and protects the user from slips, providing additional security in each step
- Uprights give additional stability by restricting motion of the leg
- Lightweight design perfectly designed for kids to keep active

### REQUIRED MEASUREMENTS

- Shoe size



## ORDERING INFORMATION

### PAEDIATRIC WALKER

| Size | Shoe Size          | Part#   |
|------|--------------------|---------|
| S    | Child 6.5 - 9.5    | PEW0400 |
| M    | Child 10 - Youth 1 | PEW0600 |
| L    | Adult 1.5 - 2.5    | PEW0800 |

### REPLACEMENT LINERS

| Size | Shoe Size        | Part#   |
|------|------------------|---------|
| S    | Child 6.5 - 9.5  | PEW0001 |
| M    | Child 10 - Youth | PEW0002 |
| L    | Adult 1.5 - 2.5  | PEW0003 |

## AFO DYNAMIC

The AFO Dynamic provides a strong, lightweight solution for patients with varying degrees of drop foot secondary to a variety of neurological conditions. Össur has taken the carbon fibre technology from its Flex-Foot® product and incorporated it into the new AFO Dynamic to provide increased support over traditional ankle foot orthoses. The energy storing properties of carbon fibre make it an ideal material for this dynamic ankle foot orthosis because it provides real energy return. For those people leading busy lives, the AFO Dynamic will ensure a more efficient walk for the user.



### INDICATIONS FOR USE

Drop foot secondary to:

- Cerebro Vascular Accident (CVA)
- Other neurological conditions
- Multiple Sclerosis
- Mild knee instability
- Charcot Marie Tooth disease (HMSN)

Contraindications:


- Severe ankle - foot deformities
- Fluctuating oedema
- Severe spasticity
- Neuropathy

### FEATURES AND BENEFITS

- Carbon fibre design provides excellent strength-to-weight ratio and dynamic motion
- The tapered carbon heel stores energy and manages forces at heel strike
- Full-length toe lever provides full support while foot is loaded, as well as reduced shock to the contralateral limb
- Creates medio-lateral stability at mid-stance
- Reduces energy consumption by reducing hip hike
- Medial strut makes it less visible in shoe
- Open anterior tibia shell provides superior anatomical fit, and reduces heat and perspiration
- Weight limit: 120 kg

### REQUIRED MEASUREMENTS

- Shoe size
- Calf height at widest part of the calf



# AFO DYNAMIC

## ORDERING INFORMATION

### AFO DYNAMIC

|      |             | Shoe Size UK |           | Shoe Size EU |           |          |          |
|------|-------------|--------------|-----------|--------------|-----------|----------|----------|
| Size | Calf Height | Male         | Female    | Male         | Female    | Left     | Right    |
| XS   | 30 - 34cm   | 2.5 - 4.5    | < 4       | 34 - 37      | < 37      | O-102021 | O-102011 |
| S    | 32 - 36cm   | 5 - 7        | 4.5 - 6.5 | 37.5 - 40    | 37.5 - 40 | O-102022 | O-102012 |
| M    | 34 - 38cm   | 7.5 - 10     | 7 - 9.5   | 40.5 - 44.5  | 40.5 - 44 | O-102024 | O-102014 |
| L    | 36 - 40cm   | 10.5 - 13.5  |           | 45 - 48.5    |           | O-102026 | O-102016 |

## ACCESSORIES

### ACCESSORY KITS

| Size | Part#    |
|------|----------|
| XS   | O-102101 |
| S    | O-102102 |
| M    | O-102104 |
| L    | O-102106 |

Note: Accessory kits come with foam padding, upper and lower Hook & Loop straps, and Hook & Loop dots.

## AFO LIGHT

AFO Light provides support only where it is needed, eliminating unnecessary bulk and providing a cosmetically-appealing look. Using the best of Flex-Foot® technology and the energy storing properties of carbon fibre, AFO Light offers a strong and lightweight solution for people with mild to moderate drop foot.



### INDICATIONS FOR USE

- |                                   |                                      |                                   |
|-----------------------------------|--------------------------------------|-----------------------------------|
| Drop foot secondary to:           | • Charcot Marie Tooth disease (HMSN) | Contraindications:                |
| • Cerebro Vascular Accident (CVA) | • Lower extremity overuse injuries   | • Severe ankle - foot deformities |
| • Other neurological conditions   |                                      | • Fluctuating edema               |
| • Multiple Sclerosis              |                                      | • Severe spasticity               |
| • Mild knee instability           |                                      | • Neuropathy                      |

### FEATURES AND BENEFITS

- Energy return and propulsion assistance of the limb at toe off
- Supports the foot for ground clearance during swing phase
- Carbon fibre design provides excellent strength-to-weight ratio and dynamic motion
- Gradual plantar flexion at initial contact
- Stiffness gradient at heel and the toe correspond to ground reaction forces as they are applied to the foot
- Easy to don and doff
- Weight limit: 100kg

### REQUIRED MEASUREMENTS

- Shoe size
- Calf height at widest part of the calf

### ORDERING INFORMATION

#### AFO LIGHT

|      |             | Shoe Size UK |           | Shoe Size EU |           |          |          |
|------|-------------|--------------|-----------|--------------|-----------|----------|----------|
| Size | Calf height | Male         | Female    | Male         | Female    | Left     | Right    |
| S    | 32 - 36cm   | < 7          | 4.5 - 6.5 | < 40         | 37.5 - 40 | O-121022 | O-121012 |
| M    | 34 - 38cm   | 7.5 - 10     | 7 - 9.5   | 40.5 - 44.5  | 40.5 - 44 | O-121024 | O-121014 |
| L    | 36 - 40cm   | 10.5 - 13.5  |           | 45 - 48.5    |           | O-121026 | O-121016 |
| XL   | 38 - 42cm   | 11.5 - 14.5  |           | 46 - 50      |           | O-121028 | O-121018 |

### ACCESSORIES

| Description | Part#    |
|-------------|----------|
| Refurb Kit  | O-122101 |

## PROLITE AFO

Comfortable, low profile prefabricated leaf spring design to correct flaccid foot drop, especially on post-stroke patients. Injection moulded polypropylene allows for different thicknesses throughout the orthosis — thicker on the vertical aspect for rigidity, thinner footplate for ease of trimming, comfort and flexibility.

### INDICATIONS FOR USE

- Drop foot
- Mild drop foot secondary to other neurological pathologies

### FEATURES AND BENEFITS

- Injection-moulded for a high quality finish
- Easily trimmed footplate
- May be heat modified for minor fitting adjustments

### REQUIRED MEASUREMENTS

- Foot length



## ORDERING INFORMATION

| Size | Foot length | Brace height | Left   | Right  |
|------|-------------|--------------|--------|--------|
| S    | 24cm        | 26.5cm       | 3500OP | 3505OP |
| M    | 28cm        | 28.5cm       | 3501OP | 3506OP |
| L    | 28cm        | 33cm         | 3502OP | 3507OP |
| XL   | 30cm        | 39cm         | 3503OP | 3508OP |

## EXTRA STRONG PROLITE

### INDICATIONS FOR USE

- Drop foot
- Mild drop foot secondary to other neurological pathologies

### FEATURES AND BENEFITS

- Made from a specific acetal resin, polyoxymethylene, ensuring no moisture absorption, outstanding fatigue strength, impact resistance
- Increased rigidity in the sagittal plane
- May be heat modified for minor fitting adjustments

### REQUIRED MEASUREMENTS

- Foot length



## ORDERING INFORMATION

| Size | Foot length | Brace height | Left    | Right   |
|------|-------------|--------------|---------|---------|
| S    | 24cm        | 26.5cm       | 3520BLU | 3525BLU |
| M    | 28cm        | 28.5cm       | 3521BLU | 3526BLU |
| L    | 28cm        | 33cm         | 3522BLU | 3527BLU |
| XL   | 30cm        | 39cm         | 3523BLU | 3528BLU |

## REBOUND FOOT-UP®

Rebound Foot-Up is a lightweight ankle/foot orthosis (AFO) that offers dynamic support for drop-foot or related disorders that require dorsiflexion support. The Rebound Foot-Up supports the ankle dorsiflexion during swing phase. Dynamic support is provided by an elastic strap connected to the shoe by either a transparent plastic inlay placed between shoelaces, or by a nut and bolt fastened through the eyelets. The positioning of the fastening allows for individualising the support by providing the option of positioning the dorsiflexion support on the medial or lateral side. When using the plastic inlay, the dorsiflexion support becomes neutral.

The Rebound Foot-Up Ankle Cuff has a smart clamp that holds the cuff in place and facilitates one-handed donning and doffing.

Optional: Rebound Foot-Up Foot Wrap can be used when barefoot or in open shoes or sandals. The adjustable elastic strap is attached to the Foot Wrap using a nut and bolt fastened through the eyelets.



### INDICATIONS FOR USE

- Drop-foot
- Disorders that require dorsiflexion support

### FEATURES AND BENEFITS

- Dynamic and discreet support for drop-foot or related disorders that require dorsiflexion support
- Improved function, ease of use, comfort, and quality
- Allows for voluntary plantar flexion
- The elastic strap is adjustable and easily replaceable
- The Rebound Foot-Up Ankle Cuff and Rebound Foot-Up Foot Wrap are made from breathable material that can be worn comfortably for long periods of time
- The Rebound Foot-Up is hand-washable
- Positioned medial or lateral support
- Contoured to the shape of leg
- Adjustable strap length
- Smart clamp facilitates one-handed donning / doffing
- Optional Rebound Foot-Up Foot Wrap can be used without shoes or in open shoes or sandals

### REQUIRED MEASUREMENTS

- Circumference of ankle, measured around the malleoli



### ORDERING INFORMATION

#### ANKLE CUFF

| Description                | Size | Ankle Circumference | Part #   |
|----------------------------|------|---------------------|----------|
| Rebound Foot-Up Ankle Cuff | S/M  | 18 - 25cm           | RFU50001 |
| Rebound Foot-Up Ankle Cuff | L/LX | 26 - 33cm           | RFU50002 |

#### FOOT WRAP

| Description               | Size     | Foot Circumference | Part #   |
|---------------------------|----------|--------------------|----------|
| Rebound Foot-Up Foot Wrap | One Size | 20 - 28cm          | RFU50101 |

### ACCESSORIES

| Description                                                                                   | Part #   |
|-----------------------------------------------------------------------------------------------|----------|
| Rebound Foot-Up Attachment Kit (includes plastic inlay, elastic strap, and 1 x nut and screw) | RFU50201 |
| Rebound Foot-Up Replacement Nuts And Bolts (5 nuts included: 2 short and 3 long screws)       | RFU50202 |

| Description            | Size | Foot Circumference | Part #    |
|------------------------|------|--------------------|-----------|
| Rebound Foot-Up Bundle | S/M  | 18 - 25cm          | RFU500011 |
| Rebound Foot-Up Bundle | L/LX | 26 - 33cm          | RFU500021 |

## FOOT-UP®

Foot-Up is a lightweight AFO that offers dynamic support for drop foot or similar complaints. The orthosis provides visible improvement in gait by supporting the foot the moment it is raised.



### INDICATIONS FOR USE

- Drop foot
- Support of dorsiflexion

### FEATURES AND BENEFITS

- Streamlined, lightweight and invisible when worn
- Breathable material, padded and comfortable to wear
- Easy to fit
- Available in black and beige
- Attached to the shoe, or available as shoeless model

### REQUIRED MEASUREMENTS

- Basic version: Circumference of ankle, measured around the malleoli
- Shoeless version: Circumference of foot immediately behind the ball of the foot

### ORDERING INFORMATION

#### FOOT-UP

| Size | Ankle Circumference | Beige   | Black   |
|------|---------------------|---------|---------|
| M    | 18 - 21cm           | 07811-1 | 07810-1 |
| L    | 22 - 26cm           | 07811-2 | 07810-2 |
| XL   | 27 - 33cm           | 07811-3 | 07810-3 |

#### SHOELESS VERSION\*

| Size | Foot Circumference | Beige   | Black   |
|------|--------------------|---------|---------|
| S    | < 20cm             | 27811-1 | 27810-1 |
| M    | 20 - 24cm          | 27811-2 | 27810-2 |
| L    | 24 - 28cm          | 27811-3 | 27810-3 |
| XL   | > 28cm             | 27811-4 | 27810-4 |

\*Shoeless version item comprises of foot wrap only and must be ordered with the Foot-Up unless the patient already has one.

### ACCESSORIES

| Description            | Beige | Black |
|------------------------|-------|-------|
| Foot-Up plastic insert | 17811 | 17810 |

| Description                         | Beige   | Black   |
|-------------------------------------|---------|---------|
| Foot-Up with Shoeless add-on Bundle | 07811xy | 07810xy |

x = size of main foot-up part

y = size of shoeless add-on

e.g. 0781012 = size 1 Foot-up with a size 2 Shoeless add-on in Black



## FORMFIT® PRO ANKLE

Össur Formfit® Pro Ankle is a compressive ankle sleeve with 2 integrated malleoli pads, designed to be worn after mild ankle injuries. The dynamic 3D MotionTech™ knit is supportive and meets certified medical grade compression requirements. Proven superior to leading brands for moisture wicking and breathability, your ankle stays cool during wear. CoolVent™ technology in the area in front of the ankle further enhances breathability and comfort.



MotionTech™



CoolVent™



SmartPad™



FoldLock™



### INDICATIONS FOR USE

- For ankle conditions that may benefit from compression, such as:
- Ankle injury rehabilitation
- Mild sprains and strains
- Mild instability


### CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

### FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit with superior moisture wicking and breathability, keeping your ankle cool during wear
- Certified medical grade compression (20-30 mmHg)
- CoolVent™ technology in the area in front of the ankle for added breathability and comfort
- 2 integrated malleoli SmartPads™ for targeted support and cushioning
- Patented FoldLock™ system for easier and more secure application

### REQUIRED MEASUREMENTS



Measure ankle circumference diagonally over the base of the heel and across the ankle bone. Ankle in 90 degrees.

| Size | Ankle circumference | Left     | Right    |
|------|---------------------|----------|----------|
| XS   | 18 - 20cm           | FP200201 | FP200101 |
| S    | 20 - 23cm           | FP200202 | FP200102 |
| M    | 23 - 26cm           | FP200203 | FP200103 |
| L    | 26 - 29cm           | FP200204 | FP200104 |
| XL   | 29 - 32cm           | FP200205 | FP200105 |

## FORMFIT® PRO ANKLE

### ORDERING INFORMATION

| Size | Left     | Right    |
|------|----------|----------|
| XS   | FP200201 | FP200101 |
| S    | FP200202 | FP200102 |
| M    | FP200203 | FP200103 |
| L    | FP200204 | FP200104 |
| XL   | FP200205 | FP200105 |



## FORMFIT® ANKLE BRACE

The Formfit Ankle brace provides stability and protection for the ankle – whether used prophylactically, for chronic instability or following an injury. A hybrid between a rigid ankle stirrup and a soft ankle support, the Formfit Ankle ensures a close fit which limits inversion, eversion, flexion and extension movements. The low profile design fits easily into shoes and the quick lace up feature makes it easy to don and doff. Breathable material ensures it can be comfortably worn all day. The Formfit Ankle with Figure 8 strap version, offers an even greater level of support, by further limiting the range of motion of the ankle.

### INDICATIONS FOR USE


- Ankle injury rehabilitation/immobilisation (Grade I/II sprains)
- Support chronic instability
- Prophylactic use

### FEATURES AND BENEFITS

- Breathable material
- Quick lace Up
- Inversion/Eversion Control
- Figure 8 strap version

### REQUIRED MEASUREMENTS

- Measure ankle circumference diagonally over the base of the heel and across the ankle bone.



## ORDERING INFORMATION

| Size | Ankle Circ. | Form Fit Ankle Brace | Form Fit Ankle with Figure 8 |
|------|-------------|----------------------|------------------------------|
| XS   | 28 - 30.5cm | B-212010001          | B-212000001                  |
| S    | 30.5 - 33cm | B-212010002          | B-212000002                  |
| M    | 33 - 35.5cm | B-212010003          | B-212000003                  |
| L    | 35.5 - 38cm | B-212010004          | B-212000004                  |
| XL   | 38 - 40.5cm | B-212010005          | B-212000005                  |

## FORMFIT® ANKLE WITH SPEEDLACE

### INDICATIONS FOR USE

- Ankle injury rehabilitation/immobilisation (Grade I/II sprains)
- Support chronic instability
- Prophylactic use


### FEATURES AND BENEFITS

- Durable materials
- Quick-Pull closure straps
- Figure 8 straps
- Removable stays
- Circumferential cuff



### REQUIRED MEASUREMENTS

- Measure ankle circumference diagonally over the base of the heel and across the ankle bone. Ankle in 90 degrees.



### ORDERING INFORMATION

| Size | Ankle Circ. | Part#   |
|------|-------------|---------|
| XXS  | 23 - 25.5cm | W-10621 |
| XS   | 25.5 - 28cm | W-10622 |
| S    | 28 - 30.5   | W-10623 |
| M    | 30.5 - 33cm | W-10625 |
| L    | 33 - 35.5cm | W-10627 |
| XL   | 35.5 - 38cm | W-10628 |
| XXL  | 38 - 40.5cm | W-10629 |

## FORMFIT® ANKLE STIRRUPS

Formfit Ankle Stirrups are lightweight and sturdy. Whether worn to accommodate an acute injury, chronic instability, or post-cast healing, Formfit Ankle Stirrups stabilise injured ankles by restricting the inversion/eversion motion. The universal design has been carefully engineered to accommodate the malleoli on both sides and ensures a great, comfortable fit. The slip resistant heel strap can be adjusted to accommodate most ankles. Our signature Flex-Edge™ overmould surrounds the stirrup's outer shell protecting patient's skin from any hard edges, further adding to a comfortable fit. The universal stirrup design is available with three different liner options: inflatable air, gel, or foam, offering each patient a solution that best fits their needs.

### INDICATIONS FOR USE

Situations requiring inversion/eversion restriction. These may include:

- Mild to severe ankle sprains
- Chronic ankle instability
- Post-cast and post-walker support
- Ankle injury rehabilitation and/or post-surgical use

### FEATURES AND BENEFITS

- Universal shell, carefully crafted to fit both right and left ankles
- Soft Flex-Edge overmould for increased comfort
- Adjustable heel strap to optimise patient fit and comfort
- Slip-resistant silicone dots to reduce slipping and stabilise the foot inside of shoe
- 2 bidirectional hook and loop straps ensure a secure and stable fit
- Choice of inflatable, foam, and gel liners to suit your needs
- Available in Adult and Youth sizes

### REQUIRED MEASUREMENTS

- Patient's height



| Size  | Patient's Height | Formfit Air Ankle Stirrup | Formfit Gel Ankle Stirrup | Formfit Foam Ankle Stirrup |
|-------|------------------|---------------------------|---------------------------|----------------------------|
| Adult | Above 160 cm     | 70350                     | 70351                     | 70352                      |
| Youth | Below 160 cm     | 71350                     | 71351                     | 71352                      |

## EXOFORM® DORSAL NIGHT SPLINT

### INDICATIONS FOR USE

- Plantar fasciitis
- Post-static pain
- General heel pain
- Drop foot
- Achilles tendonitis
- Lower extremity overuse injuries

### FEATURES AND BENEFITS

- Adjustable aluminum stay inside dorsal shell allows modification in the angle of dorsiflexion
- Soft, overmolded edges enhance comfort
- Hinged shell at the calf area prevents brace rotation
- Semi-rigid plantar pad reduces discomfort along the sides of the foot while providing an improved angle of dorsiflexion
- Non-slip surface on plantar pad prevents slipping during normal walking
- Circumferential heel and posterior strapping prevents migration while increasing overall comfort and effectiveness
- Strapping and shell perforations improve comfort and breathability



### REQUIRED MEASUREMENTS

- Shoe size

### ORDERING INFORMATION

| Size | Shoe Size UK | Shoe Size EU | Part#   |
|------|--------------|--------------|---------|
| M    | 4.5 - 8      | 37.5 - 42.5  | W-50085 |
| L    | 8.5 - 13.5   | 43 - 48.5    | W-50087 |

## FORMFIT® NIGHT SPLINT

The Formfit Night Splint, designed to help relieve the pain associated with plantar fasciitis by providing a gentle overnight stretch. Straps can be adjusted from 0° up to 10° degrees as needed. The soft breathable liner and padded straps provide comfort and help keep patient cool. Straps adjust from 0°-10° dorsiflexion

### INDICATIONS FOR USE

Conditions that may benefit from resting with the ankle in a neutral position, such as:

- Plantar fasciitis
- Achilles tendinopathy

### FEATURES AND BENEFITS

- Straps adjust from 0°-10° of dorsiflexion
- Contoured for a comfortable fit
- Soft foam interior is washable and breathable quick release buckles for easy application
- Slip resistant tread

### REQUIRED MEASUREMENTS

- Shoe size



## ORDERING INFORMATION

| Size | Shoe Size UK |            | Shoe Size EU |             | Standard   |
|------|--------------|------------|--------------|-------------|------------|
|      | Male         | Female     | Male         | Female      | with tread |
| S    | Up to 6.5    | Up to 5    | Up to 39     | Up to 38    | 50023      |
| M    | 7 - 9.5      | 5.5 - 8    | 39.5 - 44    | 38.5 - 43   | 50025      |
| L    | 10 and up    | 8.5 and up | 44.5 and up  | 43.5 and up | 50027      |

## ACCESSORIES

| Description         | Part# |
|---------------------|-------|
| Foam Wedge          | 10700 |
| Replacement Liner S | 11760 |
| Replacement Liner M | 11761 |
| Replacement Liner L | 11762 |

## DH OFFLOADING POST-OP STYLE SHOE

Designed for use in the treatment of plantar wounds, this shoe can be an excellent alternative to post-op style shoes with static pads (wool, felt, foam, etc.). It accommodates the ambulatory and compliant patient for short-term treatment (1 to 4 months).

### INDICATIONS FOR USE

- Plantar wounds
- Post-operative use
- Plantar foot ulcers

### FEATURES AND BENEFITS

- Designed to accommodate either left or right foot
- Designed around Össur's unique Active Offloading Insole
- Designed to accommodate bulky dressings
- Super soft nylon upper with reinforced heel counter provides extra comfort and stability
- Patented removable toe pieces aides in wound accessibility for daily wound site inspection



### REQUIRED MEASUREMENTS

- Shoe size

## ORDERING INFORMATION

|      | Shoe Size UK |           | Shoe Size EU |             |       |
|------|--------------|-----------|--------------|-------------|-------|
| Size | Male         | Female    | Male         | Female      | Part# |
| S    | 3.5 - 5.5    | 4 - 6     | 35.5 - 38    | 37 - 39     | 10340 |
| M    | 6 - 8        | 6.5 - 8.5 | 38.5 - 42    | 39.5 - 43.5 | 10341 |
| L    | 8.5 - 10.5   | 9 - 10    | 42.5 - 45    | 44 - 45     | 10342 |
| XL   | 11 - 13.5    | 10.5+     | 45.5 - 48.5  | 45.5-48.5   | 10343 |

## ACCESSORIES

### REPLACEMENT PARTS

| Size | Insoles |
|------|---------|
| S    | 12365   |
| M    | 12366   |
| L    | 12367   |
| XL   | 12368   |



## MESH TOP POST-OP SHOE

### INDICATIONS FOR USE

- For recovery after soft tissue procedures and post trauma

### FEATURES AND BENEFITS

- Breathable mesh upper with reinforced heel counter and padded heel collar provides increased stability overall to the rear of the foot while reducing slippage
- Semi-rigid sole provides superior protection and control after osteotomies
- Extended strap and dual side closure system accommodates extra bulky dressings

### REQUIRED MEASUREMENTS

- Shoe size



### ORDERING INFORMATION

|      | Shoe Size UK |           | Shoe Size EU |             |       |        |
|------|--------------|-----------|--------------|-------------|-------|--------|
| Size | Male         | Female    | Male         | Female      | Male  | Female |
| S    | 3.5 - 5.5    | 2.5 - 3.5 | 35.5 - 38    | 35 - 36     | 17013 | 17003  |
| M    | 6 - 8        | 4 - 5.5   | 38.5 - 42    | 36.5 - 38.5 | 17015 | 17005  |
| L    | 8.5 - 10.5   | 6 +       | 42.5 - 45    | 39 +        | 17017 | 17007  |
| XL   | 11 - 13.5    |           | 45.5 - 48.5  |             | 17018 |        |

## SOFT TOP POST-OP SHOE

### INDICATIONS FOR USE

- Post-operative use or minor trauma wounds
- Geriatric or insensate foot

### FEATURES AND BENEFITS

- The cotton-poly blend upper allows the patient's foot to breathe and provides added comfort
- Designed specifically for recovery after soft tissue procedures post trauma
- The contact closure system provides even torsion and is designed for ease of use

### REQUIRED MEASUREMENTS

- Shoe size



### ORDERING INFORMATION

|      | Shoe Size UK |           | Shoe Size EU |             |       |        |
|------|--------------|-----------|--------------|-------------|-------|--------|
| Size | Male         | Female    | Male         | Female      | Male  | Female |
| S    | 3.5 - 5.5    | 2.5 - 3.5 | 35.5 - 38    | 35 - 36     | 18013 | 18003  |
| M    | 6 - 8        | 4 - 5.5   | 38.5 - 42    | 36.5 - 38.5 | 18015 | 18005  |
| L    | 8.5 - 10.5   | 6 +       | 42.5 - 45    | 39 +        | 18017 | 18007  |
| XL   | 11 - 13.5    |           | 45.5 - 48.5  |             | 18018 |        |

# CANVAS ROCKER BOTTOM CAST SHOE

## INDICATIONS FOR USE

- For use with a short leg cast

## FEATURES AND BENEFITS

- Cushioned rocker bottom cast sandal
- Contact closure, open toe

## REQUIRED MEASUREMENTS

- Insole length



## ORDERING INFORMATION

| Size | Insole length | Blue      | Black   |
|------|---------------|-----------|---------|
| XS   | 18cm          | 308XS     |         |
| S    | 19 - 23cm     | 308SM     | B308SM  |
| S/M  | 23 - 25cm     | 308SM/MED | B308S/M |
| M    | 25.5 - 26.5cm | 308MED    | B308MED |
| M/L  | 26.5 - 28cm   | 308MLG    | B308MLG |
| L    | 28 - 30.5cm   | 308LG     | B308LG  |
| XL   | 30.5 - 35cm   | 308XLG    | B308XLG |



## TABLE OF CONTENTS

---

|                          |     |
|--------------------------|-----|
| CASTING & TRACTION ..... | 157 |
| TECHFORM® PREMIUM .....  | 159 |

## TECHFORM® PREMIUM

Techform Casting Tape was designed by and for orthopaedic casting professionals who look for a superior handling experience and a smooth, strong finished cast. It incorporates a patented rub-responsive formula that works with our unique fabric to produce the finished cast you'll ever see.

### INDICATIONS FOR USE

- Immobilisation

### FEATURES AND BENEFITS

- Glides off the roll for easier application
- Soft, conformable rolls allow for easy handling
- "Rub-responsive" resin formula adapts to the technician's technique - the more it is rubbed, the smoother the cast gets
- Excellent lamination ensures a strong and durable cast
- Smoother finished cast is non-abrasive and will not snag on clothing
- Precision welded ends do not fray and lay down easily
- Finger holes in the pouch allow for easy opening



## ORDERING INFORMATION

| WIDTH | WH    | BK    | OR    | RD    | PK    | HP    | LB    | DB    | PR    | GN    | RB    |
|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| 1.5"  | 154WH |       |       |       |       |       |       |       |       |       |       |
| 2"    | 204WH | 204BK | 204OR | 204RD | 204PK | 204HP | 204LB | 204DB | 204PR | 204GN | 204RB |
| 3"    | 304WH | 304BK | 304OR | 304RD | 304PK | 304HP | 304LB | 304DB | 304PR | 304GN | 304RB |
| 4"    | 404WH | 404BK | 404OR | 404RD | 404PK | 404HP | 404LB | 404DB | 404PR | 404GN | 404RB |
| 5"    | 504WH |       |       |       |       |       |       |       |       |       |       |

Note: Packed in boxes of 10

## NOTES

---

## NOTES


---

## NOTES

---







P-970278EN-EN

FOLLOW ÖSSUR ON



TEL 03450 065 065  
ukinfo@ossur.com

WWW.OSSUR.COM

