



Össur Legs with Dynamic TT Solutions

Digital orderform



Liner

Socket

Adapter

Tube

Foot



Össur Legs - our aim is to make doing business easy

We make it easy to order a prosthetic leg, configured by you to fit your patients’ needs. You are in control of every step of the way as Össur Legs can be tailored to fit your practice needs. Let us, help you improve:



Outcomes

In today’s O&P environment, your greatest value is in strengthening the relationships essential to growing your business. Whether providing next-level care to your patients, networking with referral sources, or collecting outcomes data to support claims, there are more demands on your time than ever.

With Össur Legs, you’ll spend less time purchasing and collecting parts, and lowers the administrative burden to comply to quality and regulatory matters. Instead, you can spend your time just where you want to, with your patients. Our high-quality products, are vigorously tested and verified to work together. On the following pages you



Efficiency

will find our product range to fulfill the specific needs of active, transfemoral users. Please contact us for our other product ranges.

7 reasons to use Össur Legs with Proprio Foot & Pro-Flex Feet:

- Components designed to be used together
- One place to select and order
- One place of contact
- All products comply to EU MDR
- Use of safe and compatible products
- Attractive package offering
- Clinical support



Patient Satisfaction

Follow the next steps to order an Össur Leg with Rheo Knee

- 1

Select the configuration.
- 2

Fill out this form > select amount, size and/or other information.
- 3

Send this form by using the ‘send button’ at the bottom of the last page.

Select your solution

Solution	Foot	Adapter(s)*	Liner	Socket
TT Solution + liner + Direct Socket	✓	✓	✓	✓
TT Solution + liner	✓	✓	✓	
TT Solution	✓	✓		

*A solution must at least contain 3 adapters.

Select the liner

Please select when not applicable

Sizing TT liners:
Standard: 16,18, 20, 22, 23.5, 25, 26.5, 28, 30, 32, 34, 36, 40

Sizing Seal Rings:
20, 22, 24, 26, 28, 30, 32, 35, 38, 41, 44, 47, 51, 55, 60, 65

TRANSTIBIAL LOCKING LINERS



Iceross Dermo Locking	Size	Pieces
Dermo Locking 3mm TT Conical		
Dermo Locking 3mm TT Standard		
Dermo Locking 3mm TT Wave		
Dermo Locking 6mm TT Std		



Iceross Synergy Locking	Size	Pieces
Synergy Locking 3mm TT Standard		
Synergy Locking 3mm TT Wave		



Iceross Comfort Locking	Size	Pieces
Comfort Locking 3mm TT Standard		
Comfort Locking 3mm TT Wave		
Comfort Locking 6mm TT Standard		

TRANSTIBIAL CUSHION LINERS



Iceross Comfort Cushion	Size	Pieces
Comfort Cushion 3mm TT Standard		
Comfort Cushion 6mm TT Standard		



Iceross Dermo Cushion	Size	Pieces
Dermo Cushion 3mm TT Standard		
Dermo Cushion 3mm TT Wave		
Dermo Cushion 6mm TT Standard		



Iceross Synergy Cushion	Size	Pieces
Synergy Cushion 3mm TT Standard		
Synergy Cushion 3mm TT Wave		

Select the socket solution

DIRECT SOCKET TT



TRANSTIBIAL SEAL-IN LINERS



Iceross Dermo Seal-In (HSM)	Size	Pieces
Dermo 3mm TT Standard		



Iceross Seal-In X5	Size	Pieces
Seal-In X5 3mm TT Standard		
Seal-In X5 3mm TT Wave		



Iceross Seal-In V	Size	Pieces
Seal-In V 3mm TT Standard		
Seal-In V 3mm TT HP		



Iceross Seal-In X*	Size	Pieces
Seal-In X 3mm TT		
Seal-In X 6mm TT		

*Order in combination with Seal-In Seal.

TRANSTIBIAL SEAL-IN X LOCKING



Iceross Seal-In X Locking*	Size	Pieces
Seal-In X Locking 3mm TT		
Seal-In X Locking 6mm TT		

*Order in combination with Seal-In Seal.

SEAL-IN SEALS



Iceross Seal-In Seals	Size	Pieces
Iceross Seal-In X-Classic S20-S65		
Iceross Seal-In X-Grip S20-S65		
Iceross Seal-In X-Volume S20-S65		

Please select when not applicable

Direct Socket Kits Basalt	Pieces
Direct Socket, 4” Basalt Fibers with standard distal connector, 200ml Black Resin	
Direct Socket, 5” Basalt Fibres with standard distal connector, 200 ml Black Resin	
Direct Socket, 7” Basalt Fibers with standard distal connector, 400ml Black Resin	

Direct Socket can only be ordered in combination with a complete leg configuration as displayed in the table 'select your solution' at the bottom left page.

Select the adapters (min. 3)



4-Hole socket adapter		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		



4-Prong socket adapter	
Weight limit	166kg
Material	Stainless Steel
Pieces	



3-Prong socket adapter	
Weight limit	166kg
Material	Stainless Steel
Pieces	



Male insert for prong		
Weight limit	100kg	166kg
Material	Stainless Steel	Titanium
Pieces		



Female insert for prong		
Weight limit	100kg	166kg
Material	Stainless Steel/Titanium	
Pieces		



Height adjustable pylon		
Weight limit	100kg	
Material	Aluminium/Stainless Steel	
Length	Standard 7,7-9,5cm	Long 9,4-12,8cm
Pieces		



4-Hole male pyramid		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		



14mm offset tube female adapter	
Weight limit	100kg
Material	Aluminium
Pieces	



4-Hole female pyramid		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		



4-Hole male pyramid with 5/16" / 8mm offset		
Weight limit	166kg	
Material	Aluminium	
Pieces		



Female pyramid tube clamp		
Weight limit	100kg	166kg
Material	Aluminium	Titanium
Pieces		



Female pylon		
Weight limit	100kg	
Material	Aluminium	
Length	Standard 250mm	Long 400mm
Pieces		



Female pylon		
Weight limit	166kg	
Material	Aluminium/ Titanium	
Length	Standard 250mm	Long 400mm
Pieces		



Icelock 562 Hybrid for Seal-In X Locking	Weight limit	Pieces
Icelock 562 Hybrid for unity	166kg	
Icelock 562 Hybrid for passive vacuum	166kg	



Icelock 200	Weight limit	Pieces
Icelock 211 clutch	n.a.	
Icelock 214 clutch 4-Hole	120kg	



Icelock 500	Pieces
Icelock Expulsion Valve 551	
Icelock Expulsion Valve 552	
Icelock Expulsion Kit 544	
Icelock 544 Direct Socket Adapter Plate	
Icelock 544 Plate Expulsion	
Icelock 544 Plate Unity	



Icelock 600 without adapter	Pieces
Icelock Ratchet 621	
Icelock Lanyard 631	
Icelock Valve 641	
Icelock Smooth 651	



Proprio Foot	
Amputation level	Transtibial and Transfemoral
Impact level	Low to Moderate
Maximum patient Weight	125kg
Build height	180mm
Heel height Adjustability	Up to 50mm
Foot size 22-30	
Category 1-8	
Foot cover beige or brown	
Side: left or right	

Order with Unity

Category selection chart								
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-125
Low Impact Level	1	1	2	3	4	5	6	7
Moderate Impact Level	1	2	3	4	5	6	7	8

Unity can be ordered with Proprio Foot in sizes 25 and above.



Pro-Flex LP Align	
Build height	147mm
Weight limit	116kg Low impact 100kg Moderate impact
Foot size 22-30	
Category 1-9	
Foot cover beige or brown	
Side: left or right	

Order with Unity

Category selection chart							
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116
Low Impact Category							
Size 22-24	1	1	2	3	4	5	N/A
Sizes 25-27	1	1	2	3	4	5	6
Sizes 28	3	3	3	3	4	5	6
Moderate Impact Category							
Size 22-24	1	2	3	4	5	N/A	N/A
Sizes 25-27	1	2	3	4	5	6	N/A
Sizes 28	3	3	3	4	5	6	N/A

PLEASE SEE NEXT PAGE FOR MORE FOOT OPTIONS

Select the foot



Pro-Flex Terra	
Build height	179mm
Weight limit	166kg
Foot size 22-30	
Category 1-9	
Foot cover beige or brown	
Side: left or right	



Pro-Flex XC	
Build height	155mm
Weight limit	166kg
Foot size 22-30	
Category 1-9	
Foot cover beige or brown	
Side: left or right	

Order with Unity for
TT solution

Order with Unity for
TF solution



Pro-Flex XC Torsion	
Build height	222mm
Weight limit	147kg
Foot size 22-30	
Category 1-9	
Foot cover beige or brown	
Side: left or right	

Order with Unity for
TT solution

Order with Unity for
TF solution



Pro-Flex Pivot	
Build height	155mm
Weight limit	125kg
Foot size 22-30	
Category 1-7	
Foot cover beige or brown	
Side: left or right	

Order with Unity for
TT solution

Order with Unity for
TF solution

Category selection chart											
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166	
Low Impact Category											
Sizes 22-24	1	1	2	3	4	5	6	N/A	N/A	N/A	
Sizes 25-27	1	1	2	3	4	5	6	7	8	9	
Sizes 28-30	3	3	3	3	4	5	6	7	8	9	
Moderate Impact Category											
Sizes 22-24	1	2	3	4	5	6	N/A	N/A	N/A	N/A	
Sizes 25-27	1	2	3	4	5	6	7	8	9	N/A	
Sizes 28-30	3	3	3	4	5	6	7	8	9	N/A	
High Impact Category											
Sizes 22-24	2	3	4	5	6	N/A	N/A	N/A	N/A	N/A	
Sizes 25-27	2	3	4	5	6	7	8	9	N/A	N/A	
Sizes 28-30	3	3	4	5	6	7	8	9	N/A	N/A	

Category selection chart											
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166	
Low Impact Category											
Sizes 22-24	1	1	2	3	4	5	6	N/A	N/A	N/A	
Sizes 25-27	1	1	2	3	4	5	6	7	8	9	
Sizes 28-30	3	3	3	3	4	5	6	7	8	9	
Moderate Impact Category											
Sizes 22-24	1	2	3	4	5	6	N/A	N/A	N/A	N/A	
Sizes 25-27	1	2	3	4	5	6	7	8	9	N/A	
Sizes 28-30	3	3	3	4	5	6	7	8	9	N/A	
High Impact Category											
Sizes 22-24	2	3	4	5	6	N/A	N/A	N/A	N/A	N/A	
Sizes 25-27	2	3	4	5	6	7	8	9	N/A	N/A	
Sizes 28-30	3	3	4	5	6	7	8	9	N/A	N/A	

Category selection chart										
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	
Low Impact Category										
Sizes 22-24	1	1	2	3	4	5	6	N/A	N/A	
Sizes 25-27	1	1	2	3	4	5	6	7	8	
Sizes 28-30	3	3	3	3	4	5	6	7	8	
Moderate Impact Category										
Sizes 22-24	1	2	3	4	5	6	N/A	N/A	N/A	
Sizes 25-27	1	2	3	4	5	6	7	8	N/A	
Sizes 28-30	3	3	3	4	5	6	7	8	N/A	
High Impact Category										
Sizes 22-24	2	3	4	5	6	N/A	N/A	N/A	N/A	
Sizes 25-27	2	3	4	5	6	7	8	N/A	N/A	
Sizes 28-30	3	3	4	5	6	7	8	N/A	N/A	

Category selection chart									
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-125	
Low Impact Category									
Sizes 22	1	1	2	3	4	N/A	N/A	N/A	
Sizes 23-24	1	1	2	3	4	5	N/A	N/A	
Sizes 25-26	2	2	2	3	4	5	6	7	
Sizes 27-28	N/A	N/A	N/A	3	4	5	6	7	
Sizes 29-30	N/A	N/A	N/A	N/A	4	5	6	7	
Moderate Impact Category									
Sizes 22	1	2	3	4	N/A	N/A	N/A	N/A	
Sizes 23-24	1	2	3	4	5	N/A	N/A	N/A	
Sizes 25-26	3	3	3	4	5	6	7	N/A	
Sizes 27-28	N/A	N/A	4	5	6	6	7	N/A	
Sizes 29-30	N/A	N/A	N/A	5	6	7	8	N/A	



Pro-Flex LP	
Build height	73mm
Weight limit	166kg
Foot size 22-30	
Category 1-9	
Foot cover beige or brown	
Side: left or right	

Order with Unity for
TT solution

Order with Unity for
TF solution

Category selection chart											
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166	
Low Impact Category											
Sizes 22-24	1	1	2	3	4	5	N/A	N/A	N/A	N/A	
Sizes 25-27	1	1	2	3	4	5	6	7	8	9	
Sizes 28-30	3	3	3	3	4	5	6	7	8	9	
Moderate Impact Category											
Sizes 22-24	1	2	3	4	5	N/A	N/A	N/A	N/A	N/A	
Sizes 25-27	1	2	3	4	5	6	7	8	9	N/A	
Sizes 28-30	3	3	3	4	5	6	7	8	9	N/A	
High Impact Category											
Sizes 22-24	2	3	4	5	N/A	N/A	N/A	N/A	N/A	N/A	
Sizes 25-27	2	3	4	5	6	7	8	9	N/A	N/A	
Sizes 28-30	3	3	4	5	6	7	8	9	N/A	N/A	



Pro-Flex LP Torsion	
Build height	147mm
Weight limit	147kg
Foot size 22-30	
Category 1-9	
Foot cover beige or brown	
Side: left or right	

Order with Unity for
TT solution

Order with Unity for
TF solution

Category selection chart									
Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
Low Impact Category									
Sizes 22-24	1	1	2	3	4	5	N/A	N/A	N/A
Sizes 25-27	1	1	2	3	4	5	6	7	8
Sizes 28-30	3	3	3	3	4	5	6	7	8
Moderate Impact Category									
Sizes 22-24	1	2	3	4	5	N/A	N/A	N/A	N/A
Sizes 25-27	1	2	3	4	5	6	7	8	N/A
Sizes 28-30	3	3	3	4	5	6	7	8	N/A
High Impact Category									
Sizes 22-24	2	3	4	5	N/A	N/A	N/A	N/A	N/A
Sizes 25-27	2	3	4	5	6	7	8	N/A	N/A
Sizes 28-30	3	3	4	5	6	7	8	N/A	N/A



You are **almost ready**

Please follow the next steps to finalise the order:

- Check the specific conditions that apply to your Össur Leg configuration.
- Please fill out the details down below and click the send button to submit the order or use the reference code when placing the order.

Conditions

The discount only applies to a configuration for one individual user/patient. However it is possible to order multiple products for trial purposes.

The invoice discount only relates to the purchase of the selected product configuration. In the situation where returning products results in incomplete configurations, the standard conditions agreed with your organisation apply. We will charge the difference of the discount.

CONTACT DETAILS

Name: _____

Contact person: _____

Phone: _____

E-mail: _____

Customer number: _____

Reference: _____

Date: _____

Please first send me a quote:

Comments:

Clicking on the yellow button should open an external window to send the e-mail. If nothing happens, please attach this orderform to an e-mail and send it to **info@ossur.com**.

Thank you very much and please contact our customer care team if you have questions.

FOLLOW ÖSSUR ON



TEL + 46 18 18 22 00
info@ossur.com

WWW.OSSUR.SE

