

Service Fabrication Order Form

1. Customer / Patient Information

Company name: _____ Fax: _____
 Customer number: _____ Össur order no.: _____
 Address: _____ *(to be completed by Össur)*
 _____ User reference: _____
 Country: _____ Date: _____
 Contact person: _____ Anatomic gender: Male Female
 Phone: _____ Weight: _____ kg
 E-mail: _____ Side: Left Right

2. Basis for Definitive Socket Production

Existing Check Socket

[New Prolutions check socket \(Please provide separate Prolutions order form\).](#)

3. Socket Type

Fiber Socket TT (PR-DEF-TT01)	Fiber Socket TF (PR-DEF-TF01)	Fiber Socket Flex KD (PR-DEF-KD02)
Fiber Socket Flex TT (PR-DEF-TT02)	Fiber Socket Flex TF (PR-DEF-TF02)	Fiber Socket Comfort KD (PR-DEF-KD03)
Fiber Socket Comfort TT (PR-DEF-TT03)	Fiber Socket Comfort TF (PR-DEF-TF03)	

4. Select Reinforcement Options For Definitive Socket

1. PETG innerlayer
2. Flex innerlayer (EVA Material)
3. Carbonlayer:

1. Laminate below adapter:

_____ Layers carbon hose;
 _____ Unidirectional carbon layer(s) (vertical reinforcement);
 _____ Unidirectional carbon layer(s) circular proximal

Recommended minimum lamination layers:
 Below adapter: 1,5 carbon layers hose, 1 unidirectional carbon layer, 1 unidirectional carbon layer circular proximal.
 Above adapter: 1 unidirectional carbon layer vertical and 2 horizontal over each prong, 1 unidirectional carbon layer vertical, 2 unidirectional carbon layers horizontal, 2 layers carbon hose - only for TF.

2. Laminate above adapter: _____ Unidirectional layer(s) vertical and _____ horizontal over each prong.
 _____ Unidirectional carbon layer(s) (vertical reinforcement);
 _____ Unidirectional carbon layers (horizontal reinforcement);
 _____ Layers carbon hose

Please define the trimline of the carbonfibre lamination on the surface of the check socket:

Posterior Cut-out height from proximal _____ cm and width _____ cm (please also mark on check socket)

5. Select Finish of Definitive Socket

Carbon layer optic Black (no additional costs)

Designs and Colours
(PR-DEF-DSGN; charged separately)

6. Select Össur-Leg Parts

Info: To make sure that your order will be classified as an Össur Leg order please make sure to select correct number of primary / secondary products according to the following table for TT/TF legs:

Össur leg definition:				
Primary products	Liner	Knee	Foot	
Secondary products		Adapter(s)	Lock	Pylon
TT Leg Minimum	2 Primary Products		1 or more Secondary Products	
TF Leg Minimum	2 or more Primary Products		1 or more Secondary Products	

6.1 Liner Selection

Send liners first for casting/measuring		No liner	
TF Locking Liner	TF Seal-In Liner	TT Locking Liner	
Iceross TF Locking stnd profile (I-7032XX) conical profile (I-7132XX)	Iceross Seal-In X5 TF stnd profile (I-TF673XX) conical profile (I-TF678XX)	Iceross Synergy Locking Wave stnd 10mm matrix (I-6303XX) special order matrix ____ cm (I-6206XX)	Iceross Comfort Locking Wave (I-5303XX)
AKOS TF stnd profile (M8955XX) conical profile (M8956XX)	Iceross Seal-In TF stnd profile (I-7532X) conical profile (I-7632XX)	Iceross Synergy Locking stnd 10mm matrix (I-6003XX) special order matrix ____ cm (I-6203XX)	Iceross Comfort Locking stnd 10mm matrix, 3mm profile (I-5406XX) special order matrix ____ cm, 3mm profile (I-5506XX)
Relax TF stnd profile (M8957XX) conical profile (M8958XX)	Iceross Seal-In X TF stnd profile (I-8532XX) conical profile (I-8632XX)	Iceross Dermo Locking Wave stnd 10mm matrix (I-4313XX) special order matrix ____ cm (I-4223XX)	Iceross Comfort Locking stnd 10mm matrix, 6mm profile (I-5006XX) special order matrix ____ cm, 6mm profile (I-5106XX)
	Iceross Seal-In X Locking TF stnd profile (I-SXLSXX) conical profile (I-SXLCXX)	Iceross Dermo Locking Wave stnd 10mm matrix (I-4313XX) special order matrix ____ cm (I-4223XX)	Iceross Original stnd 10mm matrix (I-0124XX) special order matrix ____ cm (I-2024XX)
	4 Seal Classic stnd profile (M8917XX) conical profile (M8918XX)	Iceross Dermo Locking stnd 10mm matrix, 3mm profile (I-4013XX) special order matrix ____ cm, 3mm profile (I-4213XX)	Iceross Sport stnd 10mm matrix (I-1033XX) special order matrix ____ cm (I-1233XX)
	4 Seal Senso stnd profile (M8963XX) conical profile (M8964XX)	stnd 10mm matrix, 6mm profile (I-4016XX) special order matrix ____ cm, 6mm profile (I-4216XX)	Relax Locking TT 3mm profile (M8994XX) 6mm profile (M8996XX)
		Iceross Dermo Uniform Locking Medium (I-410618) Medium+ (I-410620) Large (I-410625) Large+ (I-410628) XLarge (I-410633)	Iceform Locking 3mm profile (T-0103XX) 6mm profile (T-0106XX) 9mm profile (T-0109XX)

TT Seal-In Liner		TT Cushion Liner	
Iceross Seal-In X Locking 3mm profile (I-SXL3XX) 6mm profile (I-SXL6XX)	Iceross Seal-In X5 TT Wave (I-3664XX) Iceross Seal-In X5 (I-3663XX) Iceross Dermo Seal-In (I-4613XX)	Iceross Synergy Cushion Wave (I-CW63XX) Iceross Synergy Cushion (I-CL63XX) Iceross Dermo Cushion Wave (I-CW43XX) Iceross Dermo Cushion 3mm profile (I-CL43XX) 6mm profile (I-CL46XX) Iceross Dermo Uniform Cushion Medium (I-440618) Medium+ (I-440620) Large (I-440625) Large+ (I-440628) XLarge (I-440633)	Iceross Comfort Cushion 3mm profile (I-CL53XX) 6mm profile (I-CL56XX) Iceross Activa Cushion (I-CT62XX) Relax Cushion 3mm profile (M8997XX) 6mm profile (M8998XX) Iceform Cushion 3mm profile (T-CL03XX) 6mm profile (T-CL06XX) 9mm profile (T-CL09XX)

















Select Liner size	
Liner size	<input type="text"/>





























Seal-Rings		
Seal type	Classic (I-SXC0XX)	
	Grip (I-SXC0XX)	
	Volume (I-SXV0XX)	
Seal-Ring size	<input type="text"/>	

Sleeve		
Iceross Sleeve (F-2011XX)	Size:	<input type="text"/>
Iceflex Balance Sleeve (F-1010XX)	Size:	<input type="text"/>
Genu Sleeve (MK0020X)	Size:	<input type="text"/>
Iceform Sleeve (S-2011XX)	Size:	<input type="text"/>











6.2 Adapter Selection




(internal note: add parts from this page to socket order)



Locking Options			
Icelock 621 Ratchet & Extra Pin (L-621000) 	Icelock 651 Smooth & Extra Pin (L-651000) 	Icelock 211 Clutch & Extra Pin (L-211000) 	Icelock 234 Lanyard (L-234000) 
WORKS WITH			
Icelock 672 Male Pyramid Ti (L-672000) 	Icelock 674 4-Hole Adapter Al (L-674000) 	Icelock 214 4-Hole Clutch & Extra Pin (L-214000) 	
Icelock 673 Male Pyramid Al (L-673000) 	Icelock 675 4-Hole Adapter Ti (L-675000) 		
Icelock Clutch Options For TF			
214 AK for left leg with release button on medial side (L-214311) 	214 AK for left leg with release button on lateral side (L-214312) 	214 AK for right leg with release button on lateral side (L-214322) 	214 AK for right leg with release button on medial side (L-214321) 
Expulsion Options		Unity Elevated Vacuum Options	
Icelock 544 Expulsion Kit for Seal-In (L-544621) 	Icelock 544 Expulsion only (L-544620) 	Icelock 544 Unity Plate (L-544611) 	







Socket Adapters			
<p>Icelock 544 Socket Adapter (L-544400)</p>  	<p>3-Prong Socket Adapter (A-114030)</p>  	<p>4-Prong Socket Adapter (A-114040)</p>  	
<p>4-Hole Socket Adapter (A-122100)</p>  	<p>4-Hole Connection Plate (1651014)</p> 	<p>Loop Adapter (1970862)</p> 	
<p>Male Pyramid Socket Adapter (A-135100)</p>  	<p>4-Hole Socket Adapter (A-125100)</p>  		
Hybrid Locks			
<p>Icelock 850 Hybrid Lamination (L-850003)</p>  	<p>Icelock 850 Hybrid 4-hole (L-850002)</p>  	<p>Icelock 562 Hybrid for Unity (L-562000)</p>  	<p>Icelock 562 passive vacuum (L-562500)</p>  
Expulsion		Unity	
<p>Expulsion Valve 551 (L-551002)</p> <p>Expulsion Valve 552 (L-552000)</p>   	<p>Unity TT Valve (UTT0001)</p> <p>Unity TF Valve (UTF0001)</p>   		

(internal note: add parts from this selection to leg part order)



























Prong Adapters		
	Male Pyramid Insert For Prong (A-134300)	
	Female Pyramid Insert For Prong (A-144300)	
	Male Pyramid Insert for Prong (A-135300)	
	Female Pyramid Insert for Prong (A-145300)	
	Female Pyramid Insert for Prong LP (A-145310)	













Height Adjustable Pylons		
	Standard (A-542350)	
	Long (A-542360)	


Rotation Adapter		
	Knee Rotation Adapter (A-554300)	













Flexion Plates		
	45mm offset 10 degrees flexion (A-282440)	
	60mm offset 10 degrees flexion (A-282460)	
	75mm offset 15 degrees flexion (A-282470)	


(internal note: add parts from this selection to leg part order)

4-Hole Adapters		
	4-Hole Male (A-235100)	
	4-Hole Female (A-245100)	
	4-Hole Male w/ Rotation (A-235300)	
	4-Hole Male Pyramid, Aluminum (A-233100)	
	4-Hole Female Pyramid, Aluminum (A-242100)	
	Euro 4-Hole Adapter (A-122300)	
	4-Hole Male Pyramid With Rotation (A-234320)	
	4-Hole Pyramid Receiver With Rotation (A-244320)	
	4-Hole Male w/ 8mm Ax Offset (A-235400)	
	4-Hole Female w/ Rotation (A-245300)	
	4-Hole Female Pyramid 20/20 Offset 10° Right (A-242510)	
	4-Hole Female Pyramid 20/20 Offset 10° Left (A-242520)	
	4-Hole Female Pyramid With 20/20 Offset (A-242500)	

Single Adapters		
	Male Single Adapter Short (A-433110)	
	Male Single Adapter Long (A-433120)	
	Male Single Adapter Short (A-435110)	
	Male Single Adapter Long (A-435120)	
	Female Single Adapter Short (A-445100)	
	Female Single Adapter Long (A-442110)	






Double Adapters		
	Female Double Adapter 32mm (A-542110)	
	Female Double Adapter 45mm (A-542120)	
	Female Double Adapter 60mm (A-542130)	
	Female Double Adapter 75mm (A-542140)	
	Female Double Adapter 32mm (A-545110)	
	Female Double Adapter 45mm (A-545120)	
	Female Double Adapter 60mm (A-545130)	
	Female Double Adapter 75mm (A-545140)	

Pylons		
	Female Short Aluminum (A-742010)	
	Female Long Aluminum (A-742020)	
	Female Short Aluminum/Titan (A-748010)	
	Female Long Aluminum/Titan (A-748020)	
	Female Short Aluminum Tube (A-748030)	
	Female Long Aluminum Tube (A-748040)	

Tube Clamps		
	Male Pyramid Tube Clamp (A-335100)	
	Female Pyramid Tube Clamp (A-345100)	
	4-Hole Tube Clamp (A-325100)	
	Male Tube Clamp (A-333100)	
	Female Tube Clamp (A-342100)	
	4-Hole Tube Clamp (A-322100)	







6.3 Knee Selection

No knee

	OP4 knee (1720185)		Balance Knee OFM1 w/ IKF Adapter (1721127) w/ Male Pyramid (1721126) w/ Loop Adapter (1721128)		Balance Knee OFM2 (1721120)
	Total Knee 2000 (2000) TK2000 incl Male Pyramid (2000-OS) TK2000 incl- 4 Arm Adapter (2000-KE)		Total Knee 1900 (1900) TK 1900 incl. Male Pyramid (1900-KE) TK1900 incl. 4 Arm Adapter (1900-KE)		Total Knee 2100 (2100) TK2100 incl. Pyramid Adapter (2100-OS) TK21000 incl. 4 Arm Adapter (2100-KE)
	OHP3 Knee w/ IKF Adapter (1728161) w/ Male Pyramid (1728160) w/ Loop Adapter (1728162)		Paso Knee w/ IKF Adapter (1728181) w/ Male Pyramid (1728180) w/ Loop Adapter (1728182)		OP5 knee w/ IKF Adapter (1720182) w/ Male Pyramid (1720180) w/ Loop Adapter (1720181)
	Mauch Knee (MKN01360) Mauch Knee Plus(MKN01660)		Locking Knee (LKN12500)		
	Rheo Knee incl. 2 Years Warranty (RKN130002) Rheo Knee incl. 3 years Warranty (RKN130003) Rheo Knee incl. 5 years Warranty and additional 40 month Service Check (RKN130005)		Rheo Knee XC incl. 2 Years Warranty (RKNXC0002) Rheo Knee XC incl. 3 years Warranty (RKN0003) Rheo Knee XC incl. 5 years Warranty and additional 40 month Service Check (RKN0005)		Navii knee incl. 2 Years Warranty (NAV01002) Navii knee incl. 3 Years Warranty (NAV01003) Navii knee incl. 5 Years Warranty and additional 40 month Service Check (NAV01005)
	Power Knee incl. 2 years warranty (PKA10001) One Year Warranty Extension (PKA01GE1Y) Two years Warranty Extension incl. 40 Month Service Check (PKA01GE2Y) Three years Warranty Extension incl. 40 Month Service Check (PKA01GE3Y)				Navii knee incl. 6 years warranty and additional 40 month Service Check (NAV01006)

6.4 Foot Selection

No Foot

	Pro-Flex Terra (PTA0XYYZ)		Pro-Flex XC (PXC0XYYZ) with Unity (PXCUXYYZ)		Pro-Flex XC Torsion (PXT0XYYZ) with Unity (PXTUXYYZ)
	Pro-Flex LP (PLP0XYYZ) with Unity (PLPUXYYZ)		Pro-Flex LP Torsion (PLT0XYYZ) with Unity (PLTUXYYZ)		Pro-Flex LP Align (PLA0XYYZ)
	Pro-Flex Pivot (PFP0XYYZ) with Unity (PFPUXYYZ)		Pro-Flex ST (PST0XYYZ)		Pro-Flex Modular (PMD0XYYZ)
	ReFlex Shock (SSPEXY0) with Unity (SSPUXYY0)		ReFlex Rotate (RSPEXY0) with Unity (RSPUXYYZ)		Vari-Flex (VFPEXY0) With Unity (VFPUXYY0)
	Flex Foot Balance (FBP0XYYZ) with Unity (FBPUXYYZ)		Flex Foot Balance with D/P Flexion (FBD0XYYZ) with Unity (FBDUXYYZ)		Flex-Foot Assure (FAPEXYYZ) with Unity (FAPUXYYZ)
	Balance Foot S (BSP0XYYZ) with Unity (BSPUXYYZ)		Balance Foot S Torsion (BST0XYYZ) with Unity (BSTUXYYZ)		Balance Foot J (JBPEXYYZ) with Unity (JBPUXYYZ)
	Proprio Foot (PSX01XYYZ) with Unity (PSXU1XYYZ) Warranty Extension to 3 Years (PSXGE1Y) Warranty Extension to 4 Years incl additional 40 Month Service Check (PSXGE2Y) Warranty Extension to 5 Years incl additional 40 Month Service Check (PSXGE3Y)				

Foot Options: Foot Size: _____ cm. Category (1-9) _____ Left Right Foot Shell: Beige Brown