

Össur Leg:

K-1 BK Order Form

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Step 1 General Information


Prosthetist Name: _____ Acc Nr: _____
 Patient Identifier: _____ Weight: _____ kg Date of Birth: ____ / ____ / _____

Step 2 Select Suspension System



Seal-In®	Cushion	Locking	PLUS Suspension System
ICEROSS SEAL-IN® X 3mm 6mm Size: _____ Ring size: _____ Classic Volume Grip	ICEROSS® SYNERGY™ CUSHION Wave Size: _____ ICEROSS DERMO® CUSHION 3mm 6mm 3mm Wave Size: _____ *ICEROSS COMFORT® CUSHION 3mm 6mm Size: _____	ICEROSS® SYNERGY™ LOCKING Wave Size: _____ ICEROSS DERMO® LOCKING 3mm 6mm 3mm Wave 3mm Conical Wave Size: _____ *ICEROSS COMFORT® LOCKING 3mm 6mm 3mm Wave Size: _____	ICELOCK ICELOCK® CLUTCH 211 L-211000 OR ICELOCK® CLUTCH 4H 214 L-214000 OR *ICELOCK® 621 RATCHET L-621000 OR ICELOCK® LANYARD 631 L-631000 OR ICELOCK® SMOOTH 651 L-651000 OR ICELOCK® 562 HYBRID L-562000 OR ICELOCK® 562 HYBRID DS L-562400 OR ICELOCK® 850 HYBRID DS L-850001 OR ICELOCK® 850 HYBRID 4H L-850002 OR ICELOCK® 850 HYBRID LS L-850003 PLUS Expulsion VALVE SYSTEM *ICELOCK® EXPULSION VALVE 551 L-551002 OR PLUS UNITY SYSTEM UNITY TT VALVE KIT UTT001 UNITY PUMP
ICEROSS SEAL-IN® V Standard Profile High Profile Size: _____	ICEROSS SEAL-IN® X5 Wave Size: _____ *ICEROSS SEAL-IN DERMO® Size: _____	ICEROSS SEAL-IN X LOCKING 3mm 6mm Classic Volume Grip Ring size: _____	

Step 3 Select Socket - Direct Socket TT



Braid Selection

Liner Size	Braid Size
16 - 22	4"
23.5 - 30	5"
32 - 34	7"

ICEROSS® CASTING LINER TT

Size: _____

Socket Components Selection

600 MODULAR SPACER FOR ICELOCK 600 SERIES	L-661000
BOTTOM SEAL FOR DIRECT SOCKET	M-400206

Direct Socket Material Kit Selection

4" Glass Fibers with small distal connector, Neutral Resin	M-500104
5" Glass Fibers with medium distal connector, Neutral Resin	M-500105
7" Glass Fibers with medium distal connector, Neutral Resin	M-500107
7" Glass Fibers with large distal connector, Neutral Resin	M-501107
4" Glass Fibers with small distal connector, Black Resin	M-500114
5" Glass Fibers with medium distal connector, Black Resin	M-500115
7" Glass Fibers with medium distal connector, Black Resin	M-500117
7" Glass Fibers with large distal connector, Black Resin	M-501117
5" Carbon Fiber with medium distal connector, Black Resin	M-500205
7" Carbon Fiber with medium distal connector, Black Resin	M-500207
4" Basalt Fibers with small distal connector, Black Resin	M-500314
5" Basalt Fibers with medium distal connector, Black Resin	M-500315
7" Basalt Fibers with medium distal connector, Black Resin	M-500317
7" Basalt Fibers with large distal connector, Black Resin	M-501317

Step 4 Select Foot + Cover



*FLEX-FOOT® BALANCE



BALANCE™ FOOT S



BALANCE™ FOOT J



BREEZE™ FOOT



*FOOT COVER Beige Brown

Foot Size: _____

Left/Right: _____

Category: _____

Step 5 Select Adapter Combinations

SOCKET ADAPTERS

4-HOLE SOCKET ADAPTER AL 100kg	A-122100
4-HOLE SOCKET ADAPTER TI 166kg	A-125100
4-HOLE CONNECTION PLATE AL 136kg	1651014
4-PRONG SOCKET ADAPTER SS 166kg	A-114040
3-PRONG SOCKET ADAPTER SS 166kg	A-114030
MALE PYRAMID SOCKET ADAPTER TI 166kg	A-135100

PYRAMID ADAPTERS

PYRAMID INSERTS

MALE PYRAMID INSERT FOR PRONG TI 166kg	A-135300
FEMALE PYRAMID INSERT FOR PRONG TI 166kg	A-145300
FEMALE PYRAMID INSERT FOR PRONG Low Profile TI 166kg	A-145310
MALE PYRAMID INSERT FOR PRONG SS Insert 100kg	A-134300
FEMALE PYRAMID INSERT FOR PRONG SS 100kg	A-144300

4 HOLE ADAPTERS

4-HOLE MALE PYRAMID AL with SS Insert 100kg	A-233100
4-HOLE FEMALE PYRAMID AL 100kg	A-242100
4-HOLE TUBE CLAMP AL 100kg	A-322100
4-HOLE MALE PYRAMID WITH ROTATION SS 136kg	A-234320
4-HOLE PYRAMID RECEIVER WITH ROTATION SS 136kg	A-244320
4-HOLE MALE PYRAMID TI 166kg	A-235100
4-HOLE MALE PYRAMID WITH ROTATION TI 166kg	A-244320
4-HOLE FEMALE PYRAMID TI 166kg	A-235300
4-HOLE FEMALE PYRAMID WITH ROTATION TI 166kg	A-245100
4-HOLE TUBE CLAMP TI 166kg	A-325100
EURO 4-HOLE ADAPTER AL 125kg	A-122300

SHIFT ADAPTERS

FEMALE PYRAMID TUBE CLAMP WITH 14MM OFFSET AL 100kg	A-342430
4-HOLE FEMALE PYRAMID WITH 20/20 OFFSET AL 125kg	A-242500
4-HOLE FEMALE PYRAMID 20/20 OFFSET 10° RIGHT AL 125kg	A-242510
4-HOLE FEMALE PYRAMID 20/20 OFFSET 10° LEFT AL 125kg	A-242520
MALE-FEMALE DOUBLE ADAPTER WITH 8MM OFFSET & 10° ANGLE SS 125kg	A-554700
4-HOLE MALE PYRAMID WITH 5/16" / 8MM OFFSET TI 166kg	A-235400
SHIFT ADAPTER WITH TITANIUM TI 166kg	A-645210

PYLONS AND TUBES

100kg

FEMALE PYLON KIT SHORT A-762111 AL 100kg	A-762111
FEMALE PYLON KIT LONG A-762121 AL 100kg	A-762121
HEIGHT ADJUSTABLE PYLON SHORT 77-95MM AL 100kg	A-542350
HEIGHT ADJUSTABLE PYLON LONG 94-128MM AL 100kg	A-542350

FEMALE PYLON SHORT AL 100kg	
FEMALE DOUBLE ADAPTER 32MM AL 100kg	A-135100
FEMALE DOUBLE ADAPTER 45MM AL 100kg	A-542110
FEMALE DOUBLE ADAPTER 60MM AL 100kg	A-542120
FEMALE DOUBLE ADAPTER 75 MM AL 100kg	A-542140

MALE DOUBLE ADAPTER AL with SS Pyramid Insert 100kg	A-533100
---	----------

MALE PYRAMID TUBE CLAMP AL 100kg	A-333100
FEMALE PYRAMID TUBE CLAMP AL 100kg	A-342100

MALE SINGLE ADAPTER - SHORT AL 100kg	A-433110
MALE SINGLE ADAPTER - LONGAL 100kg	A-433120*
FEMALE SINGLE ADAPTER - SHORT AL 100kg	A-442100*
FEMALE SINGLE ADAPTER - LONG AL 100kg	A-442110*

166 KG

FEMALE PYLON SHORT TI 166kg	A-748030
FEMALE DOUBLE ADAPTER 32MM TI 166kg	A-545110
FEMALE DOUBLE ADAPTER 45MM TI 166kg	A-545120
FEMALE DOUBLE ADAPTER 60MM TI 166kg	A-545130
FEMALE DOUBLE ADAPTER 75 MM TI 166kg	A-545140

MALE DOUBLE ADAPTER TI 166kg	A-535100
------------------------------	----------

MALE PYRAMID TUBE CLAMP TI 166kg	A-335100
FEMALE PYRAMID TUBE CLAMP TI 166kg	A-345100

MALE SINGLE ADAPTER - SHORT TI 166kg	A-435110
MALE SINGLE ADAPTER - LONG TI 166kg	A-435120
FEMALE SINGLE ADAPTER TI 166kg	A-445100

* Reimbursed for PMB



WWW.OSSUR.COM/EN-ZA

Össur South Africa
Unit 1, Stonewood Business Park
9 Tanzanite Crescent
Kraaifontein
Cape Town
7550

T 0800 996 220
F 0800 981 744
E ossursa@ossur.com

