

# Össur Leg:

## K-3 BK Order Form

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilisation beyond simple locomotion.

### Step 1 General Information


Prosthetist Name: \_\_\_\_\_ Acc Nr: \_\_\_\_\_  
 Patient Identifier: \_\_\_\_\_ Weight: \_\_\_\_\_ kg Date of Birth: \_\_\_\_ / \_\_\_\_ / \_\_\_\_\_

### Step 2 Select Suspension System



Seal-In®	Cushion	Locking	PLUS Suspension System
<b>ICEROSS SEAL-IN® X</b> 3mm 6mm Size: _____ Ring size: _____ Classic Volume Grip	<b>ICEROSS® SYNERGY™ CUSHION</b> Wave Size: _____ <b>ICEROSS DERMO® CUSHION</b> 3mm 6mm 3mm Wave Size: _____	<b>ICEROSS® SYNERGY™ LOCKING</b> Wave Size: _____ <b>ICEROSS DERMO® LOCKING</b> 3mm 6mm 3mm Wave 3mm Conical Wave Size: _____	<b>ICELOCK</b> ICELOCK® CLUTCH 211 L-211000 OR ICELOCK® CLUTCH 4H 214 L-214000 OR *ICELOCK® 621 RATCHET L-621000 OR ICELOCK® LANYARD 631 L-631000 OR ICELOCK® SMOOTH 651 L-651000 OR ICELOCK® 562 HYBRID L-562000 OR ICELOCK® 562 HYBRID DS L-562400 OR ICELOCK® 850 HYBRID DS L-850001 OR ICELOCK® 850 HYBRID 4H L-850002 OR ICELOCK® 850 HYBRID LS L-850003
<b>ICEROSS SEAL-IN® V</b> Standard Profile High Profile Size: _____	<b>*ICEROSS COMFORT® CUSHION</b> 3mm 6mm Size: _____	<b>*ICEROSS COMFORT® LOCKING</b> 3mm 6mm 3mm Wave Size: _____	<b>PLUS Expulsion VALVE SYSTEM</b> *ICELOCK® EXPULSION VALVE 551 L-551002
<b>ICEROSS SEAL-IN® X5</b> Wave Size: _____	<b>PLUS SLEEVE</b> ICEROSS® SLEEVE OR ICEFLEX® BALANCE OR ICEFORM® SLEEVE OR *GENU SLEEVE Sleeve size: _____		OR <b>PLUS UNITY SYSTEM</b> UNITY TT VALVE KIT UTT0001 UNITY PUMP
<b>ICEROSS SEAL-IN X LOCKING</b> 3mm 6mm Classic Volume Grip Ring size: _____			

### Step 3 Select Socket - Direct Socket TT



#### Braid Selection

Liner Size	Braid Size
16 - 22	4"
23.5 - 30	5"
32 - 34	7"

#### ICEROSS® CASTING LINER TT

Size: \_\_\_\_\_

#### Socket Components Selection

600 MODULAR SPACER FOR ICELOCK 600 SERIES	L-661000
BOTTOM SEAL FOR DIRECT SOCKET	M-400206

#### Direct Socket Material Kit Selection

4" Glass Fibers with small distal connector, Neutral Resin	M-500104
5" Glass Fibers with medium distal connector, Neutral Resin	M-500105
7" Glass Fibers with medium distal connector, Neutral Resin	M-500107
7" Glass Fibers with large distal connector, Neutral Resin	M-501107
4" Glass Fibers with small distal connector, Black Resin	M-500114
5" Glass Fibers with medium distal connector, Black Resin	M-500115
7" Glass Fibers with medium distal connector, Black Resin	M-500117
7" Glass Fibers with large distal connector, Black Resin	M-501117
5" Carbon Fiber with medium distal connector, Black Resin	M-500205
7" Carbon Fiber with medium distal connector, Black Resin	M-500207
4" Basalt Fibers with small distal connector, Black Resin	M-500314
5" Basalt Fibers with medium distal connector, Black Resin	M-500315
7" Basalt Fibers with medium distal connector, Black Resin	M-500317
7" Basalt Fibers with large distal connector, Black Resin	M-501317

### Step 4 Select Foot + Cover



\*VARI-FLEX®

\*▲VARI-FLEX® LP

▲PRO-FLEX® PIVOT

▲PRO-FLEX® LP

▲PRO-FLEX® LP TORSION

▲PRO-FLEX® XC

▲PRO-FLEX® XC TORSION

▲PRO-FLEX® LP ALIGN

▲PRO-FLEX® TERRA

▲PRO-FLEX® ST

▲PROPRIO FOOT®

▲RE-FLEX SHOCK™ STD

▲RE-FLEX SHOCK™ TALL

\*FOOT COVER Beige Brown

Foot Size: \_\_\_\_\_

Left/Right: \_\_\_\_\_

Category: \_\_\_\_\_

▲ Please note delivery period 14 working days for this product

## Step 5 Select Adapter Combinations

### SOCKET ADAPTERS

4-HOLE SOCKET ADAPTER AL 100kg	A-122100
4-HOLE SOCKET ADAPTER TI 166kg	A-125100
4-HOLE CONNECTION PLATE AL 136kg	1651014
4-PRONG SOCKET ADAPTER SS 166kg	A-114040
3-PRONG SOCKET ADAPTER SS 166kg	A-114030
MALE PYRAMID SOCKET ADAPTER TI 166kg	A-135100

### PYRAMID ADAPTERS

#### PYRAMID INSERTS

MALE PYRAMID INSERT FOR PRONG TI 166kg	A-135300
FEMALE PYRAMID INSERT FOR PRONG TI 166kg	A-145300
FEMALE PYRAMID INSERT FOR PRONG Low Profile TI 166kg	A-145310
MALE PYRAMID INSERT FOR PRONG SS Insert 100kg	A-134300
FEMALE PYRAMID INSERT FOR PRONG SS 100kg	A-144300

### 4 HOLE ADAPTERS

4-HOLE MALE PYRAMID AL with SS Insert 100kg	A-233100
4-HOLE FEMALE PYRAMID AL 100kg	A-242100
4-HOLE TUBE CLAMP AL 100kg	A-322100
4-HOLE MALE PYRAMID WITH ROTATION SS 136kg	A-234320
4-HOLE PYRAMID RECEIVER WITH ROTATION SS 136kg	A-244320
4-HOLE MALE PYRAMID TI 166kg	A-235100
4-HOLE MALE PYRAMID WITH ROTATION TI 166kg	A-244320
4-HOLE FEMALE PYRAMID TI 166kg	A-235300
4-HOLE FEMALE PYRAMID WITH ROTATION TI 166kg	A-245100
4-HOLE TUBE CLAMP TI 166kg	A-325100
EURO 4-HOLE ADAPTER AL 125kg	A-122300

### SHIFT ADAPTERS

FEMALE PYRAMID TUBE CLAMP WITH 14MM OFFSET AL 100kg	A-342430
4-HOLE FEMALE PYRAMID WITH 20/20 OFFSET AL 125kg	A-242500
4-HOLE FEMALE PYRAMID 20/20 OFFSET 10° RIGHT AL 125kg	A-242510
4-HOLE FEMALE PYRAMID 20/20 OFFSET 10° LEFT AL 125kg	A-242520
MALE-FEMALE DOUBLE ADAPTER WITH 8MM OFFSET & 10° ANGLE SS 125kg	A-554700
4-HOLE MALE PYRAMID WITH 5/16" / 8MM OFFSET TI 166kg	A-235400
SHIFT ADAPTER WITH TITANIUM TI 166kg	A-645210

### PYLONS AND TUBES

#### 100kg

FEMALE PYLON KIT SHORT A-762111 AL 100kg	A-762111
FEMALE PYLON KIT LONG A-762121 AL 100kg	A-762121
HEIGHT ADJUSTABLE PYLON SHORT 77-95MM AL 100kg	A-542350
HEIGHT ADJUSTABLE PYLON LONG 94-128MM AL 100kg	A-542350

FEMALE PYLON SHORT AL 100kg	
FEMALE DOUBLE ADAPTER 32MM AL 100kg	A-135100
FEMALE DOUBLE ADAPTER 45MM AL 100kg	A-542110
FEMALE DOUBLE ADAPTER 60MM AL 100kg	A-542120
FEMALE DOUBLE ADAPTER 75 MM AL 100kg	A-542140

MALE DOUBLE ADAPTER AL with SS Pyramid Insert 100kg	A-533100
---	----------

MALE PYRAMID TUBE CLAMP AL 100kg	A-333100
FEMALE PYRAMID TUBE CLAMP AL 100kg	A-342100

MALE SINGLE ADAPTER - SHORT AL 100kg	A-433110
MALE SINGLE ADAPTER - LONGAL 100kg	A-433120*
FEMALE SINGLE ADAPTER - SHORT AL 100kg	A-442100*
FEMALE SINGLE ADAPTER - LONG AL 100kg	A-442110*

#### 166 KG

FEMALE PYLON SHORT TI 166kg	A-748030
FEMALE DOUBLE ADAPTER 32MM TI 166kg	A-545110
FEMALE DOUBLE ADAPTER 45MM TI 166kg	A-545120
FEMALE DOUBLE ADAPTER 60MM TI 166kg	A-545130
FEMALE DOUBLE ADAPTER 75 MM TI 166kg	A-545140

MALE DOUBLE ADAPTER TI 166kg	A-535100
------------------------------	----------

MALE PYRAMID TUBE CLAMP TI 166kg	A-335100
FEMALE PYRAMID TUBE CLAMP TI 166kg	A-345100

MALE SINGLE ADAPTER - SHORT TI 166kg	A-435110
MALE SINGLE ADAPTER - LONG TI 166kg	A-435120
FEMALE SINGLE ADAPTER TI 166kg	A-445100

\* Reimbursed for PMB



WWW.OSSUR.COM/EN-ZA

Össur South Africa  
Unit 1, Stonewood Business Park  
9 Tanzanite Crescent  
Kraaifontein  
Cape Town  
7550

T 0800 996 220  
F 0800 981 744  
E ossursa@ossur.com

