

**rössur**<sup>®</sup>  
LIFE WITHOUT LIMITATIONS

**Navii**<sup>®</sup>

Expand Your World



# Meet Navii®

**An Evolution** in Microprocessor  
Knee Technology



## Next Gen Actuator

Powerful stability,  
effortless swing



## Instantaneous Stance Control

Supporting users of all  
functional levels



## Waterproof In Any Environment

Navii is IP68 + More



## High Activity Features

Seamless walk to run, cycling and  
stair ascent



## Variable Position Manual Lock

The first of its kind



## Removable Protective Cover


Designed for form and function



## Össur Logic

Training exercises and  
advanced settings





## A Trusted Partner

From day one, Navii was created to enable a life without limitations. We designed a device that would keep moving in any environment, and expand the envelope of what's possible for the user.

### Built for Impact

Robust load cells and actuator

### Taking Waterproofing to a New Level

Corrosion resistance and easy drainage and cleaning

### High Capacity, High Durability Battery

Up to 72 hours continuous battery autonomy – battery life can be extended by powering off when not in use

## Navii is IP68 + more

**6** Ingress protection from solids:  
**Dust-tight**


**8** Ingress protection from water:  
**Continuous immersion**

**+**  
More Navii is protected in any environment:  
**saltwater, chlorine, sand and mud.**  
Just rinse and repeat!



## The Freedom to Move Confidently

Navii instills confidence with every step. Its revolutionary actuator and gait algorithms promote stance stability, consistent transitions, and unmatched freedom of movement.



### Standing Stability

High baseline torque provides increased resistance to support the user



### Gait Initiation

A consistent first step is critical for user safety. Navii initiates gait easier than other MPK's<sup>1-3</sup>

### Stumble Recovery

In the event of a stumble, Navii allows users to recover by providing strong instantaneous stance flexion resistance upon loading

### Effortless Swing

Low friction swing phase with lowered metabolic cost of walking<sup>1</sup>

### Enhanced Descent

New yielding control allows easier initiation of stance flexion and consistent resistance throughout high flexion angles

## Seamless Lifestyle Features



Walk to Run



Automatic Cycling Detection



Step-Over-Step Stair Ascent



Variable Position Manual Lock

## Elevate Outcomes with Össur Logic

With Össur Logic at your fingertips, you can see progress in motion on any iOS device. From easy set-up to advanced configuration, prosthetists can customise Navii settings to meet the changing needs of every individual.

### Biofeedback Training Exercises

Provide prosthetists, physiotherapists, and users with a tool that supports ongoing improvement

### Setup is a Breeze

Navii gets to know the user's gait in seconds

### Advanced Activity Reports

Help track progress and support reimbursement

### Enhanced Stability Control

Additional features available to support K2 functional level users and those undergoing rehabilitation



## User Removable Protective Covers



Navii protective covers, inspired by the landscapes of Iceland, add value in function and form.

- Anatomical shape and high-quality finish give natural look and feel
- Adds comfort and stability when kneeling
- Available in 3 lengths

## Upgrade the Navii Experience with Össur Foot Technology

Combining Navii in a customised Össur Leg Solution offers unique benefits across a wide range of user profiles.

”

Navii with Proprio Foot® gives me safety, stability, and adaptability on uneven terrains. I can essentially walk anywhere I want, no more planning or changing routes.\*

**Chris Olivier**

*Capetown, South Africa*



\*The statements made depict only one person's experience and is not a guarantee, promise, or reflection of the experience of others.



”

Navii allows me to go anywhere and do anything without a second thought, and with the Pro-Flex® Terra it feels like I have an ankle again.\*

**Noelle Lambert**  
Manchester, New Hampshire

”

No matter what the situation is, no matter where I am, Navii is able to sense what's going on and keep me from tripping or falling. After 62 years of being an amputee, all of a sudden my leg takes care of me.\*

**Gil Moncrief**  
Coeur d'Alene, Idaho





**Navii with Proprio Foot®**  
Synchronous knee and ankle function provides additional safety<sup>4,5</sup> and improved symmetry on varied terrain<sup>6-8</sup>



**Navii with Pro-Flex® Terra**  
Össur's premiere transfemoral solution for K3/K4 users living an active lifestyle and looking for maximum dynamics

Activity Level

Proprio Foot and Pro-Flex Family Feet are compatible with Navii, complementing its functionality for a range of users.



**K-Level**

	K1	K2	K3	K4
<b>Impact Level</b>	Low			
	Moderate			High

**User Information**

Amputation Level	Transfemoral, Knee Disarticulation, Hip Disarticulation	
Osseointegration	Yes	
Bilateral	Yes	
Impact Level	Low to Moderate	High
Maximum User Weight	136 kg	110 kg



Learn More about Navii  
[go.ossur.com/NaviiANZ](http://go.ossur.com/NaviiANZ)

References available via QR Code above

[WWW.OSSUR.COM.AU](http://WWW.OSSUR.COM.AU)  
[WWW.OSSUR.CO.NZ](http://WWW.OSSUR.CO.NZ)

Össur ANZ  
 TEL AU 1300 123 268  
 TEL NZ 0800 369 524  
[anzprosthetics@ossur.com](mailto:anzprosthetics@ossur.com)

FAX AU +61 2 9475 1114  
 FAX NZ 0800 448 265

