



## Mobility Clinic Gold Coast 2024

We recently hosted the Gold Coast Mobility Clinic at the beginning of April! It was an awesome and inspiring event jam packed with activities for amputees of all ages and abilities.

Our Mobility Clinic allowed athletes the chance to participate in running drills, TABATA workouts and even a new refreshing pool session followed by a soothing yoga workshop.

It was also inspiring to watch healthcare professionals learn and guide our athletes through mobility techniques and drills.

Visit our website to uncover all that our Mobility Clinic has to offer and stay connected with Össur for future events and updates.

[Explore Mobility Clinics](#)

## New Product Updates

## Male and female 4-hole rotation adapter updates



The adapters have changed vendor and undergone design changes. Vendor change will improve both product quality and delivery, and the design change will make these parts compatible with our Knee Rotation Adapter (A-554300). The new adapters have a 1mm higher build height over the old versions.



Please note: New item numbers have been created for the updated versions:

- A-244320 replaces M626400 (Female)
- A-234320 replaces M625400 (Male)

## Waterproof Unity: Balance Foot S & Balance Foot S Torsion



The Unity module for both Balance Foot S and Balance Foot S Torsion is now officially rated waterproof. This certification extends to both the individual feet and the combined Unity module with the foot.



It's important to note that there have been no changes made to the product itself. Instead, we have conducted rigorous testing according to our waterproof standards, which the Unity module has successfully passed.

If you have any questions regarding our latest updates, please contact us at anytime.

[Contact Us](#)

---

## Customer Care

Contact the Customer Care team today!

Email

Meet the team



### Össur ANZ

**Sydney:** 26 Ross St, North Parramatta, NSW 2151 | 1300 123 268

**Melbourne:** 9 Market Dr, Bayswater North, VIC 3153 | 1300 123 268

**Auckland:** 6/80 Westpoint Dr, Hobsonville, Auckland, 0618 | 0800 369 524

[www.ossur.com.au](http://www.ossur.com.au) | [www.ossur.co.nz](http://www.ossur.co.nz)

© 2025 Össur

[Unsubscribe](#)

Want to change the types of emails you receive from

Össur?

Selecting your preferences helps us tailor content to your interests, ensuring emails are relevant and engaging to you.

[Update your preferences](#)