



## The Future of Prosthetics is Here with Navii®



We're excited to introduce Navii, a groundbreaking Microprocessor-controlled knee set to transform the prosthetics industry with its advanced features and performance.

From day one, Navii has been engineered with waterproof capabilities in mind, making it fully waterproof and sandproof. This means it is not only perfect for daily activities like showering but also ideal for water adventures such as river fishing or enjoying the beach, pool or ocean. Navii's corrosion-resistant design ensures durability, while its easy drainage and cleaning features make maintenance a breeze. Say goodbye to limitations and hello to a prosthetic solution that seamlessly integrates into any water activity.

Please note that the Unity TF Kits have been updated to include an extra clip, enabling the connection of the Unity tube to the Navii knee cover.

[Download Brochure](#)

## Swing Into The AXIS®

### New Locking Shoulder Joint

Formerly the LTI Locking Shoulder Joint, the AXIS® has been re-engineered by College Park.

With a patented friction hinge, the shoulder joint is operated with simple switches that simplifies harnessing by eliminating the need for gross body movements. The AXIS is able to be integrated with all upper limb control systems with a special channel to conceal and protect power cables. The AXIS provides better task control, a natural range of motion, and ease of use for the user's daily life.



[Learn More](#)

---

## Total Knee 1900 Housing & Bumper Cover Updates



The Total Knee 1900 has been redesigned for increased size, strength, and durability. This update includes changes to the housing, knee caps and bumper covers.

- **Compatibility:** New knee caps are not compatible with current model Total Knee 1900, and vice versa.
- **Part Numbers:** New part numbers with an "N" suffix are available for the updated knee caps.
- **Identification:** Models with grey bumper covers and serial number MXB045502 or lower are the current version.
- **Bumper Covers:** Total Knee 1900 now feature black bumper covers for improved aesthetics.



## Össur Academy Courses

Össur Academy offers a range of courses year-round at our Melbourne facility, designed to provide in-depth education for clinicians in a distraction-free environment. In-clinic training is also available through your Össur Clinical Specialist. Our upcoming courses include:

- [Direct Socket TF Course \(Mar 24 - 25\)](#)
- [Connect TF and Mechanical Knees \(Apr 1\)](#)
- [Pro-Flex Family - Introducing the New Pro-Flex Terra \(Apr 8\)](#)
- [Össur TT Solutions \(May 5-6\)](#)

[View Full Calendar](#)

## 2025 Product Catalogue

Our updated 2025 catalogue is now available online.  
To access this, please click the link below.

If you prefer the hard copy, please reach out to your local representative.



[Download Here](#)

## Customer Care

Contact the Customer Care team today!

Email

Meet the team



### Össur ANZ

**Sydney:** 26 Ross St, North Parramatta, NSW 2151 | 1300 123 268

**Melbourne:** 9 Market Dr, Bayswater North, VIC 3153 | 1300 123 268

**Auckland:** 6/80 Westpoint Dr, Hobsonville, Auckland, 0618 | 0800 369 524

[www.ossur.com.au](http://www.ossur.com.au) | [www.ossur.co.nz](http://www.ossur.co.nz)

© 2025 Össur

[Unsubscribe](#)

Want to change the types of emails you receive from Össur?

Selecting your preferences helps us tailor content to your interests,  
ensuring emails are relevant and engaging to you.

[Update your preferences](#)