




Pro-Flex Terra

[Explore the Innovations](#)

Innovation Rooted in Technology, Inspired by People

Since the launch of the first carbon fiber prosthetic foot, we've been committed to solving real-world challenges for amputees. Our latest Pro-Flex advancements continue to redefine comfort, mobility, and energy return - helping users move more naturally and with less strain. Discover how our cutting-edge designs, including the Pro-Flex Terra, are transforming lives and setting new industry standards.

[Discover the Innovations](#)



College Park Torsion Adapters Now Available!

We're excited to share that 5 College Park Torsion Adapters are now available for customer orders **effective immediately**.

					
	<i>Torsion Adapter, Pyramid Adapter to Pyramid Receiver, Ti</i>	<i>Torsion Adapter, 30mm Tube Adapter, Ti, 275mm</i>	<i>Torsion Adapter, 30mm Tube Adapter, Ti, 455mm</i>	<i>Torsion Adapter, 34mm Tube Adapter, Ti, 275mm</i>	<i>Torsion Adapter, 34mm Tube Adapter, Ti, 455mm</i>
Part No.	TOR PAPR T	TOR T275 30	TOR T455 30	TOR T275 34	TOR T455 34
Material	Titanium	Titanium (receiver), Aluminum (tube)	Titanium (receiver), Aluminum (tube)	Titanium (receiver), Aluminum (tube)	Titanium (receiver), Aluminum (tube)
Build Height	43 mm (1.7 in)	275 mm (10.8 in)	455 mm (18.0 in)	275 mm (10.8 in)	455 mm (18.0 in)
Cut Length (MIN)	n/a	83 mm (3.3 in)	83 mm (3.3 in)	79 mm (3.1 in)	79 mm (3.1 in)
Assembly Weight	393 g	567 g	670 g	569 g	685 g
Patient Weight Limit	150 kg (330 lbs)	136 kg (300 lbs)	136 kg (300 lbs)	150 kg (330 lbs)	150 kg (330 lbs)
Axial Rotation	+/- 20°	+/- 20°	+/- 20°	+/- 20°	+/- 20°
Warranty	2 Year, limited	2 Year, limited	2 Year, limited	2 Year, limited	2 Year, limited


Explore College Park

i-Limb Hand Update New Colour Options Available!

We're pleased to announce that new colour options for the **i-Limb hand** (all sizes and variants) are now available.

In addition to the standard black option, the dorsal fairing is also available in the following colours: **Titanium, Carbon and White**.

The introduction of these additional colours will give the user more choice in the appearance and personalisation of their prosthetic hand.



New Transparent



Original Clear

Availability and Ordering

New orders: These colour options can be selected at no additional cost when included in a new i-Limb order.

Existing Devices: For existing hands under warranty, a cover colour update can be requested for a small fee when returned for service.

Transparent Glove Option

A new **transparent i-Limb Skin Active glove** is also available, designed to showcase the chosen fairing colour underneath. This is recommended for use with the new fairings to enhance visibility.

Please note, TS tip option for the transparent i-Limb Skin Active is not available.

[Learn More](#)



Mobility Clinics 2025

Help Your Clients Plan Ahead

Össur Mobility Clinics are in full swing with new costings, fresh coaches, and a revised agenda. These events bring amputees and clinicians together to move, connect, and

grow. Join us, take part in hands-on sessions, and support your clients in a truly powerful way.

Event	Date	Location
 On the Snow	1–4 Aug	Thredbo
 Mobility Clinic	20–21 Sept	Melbourne
 Mobility Clinic	1–2 Nov	Adelaide
 On the Waves	22 Nov	Newcastle
 On the Green	15 Dec	Melbourne

[Register Now!](#)

Flyers and Posters

Help us spread the word. You can request physical copies or download them directly.

[Contact Us for Physical Copies](#)

[Download Digital Poster](#)

Customer Care

Contact the Customer Care team today!

Email

[Meet the team](#)



Össur ANZ

Sydney: 26 Ross St, North Parramatta, NSW 2151 | 1300 123 268

Melbourne: 9 Market Dr, Bayswater North, VIC 3153 | 1300 123 268

Auckland: 6/80 Westpoint Dr, Hobsonville, Auckland, 0618 | 0800 369 524

www.ossur.com.au | www.ossur.co.nz

© 2025 Össur

[Unsubscribe](#)

Want to change the types of emails you receive from Össur?

Selecting your preferences helps us tailor content to your interests,
ensuring emails are relevant and engaging to you.

[Update your preferences](#)