



Navii Empowers Users to Expand Their World

- **Waterproof:** Navii is rated waterproof to IP68. Fresh, salt, chlorinated water, plus sand and mud are no problem for Navii. Just rinse and repeat!
- **Actuator:** The new, stronger actuator provides powerful stability and freedom of movement. It allows for improved consistency and predictability in stair and ramp descent and variable terrains and gives a smooth swing for natural gait dynamics.
- **Lock:** The three-stage lock mechanism – 20°, 10°, and fully extended at 0° - gives total security and stability, and assists in challenging environments like standing on public transport, leaning, sitting, or playing a variety of sports
- **Covers:** Navii's magnetic protective covers are available in five colours - **Adobe Red**, **Wave Blue**, **Mosi Green**, **Luna Grey**, and **Onyx Black**



Explore the Innovations

Prosthetics Academy Courses: Össur MPK Solutions

Dates

Wednesday 21 May, 2025


Wednesday 17 September, 2025

Location

Össur Academy facility Melbourne

Join us for a hands-on session to explore Össur's latest waterproof MPK technology and advancements. This expert-led training will equip you with the skills to optimise client outcomes through effective prescription, alignment, and gait training. You'll gain insight into:

- Key features and benefits of Össur's waterproof MPK
- Prescription considerations for optimal results
- Alignment and programming techniques
- Gait training methods
- A demonstration of the Össur Power Knee



AOPA Endorsed Education

2 CPD Points

[Register Now](#)

Need Support?

Contact our Customer Care team today!

Email

Meet the team



Össur ANZ

Sydney: 26 Ross St, North Parramatta, NSW 2151 | 1300 123 268

Melbourne: 9 Market Dr, Bayswater North, VIC 3153 | 1300 123 268

Auckland: 6/80 Westpoint Dr, Hobsonville, Auckland, 0618 | 0800 369 524

www.ossur.com.au | www.ossur.co.nz

© 2025 Össur

[Unsubscribe](#)

Want to change the types of emails you receive from Össur?

Selecting your preferences helps us tailor content to your interests, ensuring emails are relevant and engaging to you.

[Update your preferences](#)